The image displays three distinctively colored EZ-IO needles arranged vertically. Each needle features a cylindrical hub with a ribbed grip and a sharp, straight needle tip extending from the side. The top needle is pink, the middle one is blue, and the bottom one is yellow.

The image displays the Arrow EZ-IO® product line. It features three separate cannula components: a pink cannula with a white hub, a blue cannula with a white hub, and a yellow cannula with a white hub. Each cannula has a small white tube attached to its side. Below these components is a larger hub assembly consisting of a grey base with a central tube and two smaller tubes extending from the sides.

**INDICATIONS FOR USE:**  
For intravenous access anytime in which vascular access is difficult to obtain in emergent, urgent or medically necessary cases.

<b>ADULTS</b> <ul style="list-style-type: none"> <li>• Proximal humerus</li> <li>• Proximal tibia</li> <li>• Distal tibia</li> </ul>	<b>PEDIATRICS</b> <ul style="list-style-type: none"> <li>• Proximal humerus</li> <li>• Proximal tibia</li> <li>• Distal tibia</li> <li>• Distal femur</li> </ul>
--	--

**CONTRAINdicATIONS FOR USE:**

- Fracture in target bone
- Previous significant orthopedic procedure at the site, prosthetic limb or joint.
- 10° catheter use in >48 hours of the target bone.
- Infection in the area of insertion.
- Excessive tissue (severe obesity) and/or absence of adequate anatomical landmarks.

**WARNINGS AND PRECAUTIONS FOR EZ-10® INTRAVASCULAR ACCESS SYSTEM:** 

**CAUTION:** Use aseptic technique.

**CAUTION:** Check skin, adipose and muscle thickness before insertion.

**CAUTION:** Extra care should be taken during insertion and site monitoring when used in patients with bone diseases that increase the likelihood of fracture, extravasation and dislodgement.

**CAUTION:** Do not recap Needle Sets or reconnect separated components. Use biohazard and sharps disposal precautions. Re-use of content may cause transmission of disease, leading to patient risk at time of application(s).

**CAUTION:** Before administering heparin, insulin, or highly-concentrated drugs, check the 10° Catheter again for placement and integrity.

**CAUTION:** Use caution with chemotherapeutic agents.

**CAUTION:** Monitor insertion site frequently for extravasation.

**CAUTION:** Stylet and Catheter are **NOT MRI compatible**.

**CAUTION:** Do not leave the Catheter inserted for longer than 72 hours.

**CAUTION:** Needle Sets are single use only; serious medical consequences (e.g. life-threatening infection) and reduced performance (e.g. blunted needles) may occur if compliance to this warning is not followed. (For a complete listing of these serious medical and performance consequences, please contact Teleflex.)


**STORAGE:** Store in a cool dry place.

**ARROW EZ-IO**  
INTRAVASCULAR VASCULAR ACCESS

**EZ-10® NEEDLE SETS: DESCRIPTION**

- Comprised of Catheter with Luer-lock connection, Stylet, Safety Cap.
- 15 gauge, 304 stainless steel in 15 mm, 25 mm and 45 mm lengths.
- Sterile, non-pyrogenic, in protective packaging.
- Intended for use with EZ-10® Power Driver.

**EZ-10® Power Driver and Needle Sets: Description**

**Insertion Sites**

For additional clinical educational resources please visit [Teleflex.com/EZ10education](http://Teleflex.com/EZ10education)

**1. Clean insertion site**

**2. Prepare supplies.**

- a. Prime EZ-Connex
- Unlock clamp
- Prime set and
- b. Open EZ-Stabilizer

**3. Attach Needle Set**

**IMPORTANT:** Only handle the hub.

**IMPORTANT:** Control the needle.

**4. Push Needle Set into Hub**

**5. 5 mm of the Catheter**

**IMPORTANT:** The use of depth markings measuring guides.

**Neonates and Infants**

General recommendations for marking verification:

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The diagram illustrates the EZ-10 Power Driver, which is a red, ergonomic device with a handle and a mounting slot for the needles. Four different needle sizes are shown: EZ-10® 15 mm (3-15), EZ-10® 25 mm (3-25), EZ-10® 40 mm (3-40), and EZ-10® 45 mm (3-45). Each needle is shown with its corresponding cap (Turvalasitus) and a catheter (Catheter) attached. Labels include: Neulasjori, Turvalasitus, EZ-10® 15 mm (3-15), EZ-10® 25 mm (3-25), EZ-10® 40 mm (3-40), EZ-10® 45 mm (3-45) (note: the 45 mm label is rotated), Mandriini, Napa, and Catheter.

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The diagram shows three anatomical sites for intraosseous insertion: 
 

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**constitutional protocol/policy.**

**0° Power Driver and remove Safety Cap from Catheter.**

Stylet Set by the plastic Hub.

movement prior to and during procedure.

skin until tip touches bone.

(at least one black line) must be visible outside the skin.

**determine correct needle selection is**  
**depth markings on each catheter function as depth**  
**marker soft tissue depth overlying bone (see above),**

):

for needle set selection include the following (depth  
 will be done prior to insertion):

- infants proximal and distal tibia.
- infants in distal femur, proximal and distal tibia.

**6. Squeeze trigger and apply moderate steady pressure.**

**IMPORTANT: DO NOT USE EXCESSIVE FORCE.** Use moderate steady downward pressure and allow Needle Set rotation to penetrate the bone.

**Note:** If Driver stalls and Needle Set will not penetrate the bone, operator may be applying much downward pressure to penetrate bone.

**Note:** In the event of a Driver failure, disconnect the Power Driver, grasp the Needle Set handle and advance into the medullary space while twisting.

**7. Advance Needle Set and release Trigger.**

**Pediatrics:** Release Trigger when sudden "give" or "pop" is felt, indicating entry into medullary space.

**Adult:** Advance Needle Set approximately 1-2 cm after entry into medullary space; Proximal humerus for most adults needle should be advanced 2 cm or until Hub is flush or against the skin.

**8. Stabilize Needle Set Hub, disconnect Driver, and remove Stylet.**

**9. Place Stylet into NeedleVISE® for sharp containment.**

**Note:** Place the NeedleVISE® on a flat stable surface. Immediately following use of a sharp, still holding it with one hand away from the sharp end, firmly snap the sharp point straight down into the opening in the NeedleVISE® until it stops, making sure to KEEP FREE HAND AWAY FROM THE SHARP SECURING DEVICE DURING INSERTION. DO NOT USE NEEDLEVISE® WITH FREE HAND WHILE INSERTING NEEDLE. ALWAYS USE ONE-HAND TECHNIQUE WHEN INSERTING SHARP into NeedleVISE®. Dispose of opened sharp into NeedleVISE® whether or not it has been used.

**10. Obtain samples for lab analysis, if needed.**

**Note:** Only attach a Syringe directly to the EZ-10® Catheter Hub when drawing blood for laboratory analysis (stabilize Catheter) or removal.

**11. Use of the EZ-Stabilizer® is strongly recommended for all EZ-10® insertions.**

- Place Stabilizer over Catheter Hub.

**12. Attach a primed EZ-Connect® extension set to the Hub, firmly secure by twisting clockwise.**

**Note:** DO NOT use any instruments to tighten connections.

**Note:** To prevent valve damage, DO NOT use needles or blunt cannula to access the swabbable valve. Non-standard syringes or connectors can damage the swabbable valve.

**Note:** Operator may use a sterile alcohol wipe to swab the EZ-Connect® valve and let it dry before connecting to the hub.

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
©2014, kaikki oikeudet  
EZ-Stabilizer™ |


The image shows a white medical device package. At the top left is the 'teleflex' logo. To the right is the text 'Customer Service: 1.866.479.8500'. Below this is a small icon of a person in a hospital bed with the text 'Manufactured for: Teleflex Medical'. The main label features a central illustration of a medical device with a red handle and a black probe. To the left of the illustration is a circular 'CE' mark. To the right is a circular 'Single Use' label. Below the illustration are three text labels: 'Do Not Use If Package is Damaged' (with a crossed-out trash bin icon), 'Keep Dry' (with a raindrop icon), and 'Caution' (with a warning triangle icon). At the bottom of the label, it says 'For The Medical Device Directive (EU) 93/42/EEC'. The bottom edge of the package has the text 'Arrow™, EZ-10 Intraoperative Vascular System, and EZ-Connect™ are trademarks of Teleflex Inc.'

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