ARROW®

INTERNATIONAL

EN	Arterial Catheterization Product	1
FR	Produit de cathétérisme artériel	3
DE	Besteck zur Katheterisierung von Arterien	5
IT	Prodotto per cateterismo arterioso	7
PL	Przyrząd do cewnikowania tętnic	9
PT	Produto de Cateterização Arterial	11
RU	Устройство для артериальной катетеризации	13
SL	Pripomoček za arterijsko kateterizacijo	15
ES	Dispositivo para cateterismo arterial	17
SV	Produkt för artärkatetrering	19
TR	Arteriyel Kateterizasyon Ürünü	21

Arterial Catheterization Product

Safety and Efficacy Considerations:

Do not use if package has been previously opened or damaged.

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Do not alter catheter or any other kit/set component during insertion, use, or removal.

Procedure must be performed by trained personnel well versed in anatomical landmarks, safe technique, and potential complications.

Indications for Use:

The Arrow® arterial catheterization device permits access to the peripheral arterial circulation.

Contraindications:

Inadequate collateral arterial circulation.

Warnings and Precautions:*

- Warning: Sterile, Single use: Do not reuse, reprocess or resterilize. Reuse of device creates a potential risk of serious injury and/or infection which may lead to death.**
- Warning: Practitioners must be aware of complications associated with arterial procedures^{1,7,8,10} including bacteremia,² septicemia, vessel wall perforation, intravascular clotting and embolization, hematoma, arterial spasm, tissue necrosis, hemorrhage, thrombosis, peripheral ischemia and infarction, peripheral nerve damage⁹ and air embolism.³**
- Warning: In brachial procedures collateral flow cannot be guaranteed and therefore intravascular clotting can result in tissue necrosis.¹²**
- Warning: In radial artery procedures, practitioners must ascertain that definite evidence of collateral ulnar flow exists.**
- Warning: Care should be exercised that indwelling catheter is not inadvertently kinked at hub area when securing catheter to patient. Kinking may weaken the wall of the catheter and cause a fraying or fatigue of material, leading to possible separation of catheter.**
- Precaution: Do not suture directly to outside diameter of catheter body to minimize the risk of damaging the catheter or adversely affecting monitoring capabilities.**
- Warning: Due to the risk of exposure to HIV (Human Immunodeficiency Virus) or other blood borne pathogens, health care workers should routinely use universal blood and body-fluid precautions in the care of all patients.**
- Precaution: To minimize the risk of disconnects, it is recommended that only Luer-Lock connecting tubing be used with this device.⁵**
- Precaution: Indwelling catheter should be routinely inspected for desired flow rate, security of dressing, and possible migration. Do not use scissors to remove dressing to minimize the risk of cutting catheter.**

A Suggested Procedure:

Use sterile technique.

- Assess for adequate collateral arterial circulation.
- Prep and drape anticipated puncture site per hospital protocol.
- Utilize local anesthetic as required.
- Remove protective shield from needle/catheter assembly. Use of catheter hub wing clip is optional. It may be removed from suture ring groove, if desired, and the suture ring on catheter may be used for suturing catheter to skin (refer to Fig. 1).


Fig. 1

- Locate and puncture artery with needle/catheter assembly using a continuous, controlled, slow, forward motion (refer to Fig. 2). Avoid transfixing both vessel walls. **Precaution: If both vessel walls are punctured, subsequent advancement of spring-wire guide could result in inadvertent sub-arterial placement.** Blood flashback into hub of introducer needle indicates successful entry into vessel.


Fig. 2

- Holding introducer needle hub in place, carefully advance only the catheter 1-2 mm into vessel.
- Holding catheter in place, remove introducer needle. Pulsatile blood flow indicates successful entry into artery. **Warning: Exposure of arterial circulation to atmospheric pressure may result in entry of air into circulation.** Apply pressure over artery to minimize blood loss. **Precaution: Do not reinsert needle into catheter to minimize the risk of catheter damage.**
- Insert feed tube with spring-wire guide into catheter hub. Advance guide wire slowly to required position in artery using actuating lever (refer to Fig. 3). **Precaution: Do not advance guide wire unless there is free blood flashback.**


Fig. 3

When reference mark on clear feed tube coincides with edge of internal cylinder of actuating lever, the tip of spring-wire guide is located at catheter tip. **Precaution: If resistance is encountered while advancing spring-wire guide, do not force feed. Withdraw entire unit and attempt new puncture.**

9. Holding spring-wire guide and feed tube in position, advance catheter forward with a slight rotating motion over spring-wire guide into vessel (refer to Fig. 4).


Fig. 4

10. Hold catheter in place and remove spring-wire guide assembly. Pulsatile blood flow indicates positive arterial placement. **Warning: Exposure of arterial circulation to atmospheric pressure may result in entry of air into circulation.**
11. Attach desired stopcock, injection cap or connecting tubing to catheter hub. Secure catheter to patient in preferred manner using suture ring or wing clip. **Warning: Care should be exercised that indwelling catheter is not inadvertently kinked at hub area when securing catheter to patient. Kinking may weaken the wall of the catheter and cause a fraying or fatigue of material, leading to possible separation of catheter. Precaution: Do not suture directly to outside diameter of catheter body to minimize the risk of damaging catheter or adversely affecting monitoring capabilities.**
12. Dress puncture site according to hospital protocol.
13. Document catheter insertion procedure.

Catheter Removal Procedure:

1. Remove dressing. **Precaution: To minimize the risk of cutting catheter do not use scissors to remove dressing.**
2. **Warning: Exposure of arterial circulation to atmospheric pressure may result in entry of air into circulation.**

Remove suture(s). Be careful not to cut catheter. Remove catheter slowly, pulling parallel to skin. As catheter exits site, apply pressure with a dressing until hemostasis occurs. Apply light pressure dressing to site.

3. Upon removal of catheter, inspect to make sure entire length has been withdrawn.
4. Document removal procedure.

References:

1. Abadir AR, Ung KA. Complications of radial artery cannulation. *Anesthesiology Rev*. 1980;7: 11-16.
2. Band JD, Maki DG. Infections caused by arterial catheters used for hemodynamic monitoring. *Am J Med*. 1979;67:735-741.
3. Chang C, Dughi J, Shitabata P, Johnson G, Coel M, McNamara JJ. Air embolism and the radial arterial line. *Crit Care Med*. 1988;16:141-143.
4. Clark CA, Harman EM. Hemodynamic monitoring: arterial catheters. In: Taylor RW, Civetta JM, Dirby RR, eds. *Techniques and Procedures in Critical Care*. Philadelphia, PA: JB Lippincott; 1990:218-230.
5. Daily EK, Schroeder JS. *Techniques in Bedside Hemodynamic Monitoring*. 5th ed. St. Louis, MO: Mosby; 1994:71-77.
6. Falk PS, Scuderri PE, Sherertz RJ, Motsinger SM. Infected radial artery pseudoaneurysms occurring after percutaneous cannulation. *Chest*. 1992;101: 490-495.
7. Kaye W. Invasive monitoring techniques: arterial cannulation, bedside pulmonary artery catheterization, and arterial puncture. *Heart Lung*. 1983;12:395-427.
8. Kaye WE, Dubin HG. Vascular cannulation. In: Taylor RW, Civetta JM, Dirby RR, eds. *Techniques and Procedures in Critical Care*. Philadelphia, PA: JB Lippincott; 1990:204-208.
9. Marshall G, Edelstein G, Hirshman CA. Median nerve compression following radial arterial puncture. *Anesth Analg*. 1980;59:953-954.
10. Puri VK, Carlson RW, Bander JJ, Weil MH. Complications of vascular catheterization in the critically ill. *Crit Care Med*. 1980;8:495-499.
11. Norwood SH, Cormier B, McMahon NG, Moss A, Moore V. Prospective study of catheter-related infection during prolonged arterial catheterization. *Crit Care Med*. 1988;16:836-839.
12. Venus B, Satish P. Vascular Cannulation. In: Civetta JM, Taylor RW, Kirby RR, eds. *Critical Care*. 3rd ed. Philadelphia, PA: Lippincott-Raven Publishers; 1997:539-542.

Arrow International, Inc. recommends that the user be acquainted with the reference literature.

*If you have any questions or would like additional reference information, please contact Arrow International, Inc.

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Arrow International, Inc. kullanıcının referans literatüre aşına olmasını önerir.
*Sorularınız varsa veya ek referans bilgisi isterseniz lütfen Arrow International, Inc. ile irtibat kurun.