
INSTRUCTIONS

EN

English

1. **Select** the area for the placement of the Grip-Lok.
Note: The catheter should be already inserted in the patient with Luer connection attached.
2. **Prepare** the skin according to the standard hospital protocol for dressing application. Hair removal may be required on some patients for better adhesion.
3. **Open** the Grip-Lok and fold back the top liner to expose center section of adhesive.
4. **Slide** the Grip-Lok under the catheter and the Luer connector centering in the exposed adhesive area.
5. **Hold** the Grip-Lok and catheter in position and remove the folded liner.
6. **Hold** the Grip-Lok in position and remove one side of the bottom liner then remove the other side of the bottom liner and gently press adhesive to skin.
7. **Wrap** top adhesive tab around the catheter pressing to conform to the shape of the catheter.
8. **Wrap** the remaining top adhesive tab around the Luer connector pressing to conform to the shape of the connector.

Direction: for added skin protection, skin prep may be used prior to application of securement device.
Note: Re-use of this device may change its mechanical or biological features and may cause device failure, allergic reactions or bacterial infections.
Note: replace securement device if soiled or saturated in fluid.
Note: use of an alcohol swab will ease removal of securement device from the skin.

This product does not contain any natural rubber latex.

This product is indicated for single patient use.

Grip-Lok is a registered trademark of TIDI Products, LLC.
Patent Pending.

G-04020-100A (1/18)

REF

Description

3300MIV	Peripheral IV & Arterial Catheter Securement
3300MART	Arterial Catheter & Peripheral IV Securement

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