Right Atrial Electrocardiography (RAECG) Adapter

Safety and Efficacy Considerations:

Do not use if package has been previously opened or damaged. **Warning: Prior to use read all package insert warnings, precautions, and instructions. Failure to do so may result in severe patient injury or death.**

Do not alter the adapter during use.

Procedure must be performed by trained personnel well versed in anatomical landmarks, safe technique and potential complications.

Indications for Use:

The RAECG Adapter permits verification of catheter tip position using electrocardiography.

Contraindications:

None known.

Warnings and Precautions:*

- Warning: Sterile, Single use: Do not reuse, reprocess or resterilize. Reuse of device creates a potential risk of serious injury and/or infection which may lead to death.**
- Warning: Due to the risk of exposure to HIV (Human Immunodeficiency Virus) or other blood borne pathogens, health care workers should routinely use universal blood and body-fluid precautions in the care of all patients.**

Carefully read all warnings and precautions throughout procedure instructions.

RAECG Adapter

Electro-Shock Warnings:

- Warning: Strict electrical isolation of monitoring equipment and good grounding techniques are mandatory when using a saline-filled catheter as a transvenous ECG electrode. All equipment should comply with American Heart Association recommendation of being incapable of delivering more than 10 microamp leakage current to the patient. An electrocautery protection filter should be inserted between ECG cable and ECG machine. Rubber cap must be over stainless steel male**

RAECG button when adapter is in-line but not in use.

- Warning: Physician must be aware of complications associated with right atrial and inadvertent right ventricular catheterization whenever central venous catheter being introduced is intentionally passed beyond junction of superior vena cava and right atrium.**
- Warning: Any premature heart beat might indicate catheter touching heart; withdraw catheter and check position.**

A Suggested Procedure: Use sterile technique.

- Aseptically place catheter intravenously.
- Sterile Operator – Pass distal or female end of a sterile length of IV tubing to non-sterile assistant.
- Non-sterile Assistant – Luer-Lock adapter to IV tubing and flush system with an ionic salt solution, e.g. 0.9% saline solution, 4% sodium chloride or 8.4% sodium bicarbonate solution. Meticulously remove all air bubbles and blood clots.
- Sterile Operator – Luer-Lock male or proximal end of sterile IV tubing to catheter.
- Non-sterile Assistant – Flush entire assembly.
- Sterile Operator – Remove rubber protective cap from stainless steel male button on adapter.
- Non-sterile Assistant – For ECG machines without lead selectors, connect negative terminal lead of ECG to adapter, the positive terminal to left leg, and the ground terminal lead to left arm. For ECG machines with lead selectors, choose lead II on ECG machine and place right arm lead to adapter.
- Sterile Operator – Advance or withdraw saline filled catheter until size of the P wave of ECG trace just begins to progressively get larger. At this point, catheter tip should be at cavoatrial junction. If tip is advanced further into right atrium, P wave will either be as large as R wave or become biphasic. If catheter enters right ventricle, P wave will disappear and a wider QRS will appear.
- Continue with procedure.

Trouble Shooting Tips:

- Salt bridge solution – The solution in the catheter must be highly ionic to conduct the electrical impulse. Do not use dextrose in water since it is non-ionic.
- Flush fluids – Discontinue all flush fluids to minimize the risk of electrical deflection.
- Grounding – Ensure good skin connection with ECG electrodes to optimize ECG tracing.

- Internal catheter movement – Movement of chest during respiration or movement of heart during contraction can cause catheter to move and produce an undulating trace. Repositioning chest leads may alleviate problem.

*If you have any questions or would like additional reference information, please contact Arrow International, Inc.

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* Om du har frågor eller önskar ytterligare referensinformation, kontakta Arrow International, Inc.