## Emergency Infusion Device - EID™

### Safety and Efficacy Considerations:

Do not use if package has been previously opened or damaged.

**Warning:** Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Do not alter EID™ or any kit/set component during insertion, use or removal.

Procedure must be performed by trained personnel well versed in anatomical landmarks, safe technique and potential complications.

### Indications for Use:

The Arrow Emergency Infusion Device (EID) permits placement of a sheath for rapid volume infusion.

### Contraindications:

None known.

### Warnings and Precautions:\*

1. **Warning: Sterile, Single use: Do not reuse, reprocess or resterilize.** Reuse of device creates a potential risk of serious injury and/or infection which may lead to death.
2. **Warning: Practitioners must be aware of complications associated with percutaneous sheath introduction including vessel wall perforation, pleural and mediastinal injuries, air embolism, sheath embolism, thoracic duct laceration, bacteremia, septicemia, thrombosis, inadvertent arterial puncture, nerve damage, hematoma formation, and hemorrhage.**
3. **Warning: Do not apply excessive force in removing guide wire, dilator or sheath. If withdrawal cannot be easily accomplished a chest x-ray should be obtained and further consultation requested.**
4. **Warning: The practitioner must be aware of potential air embolism problems associated with leaving open needles, sheaths, or catheters in venous puncture sites or as a consequence of inadvertent disconnects. To lessen the risk of disconnects, only securely tightened Luer-Lock connections should be used with this device. Follow hospital protocol for all sheath maintenance to guard against air embolism.**
5. **Warning: Due to the risk of exposure to Human Immunodeficiency Virus (HIV) or other blood borne pathogens, health care workers should routinely use universal blood and body-fluid precautions in the care of all patients.**
6. **Precaution: Do not suture directly to the outside diameter of the sheath to minimize the risk of cutting or damaging the sheath or impeding sheath flow.**

7. **Precaution: Indwelling sheaths should be routinely inspected for desired flow rate, security of dressing, correct position and for proper Luer-Lock connection.**
8. **Precaution: Maintain the insertion site with regular meticulous redressing using aseptic technique.**

### A Suggested Procedure: Use sterile technique.

1. **Precaution: Place patient in slight Trendelenburg position as tolerated to reduce the risk of air embolism. If femoral approach is used, place patient in supine position.**
2. Prep area of anticipated venipuncture.
3. Drape puncture site as required.
4. Perform skin wheal using desired needle.
5. Tighten Luer-Lock connection between syringe and Y-connector.
6. Use scalpel to nick skin at puncture site.
7. Do not remove guard from sheath. Using guard tab, pull guard back toward Y-connector to expose approximately one-half of sheath.
8. Puncture vessel using a continuous, controlled, slow, forward motion. **Precaution: Avoid transfixing both vessel walls. If both vessel walls are punctured, advancement of spring-wire guide could result in inadvertent subvessel placement.**  
The guard can be used to support the sheath during insertion. To remove guard continue to pull back. Guard will 'peel off' sheath.  
Retracting syringe will hasten blood flashback that can be noted in the clear hub of introducer needle and syringe, indicating successful entry into vessel. **Precaution: Do not advance guide wire unless venous blood can be aspirated by retracting syringe plunger.**
9. Stabilize position of introducer needle and carefully advance spring-wire guide as far as required into vessel.
10. Advance spring-wire by firmly gripping wire and sterility shield approximately 5 cm away from rubber seal and pushing forward. The "J" tip of spring-wire is positioned directly behind the heel of introducer needle. Approximately 1.5 cm of forward movement will fully expose the "J" tip from the lumen of the needle. **Precaution: Grasping wire too close to proximal end of wire will result in bending the wire. Precaution: If resistance is encountered while advancing spring-wire guide, do not force feed and do not withdraw spring-wire guide against needle bevel to minimize the risk of possible severing or damaging of spring-wire guide. Withdraw entire unit and attempt new puncture.**

11. Advance EID 1 to 2 mm further into vessel.
  12. Hold needle hub in position and advance sheath/dilator assembly forward into vessel using slight twisting motion. **Precaution: Do not withdraw dilator until sheath is well within vessel to minimize the risk of damage to sheath tip.**
  13. Advance sheath off dilator into vessel, again grasping near skin and using slight twisting motion.
  14. Holding sheath in place, remove spring-wire guide, dilator and needle/syringe assembly. Free blood flow indicates placement in vessel. Cover sheath hub connection with sterile-gloved finger to minimize the risk of air embolism. **Warning: Do not leave tissue dilator in place as an indwelling catheter. Leaving tissue dilator in place puts patient at risk for possible vessel wall perforation. Warning: Although the incidence of spring-wire guide failure is extremely low, practitioner should be aware of the potential for breakage if undue force is applied to the wire.**
  15. Flush infusion tubing and securely attach to sheath hub.
  16. Use suture to secure sheath and/or anchor with a purse string suture around the sheath suture ring. **Precaution: Do not suture directly to the outside diameter of the sheath to minimize the risk of cutting or damaging the sheath or impeding flow.**
  17. Cover puncture site with sterile dressing.
  18. Record insertion procedure on patient's chart.
3. If applicable, remove sutures from sheath. **Precaution: Be careful not to cut sheath.**
  4. **Warning: Exposure of the central vein to atmospheric pressure may result in entry of air into the central venous system.** Remove sheath slowly, pulling it parallel to the skin. As sheath exits the site, apply pressure with a dressing impermeable to air e.g. VASELINE®\* gauze. Because the residual sheath track remains an air entry point until completely sealed, the occlusive dressing should remain in place for a least 24-72 hours dependent upon the amount of time the sheath was indwelling.<sup>1,2,3</sup>
  5. Upon removal of the sheath, inspect it to make sure that the entire length has been withdrawn.
  6. Verify that the sheath was intact upon removal.
  7. Document removal procedure.

### References:

1. Kashuk JL, Penn I. Air embolism after central venous catheterization. *Surg Gynecol Obstet.* September 1984;159:249-252.
2. Phifer TJ, Bridges M, Conrad SA. The residual central venous catheter track – an occult source of lethal air embolism: case report. *J Trauma.* 1991;31:1558-1560.
3. Thielen JB, Nyquist J. Subclavian catheter removal. *JIN.* March/April 1991;14:114-118.

Arrow International, Inc. recommends that the user be acquainted with the reference literature.

\*If you have any questions or would like additional reference information, please contact Arrow International, Inc.

†A registered trademark of Unilever Supply Chain, Inc.

### Sheath Removal Procedure:

1. **Precaution: Place the patient in a supine position.**
2. Remove dressing. **Precaution: To minimize the risk of cutting the sheath, do not use scissors to remove the dressing.**

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\*Om du har frågor eller önskar ytterligare referensinformation, kontakta Arrow International, Inc.