# Emergency Infusion Device - EID

For short term use (<30 days)

Rx only.

## Indications for Use:

The Arrow Emergency Infusion Device (EID) permits placement of a sheath for rapid volume infusion.

## Contraindications:

None known.

## Clinical Benefits to be Expected:

The ability to access into the circulation and infuse large fluid volumes rapidly into a patient for treatment of shock or trauma, as examples.

## ⚠️ General Warnings and Precautions

### Warnings:

1. Sterile, Single use: Do not reuse, reprocess or resterilize. Reuse of device creates a potential risk of serious injury and/or infection which may lead to death. Reprocessing of medical devices intended for single use only may result in degraded performance or a loss of functionality.
2. Read all package insert warnings, precautions and instructions prior to use. Failure to do so may result in severe patient injury or death.
3. Do not apply excessive force in removing guidewire, dilator or sheath. If withdrawal cannot be easily accomplished a chest x-ray should be obtained and further consultation requested.
4. The practitioner must be aware of potential air embolism problems associated with leaving open needles, sheaths, or catheters in venous puncture sites or as a consequence of inadvertent disconnects. To lessen the risk of disconnects, only securely tightened Luer-Lock connections should be used with this device. Follow hospital protocol for all sheath maintenance to guard against air embolism.
5. Clinicians must be aware of complications/undesirable side-effects associated with percutaneous sheath introduction including, but not limited to:

- vessel wall perforation
- pleural and mediastinal injuries
- air embolism
- sheath embolism
- thoracic duct laceration
- bacteremia
- septicemia
- thrombosis
- inadvertent arterial puncture
- nerve damage
- hematoma formation
- hemorrhage

### Precautions:

1. Do not alter the sheath or any other kit/set component during insertion, use, or removal.
2. Procedure must be performed by trained personnel well versed in anatomical landmarks, safe technique and potential complications.
3. Use standard precautions and follow institutional policies for all procedures including safe disposal of devices.
4. Do not suture directly to the outside diameter of the sheath to minimize the risk of cutting or damaging the sheath or impeding sheath flow.

5. Indwelling sheaths should be routinely inspected for desired flow rate, security of dressing, correct position and for proper Luer-Lock connection.
6. Maintain the insertion site with regular meticulous redressing using aseptic technique.
7. Some disinfectants used at device insertion site contain solvents which can weaken the device material. Alcohol, acetone, and polyethylene glycol can weaken the structure of polyurethane materials. These agents may also weaken the adhesive bond between stabilization device and skin.
  - Do not use acetone on device surface.
  - Do not use alcohol to soak device surface or allow alcohol to dwell in a device lumen to restore patency or as an infection prevention measure.
  - Do not use polyethylene glycol containing ointments at insertion site.
  - Take care when infusing drugs with a high concentration of alcohol.
  - Allow insertion site to dry completely prior to applying dressing.

## A Suggested Procedure: Use Sterile Technique.

⚠️ **Precaution:** Place patient in slight Trendelenburg position as tolerated to reduce the risk of air embolism. If femoral approach is used, place patient in supine position.

1. Prep area of anticipated venipuncture.
2. Drape puncture site as required.
3. Perform skin wheal using desired needle.
4. Tighten Luer-Lock connection between syringe and Y-connector.
5. Use scalpel to nick skin at puncture site.
6. Do not remove guard from sheath. Using guard tab, pull guard back toward Y-connector to expose approximately one-half of sheath.
7. Puncture vessel using a continuous, controlled, slow, forward motion.

⚠️ **Precaution:** Avoid transfixing both vessel walls. If both vessel walls are punctured, advancement of guidewire could result in inadvertent subvessel placement.

The guard can be used to support the sheath during insertion. To remove guard continue to pull back. Guard will 'peel off' sheath.

Retracting syringe will hasten blood flashback that can be noted in the clear hub of introducer needle and syringe, indicating successful entry into vessel.

⚠️ **Precaution:** Do not advance guidewire unless venous blood can be aspirated by retracting syringe plunger.

8. Stabilize position of introducer needle and carefully advance guidewire as far as required into vessel.
9. Advance guidewire by firmly gripping wire and sterility shield approximately 5 cm away from rubber seal and pushing forward. The "J" tip of spring-wire is positioned directly behind the heel of introducer needle. Approximately 1.5 cm of forward movement will fully expose the "J" tip from the lumen of the needle.

⚠️ **Precaution:** Grasping wire too close to proximal end of wire will result in bending the wire.

⚠️ **Precaution:** If resistance is encountered while advancing guidewire, do not force feed and do not withdraw guidewire against needle bevel to minimize the

risk of possible severing or damaging of guidewire. Withdraw entire unit and attempt new puncture.

10. Advance EID 1 to 2 mm further into vessel.
11. Hold needle hub in position and advance sheath/dilator assembly forward into vessel using slight twisting motion.
  - ⚠ **Precaution: Do not withdraw dilator until sheath is well within vessel to minimize the risk of damage to sheath tip.**
12. Advance sheath off dilator into vessel, again grasping near skin and using slight twisting motion.
13. Holding sheath in place, remove guidewire, dilator and needle/syringe assembly. Free blood flow indicates placement in vessel. Cover sheath hub connection with sterile-gloved finger to minimize the risk of air embolism.
  - ⚠ **Warning: Do not leave tissue dilator in place as an indwelling catheter. Leaving tissue dilator in place puts patient at risk for possible vessel wall perforation. Warning: Although the incidence of guidewire failure is extremely low, practitioner should be aware of the potential for breakage if undue force is applied to the wire.**
14. Flush infusion tubing and securely attach to sheath hub.
15. Use suture to secure sheath and/or anchor with a purse string suture around the sheath suture ring.
  - ⚠ **Precaution: Do not suture directly to the outside diameter of the sheath to minimize the risk of cutting or damaging the sheath or impeding flow.**
16. Cover puncture site with sterile dressing.
17. Record insertion procedure on patient's chart.

### Sheath Removal Procedure:

- ⚠ **Precaution: Place the patient in a supine position.**
1. Remove dressing.

- ⚠ **Precaution: To minimize the risk of cutting the sheath, do not use scissors to remove the dressing.**
2. If applicable, remove sutures from sheath.
  - ⚠ **Precaution: Be careful not to cut sheath.**
  3. Remove sheath slowly, pulling it parallel to the skin. As sheath exits the site, apply pressure with a dressing impermeable to air. Because the residual sheath track remains an air entry point until completely sealed, the occlusive dressing should remain in place for a least 24-72 hours dependent upon the amount of time the sheath was indwelling.
    - ⚠ **Warning: Exposure of the central vein to atmospheric pressure may result in entry of air into the central venous system.**
  4. Upon removal of the sheath, inspect it to make sure that the entire length has been withdrawn.
  5. Verify that the sheath was intact upon removal.
  6. Document removal procedure.

For reference literature concerning patient assessment, clinician education, insertion technique, and potential complications associated with this procedure, consult standard textbooks, medical literature, and Arrow International LLC website: [www.teleflex.com](http://www.teleflex.com)

A pdf copy of this IFU is located at [www.teleflex.com/IFU](http://www.teleflex.com/IFU)

For a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation 2017/745/EU on Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority. The contacts of national competent authorities (Vigilance Contact Points) and further information can be found on the following European Commission website: [https://ec.europa.eu/growth/sectors/medical-devices/contacts\\_en](https://ec.europa.eu/growth/sectors/medical-devices/contacts_en)

Symbol Glossary: Symbols are in compliance with ISO 15223-1.

Some symbols may not apply to this product. Refer to product labeling for symbols that apply specifically to this product.

Caution	Medical device	Consult instructions for use	Do not reuse	Do not resterilize	Sterilized by ethylene oxide	Single sterile barrier system	Keep away from sunlight	Keep dry
Do not use if package is damaged	Not made with natural rubber latex	Catalogue number	Lot number	Use by	Manufacturer	Date of manufacture	Importer	

Arrow, the Arrow logo, Teleflex and Teleflex logo are trademarks or registered trademarks of Teleflex Incorporated or its affiliates, in the U.S. and/or other countries. © 2021 Teleflex Incorporated. All rights reserved.

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