Arrow® Raulerson Spring-Wire Introduction Syringe

Rx only

Intended Purpose:

The Arrow® Raulerson Syringe is indicated for use with spring-wire guides ranging from .025" (.64 mm) to .038" (.96 mm) and spring-wire guide introduction needles sized to place these wires. The syringe allows for vessel location and passage of the guidewire through the syringe obviating the need for separating the introducer needle from the syringe.

Patient Target Group:

Patients requiring central catheter placement.

Contraindications:

None known.

Clinical Benefits to be Expected:

The ability to use as an adjunct in the placement of a spring-wire guide using the Seldinger technique to reduce the potential for air emboli, and blood exposure for the user.

Performance Characteristics:

To facilitate spring-wire guide into central anatomy allowing successful Seldinger insertion.

Contains Hazardous Substance:

Components manufactured using Stainless Steel can contain > 0.1% weight by weight of Cobalt (CAS # 7440-48-4) which is considered a category 1B CMR (Carcinogenic, mutagenic or toxic to reproduction) substance. The amount of Cobalt in the Stainless Steel components has been evaluated and considering the intended purpose and toxicological profile of the devices there is no biological safety risk to patients when using the devices as instructed within this IFU.

⚠ General Warnings and Precautions

Warnings:

1. Sterile, Single use: Do not reuse, reprocess or resterilize. Reuse of device creates a potential risk of serious injury and/or infection which may lead to death. Reprocessing of medical devices intended for single use only may result in degraded performance or a loss of functionality.
2. Read all package insert warnings, precautions and instructions prior to use. Failure to do so may result in severe patient injury or death.
3. Passage of guidewire into the right heart can cause dysrhythmias, right bundle branch block, and a perforation of vessel, atrial or ventricular wall.
4. Air embolism can occur if air is allowed to enter a vascular access device or vein. Do not leave open needles in venous puncture site. Use only securely tightened Luer-Lock connections with any vascular access device to guard against inadvertent disconnection.

5. Clinicians must be aware of complications/undesirable side-effects associated with vessel access and guidewire placement including, but not limited to:

- vessel wall perforation
- air embolism
- guidewire embolism
- bacteremia
- septicemia
- inadvertent arterial puncture
- guidewire malposition
- hematoma
- dysrhythmias
- nerve damage
- anaphylaxis
- tissue damage/trauma
- local tissue inflammation
- peripheral trauma

Precautions:

1. Do not alter any kit/set component during insertion, use or removal.
2. Procedure must be performed by trained personnel well versed in anatomical landmarks, safe technique and potential complications.
3. Use standard precautions and follow institutional policies for all procedures including safe disposal of devices.
4. If the package is damaged or unintentionally opened before use do not use the device. Dispose of the device.
5. Storage conditions for these devices require that they are kept dry and out of direct sunlight.

A Suggested Procedure: Use Sterile Technique.

Prep Puncture Site:

1. Position patient as appropriate for insertion site.
2. Prepare clean skin with an appropriate antiseptic agent and allow to dry.
3. Drape puncture site.
4. Administer local anesthetic per institutional policies and procedures.
5. Dispose of needle.
6. Insert introducer needle or catheter/needle with attached Arrow Raulerson Syringe into vein and aspirate.

⚠ Warning: Do not leave open needles or uncapped, unclamped catheters in venous puncture site. Air embolism can occur if air is allowed to enter a vascular access device or vein.

⚠ Precaution: Do not reinsert needle into introducer catheter (where applicable) to reduce risk of catheter embolus.

Verify Venous Access:

Utilize one of the following techniques to verify venous access because of the potential for inadvertent arterial placement:

- Central Venous Waveform:
 - Insert fluid primed blunt tip pressure transduction probe into rear of plunger and through valves of Arrow Raulerson Syringe and observe for central venous pressure waveform.
 - Remove transduction probe if using Arrow Raulerson Syringe.
- Pulsatile Flow (if hemodynamic monitoring equipment is not available):
 - Use transduction probe to open syringe valving system of Arrow Raulerson Syringe and observe for pulsatile flow.
 - Disconnect syringe from needle and observe for pulsatile flow.

⚠ Warning: Pulsatile flow is usually an indicator of inadvertent arterial puncture.

⚠ Precaution: Do not rely on blood aspirate color to indicate venous access.

7. Straighten "J" Tip of guidewire for introduction of the guidewire into Arrow Raulerson Syringe.
8. Advance guidewire into Arrow Raulerson Syringe approximately 10 cm until it passes through syringe valves or into introducer needle.
 - Advancement of guidewire through Arrow Raulerson Syringe may require a gentle twisting motion.
9. Use centimeter markings (where provided) on guidewire as a reference to assist in determining how much guidewire has been inserted.
 - ⚠ **Precaution:** Maintain firm grip on guidewire at all times. Keep sufficient guidewire length exposed for handling purposes. A non-controlled guidewire can lead to wire embolus.
 - ⚠ **Warning:** Do not aspirate Arrow Raulerson Syringe while guidewire is in place; air may enter syringe through rear valve.
 - ⚠ **Precaution:** Do not reinfuse blood to reduce risk of blood leakage from rear (cap) of syringe.
 - ⚠ **Warning:** Do not withdraw guidewire against needle bevel to reduce risk of possible severing or damaging of guidewire.
10. Remove introducer needle and Arrow Raulerson Syringe (or catheter) while holding guidewire in place.
11. Proceed with procedure per manufacturer's instructions or institutional policies and procedures.

For reference literature concerning patient assessment, clinician education, insertion technique, and potential complications associated with this procedure, consult standard textbooks, medical literature, and Arrow International LLC website: www.teleflex.com

A pdf copy of this IFU is located at www.teleflex.com/IFU

For a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation 2017/745/EU on Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority. The contacts of national competent authorities (Vigilance Contact Points) and further information can be found on the following European Commission website: https://ec.europa.eu/growth/sectors/medical-devices/contacts_en

en

Symbol Glossary: Symbols are in compliance with ISO 15223-1.

Some symbols may not apply to this product. Refer to product labeling for symbols that apply specifically to this product.

Caution	Medical device	Consult instructions for use	Contains hazardous substances	Do not reuse	Do not sterilize	Sterilized by ethylene oxide	Single sterile barrier system	Keep away from sunlight
Keep dry	Do not use if package is damaged	Not made with natural rubber latex	Catalogue number	Lot number	Use by	Manufacturer	Date of manufacture	Importer

Arrow, the Arrow logo, Teleflex and the Teleflex logo are trademarks or registered trademarks of Teleflex Incorporated or its affiliates, in the U.S. and/or other countries. © 2024 Teleflex Incorporated. All rights reserved.

"Rx only" is used within this labeling to communicate the following statement as presented in the FDA CFR : Caution: Federal law restricts this device to sale by or on the order of a licensed healthcare practitioner.

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