Arrow® Tissue Dilator

Rx only

Intended Purpose:

The Arrow® Dilator is used during vascular or cavity drainage device insertion procedure to enlarge the subcutaneous tissue, to dilate or distend an opening to allow the catheter to pass with ease.

Patient Target Group:

Intended to be used in patients with anatomy suitable for use with the device.

Contraindications:

None known.

Clinical Benefits to be Expected:

Enlarges the subcutaneous tissue and vessel site to dilate or distend an opening to allow a catheter or sheath device to pass with ease.

Performance Characteristics:

To facilitate tissue dilation for advancement of vascular access related devices.

⚠️ General Warnings and Precautions

Warnings:

1. Sterile, Single use: Do not reuse, reprocess or resterilize. Reuse of device creates a potential risk of serious injury and/or infection which may lead to death. Reprocessing of medical devices intended for single use only may result in degraded performance or a loss of functionality.
2. Read all package insert warnings, precautions and instructions prior to use. Failure to do so may result in severe patient injury or death.
3. Do not use excessive force when introducing tissue dilator as this can lead to vessel perforation, bleeding or component damage.
4. Do not leave tissue dilator in place as an indwelling catheter. Leaving tissue dilator in place puts patient at risk for possible vessel wall perforation.
5. Do not leave dilators uncapped in venous puncture site. Air embolism can occur with these practices.
6. Clinicians must be aware of complications/undesirable side-effects associated with vascular access including but not limited to:

• vessel perforation	• thrombosis
• hemorrhage	• inadvertent arterial puncture
• infection	• nerve damage
• air embolism	• hematoma
• dilator embolism	
• mediastinal injury	

Precautions:

1. Do not alter the dilator during insertion, use or removal.
2. Procedure must be performed by trained personnel well versed in anatomical landmarks, safe technique and potential complications.
3. Use standard precautions and follow institutional policies for all procedures including safe disposal of devices.
4. If the package is damaged or unintentionally opened before use do not use the device. Dispose of the device.

5. Storage conditions for these devices require that they are kept dry and out of direct sunlight.

A Suggested Procedure: Use Sterile Technique

Dilator use for Sheath Placement:

1. Insert entire length of dilator through hemostasis valve into device pressing hub of dilator firmly into hub of hemostasis valve assembly. Place assembly on sterile field awaiting final placement.
2. Insert guidewire per manufacturer's instructions.
3. Enlarge cutaneous puncture site with cutting edge of scalpel, positioned away from guidewire.
 - ⚠️ Warning: Do not cut guidewire to alter length.
 - ⚠️ Warning: Do not cut guidewire with scalpel.
4. Thread tapered tip of dilator assembly over guidewire.
 - ⚠️ Precaution: Sufficient guidewire length must remain exposed at hub end of device to maintain a firm grip on guidewire.
5. Grasping near skin, advance assembly with slight twisting motion to a depth sufficient to enter vessel. Dilator may be partially withdrawn to facilitate advancement of device through tortuous vessel.
 - ⚠️ Precaution: Do not withdraw dilator until the device is well within the vessel to reduce the risk of damaging tip.
6. Advance access device assembly off dilator into vessel, again grasping near skin and using slight twisting motion.
7. To check for proper device placement within the vessel, attach syringe to distal side port for aspiration. Hold assembly in place and withdraw guidewire and dilator sufficiently to allow venous blood flow to be aspirated into distal side port.
 - ⚠️ Precaution: Maintain firm grip on guidewire at all times.
8. Holding assembly in place, remove guidewire and dilator as a unit. Place sterile-gloved finger over hemostasis valve until an obturator or catheter is inserted.
 - ⚠️ Warning: Hemostasis valve must be occluded at all times to reduce the risk of air embolism or hemorrhage. If catheter introduction is delayed, temporarily cover valve opening with sterile gloved finger until obturator is inserted. Use Arrow obturator, sold separately, to occlude hemostasis valve assembly. This will ensure that leakage does not occur and inner seal is protected from contamination.

Tissue Dilator Insertion only:

1. Insert guidewire per manufacturer's instructions.
2. If necessary, enlarge cutaneous puncture site with cutting edge of scalpel, positioned away from guidewire.
 - ⚠️ Warning: Do not cut guidewire to alter length.
 - ⚠️ Warning: Do not cut guidewire with scalpel.
3. Use tissue dilator to enlarge tissue tract to the vein as required. Follow the angle of the guidewire slowly through the skin.
4. Holding wire in place, remove dilator, place sterile gloved finger over puncture site.
5. Continue with catheter insertion per manufacturer's instructions.

For reference literature concerning patient assessment, clinician education, insertion technique, and potential complications associated with this procedure, consult standard textbooks, medical literature, and Arrow International LLC website: www.teleflex.com

A pdf copy of this IFU is located at www.teleflex.com/IFU

For a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation 2017/745/EU on Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority. The contacts of national competent authorities (Vigilance Contact Points) and further information can be found on the following European Commission website: https://ec.europa.eu/growth/sectors/medical-devices/contacts_en

en

Symbol Glossary: Symbols are in compliance with ISO 15223-1. Some symbols may not apply to this product. Refer to product labeling for symbols that apply specifically to this product.

								
Caution	Medical device	Consult instructions for use	Do not reuse	Do not resterilize	Sterilized by ethylene oxide	Single sterile barrier system	Keep away from sunlight	Keep dry
								
Do not use if package is damaged	Not made with natural rubber latex	Catalogue number	Lot number	Use by	Manufacturer	Date of manufacture	Importer	

Arrow, the Arrow logo, Teleflex, and the Teleflex logo are trademarks or registered trademarks of Teleflex Incorporated or its affiliates, in the U.S. and/or other countries. © 2023 Teleflex Incorporated. All rights reserved.

"Rx only" is used within this labeling to communicate the following statement as presented in the FDA CFR : Caution: Federal law restricts this device to sale by or on the order of a licensed healthcare practitioner.

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