3M™ Tegaderm™ Transparent Film Dressing

Indications for Use (USA):

Tegaderm™ Film can be used to cover and protect catheter sites and wounds, to maintain a moist environment for wound healing or to facilitate autolytic debridement, as a secondary dressing, as a protective cover over at-risk skin, to secure devices to the skin, to cover first and second degree burns, and as a protective eye covering. Do not use the dressing as a replacement for sutures and other primary wound closure methods.

Intended Purpose:

The transparent film dressing is intended to be used in Arrow kits as a mechanical barrier, to cover and protect a device insertion site.

Contraindications:

None known.

Clinical Benefits Statement:

Transparent film dressings are used to protect IV sites.

Warnings and Precautions

Warnings:

1. **Sterile, Single use:** Do not reuse, reprocess or resterilize. Reuse of device creates a potential risk of serious injury and/or infection which may lead to death. Reprocessing of medical devices intended for single use only may result in degraded performance or a loss of functionality.
2. Read all package insert warnings, precautions and instructions prior to use. Failure to do so may result in severe patient injury or death.

Precautions:

1. Stop any bleeding at the site before applying the dressing.
2. Do not stretch the dressing during application as tension can cause skin trauma.

3. Make sure the skin is clean, free of soap residue and lotion and allowed to dry thoroughly before applying the dressing to prevent skin irritation and to ensure good adhesion.
4. The dressing may be used on an infected site, only when under the care of a health care professional.
5. Antimicrobial ointments containing polyethylene glycols may compromise the strength of the Tegaderm™ Transparent Film Dressings.
6. Tegaderm™ Transparent Dressings should not be re-sterilized by gamma, E-beam or steam methods.

Application:

1. Ensure at least a one inch margin of dressing adheres to healthy, dry skin.
2. Prepare the catheter insertion site according to your institution's approved protocol.
3. To ensure good adhesion, clip excess hair where the dressing will be placed. Do not shave the skin because of the potential for microabrasions.
4. Make sure skin is free of soaps, detergents, and lotions. Allow all preps and protectants to dry thoroughly before applying the dressing. Wet preps and soap residues can cause irritation if trapped under the dressing. Additionally, adhesive products do not adhere well to wet or oily surfaces.
5. Do not stretch the Tegaderm™ dressing during application. Applying an adhesive product with tension can produce mechanical trauma to the skin. Stretching can also cause adhesion failure.
6. The adhesive of Tegaderm™ dressing is pressure-sensitive. To ensure best adhesion, always apply firm pressure to the dressing from the center out to the edges.

A pdf copy of this IFU is located at www.teleflex.com/IFU

3M and Tegaderm are trademarks of 3M.

Arrow is a registered trademark of Teleflex Incorporated or its affiliates.

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**Arrow International LLC**