

ANY CIRCUMSTANCES.  
TYPE ARE NOT APPROVED FOR SALE IN THE U.S. UNDER  
CHLORHEXIDINE GLUCONATE, PREPARATIONS OF THIS  
CHLORHEXIDINE GLUCONATE, WHICH WERE USED DURING  
IN PATIENTS TREATED WITH LUBRICANTS CONTAINING  
REPORTED IN SEVERAL COUNTRIES. THE MOST SERIOUS  
TOPICAL USE OF CHLORHEXIDINE GLUCONATE HAVE BEEN  
HYPERSENSITIVITY REACTIONS ASSOCIATED WITH THE  
OCCUR, DISCONTINUE USE OF THE DRESSING IMMEDIATELY.  
REACTIONS ARE VERY RARE, BUT IF ANY SUCH REACTIONS  
HYPERSENSITIVITY, AND GENERALIZED ALLERGIC  
CHLORHEXIDINE GLUCONATE, ADVERSE REACTIONS TO  
ON PATIENTS WITH A KNOWN SENSITIVITY TO  
DO NOT USE BIOPATCH® DIRECTLY OVER BURN INJURY  
BEEN ESTABLISHED IN CHILDREN UNDER 16 YEARS OF AGE.  
THE SAFETY AND EFFECTIVENESS OF BIOPATCH® HAS NOT  
MEMBRANES.  
FOR EXTERNAL USE ONLY. DO NOT ALLOW THIS PRODUCT  
TO CONTACT THE EYES, EARS, MOUTH, OR MUCOUS  
MEMBRANES.  
NECROSIS OF THE SKIN.  
HAS RESULTED IN HYPERSENSITIVITY REACTIONS AND  
INFANTS. USE OF THIS PRODUCT ON PREMATURE  
WARNING: DO NOT USE BIOPATCH® ON PREMATURE

#### Warnings

#### Labeling Symbols

	Do not re-use		Do not resterilize
	Caution		Do not use if package is damaged
	Batch code		Use-by date
	Catalogue number		Manufacturer
	Temperature limit		
	Sterilized using ethylene oxide		
	Caution: Federal law restricts this device to sale by or on the order of a licensed healthcare practitioner.		

device-related infections.  
BIOPATCH® should not be placed over infected wounds. It is not intended to be used as a treatment of percutaneous  
**Precautions**  
CRBSI, in patients with central venous or arterial catheters.  
skin colonization of microorganisms commonly related to  
catheter-related blood stream infections (CRBSI), and  
catheters. It is also intended to reduce local infections,  
tubes, externally placed orthopedic pins, and epidural  
coronary catheters, mid-line catheters, drains, chest  
arterial catheters, dialysis catheters, peripherally inserted  
medical devices such as: IV catheters, central venous lines,  
by the use of vascular and non-vascular percutaneous  
intended for use as a hydrophilic wound dressing that  
BIOPATCH® containing Chlorhexidine Gluconate is  
**Indication For Use**  
Chlorhexidine Gluconate is a well-known antiseptic agent  
with broad-spectrum antimicrobial and antifungal activity,  
under the dressing.  
incorporated into the dressing inhibits bacterial growth  
eight times its own weight in fluid, while the CHG  
Gluconate (CHG). The foam material absorbs up to  
polyurethane absorptive foam with CHG is a hydrophilic  
BIOPATCH® Protective Disk with CHG is a hydrophilic  
**Product Description**  
(Please Read Carefully Before Using)  
**Instructions For Use**

Protective Disk with CHG

Manufactured for  
ETHICON, INC.  
Raritan, NJ 08869 USA  
© Ethicon, Inc. 2012

U.S. customers: to order product call  
1-800-255-2500;  
for product quality and  
technical questions call  
1-877-384-4266.

LAB0010999v4 09/2020 10558-731-04-Rev. A

CAUTION SHOULD BE USED WHEN USING CHLORHEXIDINE-CONTAINING PREPARATIONS AND THE PATIENT SHOULD BE OBSERVED FOR THE POSSIBILITY OF HYPERSENSITIVITY REACTIONS. THE GOVERNMENT OF JAPAN HAS REPORTED ANAPHYLACTOID-TYPE ADVERSE EVENTS IN 13 PATIENTS WHILE USING CENTRAL VENOUS CATHETERS IMPREGNATED WITH CHLORHEXIDINE.

Clinical Trial Results

A controlled, randomized, clinical trial consisting of 687 subjects with 1699 central venous or arterial catheter insertion sites was conducted at two centers.<sup>1</sup> Results showed that the use of BIOPATCH® resulted in a statistically significant 44% reduction in the incidence of local infection (p≤0.0001).

Table 1: Summary of local infections in 1401 evaluable lines

	No Local Infection # of lines (%)	Local Infection # of lines (%)	Total
BIOPATCH®	556 (83.6%)	109 (16.4%)	665
Control	520 (70.7%)	216 (29.3%)	736
Total	1076	325	1401

Results also showed that the use of BIOPATCH® resulted in a statistically significant 60% reduction in the incidence of catheter-related blood stream infections (p=0.026).

Table 2: Summary of catheter-related blood stream infections (CRBSI) in 589 evaluable subjects

	No CRBSI Frequency (%)	CRBSI <sup>1</sup> Frequency (%)	Total
BIOPATCH®	288 (97.6%)	7 (2.4%)	295
Control	276 (93.9%)	18 (6.1%)	294
Total	564	25	589

<sup>1</sup>Clinical diagnosis based on positive blood cultures and DNA typing.

Results of this study also showed that use of BIOPATCH® resulted in statistically significant reduction in skin colonization of microorganisms commonly associated with CRBSI (p≤0.05). Patients randomized to the BIOPATCH® Treatment Group experienced no serious device-related adverse events.

Information regarding the use of BIOPATCH® on patients <16 years of age is limited. A study performed on 16 patients, ages 3 days to 15 years, was performed to evaluate the effectiveness of BIOPATCH® in the management of insertion or exit sites of indwelling CVCs. No cases of catheter-related infections were reported during the course of the trial. Compared to the institution's standard therapy, BIOPATCH® resulted in better appearance of entrance/exit sites in 56% of cases (p=0.002); less irritation of entrance/exit sites in 50% of cases (p=0.011); better entrance/exit site protection in 94% of cases (p<0.001). BIOPATCH® was the preference of the investigators over standard therapy in 81% of cases (p<0.001).

<sup>1</sup>Maki DG, Mermel L, Genthner D, Hua S, Chiacchierini RP: An evaluation of BIOPATCH® Antimicrobial Dressing compared to routine standard of care in the prevention of catheter-related blood stream infection. Ethicon, Inc. 2000.

Directions For Use

1. Prepare the skin surrounding the percutaneous device according to hospital protocol.
2. Remove BIOPATCH® from the sterile package using aseptic technique.
3. Place BIOPATCH® around the device, making sure the BLUE PRINTED side is facing upward. The WHITE foam side releases the Chlorhexidine Gluconate (CHG) and should be in contact with the patient's skin.
4. In order to ensure easy removal when used with a film dressing, place BIOPATCH® around the device site in such a way that the device rests upon the slit portion of the BIOPATCH®. The edges of the radial slit must be pushed together and remain in contact to maximize efficacy.
5. Secure the device and BIOPATCH® to the skin. Ensure complete contact between the skin and BIOPATCH®.
6. Change the patch as necessary, in accordance with facility protocol; dressing changes should occur at a minimum of every 7 days. Dressing changes will be needed more frequently with highly exuding wounds.
7. To remove the transparent film dressing, pick up the corner of the dressing and stretch the dressing away from the device, holding the device in place. (Dressing will partially lift.) Peel back until resistance is felt. Repeatedly stretch and peel as necessary until the dressing is removed.
8. BIOPATCH® will remain attached to the transparent film dressing, so removal will be simultaneous.

Storage Information

- Store between 15°C and 30°C (59°F and 86°F).
- It is to be stored in its original packaging.
- Expiration date of the product is indicated as year (4 digits), month (2 digits) and day (2 digits).
- Do not resterilize. Discard all open and unused portions of the device.
- Do not use if the package is opened or damaged. Do not use if seal is broken or compromised.
- After use, handle and dispose of all unused product and packaging in accordance with accepted medical practice and applicable local, state, and federal laws and regulations.

NOTE: Over time, the BIOPATCH® may turn yellow in color. This coloration does not reduce the antimicrobial efficacy of the dressing.

How Supplied

BIOPATCH® is supplied sterile. Each package contains a single disk. BIOPATCH® is intended for single use only. Do not resterilize.

Product Code Ending	BIOPATCH® Size	Maximum Amount of CHG per Dressing
-150	1" DISK (2.5 cm) 4.0 mm center hole	92 mg
-151	3/4" DISK (1.9 cm) 1.5 mm center hole	53 mg
-152	1" DISK (2.5 cm) 7.0 mm center hole	86.8 mg