# Pressure Injectable PICC Product

---

Venous Access | Critical Care

EN

NL

FR

DE

EL

IT

PL

PT

RU

ES

SV

TR

ARROW®

Caution

Do not reuse

Do not resterilize

Sterilized  
by ethylene oxide

Do not use  
if package is damaged

Does not contain  
natural rubber latex

Consult  
instructions for use

**Arrow International, Inc.**

© 2010 Arrow International, Inc. All rights reserved. Printed in the USA.

An issued or revision date for these instructions is included for user information. If two years have elapsed between this date and product use, the user should contact Arrow International, Inc. to see if additional product information is available.

Issued Date: August 2015

**Rx only.**

# Table of Contents

<i>Product Description</i>	1	<i>Pre-PICC Insertion &amp; Patient Assessment Activities</i>	6
Indications / Contraindications . . . . .	.1	Procedural Pause . . . . .	.6
<hr/>		Preparing for PICC Insertion . . . . .	.7
<i>Peripherally Inserted Central Catheter Warnings and Precautions</i>	1	Prep Puncture Site . . . . .	.7
General . . . . .	.1	Prepare All Equipment . . . . .	.7
Catheter . . . . .	.2	Trim Catheter . . . . .	.7
Placement Wire & Guidewire / SWG . . . . .	.2	Flush Catheter. . . . .	.7
Tissue Dilator . . . . .	.2	<hr/>	
Peel-Away Sheath over Tissue Dilator . . . . .	.2	<i>Catheter Insertion Instructions</i>	7
Pressure Injection . . . . .	.2	Gain Initial Venous Access . . . . .	.7
Possible Complications . . . . .	.3	Place Sheath . . . . .	.8
<hr/>		Advance Catheter . . . . .	.8
<i>Accessory Component Instructions</i>	3	Placement Wire Use . . . . .	.8
Catheter Stabilization Device . . . . .	.3	Verify Catheter Tip Placement . . . . .	.8
Catheter Trimmer . . . . .	.3	Complete Catheter Insertion . . . . .	.8
Dressing . . . . .	.4	<hr/>	
Echogenic Needle . . . . .	.4	<i>Documentation</i>	9
Filter Straw / Filter Needle . . . . .	.4	<hr/>	
Guidewire / SWG . . . . .	.4	<i>Pressure Injection Instructions</i>	9
80 cm . . . . .	.4	<hr/>	
130 cm . . . . .	.4	<i>Care and Maintenance</i>	9
Luer-Activated Connector . . . . .	.5	Dressing . . . . .	.9
Maximal Barrier Drape. . . . .	.5	Catheter Patency . . . . .	.9
Protected Needle . . . . .	.6	<hr/>	
SharpsAway II™ Locking Disposal Cup . . . . .	.6	<i>Catheter Removal Instructions</i>	9

For convenience, procedural and general Warnings and Precautions are listed at the beginning of the instructions. Please review all content before performing the procedure.

For reference literature concerning patient assessment, clinician education, insertion techniques and potential complications associated with this procedure refer to Arrow International, Inc. website: [www.arrowintl.com](http://www.arrowintl.com)

# Pressure Injectable Peripherally Inserted Central Catheter (PICC) Product

EN

## Product Description

The Arrow® Pressure Injectable PICC is a peripherally inserted central venous catheter manufactured with medical grade, flexible polyurethane. The Arrow® PICC has a non-tapered catheter body with either a blunt tip or a Blue FlexTip® that is softer than a cut tip with a contour design to enhance maneuverability. The Blue FlexTip® also provides visual confirmation of an intact catheter upon removal. The kit components assist the clinician in maintaining maximal sterile barrier precautions.

## Indications:

The Arrow® Pressure Injectable PICC is indicated for short or long term peripheral access to the central venous system for intravenous therapy, blood sampling, infusion, and power injection of contrast media. The maximum pressure of power injector equipment used with the Pressure Injectable PICC may not exceed 300 psi.

## Contraindications:

This device is contraindicated wherever there is presence of device related infections, previous or current thrombosis. Clinical assessment of patient must be completed to ensure no contraindications exists.

## Peripherally Inserted Central Catheter

### ⚠ Warnings and Precautions:

Do not place the catheter into or allow it to remain in the right atrium or right ventricle (refer to Figure 1).

Figure 1

## General Warnings and Precautions

### Warnings:

1. Sterile, Single use: Do not reuse, reprocess or resterilize. Reuse of device creates a potential risk of serious injury and/or infection which may lead to death.
2. Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.
3. Practitioners must be aware of complications associated with central vein catheters including but not limited to: cardiac tamponade secondary to vessel wall, atrial or ventricular perforation, pleural and mediastinal injuries, air embolism, catheter embolism, catheter occlusion, thoracic duct laceration, bacteremia, septicemia, thrombosis, inadvertent arterial puncture, nerve damage, hematoma, hemorrhage, and dysrhythmias.
4. Practitioners must be aware of clinical conditions that may limit use of PICCs including but not limited to: dermatitis, cellulitis, and burns at or about the insertion site, previous ipsilateral venous thrombosis, radiation therapy at or about insertion site, contractures, mastectomy, and potential use for AV fistula.
5. Do not place Central Venous Catheter (CVC) or Peripherally Inserted Central Catheter (PICC) into or allow them to remain in the right atrium or right ventricle. X-ray exam or other method in compliance with hospital/institutional protocol must show catheter tip located in right side of mediastinum in the Superior Vena Cava (SVC) above its junction with right atrium and parallel to vessel wall and its distal tip positioned at a level above either azygos vein or carina of the trachea, whichever is better visualized. Although cardiac tamponade secondary to pericardial effusion is uncommon, there is a high mortality rate associated with it. Improper advancement of guidewire into the heart has also been implicated in causing cardiac perforation and tamponade.

6. Ensure catheter tip has not entered the heart or no longer lies parallel to vessel wall by performing an x-ray exam or other method in compliance with hospital/institutional protocol. If catheter position has changed, immediately re-evaluate.
7. Practitioners must be aware of the potential for entrapment of guidewire by any implanted device in circulatory system (i.e., vena cava filters, stents). Review patient's history before catheterization procedure to assess for possible implants. Care should be taken regarding length of guidewire inserted. It is recommended that if patient has a circulatory system implant, catheter procedure be done under direct visualization to reduce risk of guidewire entrapment.
8. Choose appropriate sized catheter for size of vessel to be cannulated.
9. Catheter tip must be located in central circulation when administering > 10% glucose solution, total parenteral nutrition, continuous vesicant therapy, infusates with pH less than 5 or greater than 9, and infusates with an osmolality above 600 mOsm/L, or any medication known to be irritating to vessels proximal to the vena cava.
10. Infusion of incompatible drugs through a non "staggered port" may cause precipitation.
11. Be aware of the risk of chemically induced thrombophlebitis when catheter is placed with distal end located in a vessel proximal to the SVC.
12. Do not leave open needles or uncapped, unclamped catheters in central venous puncture site. Air embolism can occur with these practices.
13. Use only securely tightened Luer-Lock connections with any Venous Access Device (VAD) to guard against inadvertent disconnect.
14. Use Luer-Lock connectors to help guard against air embolism and blood loss.
15. Pulsatile flow is usually an indicator of inadvertent arterial puncture.

### Cautions:

1. Do not use if package has been previously opened or damaged.
2. Do not alter the catheter, guidewire, or any other kit/set component during insertion, use, or removal (except as instructed).
3. Procedure must be performed by trained personnel well versed in anatomical landmarks, safe technique, and potential complications.
4. Assess patient for heparin sensitivity. Heparin-Induced Thrombocytopenia (HIT) has been reported with the use of heparin flush solutions.
5. Do not use topical antibiotic ointment or creams on insertion sites (except when using dialysis catheters) because of their potential to promote fungal infections and antimicrobial resistance.
6. Temporarily shut off remaining port(s) through which solutions are being infused before blood sampling.
7. Do not rely on blood aspirate color to indicate venous access.
8. Do not reinsert needle into introducer catheter to reduce risk of catheter embolism.
9. Maintain insertion site with regular meticulous redressing using aseptic technique.
10. Engage safety and/or locking feature of scalpel (where provided) when not in use to reduce risk of sharps injury.
11. Perform hand hygiene:
  - before and immediately after all clinical procedures
  - before and after donning and removal of gloves
12. Properly handle and dispose of sharps in sharps container in accordance with US OSHA or other governmental standards for blood borne pathogens and/or hospital/institutional policy.
13. Hands must remain behind the needle at all times during use and disposal.
14. Use universal blood and body-fluid precautions in the care of all patients due to the risk of exposure to Human Immunodeficiency Virus (HIV) or other blood borne pathogens.

## Catheter Warnings and Precautions

### Warnings:

1. Only utilize catheters indicated for high pressure injection applications for such applications. Utilizing catheters not indicated for high pressure applications can result in inter-lumen crossover or rupture with potential for injury.
2. Do not apply excessive force in placing or removing catheter. Excessive force can cause catheter breakage. If placement or withdrawal cannot be easily accomplished, an x-ray should be obtained and further consultation requested.
3. Do not secure, staple, and/or suture directly to outside diameter of catheter body or extension lines to reduce risk of cutting or damaging the catheter or impeding catheter flow. Secure only at indicated stabilization locations.
4. Do not cut catheter to alter catheter length unless procedure requires it.
5. Do not attach catheter clamp and fastener (where provided) until either guidewire or placement wire is removed.
6. Do not use scissors to remove dressing to reduce risk of cutting catheter.
7. Open catheter clamp prior to infusion through lumen to reduce risk of damage to extension line from excessive pressure.
8. Do not attempt to advance or reinsert placement wire (where provided) into catheter through the septum, if placement wire has been removed prior to catheter insertion. Attempting placement wire advancement or reinsertion increases risk of damaging catheter or wire.
9. Do not clamp extension line(s) when placement wire is in catheter to reduce risk of placement wire kinking.
10. Practitioners should remove slide clamp(s), where provided, when not in use. Slide clamp(s) may be inadvertently removed and aspirated by children or confused adults.
11. Residual catheter track remains an air entry point until completely sealed, occlusive dressing should remain in place for at least 24 - 72 hours dependent upon amount of time catheter was indwelling.

### Cautions:

1. Check ingredients of prep sprays and swabs before using. Some disinfectants used at catheter insertion site contain solvents which can attack the catheter material. Alcohol and acetone can weaken the structure of polyurethane materials. These agents may also weaken the adhesive bond between catheter stabilization device and skin.
  - Acetone: Do not use acetone on catheter surface.
  - Alcohol: Do not use alcohol to soak catheter surface or to restore catheter patency.

Take care when instilling drugs containing high concentration of alcohol. Allow insertion site to dry completely prior to applying dressing.
2. Ensure catheter patency prior to injection. Do not use syringes smaller than 10 mL (a fluid filled 1 mL syringe can exceed 300 psi), to reduce risk of intraluminal leakage or catheter rupture.
3. Remove catheter clamp and fastener (where provided) prior to attempting a catheter exchange procedure.
4. Do not exert excessive force while removing the catheter, to reduce risk of catheter breakage.
5. Continuously monitor indwelling catheter for:
  - desired flow rate
  - security of dressing
  - adherence of stabilization device to skin and connection to catheter
  - correct catheter position; use centimeter markings to identify if catheter position has changed
  - secure Luer-Lock connection(s)
6. Minimize catheter manipulation throughout procedure to maintain proper catheter tip position.
7. Provide strain relief of catheter under dressing to decrease catheter movement and assist in maintaining proper catheter tip position.
8. Inject a small amount of radiopaque dye to locate catheter tip if difficulty is encountered in visualizing the catheter tip.

## Placement Wire & Guidewire / SWG Warnings and Precautions

### Warnings:

1. Do not insert stiff end of guidewire into vessel as this may result in vessel damage.
2. Do not cut guidewire to alter length.
3. Do not withdraw guidewire against needle bevel to reduce risk of possible severing or damaging of guidewire.
4. Do not use excessive force when introducing guidewire or tissue dilator as this can lead to vessel perforation and bleeding.
5. Passage of guidewire into the right heart can cause dysrhythmias, right bundle branch block, and a perforation of vessel, atrial, or ventricular wall.
6. Do not apply undue force on guidewire to reduce risk of possible breakage.
7. Do not apply excessive force in removing guidewire or catheter. If withdrawal cannot be easily accomplished, a visual image should be obtained and further consultation requested.
8. Do not cut guidewire with scalpel.
  - Position cutting edge of scalpel away from guidewire.
  - Engage safety and/or locking feature of scalpel (where provided) once cutaneous puncture site is enlarged, to reduce risk of cutting the guidewire (refer to Figure 22).
9. Do not attempt to advance placement wire through septum.
10. Do not cut placement wire when trimming catheter to reduce risk of foreign embolism.
11. Remove placement wire and Luer-Lock sidearm assembly as a unit (refer to Figure 23). Failure to do so may result in wire breakage.

### Caution:

1. Maintain a firm grip on guidewire at all times. Keep sufficient guidewire length exposed at hub for handling purposes. A non-controlled guidewire can lead to wire embolism.
2. Check that there is no wire in cut catheter segment, after trimming catheter. If there is any evidence that placement wire has been cut or damaged, catheter should not be used.

## Tissue Dilator Warnings

### Warnings:

1. Do not leave tissue dilator in place as an indwelling catheter. Leaving tissue dilator in place puts patient at risk for possible vessel wall perforation.
2. Do not use excessive force when introducing guidewire or tissue dilator as this can lead to vessel perforation and bleeding.

## Peel-Away Sheath over Tissue Dilator Precaution

### Caution:

1. Do not withdraw tissue dilator until the sheath is well within the vessel to reduce risk of damage to sheath tip. Sufficient guidewire length must remain exposed at hub end of sheath to maintain a firm grip on guidewire.

## Pressure Injection Warnings and Precautions

### Warnings:

1. Assess each patient for appropriateness of a pressure injection procedure. Pressure injection procedures must be performed by trained personnel well versed in safe technique and potential complications.
2. Obtain a visual image to confirm catheter tip position prior to each pressure injection.
3. Ensure catheter patency prior to pressure injection to reduce risk of catheter failure and/or patient complications.

4. Discontinue pressure injections at first sign of extravasation or catheter deformation. Follow hospital/institutional protocol for appropriate medical intervention.
5. Use only lumen(s) labeled "Pressure Injectable" for pressure injection to reduce risk of catheter failure and/or patient complications.

**Cautions:**

1. Do not exceed the maximum pressure of 300 psi on power injector equipment to reduce risk of catheter failure and/or tip displacement.
2. Do not exceed 5 injections or catheter's maximum recommended flow rate located on product labeling and catheter luer hub to reduce risk of catheter failure and/or tip displacement.
3. Warm contrast media to body temperature prior to pressure injection to reduce risk of catheter failure.
4. Pressure limit settings on injector equipment may not prevent over pressurizing an occluded or partially occluded catheter.
5. Use appropriate administration set tubing between catheter and pressure injector equipment to reduce risk of catheter failure.
6. Follow the contrast media manufacturer's specified instructions for use, contraindications, warnings, and precautions.

**Possible Complications (but not limited to):**

- ◆ cardiac tamponade secondary to vessel wall, atrial or ventricular perforation
- ◆ pleural injury
- ◆ air embolism
- ◆ catheter embolism
- ◆ bleeding / hemorrhage
- ◆ bacteremia
- ◆ thrombosis
- ◆ hematoma
- ◆ brachial plexus injury
- ◆ fibrin sheath formation
- ◆ vessel erosion
- ◆ mediastinal injury
- ◆ nerve injury
- ◆ thoracic duct laceration
- ◆ occlusion
- ◆ septicemia
- ◆ inadvertent arterial puncture
- ◆ dysrhythmias
- ◆ exit site infection
- ◆ phlebitis
- ◆ catheter tip malposition

**Accessory Component Instructions**

Review the list of components that will be utilized before beginning the Arrow® Pressure Injectable PICC insertion procedure. Kits / Sets may not contain all accessory components detailed in this section. Become familiar with instructions for individual component(s) before beginning the actual insertion procedure.

**Catheter Stabilization Device:**

A catheter stabilization device should be used in accordance with manufacturer's instructions for use.

**⚠ Caution:** Minimize catheter manipulation throughout procedure to maintain proper catheter tip position.

- Cleanse and prep anticipated dressing site per hospital/institutional protocol.
- Apply skin prep (where provided) to coat skin and maximize adherence.
- Allow to dry thoroughly.
- Position catheter stabilization device appropriate distance from catheter insertion site to permit ease of insertion site care and maintenance.
- Secure catheter by using the primary suture hub.
- Place suture hub wings over posts and press down (refer to Figure 2). Snap retainer wings to closed position to secure suture hub (refer to Figure 3).


Figure 2


Figure 3

- Remove paper backing from one half of catheter stabilization device pad and press onto dry, prepared skin. Repeat process for other half.
- Complete sterile insertion site dressing according to established hospital/institutional protocol.
- Document dressing application on patient's chart.
- Replace dressing per hospital/institutional protocol. Catheter stabilization device should be replaced at least every 7 days to ensure maximum adherence.

*Refer to individual manufacturer's instructions for more information and specific detailed instructions.*

**Catheter Trimmer:**

A catheter trimmer is a one time use trimming device.

**⚠ Warning:** Do not attempt to advance placement wire through septum.

**⚠ Warning:** Do not cut placement wire when trimming catheter to reduce risk of foreign embolism.

**⚠ Caution:** Check that there is no wire in cut catheter segment, after trimming catheter. If there is any evidence that placement wire has been cut or damaged, catheter should not be used.

**NOTE:** There should be very limited resistance when cutting catheter with supplied trimmer. Any greater resistance is likely to be caused by the placement wire – which has not been sufficiently retracted. If so, do not use catheter.

- Retract placement wire 1-1/2 inches minimum (4 cm) behind where catheter is to be cut. The placement wire is to be withdrawn through septum (refer to Figure 4).


Figure 4

- Kink proximal end of placement wire at connector with side-port (refer to Figure 5) to reduce risk of placement wire extending beyond distal tip of catheter during insertion.


Figure 5

- Peel back contamination guard exposing catheter portion to be trimmed. Using trimming device, cut catheter straight across (90° to catheter cross-section) to maintain a blunt tip.

## Dressing:

A transparent dressing should be used in accordance with manufacturer's instructions for use.

- Prepare site. Allow all preps to dry completely.
- Peel liner from dressing to expose adhesive.
- Adhere center of transparent window over insertion site, while holding notched portion off the skin (refer to Figure 6).


Figure 6

- Overlap softcloth tabs under catheter to form a tight seal around catheter hub and extension line(s) (refer to Figure 7).


Figure 7

- Press dressing into place.
- Slowly remove frame while smoothing down dressing edges. Smooth dressing from center toward edges, use firm pressure to enhance adhesion (refer to Figure 8).


Figure 8

- Use sterile tape strips to secure hub, extension line(s), and/or tubing (refer to Figure 9).


Figure 9

- Label dressing according to protocol.

*Refer to individual manufacturer's instructions for more information and specific detailed instructions for dressing removal (not included).*

## Echogenic Needle:

An echogenic needle is used to provide greater needle visibility under ultrasound. The needle tip is enhanced for approximately 1 cm for clinician to identify exact needle tip location when puncturing the vessel under ultrasound.

## Filter Straw / Filter Needle:

A filter straw/filter needle (5 micron) is utilized to aspirate solution from glass ampule and reduce risk of glass particulate from entering the solution.

- Open glass ampule using appropriate sterile and sharps protection technique.
- Attach filter straw/filter needle to syringe.
- Insert filter straw/filter needle into ampule.
- Aspirate contents from ampule.
- Remove and discard filter straw/filter needle.
- Attach appropriate needleless connector or cannula to syringe.
- Purge air from syringe.
- Label syringe appropriately.

## Guidewire / SWG:

Kits/Sets are available with a variety of guidewires/SWGs. Guidewires are provided in different diameters, lengths, and tip configurations for specific insertion techniques. Become familiar with the guidewire(s) to be used with the specific technique before beginning the actual PICC insertion procedure.

Image guidance may be used to gain initial venous access.

### Catheter Insertion with an 80 cm Guidewire:

Use single 45 cm guidewire for venous access and 80 cm soft tip guidewire for catheter placement. Image guidance or fluoroscopy is used to gain initial venous access; catheter placement with 80 cm guidewire is done under fluoroscopy.

- Gain venous access with 45 cm guidewire and peel-away sheath.
- Load PICC onto 80 cm guidewire until soft tip of wire extends beyond tip of catheter.
- While maintaining control of distal end of guidewire, advance soft tip/catheter tip as a unit through peel-away sheath to desired depth.
- Once catheter is in desired location, remove guidewire.

### Catheter Insertion with an 130 cm Guidewire:

Use single 45 cm guidewire for venous access and 130 cm soft tip guidewire for catheter placement. Image guidance or fluoroscopy is used to gain initial venous access; catheter placement with 130 cm guidewire is done under fluoroscopy.

- Gain venous access with 45 cm guidewire.
- Insert soft end of 130 cm guidewire through peel-away sheath to desired depth.
- Thread catheter over guidewire and advance catheter over guidewire through sheath into vessel into correct position.
- Once catheter is in desired location, remove guidewire.

**NOTE:** Some clinicians will gain access with 130 cm guidewire and thread catheter over guidewire once wire has been correctly positioned in the SVC. This technique is done under fluoroscopy.

**Luer-Activated Connector:**

A luer-activated connector should be used to reduce risk of accidental disconnects and potential for air embolism.

- ⚠ Caution:** Do not use needles in connector unless specifically indicated.
- ⚠ Caution:** Do not cap connector unless specifically indicated.
- ⚠ Caution:** Check for clamping requirement to prevent fluid movement prior to connector disconnect and use.

**General Guidelines for Connector Use:**

Use aseptic technique.

- Remove connector from packaging.
- Inspect for integrity.
- If applicable:
  - prime in accordance with hospital/institutional protocol
  - expel air
- Disinfect desired extension set or venous access device; make sure to rotate swab to ensure cleaning of threads.
- Attach connector to desired extension set or venous access device.
- Secure connector. If using a rotating luer connector, make sure that the spin collar is tightened.
- Position connector to achieve the most comfortable position on the patient’s skin.
- Access: Swab connector with desired disinfectant in accordance with hospital/institutional and connector manufacturer protocol.
- Utilize only fully primed syringe or administration set.
- Disconnect from venous access device or administration set, grasp connector/cap and twist until loose.
- Flush after each use with normal saline or in accordance with hospital/institutional protocol.
- Repeat above steps for subsequent connections.

*Refer to individual manufacturer’s instructions for specific connector details, including priming volumes, dead space, and flow rates.*

**Maximal Barrier Drape:**

Drape(s) provide a maximal sterile barrier. Follow CDC Category 1A Recommendation.

- Drape provided is either:
  - single extra-large drape with fenestration
  - two-piece drape consisting of an arm drape with fenestration and a body drape (body drape is used to appropriately drape torso and upper-lower extremities)
- Unfold the Maximal Barrier Drape:
  - Peel off fenestration backing (refer to Figure 10).


Figure 10

- Position fenestration over intended insertion site (refer to Figure 11).


Figure 11

- Unfold width (refer to Figure 12).


Figure 12

- Unfold towards head (refer to Figure 13).


Figure 13

- Unfold towards hand (refer to Figure 14).


Figure 14

- Perform sterile procedure.
- Tear along seam to remove drape (refer to Figure 15).


Figure 15

## Protected Needle:

See individual manufacturer's instructions for product use, when used as a single product and not as a kit component.

- Warning:** Hands must remain behind needle at all times during use and disposal.
- Caution:** Make sure all needles are used in accordance with OSHA and hospital/institutional safety protocols.
- Caution:** Do not attempt to override or defeat the safety locking mechanism of a protected needle.
- Caution:** Discard in an approved sharps collector in accordance with applicable regulations and hospital/institutional policy.

### General Guidelines for Protected Needle Use:

- Aspirate medication into syringe using aseptic technique.
- If necessary to transport filled syringe to point of administration, use safe, passive recapping technique to cover needle before transport to point of use. In accordance with OSHA standards, such recapping must be accomplished by a one-handed technique, i.e., do not hold needle shield during recapping process.
- Administer injection following established technique.
- Immediately activate needle protection device upon withdrawal from patient by pushing lever arm completely forward until needle tip is fully covered (refer to Figure 16).


Figure 16

- Visually confirm lever arm has fully advanced and needle tip is covered. If unable to activate, discard immediately into approved sharps collector.
- Activation of protective mechanism may cause minimal splatter of fluid that may remain on needle after injection.
- After single use, discard in an approved sharps collector in accordance with applicable regulations and hospital/institutional policy. For greatest safety, use a one-handed technique and activate away from self and others.

## SharpsAway II™ Locking Disposal Cup:

The SharpsAway II™ Locking Disposal Cup is used for disposal of needles (15 Ga. - 30 Ga.).

- Caution:** Do not attempt to remove needles that have been placed into SharpsAway II™ Locking Disposal Cup. These needles are secured in place. Damage may occur to needles if they are forced out of disposal cup.
- Caution:** Do not re-use needles after they have been placed into the foam SharpsAway® system. Particulate matter may adhere to needle tip.

- Using one-handed technique, firmly push needles into disposal cup holes (refer to Figure 17).


Figure 17

- Once placed into disposal cup, needles will be automatically secured in place so that they cannot be reused.
- Discard entire cup, at completion of procedure, into an approved sharps container.
- Where provided, a foam SharpsAway® system may be utilized by pushing needles into foam after use.

## Pre-PICC Insertion & Patient Assessment Activities

### A Suggested Procedure:

- Warning:** Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Clinical assessment of patient must be completed to ensure no contraindications exist e.g. allergies. This device is not recommended for use in the presence of device related infections or previous/current thrombosis.

- Caution:** Perform hand hygiene:
  - before and immediately after all clinical procedures
  - before and after donning and removal of gloves

A procedural checklist is included in many Arrow® products. Check individual product labeling to see if one is included.

### Procedural Pause:

1. Verify physician order:
  - Confirm correct patient.
  - Confirm correct diagnosis.
  - Confirm correct procedure.

Physician order must include post placement assessment of catheter tip placement (direct visualization technique or other method in compliance with hospital/institutional protocol).
2. Educate patient: Explain procedure to patient. Make sure information is presented with respect to patient's level of understanding, culture, and language.
3. Have informed consent signed, if required.
4. Identify insertion vein:
  - Apply tourniquet above anticipated insertion vein.
  - Identify appropriate vein for insertion. Use direct visualization technologies, e.g. ultrasound or fluoroscopy, if available; and assess vein health.

**NOTE:** PICCs are typically inserted into basilic, brachial, or cephalic veins (refer to Figure 18).


Figure 18

5. Release tourniquet and leave in place beneath the arm.
6. Measure patient to assure placement of catheter in the SVC:
  - Extend arm laterally 45 to 90 degrees from trunk.
  - Measure distance from insertion site along presumed anatomical course of vessel to be catheterized.
  - Catheter tip should lie in distal one-third of SVC above right atrium and parallel to SVC wall.
    - ◊ If a catheter stabilization device will be used, add ½ to 1 inch (1.2 to 2.5 cm) to catheter measurement (STATLOCK®); if another device is used, check manufacturer recommendations.
    - ◊ If using upper arm circumference assessment; for consistency in measurement, measure from an anatomical point and record.

7. Position patient as appropriate for insertion site:
  - Extend arm laterally 45 to 90 degrees from trunk.
8. Prepare work area.

**Preparing for PICC Insertion:**

- Clinicians should use sterile technique, maximal sterile barrier precautions throughout the procedure, and dress in protective clothing:
  - mask
  - eye protection
  - sterile gown
  - sterile gloves
  - hair cover

**⚠ Caution:** Perform hand hygiene:
 

- before and immediately after all clinical procedures
- before and after donning and removal of gloves

**⚠ Caution:** Use universal blood and body-fluid precautions in the care of all patients due to the risk of exposure to Human Immunodeficiency Virus (HIV) or other blood borne pathogens.

**⚠ Caution:** Properly handle and dispose of sharps in sharps container in accordance with US OSHA or other governmental standards for blood borne pathogens and/or hospital/institutional policy.

**Prep Puncture Site:**

1. Prep puncture site with appropriate antiseptic/agent.
2. Drape puncture site.

See unfolding instructions for Maximal Barrier Drape (where provided) under Accessory Component Instructions section.

3. Perform skin wheal using desired needle and local anesthetic.
4. Dispose of needle.
  - In kits where provided, the SharpsAway II™ Locking Disposal Cup is used for disposal of needles (15 Ga. - 30 Ga.).

See SharpsAway II™ Locking Disposal Cup under Accessory Component Instructions section.

**⚠ Caution:** Do not attempt to remove needles that have been placed into SharpsAway II™ Locking Disposal Cup. These needles are secured in place. Damage may occur to needles if they are forced out of disposal cup.

**⚠ Caution:** Do not re-use needles after they have been placed into the foam SharpsAway® system. Particulate matter may adhere to needle tip.

**Prepare All Equipment:**

Prepare Catheter with Placement Wire for Insertion, where provided (refer to Figure 19).


Figure 19

- Add desired luer-activated connector.
- Remove catheter tip protector.

**Trim Catheter:**

If necessary, review detailed instructions for Catheter Trimmer device under Accessory Component Instructions section.

1. Identify catheter type:
  - BFT (Blue FlexTip®)
  - Non-BFT
2. Peel back contamination guard exposing catheter portion to be trimmed.
3. Review catheter marking pattern below. The catheter is marked so clinician can easily identify desired amount of catheter to be trimmed.
  - BFT double numbering pattern (refer to Figure 20):


Figure 20

- ◊ First number designates centimeters from tip of catheter.
  - ◊ Second number designates centimeters from hub of catheter.
  - ◊ This double numbering pattern permits clinician to easily identify centimeters of catheter to be trimmed and also identifies centimeters of catheter remaining.
  - ◊ Record both numbers.
- Non-BFT numbering pattern (refer to Figure 21):


Figure 21

- ◊ Number designates centimeters of catheter to be trimmed and also gives amount of catheter remaining.
4. Using the trimming device, cut catheter straight across (90° to catheter cross-section) to maintain a blunt tip.

**NOTE:** There should be very limited resistance when cutting catheter with supplied trimming device. If using a catheter with a placement wire, any greater resistance is likely to be caused by the placement wire which has not been sufficiently retracted. If so, do not use catheter.

5. Inspect cut surface for clean cut and no loose material.

**⚠ Warning:** Do not cut placement wire when trimming catheter to reduce risk of foreign embolism.

**⚠ Caution:** Check that there is no wire in cut catheter segment, after trimming catheter. If there is any evidence that placement wire has been cut or damaged, catheter should not be used.

**Flush Catheter:**

1. Use filter straw/filter needle to withdraw solution from glass ampules.
2. Attach syringe to sidearm and flush distal lumen with sterile saline solution. Leave syringe in place.
3. Flush remaining lumen(s) with sterile saline solution, to establish patency and prime lumen(s).
4. Clamp or attach luer-activated connector(s) to extension line(s) to contain saline within lumen(s).

**Catheter Insertion Instructions**

1. Reapply tourniquet and replace sterile gloves.
2. Locate vein for insertion:
  - Use image guidance, if available.

An echogenic needle is included for access.
3. Insert introducer needle into vein.
4. Check for pulsatile flow.

**⚠ Warning:** Pulsatile flow is usually an indicator of inadvertent arterial puncture.

**⚠ Caution:** Do not rely on blood aspirate color to indicate venous access.

**Gain Initial Venous Access:**

See specific guidewire instructions for Guidewire/SWG Insertion Technique under Accessory Component Instructions section.

1. Insert soft tip of guidewire through introducer needle into vein. Advance guidewire to desired depth.

**⚠ Warning:** Do not insert stiff end of soft tip guidewire into vessel as this may result in vessel damage.

**⚠ Warning:** Do not cut guidewire to alter length.

**⚠ Warning:** Do not withdraw guidewire against needle bevel to reduce risk of possible severing or damaging of guidewire.

- Remove needle:

Hold guidewire in place while removing introducer needle.

**⚠ Caution:** Maintain a firm grip on guidewire at all times. Keep sufficient guidewire length exposed at hub for handling purposes. A non-controlled guidewire can lead to wire embolism.

- Enlarge cutaneous puncture site with cutting edge of scalpel, if necessary, positioned away from guidewire.

**⚠ Warning:** Do not cut guidewire to alter length.

**⚠ Warning:** Do not cut guidewire with scalpel.

- Position cutting edge of scalpel away from guidewire.
- Engage safety and/or locking feature of scalpel (where provided) once cutaneous puncture site is enlarged, to reduce risk of cutting the guidewire (refer to Figure 22).


Figure 22

### Place Sheath:

- Thread tapered tip of peel-away sheath/dilator assembly over guidewire. Grasping near skin advance assembly with slight twisting motion to a depth sufficient to enter vessel. Dilator may be partially withdrawn to further facilitate advancement of sheath into the vessel. A slight twisting motion of the peel-away sheath might help advancement.

**⚠ Caution:** Do not withdraw tissue dilator until the sheath is well within the vessel to reduce risk of damage to sheath tip. Sufficient guidewire length must remain exposed at hub end of sheath to maintain a firm grip on guidewire.

- Check sheath placement by holding sheath in place, withdraw guidewire and dilator sufficiently to allow venous blood flow. Holding sheath in place, remove guidewire and dilator as a unit.

**⚠ Warning:** Do not leave tissue dilator in place as an indwelling catheter. Leaving tissue dilator in place puts patient at risk for possible vessel wall perforation.

**⚠ Warning:** Do not apply undue force on guidewire to reduce risk of possible breakage.

### Advance Catheter:

Advance catheter according to the guidewire used. Review detailed instructions for 80 cm and 130 cm guidewire usage under Accessory Component Instructions section.

**⚠ Warning:** Do not apply excessive force in placing or removing catheter. Excessive force can cause catheter breakage. If placement or withdrawal cannot be easily accomplished, an x-ray should be obtained and further consultation requested.

- Remove catheter tip protector.
- Insert catheter through peel-away sheath:
  - Retract and/or gently flush while advancing catheter if resistance is met.
  - Contamination guard will peel away as catheter is advanced through peel-away sheath.
- Stop advancing catheter 5 inches (13 cm) before reaching pre-established insertion length.
- Withdraw peel-away sheath over catheter until free from venipuncture site.
- Grasp tabs of peel-away sheath and pull apart, away from catheter, until sheath splits down entire length.
- Advance catheter to final indwelling position.

### Placement Wire Use (where provided):

**⚠ Caution:** Do not clamp extension line(s) when placement wire is in catheter to reduce risk of placement wire kinking.

- Complete catheter insertion.
- Remove placement wire.

**⚠ Warning:** Remove placement wire and Luer-Lock sidearm assembly as a unit (refer to Figure 23). Failure to do so may result in wire breakage.


Figure 23

**⚠ Caution:** Do not attach catheter clamp and fastener (where provided) until either guidewire or placement wire is removed.

- Examine tip of placement wire after removal to ensure wire has not been altered (refer to Figure 24).


Figure 24

### Verify Catheter Tip Placement:

- Check catheter placement with syringe by aspirating through distal lumen until free flow of venous blood is observed.

**⚠ Caution:** Do not rely on blood aspirate color to indicate venous access.

### Complete Catheter Insertion:

- Flush lumen(s) to completely clear blood from catheter.
- Connect all extension line(s) to appropriate Luer-Lock line(s) as required. Unused port(s) may be "locked" through luer-activated connector(s) using standard hospital/institutional protocol.
  - Slide clamps are provided on extension lines to occlude flow through each lumen during line and luer-activated connector changes.

**⚠ Warning:** Practitioners should remove slide clamp(s), where provided, when not in use. Slide clamp(s) may be inadvertently removed and aspirated by children or confused adults.

- Cleanse insertion site per hospital/institutional protocol.
- Ensure insertion site is dry before applying dressing. Apply skin protectant as needed.

**⚠ Caution:** Do not use topical antibiotic ointment or creams on insertion sites (except when using dialysis catheters) because of their potential to promote fungal infections and antimicrobial resistance.

- Secure catheter: Use a catheter clamp, fastener, catheter stabilization device or sterile tape strips where provided.

- Assess placement of catheter tip in compliance with hospital/institutional protocol.

**⚠ Warning:** Do not place Central Venous Catheter (CVC) or Peripherally Inserted Central Catheter (PICC) into or allow them to remain in the right atrium or right ventricle. X-ray exam or other method in compliance with hospital/institutional protocol must show catheter tip located in right side of mediastinum in the Superior Vena Cava (SVC) above its junction with right atrium and parallel to vessel wall and its distal tip positioned at a level above either azygos vein or carina of the trachea, whichever is better visualized. Although cardiac tamponade secondary to pericardial effusion is uncommon, there is a high mortality rate associated with it. Improper advancement of guidewire into the heart has also been implicated in causing cardiac perforation and tamponade.

- If catheter tip is malpositioned, reposition, redress, and re-verify.

## Documentation

Hospitals/institutions must establish a permanent medical record that documents the entire procedure, based upon their policy, procedures, and Best Practices. The actual format can differ from institution to institution. Report any product defects/failures to organization risk management, manufacturers, and appropriate regulatory agencies.

Documentation generally includes (but is not limited to) the following information:

1. Device specifics:
  - type, brand and lot number
  - length and size of Vascular Access Device (VAD)
  - internal/external catheter length
  - antimicrobial or not
  - trimmed catheter or not
2. Procedure specifics:
  - time out or procedural pause
  - informed consent, as required
  - date, time of insertion, insertion site, number and site attempts, inserter's identification
  - use of visualization and guidance technologies
  - site preparation and technique
3. Patient assessment and response:
  - pertinent dx, assessment, vital signs
  - understanding of procedure, patient's response to procedure
  - complications and barriers to care
4. Therapy specifics:
  - type of therapy, drug dose, rate, time
  - route and method of administration
  - laboratory specimen collected
5. Visual confirmation:
  - verification of appropriate tip location prior to initial use

Monitor patient for post catheter insertion complications.

## Pressure Injection Instructions

Use sterile technique.

**Warning:** Obtain a visual image to confirm catheter tip position prior to each pressure injection.

1. Remove luer-activated connector from appropriate extension line of catheter.
2. Check for catheter patency:
  - Attach 10 mL syringe filled with sterile normal saline.
  - Aspirate catheter for adequate blood return.
  - Vigorously flush catheter.

**Warning:** Ensure catheter patency prior to pressure injection to reduce risk of catheter failure and/or patient complications.

3. Detach syringe.
4. Attach pressure injection administration set tubing to appropriate extension line of catheter according to manufacturer's recommendations.

**Warning:** Use only lumen(s) labeled "Pressure Injectable" for pressure injection to reduce risk of catheter failure and/or patient complications.

**Caution:** Do not exceed the maximum pressure of 300 psi on power injector equipment to reduce risk of catheter failure and/or tip displacement.

**Caution:** Do not exceed 5 injections or catheter's maximum recommended flow rate located on product labeling and catheter luer hub to reduce risk of catheter failure and/or tip displacement.

5. Inject contrast media in accordance with hospital/institutional protocol.

**Caution:** Warm contrast media to body temperature prior to pressure injection to reduce risk of catheter failure.

6. Disconnect catheter from pressure injector equipment.
7. Flush catheter using 10 mL syringe or larger filled with sterile normal saline.
8. Disconnect syringe and replace with sterile luer-activated connector on catheter extension line.

## Care and Maintenance

### Dressing:

Replace dressing according to hospital/institutional policies, procedures, and practice guidelines. Change immediately if the integrity becomes compromised e.g. dressing becomes damp, soiled, loosened, or no longer occlusive.

**Caution:** Maintain insertion site with regular meticulous redressing using aseptic technique.

- Consult manufacturer's recommendations for dressing specifics.
- Change transparent semipermeable membrane dressing every 7 days.
- Change gauze and tape every 48 hours.
- Label dressing with type, size, and length of catheter; date and time; and initials of the clinician performing dressing change.

### Catheter Patency:

Maintain catheter patency according to hospital/institutional policies, procedures, and practice guidelines. All personnel who care for patients with catheters must be knowledgeable about effective management to prolong catheter's dwell time and prevent injury.

**Caution:** Perform hand hygiene:

- before and immediately after all clinical procedures
- before and after donning and removal of gloves

- Solution and frequency of flushing a venous access catheter should be established in hospital/institutional policy.
- Establish and maintain catheter patency by:
  - flushing intermittently via syringe with heparinized saline or preservative-free 0.9% sodium chloride
  - continuous drip
  - positive pressure device
- Amount of heparin depends on:
  - physician preference,
  - hospital/institutional protocol,
  - patient condition

**Caution:** Assess patient for heparin sensitivity. Heparin-Induced Thrombocytopenia (HIT) has been reported with use of heparin flush solutions.

- Volume of flush solution should be:
  - equal to at least twice the priming volume of the catheter and any add-on devices

**NOTE:** Catheter priming volume is printed on product packaging.

- Properly flush (heparinization) using a positive-pressure flushing technique to help prevent occlusion, when using any VAD for intermittent infusion therapy.

**NOTE:** Neutral as well as positive displacement valve systems have also been shown to help prevent occlusion.

- Properly cleanse all valves with an appropriate antiseptic before being accessed.
- The SASH or SAS method of flushing will help eliminate occlusions due to incompatible solutions:
  - Saline • Administer drug • Saline • Heparin (if used)

## Catheter Removal Instructions

1. Perform catheter removal:
  - following order of authorized prescriber
  - in accordance with hospital/institutional policies, procedures, and practice guidelines
2. Remove catheter immediately upon patient assessment for:
  - suspected contamination i.e. when catheters are inserted during a medical emergency or if adherence to aseptic technique cannot be ensured
  - unresolved complication(s)
  - discontinuation of therapy
  - source of infection
3. Place patient in supine position, as clinically indicated to reduce risk of potential air embolism.
4. Remove dressing.

**⚠ Warning:** Do not use scissors to remove dressing, to reduce risk of cutting catheter.

5. Open catheter stabilization device retainer wings and remove catheter from catheter stabilization device posts.
6. Remove catheter by slowly pulling it parallel to skin. If resistance is met when removing, catheter should not be forcibly removed and physician should be notified.

**⚠ Caution:** Do not exert excessive force while removing catheter, to reduce risk of catheter breakage.

7. Upon removal of catheter:
  - measure and inspect
  - ensure entire catheter length has been removed
8. Apply direct pressure to site until hemostasis is achieved.
9. Apply alcohol swab to catheter stabilization device adhesive and gently lift pad off of skin (if applicable).

10. Dress insertion site. Apply sterile air occlusive dressing and assess site every 24 hours until site is epithelialized.

**⚠ Warning:** Residual catheter track remains an air entry point until completely sealed, occlusive dressing should remain in place for at least 24 - 72 hours dependent upon amount of time catheter was indwelling.

11. Document catheter removal procedure on patient's chart per hospital/institutional protocol.

Include:

- catheter condition
- length of catheter removed/presence of intact catheter tip
- patient's tolerance of the procedure
- any interventions needed for removal

**For reference literature concerning patient assessment, clinician education, insertion techniques and potential complications associated with this procedure refer to Arrow International, Inc. website: [www.arrowintl.com](http://www.arrowintl.com)**

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