Peripherally Inserted Central Catheter (PICC) Product

Safety and Efficacy Considerations:

Do not use if package has been previously opened or damaged. **Warning: Prior to use read all package insert warnings, precautions, and instructions. Failure to do so may result in severe patient injury or death.**

Do not alter the catheter or any other kit/set component during insertion, use or removal.

Procedure must be performed by trained personnel well versed in anatomical landmarks, safe technique, and potential complications.

Warning: Do not place the catheter into or allow it to remain in the right atrium or right ventricle (refer to Fig. 1)


Fig. 1

Cardiac Tamponade: It has been documented by many authors that placement of indwelling catheters in the right atrium is a dangerous practice^{1,2,4,5,7,16,17} that may lead to cardiac perforation and tamponade.^{1,2,4,5,16,17} Although cardiac tamponade secondary to pericardial effusion is uncommon, there is a high mortality rate associated with it.²⁰ Practitioners placing central venous catheters must be aware of this potentially fatal complication before advancing the catheter too far relative to patient size. No particular route or catheter type is exempt from this potentially fatal complication.¹⁷ The actual position of the tip of the indwelling catheter should be confirmed by x-ray after insertion.^{1,2,5,16,17,21} Central venous catheters should be placed in the superior vena cava^{1,2,4,5,7,16,22} above its junction with the right atrium and parallel to the vessel wall^{10,22} and its distal tip positioned at a level above either the azygos vein or the carina of the trachea, whichever is better visualized. Central venous catheters should not be placed in the right atrium unless specifically required for special relatively short term procedures such as aspiration of air emboli during neurosurgery. Such procedures are, nevertheless, risk prone and should be closely monitored and controlled.

Indications for Use:

A Peripherally Inserted Central Catheter permits venous access to the central circulation through a peripheral vein. It offers an alternative method of intravenous therapy for select adult and pediatric patients.

Contraindications:

None known.

Warnings and Precautions:*

1. **Warning: Sterile, Single use: Do not reuse, reprocess or resterilize. Reuse of device creates a potential risk of serious injury and/or infection which may lead to death.**
2. **Warning: Do not place the catheter into or allow it to remain in the right atrium or right ventricle. Central vein catheters should be positioned so that the distal tip of the catheter is in the superior vena cava (SVC) above the junction of the SVC and the right atrium and lies parallel to the vessel wall.**

3. **Warning:** Practitioners must be aware of complications associated with central vein catheters including cardiac tamponade secondary to vessel wall, atrial or ventricular perforation, pleural and mediastinal injuries, air embolism, catheter embolism, thoracic duct laceration, bacteremia, septicemia, thrombosis, inadvertent arterial puncture, nerve damage, hematoma, hemorrhage, dysrhythmias and catheter occlusion.
4. **Warning:** Do not apply excessive force in placing or removing catheter. If placement or withdrawal cannot be easily accomplished, an x-ray should be obtained and further consultation requested.
5. **Warning:** The practitioner must be aware of potential air embolism associated with leaving open needles or catheters in central venous puncture sites or as a consequence of inadvertent disconnects. To lessen the risk of air embolism, only securely tightened Luer-Lock connections should be used with this device. Follow hospital/agency protocol to guard against air embolism for all catheter maintenance.
6. **Warning:** The practitioner must be aware of clinical conditions that may limit the use of PICC catheters such as: dermatitis, cellulitis and burns at or about the insertion site, previous ipsilateral venous thrombosis, radiation therapy, contractures, mastectomy, and potential use for AV fistula.
7. **Warning:** Passage of the guide wire into the right heart can cause dysrhythmias, right bundle branch block,¹⁶ and a perforation of the vessel wall, atrial or ventricular.
8. **Warning:** Practitioners must be aware of the potential for entrapment of guide wire by any implanted device in the circulatory system (ie. vena cava filters, stents). Review patient's history before catheterization procedure to assess for possible implants. Care should be taken regarding the length of spring-wire guide inserted.¹ It is recommended that if patient has a circulatory system implant, catheter procedure be done under direct visualization to minimize the risk of guidewire entrapment.
9. **Warning:** Due to the risk of exposure to HIV (Human Immunodeficiency Virus) or other blood borne pathogens, health care workers should routinely use universal blood and body-fluid precautions in the care of all patients.
10. **Precaution:** Indwelling catheters should be routinely inspected for desired flow rate, security of dressing, correct catheter position and for secure Luer-Lock connection. Use centimeter markings to identify if the catheter position has changed.
11. **Precaution:** Only x-ray examination of the catheter placement can ensure that the catheter tip has not entered the heart or no longer lies parallel to the vessel wall. If catheter position has changed, immediately perform x-ray examination to confirm catheter tip position.
12. **Precaution:** For blood sampling from a two-lumen catheter, temporarily shut off remaining port through which solutions are being infused.
13. **Precaution:** Alcohol and acetone can weaken the structure of polyurethane materials. Check ingredients of prep sprays and swabs for acetone and alcohol content.
Acetone: Do not use acetone on catheter surface. Acetone may be applied to skin but must be allowed to dry completely prior to applying dressing.
Alcohol: Do not use alcohol to soak catheter surface or to restore catheter patency. Care should be taken when instilling drugs containing high concentration of alcohol. Always allow alcohol to dry completely prior to applying dressing.
These agents may also weaken the adhesive bond between the StatLock® Anchoring Device and skin.
14. **Precaution:** Some disinfectants used at the catheter insertion site contain solvents, which can attack the catheter material. Assure insertion site is dry before dressing.
15. **Precaution:** The indications for use in children are the same as adults; however, insertion techniques are often modified according to the age and size of a child. If the practitioner is inexperienced in utilizing this product in a child, appropriate consultation should be sought.
16. **Precaution:** When catheter is placed with distal end located in a vessel proximal to the SVC (midclavicular or midline), practitioners must be aware of the risk of chemically induced thrombophlebitis.¹¹ When administering > 10% glucose solution/total parenteral nutrition, any medication with an osmolarity above 600 mOsm/liter, medications with an osmolarity between 400 and 600 mOsm/liter for more than 4 weeks, or any medication known to be irritating to vessels proximal to the vena cava, the catheter tip must be located in the central circulation.¹⁹
17. **Precaution:** To minimize the risk of pressure induced damage to catheter, do not expose to pressures above 50 psi. Common sources of potentially high pressure include: syringes smaller than 10 cc used to irrigate or declot an occluded catheter (a fluid filled 1 cc syringe can exceed 300 psi⁶), certain radiographic procedures, and infusion pumps with occlusion pressure limits above 50 psi.

18. Precaution: The StatLock® Catheter Securement Device should be routinely inspected for adherence to the skin and connection to the catheter. At the same time, the catheter position should be checked using the centimeter marks on the catheter body.

**A Suggested Procedure:
Use sterile technique.**

1. Explain procedure to patient. Have informed consent signed as required.
2. Measure and record upper arm circumference.
3. Measure patient to determine length of catheter required to place distal tip in SVC.

Adult: Extend arm approximately 45 to 90 degrees from trunk. Measure distance from insertion site along presumed anatomical course of vessel to be catheterized. Catheter tip should lie in distal one-third of SVC above right atrium and parallel to SVC wall. If a StatLock® Catheter Securement Device will be used, add 1 to 1-1/2 inches (2.5 to 3.8 cm) to catheter measurement.

Pediatric: Measure distance from insertion site along presumed anatomical course of vessel to be catheterized. Catheter tip should lie in distal one-third of SVC above right atrium and parallel to SVC wall. When insertion site is saphenous vein, catheter tip should lie in distal one-third of inferior vena cava (IVC) below right atrium and parallel to IVC wall. If a StatLock® Catheter Securement Device will be used, add 1 to 1-1/2 inches (2.5 to 3.8 cm) to catheter measurement.

4. Position patient for insertion.
Adult: Position patient as appropriate for insertion site. Extend arm laterally 45 to 90 degrees from trunk.

Pediatric: Position patient as appropriate for insertion site.

5. Perform surgical scrub. Dress in protective clothing (mask, goggles, sterile gown, sterile gloves, hair cover).
6. Prepare venipuncture site utilizing appropriate antiseptics.
7. Perform a skin wheal with a local anesthetic as needed. In kits where provided, a SharpsAway® disposal cup is used for disposal of needles. Push needles into foam after use. Discard entire cup at completion of procedure. **Precaution: Do not reuse needles after they have been placed into disposal cup. Particulate matter may adhere to needle tip.**
8. Prepare all equipment. For two-lumen catheter, flush proximal lumen with sterile saline solution. Clamp or attach injection cap to proximal lumen

pigtail. Flush distal lumen and leave syringe in place.

9. Insert introducer needle into vein and aspirate. Ultrasound may be utilized to locate vessel using echogenic tip of needle. **Precaution: The color of blood is not always a reliable indicator of venous access.**¹³ Pulsatile flow is usually an indicator of inadvertent arterial puncture.

10a. Catheter Insertion Using a Single 45 cm Wire for Vessel Access:

- Insert floppy tip of spring-wire guide through introducer needle into vein. Advance spring-wire guide in routine fashion to desired depth. **Warning: Do not cut spring-wire guide to alter length. Do not withdraw spring-wire guide against needle bevel to minimize the risk of possible severing or damaging of spring-wire guide.**

- Hold spring-wire guide in place and remove introducer needle. **Precaution: Maintain firm grip on spring-wire guide at all times.**

- Thread tapered tip of peel-away sheath/dilator assembly over spring-wire guide. Grasping near skin advance assembly with slight twisting motion to a depth sufficient to enter vessel.

Precaution: Do not withdraw dilator until the sheath is well within the vessel to minimize the risk of damage to sheath tip.
Precaution: Sufficient guide-wire length must remain exposed at hub end of sheath to maintain a firm grip on guide-wire.

- Advance peel-away sheath over dilator into vessel, again grasping near skin and using slight twisting motion.

- Holding sheath in place, remove guide-wire and dilator as a unit. **Warning: Do not leave tissue dilator in place as an indwelling catheter to minimize the risk of possible vessel wall perforation.** **Warning: Although the incidence of spring-wire guide failure is extremely low, practitioner should be aware of the potential for breakage if undue force is applied.**

- Advance catheter through peel-away sheath.

10b. Catheter Insertion Using a Single 45 cm Wire for Vessel Access and an 80 cm Wire for Catheter Placement:

- Use 45 cm wire to gain vascular access (see 45 cm wire procedure - 10a.).

- Load PICC catheter onto 80 cm wire until floppy tip of wire extends beyond tip of catheter.

- While maintaining control of distal end of wire, advance floppy tip/catheter tip as a unit through peel-away sheath to desired depth.

- Remove wire.
- Once catheter is in desired location, remove spring-wire guide.

10c. **Catheter Insertion Using a Single 45 cm Wire for Vessel Access and a 130 cm Wire for Catheter Placement:**

- Use 45 cm wire to gain vascular access (see 45 cm wire procedure - 10a.).
 - Insert floppy end of 130 cm spring-wire guide through peel-away sheath to desired depth.
 - Thread catheter over guide-wire and advance catheter over wire through sheath into vessel into correct position. **Precaution: Maintain firm grip on spring-wire guide at all times.**
 - Once catheter is in desired location, remove spring-wire guide.
11. Stop advancing catheter 5 inches (13 cm) before reaching pre-established insertion length.
 12. Withdraw peel-away sheath until free from venipuncture site.
 13. Grasp tabs of peel-away sheath and pull apart, away from catheter, until sheath splits down entire length.
 14. Advance catheter to final indwelling position.
 15. Check catheter placement with attached syringe by aspirating through distal lumen until free flow of venous blood is observed. For two-lumen catheter remove proximal lumen cap and aspirate as permitted by size of catheter lumen. **Precaution: The color of blood aspirated is not always a reliable indicator of venous access.¹³**
 16. Flush lumen(s) with sufficient volume of solution to completely clear blood.
 17. Connect pigtail(s) to appropriate Luer-Lock line(s) as required. Alternately, unused port(s) may be “locked” through injection cap(s) using standard hospital/agency protocol. **Warning: This product contains slide clamps which may be inadvertently removed, and potentially aspirated, by children or confused adults. In such situations, practitioners should remove clamps when not in use.** Slide clamps are provided on pigtails to occlude flow through each lumen during line and injection cap changes. **Precaution: To minimize the risk of damage to pigtail(s) from excessive pressure, each clamp must be opened prior to infusing through that lumen.**
 18. Cleanse insertion site per hospital/agency protocol.
 19. Secure catheter. Where provided, a catheter clamp, fastener, StatLock® Catheter Securement Device or Steri-Strips®† may be used.

StatLock® Catheter Securement Device Instructions:

Initial Application:

- Cleanse and prep anticipated dressing site per hospital/agency protocol. Skin prep should be

applied to coat skin and maximize StatLock® adherence. Allow to dry thoroughly. The anchor pad will be placed so center of pad is within 1 to 1-1/2 inches (2.5 to 3.8 cm) of catheter insertion site.

- The catheter can be secured to StatLock® by using the primary suture hub. **Precaution: Minimize catheter manipulation throughout this procedure to maintain proper catheter tip position. Provide strain relief placement of catheter under dressing.**

Primary Suture Hub Securement:

From earlier measurements, after catheter is properly positioned the primary suture hub is within 1 to 1-1/2 inches (2.5 to 3.8 cm) from insertion site. Place suture hub wings over StatLock® posts and press down (refer to Fig. 2). Snap StatLock® retainer wings to closed position to secure suture hub (refer to Fig. 3).


Fig. 2


Fig. 3

Flexible Catheter Clamp and Rectangular, Rigid Fastener Securement:

If primary suture hub cannot be used on StatLock®, the catheter clamp and fastener (clamp assembly) should be used to secure catheter. Clamp assembly should be applied to area of catheter that lies over StatLock® posts.

- To apply flexible catheter clamp, spread wings of clamp and position on catheter as required to

ensure proper placement over StatLock® posts (refer to Fig. 4).


Fig. 4

- Snap rectangular, rigid fastener onto flexible catheter clamp (refer to Fig. 5).


Fig. 5

- As a unit, snap catheter clamp assembly onto StatLock® posts (refer to Fig. 6).


Fig. 6

- Snap StatLock® retainer wings to closed position to secure suture hub (see Fig. 7).


Fig. 7

- Remove paper backing from one half of StatLock® anchoring pad and press onto dry, prepared skin. Repeat process for other half of StatLock®.
- Complete sterile insertion site dressing according to established hospital/agency protocol.
- Document StatLock®/dressing application on patient's chart.
- Replace StatLock®/dressing per hospital/agency protocol. StatLock® Catheter Securement Device should be replaced at least every 7 days to ensure maximum adherence.

Alternate Technique:

Apply Steri-Strips®† and sterile dressing or apply suture and sterile dressing according to hospital/agency protocol. **Precaution: Do not suture directly to outside diameter of catheter to minimize the risk of cutting or damaging the catheter or impeding catheter flow.**

StatLock® Catheter Securement Device Removal:

- Remove dressing and sutures (if applicable).
 - Open StatLock® retainer wings and remove catheter from StatLock® posts.
 - Apply alcohol swab to StatLock® adhesive and gently lift pad off of skin.
 - Document StatLock® removal procedure on patient's chart.
20. Apply PICC label to dressing.
 21. Obtain chest x-ray immediately to verify tip placement. **Precaution: X-ray exam must show the catheter located in the right side of the mediastinum in the SVC above its junction with the right atrium and parallel to the vessel wall^{10,22} and its distal tip positioned at a level above either the azygos vein or the carina of the trachea, whichever is better visualized.** If catheter tip is malpositioned, reposition and re-verify. **Precaution: If difficulty is encountered in visualizing the catheter tip with x-ray, a small amount of radiopaque dye may be injected to assist in locating the tip.¹³**
 22. Complete adhesive chart label and affix to patient's chart. Document assessments and insertion procedure.
 23. Check catheter tip position routinely using centimeter marks on catheter.

Heparinization:

1. A variety of "locking" solution concentrations may be utilized to maintain patency of catheter. The amount of heparin used, if any, and frequency of flushing depends on physician preference, hospital/agency protocol, and patient condition.^{3,8}
2. The volume of heparin solution should be equal to, or slightly more than, the volume of the lumen that

is being locked. Catheter priming volume is printed on product packaging.

Catheter Removal Procedure:

1. Remove dressing and sutures (if applicable). **Precaution: To minimize the risk of cutting the catheter, do not use scissors to remove dressing.**
2. Open StatLock® retainer wings and remove catheter from StatLock® posts (if applicable).
3. Remove catheter by slowly pulling it parallel to skin. **Precaution: To minimize the risk of catheter breakage, do not exert excessive force if difficulty is encountered upon removal.** If resistance is met, apply heat for 20-30 minutes to area.^{15,24} Gently begin pulling catheter parallel to skin. If further difficulty is encountered, obtain an x-ray and consult physician.
4. Upon removal of catheter, measure and inspect to ensure that entire catheter length has been removed.
5. Apply alcohol swab to StatLock® adhesive and gently lift pad off of skin (if applicable).
6. Dress insertion site.
7. Document catheter removal and StatLock® removal (if applicable) removal procedures on patient's chart.

References:

1. Bar-Joseph G, Galvis AG. Perforation of the heart by central venous catheters in infants: guidelines to diagnosis and management. *J Pediatr Surg.* 1983;18:284-287.
2. Blitt CD, ed. *Monitoring in Anesthesia and Critical Care Medicine.* Central venous pressure monitoring. New York, NY: Churchill Livingstone; 1985:121-165.
3. Bonn J. Venous access: Peripherally inserted central catheters. 1994 Program, SCIVR '94, 19th Annual Scientific Meeting. San Diego, CA
4. Brandt RL, Foley WJ, Fink GH, Regan WJ. Mechanism of perforation of the heart with production of hydropericardium by a venous catheter and its prevention. *Am J Surg.* 1970;119:311-316.
5. Collier PE, Ryan JJ, Diamond DL. Cardiac tamponade from central venous catheters – report of a case and review of the English literature. *Angiology.* 1984;35:595-600.
6. Conn C. The importance of syringe size when using implanted vascular access devices. *J Vasc Access Nurs.* 1993;3:11-18.
7. Curelaru J, Linder LE, Gustavsson B. Displacement of catheters inserted through internal jugular veins with neck flexion and extension. *Intens Care Med.* 1980;6:179-183.

8. Danek GD, Noris EM. Pediatric IV catheters: Efficacy of saline flush. *Pediatric Nursing.* 1992;18:111-113.
9. Frey MA. Pediatric peripherally inserted central catheter program report – a summary of 4,536 catheter days. *JIN.* 1995;18:280-291.
10. Iberti TJ, Katz LB, Reiner MA, Brownie T, Kwun KB. Hydrothorax as a late complication of central venous indwelling catheters. *Surgery.* 1983;842-846.
11. Intravenous Nurses Society (INS). Midline and midclavicular catheters. *JIN.* July/August 1997;20:175-178.
12. James L, Bledsoe L, Hadaway LC. A retrospective look at tip location and complications of peripherally inserted central catheter lines. *JIN.* 1993;16:104-109.
13. Jobs DR, Schwartz AJ, Greenhow DE, Stephenson LW, Ellison N. Safer jugular vein cannulation: recognition of arterial punctures and preferential use of the external jugular route. *Anesthesiology.* 1983;59:353-355.
14. LaFortune S. The use of confirming x-rays to verify tip position for peripherally inserted catheters. *JIN.* 1993;16:246-250.
15. Marx M. The management of the difficult peripherally inserted central venous catheter line removal. *JIN.* 1995;18:246-249.
16. Maschke SP, Rogove HJ. Cardiac tamponade associated with a multilumen central venous catheter. *Crit Care Med.* 1984;12:611-612.
17. Peters JL, ed. *A Manual of Central Venous Catheterization and Parenteral Nutrition.* Boston, Mass: John Wright PSG; 1983:140-146,156.
18. Ryder MA. Peripherally inserted central venous catheters. *Nursing Clinics of North America.* 1993;28:937-971.
19. Ryder M. Total parenteral nutrition. In: Boggs RL, Woodridge-Kin M, eds. *AACN Procedural Manual for Critical Care Nursing.* 3rd ed. Philadelphia, PA: WB Saunders Co: 1993:771-812.
20. Sheep RE, Guiney WB Jr. Fatal cardiac tamponade. *JAMA.* 1982;248:1632-1635.
21. Sigurdsson J, Riba P, Sigurdsson S. The wandering central venous catheter. *Intensive Care Med.* 1985;11:263-264.
22. Tocino IM, Watanabe A. Impending catheter perforation of superior vena cava: radiographic recognition. *Am J Roentgenology.* 1986;146:487-490.
23. Todres ID, Cote CJ. Procedures. In: Cote CJ, Ryan JF, Todres ID, Goudsouzian NG, eds. *A Practice of Anesthesia for Infants and Children.* 2nd ed. Philadelphia, PA: WB Saunders Co; 1993:505-523.

24. Wall JL, Kierstead VL. Peripherally inserted central catheters resistance to removal: A rare complication. *JIN*. 1995;18:251-254.

Arrow International, Inc., recommends that the user be acquainted with the reference literature.

*If you have any questions or would like additional reference information, please contact Arrow International, Inc.

†Steri-Strip® is a registered trademark of 3M company.

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Arrow International, Inc. rekommenderar att användaren studerar referenslitteraturen.
*Om du har frågor eller önskar ytterligare referensinformation, kontakta Arrow International, Inc.