ARROW[®]

Pediatric Single-Lumen Central Venous Catheterization Product

Safety and Efficacy Considerations:

Do not use if package has been previously opened or damaged. **Warning: Prior to use read all package insert warnings, precautions, and instructions. Failure to do so may result in severe patient injury or death.**

Do not alter the catheter, spring-wire guide, or any other kit/set component during insertion, use, or removal.

Procedure must be performed by trained personnel well versed in anatomical landmarks, safe technique, and potential complications.

Warning: Do not place the catheter into or allow it to remain in the right atrium or right ventricle (refer to Fig. 1).


Fig. 1

Cardiac Tamponade: It has been documented by many authors that placement of indwelling catheters in the right atrium is a dangerous practice^{1,2,3,5,6,11,14,15,16,23} that may lead to cardiac perforation and tamponade.^{1,2,3,5,6,11,14,15,16,23} Although cardiac tamponade secondary to pericardial effusion is uncommon, there is a high mortality rate associated with it.²² Practitioners placing central venous catheters must be aware of this potentially fatal complication before advancing the catheter too far relative to patient size.

No particular route or catheter type is exempt from this potentially fatal complication.¹⁸ The actual position of the tip of the indwelling catheter should be confirmed by x-ray after insertion.^{1,2,3,6,11,16,23} Central venous catheters should be placed in the superior vena cava^{1,2,3,6,11,15,16,20,23} above its junction with the right atrium and parallel to the vessel wall and its distal tip positioned at a level above either the azygos vein or the carina of the trachea, whichever is better visualized.

Central venous catheters should not be placed in the right atrium unless specifically required for special relatively short term procedures, such as aspiration of air emboli during neurosurgery. Such procedures are, nevertheless, risk prone and should be closely monitored and controlled.

Indications for Use:

The single-lumen catheter permits venous access to the central circulation of pediatric patients.

Contraindications:

None known.

Warnings and Precautions:*

- Warning: Sterile, Single use: Do not reuse, reprocess or resterilize. Reuse of device creates a potential risk of serious injury and/or infection which may lead to death.**
- Warning: Do not place the catheter into or allow it to remain in the right atrium or right**

ventricle. Central vein catheters should be positioned so that the distal tip of the catheter is in the lower 1/2 to 1/3 of the superior vena cava (SVC) above the junction of the SVC and the right atrium and lies parallel to the vessel wall. For femoral vein approach, the catheter should be advanced into the vessel so that the catheter tip lies parallel to the vessel wall and does not enter the right atrium.

Due to the variety of pediatric catheter lengths available, patient size must be carefully considered relative to actual length of catheter introduced.

3. **Warning:** Practitioners must be aware of complications associated with central vein catheters including **cardiac tamponade** secondary to vessel wall, atrial or ventricular perforation, pleural and mediastinal injuries, air embolism, catheter embolism, thoracic duct laceration, bacteremia, septicemia, thrombosis, inadvertent arterial puncture, nerve damage, hematoma, hemorrhage, and dysrhythmias.
4. **Warning:** Do not apply excessive force in removing guide wire or catheters. If withdrawal cannot be easily accomplished, a chest x-ray should be obtained and further consultation requested.
5. **Warning:** The practitioner must be aware of potential air embolism associated with leaving open needles or catheters in central venous puncture sites or as a consequence of inadvertent disconnects. To lessen the risk of disconnects, only securely tightened Luer-Lock connections should be used with this device. Follow hospital protocol to guard against air embolism for all catheter maintenance.
6. **Warning:** Passage of the guide wire into the right heart can cause dysrhythmias, right bundle branch block,⁹ and vessel wall, atrial or ventricular perforation.
7. **Warning:** Due to the risk of exposure to HIV (Human Immunodeficiency Virus) or other blood borne pathogens, health care workers should routinely use universal blood and body-fluid precautions in the care of all patients.
8. **Precaution:** Indwelling catheters should be routinely inspected for desired flow rate, security of dressing, correct catheter position, and for secure Luer-Lock connection. Use centimeter markings to identify if the catheter position has changed.

9. **Precaution:** Only x-ray examination of the catheter placement can ensure that the catheter tip has not entered the heart or no longer lies parallel to the vessel wall. If catheter position has changed, immediately perform chest x-ray examination to confirm catheter tip position.

10. **Precaution:** Alcohol and acetone can weaken the structure of polyurethane materials. Check ingredients of prep sprays and swabs for acetone and alcohol content.

Acetone: Do not use acetone on catheter surface. Acetone may be applied to skin but must be allowed to dry completely prior to applying dressing.

Alcohol: Do not use alcohol to soak catheter surface or to restore catheter patency. Care should be taken when instilling drugs containing high concentration of alcohol. Always allow alcohol to dry completely prior to applying dressing.

11. **Precaution:** Some disinfectants used at the catheter insertion site contain solvents, which can attack the catheter material. Assure insertion site is dry before dressing.

12. **Precaution:** Use of a syringe smaller than 10 ml to irrigate or declot an occluded catheter may cause intraluminal leakage or catheter rupture.⁷

A Suggested Procedure: Use sterile technique.

1. **Precaution:** Place patient in slight Trendelenburg position as tolerated to reduce the risk of air embolism. If femoral approach is used, place patient in supine position.
2. Measure the patient's external anatomy to estimate the length of catheter required to place the distal tip in the SVC.
3. Prep and drape puncture site as required.
4. Perform skin wheal with desired needle.
5. **Precaution: Remove catheter stylet.** Prepare the catheter for insertion by flushing the lumen. Leave the catheter uncapped for guide wire passage. **Warning: Do not cut the catheter to alter length.**
6. Locate central vein with a needle and syringe appropriate to the size of the patient.
7. Insert introducer catheter/needle with attached syringe into vein beside locator needle and aspirate. Remove locator needle. Withdraw needle from introducer catheter. If no free flow of venous blood is observed after needle is removed, attach syringe

to the catheter and aspirate until good venous blood flow is established. **Precaution: The color of the blood aspirated is not always a reliable indicator of venous access.¹² Do not reinsert needle into introducer catheter.** Because of the potential for inadvertent arterial placement, verify venous access via a wave form obtained by a calibrated pressure transducer (refer to Fig. 2).


Fig. 2

If hemodynamic monitoring equipment is not available to permit transducing a central venous wave form, disconnect the syringe and check for pulsatile flow. Pulsatile flow is usually an indicator of inadvertent arterial puncture.

Alternate Technique:

Introducer needle may be used in the standard manner as alternative to catheter/needle assembly.

8. Using the two-piece Arrow Advancer™, advance spring-wire guide through guide wire introducer needle or catheter into vein. **Precaution: Due to the fragile nature of the spring-wire guide contained in this product, inspect “J”-tip for damage prior to insertion.** Advance spring-wire guide to a depth appropriate to patient size. Advancement of “J” tip may require a gentle rotating motion. **Warning: Do not cut spring-wire guide to alter length. Do not withdraw spring-wire guide against needle bevel to minimize the risk of possible severing or damaging of spring-wire guide.**

Arrow Two-Piece Advancer™ Instructions:

- Using your thumb, straighten the “J” by retracting the spring-wire guide into the Advancer™ (refer to Figs. 3, 4).


Fig. 3


Fig. 4

When the tip is straightened, the spring-wire guide is ready for insertion.

Alternate Technique:

If a simple straightening tube is preferred, the straightening tube portion of the Advancer™ can be disconnected from the unit and used separately.

Separate the Advancer™ tip or straightening tube from the blue Advancer™ unit. If the “J”-tip portion of the spring-wire guide is used, prepare for insertion by sliding the plastic tube over the “J” to straighten. The spring-wire guide should then be advanced in the routine fashion to the desired depth.

9. Hold spring-wire guide in place and remove introducer needle or catheter. **Precaution: Maintain firm grip on spring-wire guide at all times.**
10. Enlarge cutaneous puncture site with cutting edge of scalpel positioned away from the spring-wire guide. **Precaution: Do not cut guide wire.** Use vessel dilator to enlarge site as required. **Warning: Do not leave vessel dilator in place as an indwelling catheter to minimize the risk of possible vessel wall perforation.**
11. Thread tip of catheter over spring-wire guide. Sufficient guide wire length must remain exposed at hub end of catheter to maintain a firm grip on guide wire. Grasping near skin, advance catheter into vein with slight twisting motion. **Precaution: If provided, catheter clamp and fastener must not be attached to catheter until spring-wire guide is removed.**
12. Advance the catheter to the final indwelling position using centimeter marks (where provided) as positioning reference points.
13. Hold catheter at desired depth and remove spring-wire guide. The Arrow catheter included in this product has been designed to freely pass over the spring-wire guide. If resistance is encountered

when attempting to remove the spring-wire guide after catheter placement, the spring-wire may be kinked about the tip of the catheter within the vessel (refer to Fig. 5).


Fig. 5

In this circumstance, pulling back on the spring-wire guide may result in undue force being applied resulting in spring-wire guide breakage. If resistance is encountered, withdraw the catheter relative to the spring-wire guide about 2-3 cm and attempt to remove the spring-wire guide. If resistance is again encountered remove the spring-wire guide and catheter simultaneously. **Warning: Although the incidence of spring-wire guide failure is extremely low, practitioner should be aware of the potential for breakage if undue force is applied to the wire.**

14. Verify that the entire spring-wire guide is intact upon removal.
15. Check lumen placement by attaching a syringe to the catheter and aspirate until free flow of venous blood is observed. Connect appropriate Luer-Lock line, as required, or “lock” through injection cap using standard hospital protocol.
16. Secure and dress catheter temporarily.
17. Verify catheter tip position by chest x-ray immediately after placement. **Precaution: X-ray exam must show the catheter located in the right side of the mediastinum in the SVC with the distal end of the catheter parallel to the vena cava wall and its distal tip positioned at a level above either the azygos vein or the carina of the trachea, whichever is better visualized.** If catheter tip is malpositioned, reposition and re-verify.
18. Secure catheter to patient. Use integral suture ring and/or side wings as primary suture site. In kits where provided, the catheter clamp and fastener should be utilized as a secondary suture site as necessary. **Precaution: Do not suture directly to the outside diameter of the catheter to minimize the risk of cutting or damaging the catheter or impeding catheter flow.**

Catheter Clamp and Fastener (where provided)

Instructions for Use:

- After spring-wire guide has been removed and the line has been connected or locked, spread wings of rubber clamp and position on catheter, as required, to ensure proper tip location (refer to Fig. 6).


Fig. 6

- Snap rigid fastener onto catheter clamp (refer to Fig. 7).


Fig. 7

- Secure catheter to patient by suturing the catheter clamp and fastener to the skin, using side wings to minimize the risk of catheter migration (refer to Fig. 8).


Fig. 8

19. Dress puncture site per hospital protocol. **Precaution: Maintain the insertion site with regular meticulous redressing using aseptic technique.**

- Record on the patient's chart the indwelling catheter length using centimeter markings (where provided) on catheter where it enters the skin. Frequent visual reassessment should be made to ensure that the catheter has not moved.

Catheter Exchange Procedure:

- Use sterile technique.
- Precaution: Prior to attempting a catheter exchange procedure, remove catheter clamp and fastener, where provided.**
- Proceed per hospital protocol. Cutting the catheter is not recommended due to the potential for catheter embolism.

Catheter Removal Procedure:

- Precaution: Place patient in a supine position.**
- Remove dressing. **Precaution: To minimize the risk of cutting catheter, do not use scissors to remove dressing.**
- Warning: Exposure of central vein to atmospheric pressure may result in entry of air into central venous system.** Remove suture(s) from catheter clamp, where provided, and primary suture site. Be careful not to cut catheter. Remove catheter slowly, pulling parallel to skin. As catheter exits site, apply pressure with a dressing impermeable to air, e.g. Vaseline® gauze. Because residual catheter track remains an air entry point until completely sealed, the occlusive dressing should remain in place for at least 24-72 hours dependent upon amount of time catheter was indwelling.^{13,17,19,24}
- Upon removal of catheter, inspect it to make sure entire length has been withdrawn.
- Document removal procedure.

References:

- Agarwal KC, Khan MAA, Falla A, Amato JJ. Cardiac perforation from central venous catheters: survival after cardiac tamponade in an infant. *Pediatrics*. March 1984;73:333-338.
- Bar-Joseph G, Galvis AG. Perforation of the heart by central venous catheters in infants: guidelines to diagnosis and management. *J Pediatr Surg*. June 1983;18:284-287.
- Brandt RL, Foley WJ, Fink GH, Regan WJ. Mechanism of perforation of the heart with production of hydropericardium by a venous catheter and its prevention. *Am J Surg*. 1970;119:311-316.

- Carey BE. Major complications of central lines in neonates. *Neonatal Network*. June 1989;7:17-28.
- Casado-Flores J, Valdiviuelso-Serna A, Perez-Jurado L, et al. Subclavian vein catheterization in critically ill children: analysis of 322 cannulations. *Intensive Care Med*. 1991;17:350-354.
- Collier PE, Ryan JJ, Diamond DL. Cardiac tamponade from central venous catheters – report of a case and review of the English literature. *Angiology*. September 1984;35:595-600.
- Conn C. The importance of syringe size when using implanted vascular access devices. *J Vasc Access Nurs*. Winter 1993;3:11-18.
- Daniels SR, Hannon DW, Meyer RA, Kaplan S. Paroxysmal supraventricular tachycardia. *AJDC*. 1984;138:474-475.
- Eissa NT, Kvetan V. Guide wire as a cause of complete heart block in patients with preexisting left bundle branch block. *Anesthesiology*. 1990;73:772-774.
- Franciosi RA, Ellefson RD, Uden D, Drake RM. Sudden unexpected death during central hyperalimentation. *Pediatrics*. March 1982;69:305-307.
- Fuhrman BP. Perforation by vascular catheters. *Perspectives in Critical Care*. 1988;1:147-154.
- Jobs DR, Schwartz AJ, Greenhow DE, Stephenson LW, Ellison N. Safer jugular vein cannulation: recognition of arterial punctures and preferential use of the external jugular route. *Anesthesiology*. 1983;59:353-355.
- Kashuk JL, Penn I. Air embolism after central venous catheterization. *Surg Gynecol Obstet*. September 1984;159:249-252.
- Leibovitz E, Ashkenazi A, Levin S, Nissim F. Fatal pericardial tamponade complicating total parenteral nutrition via a silastic central vein catheter. *J Pediatr Gastroenterol Nutr*. 1988;7:306-307.
- Mupanemunda RH, Mackanjee HR. A life-threatening complication of percutaneous central venous catheters in neonates. *AJDC*. 1992;146:1414-1415. Letter.
- Otto CW. Central venous pressure monitoring. In: Blitt CD, ed. *Monitoring in Anesthesia and Critical Care Medicine*. New York, NY: Churchill Livingstone, 1985:139.
- Paskin DL, Hoffman WS, Tuddenham WJ. A new complication of subclavian vein catheterization. *Ann Surg*. March 1974;179:266-268.

18. Peters JL, ed. *A Manual of Central Venous Catheterization and Parenteral Nutrition*. Boston, MA: John Wright PSG, 1983:58-61, 155-157.
 19. Phifer TJ, Bridges M, Conrad SA. The residual central venous catheter track – an occult source of lethal air embolism: case report. *J Trauma*. 1991;31:1558-1560.
 20. Rogers BB, Berns SD, Maynard EC, Hansen TWR. Pericardial tamponade secondary to central venous catheterization and hyperalimentation in a very low birth weight infant. *Pediatr Pathology*. 1990;10:819-823.
 21. Ross P Jr, Seashore JH. Bilateral hydrothorax complicating central venous catheterization in a child: case report. *J Pediatr Surg*. 1989;24:263-264.
 22. Sheep RE, Guiney WB Jr. Fatal cardiac tamponade. *JAMA*. 1982;248:1632-1635.
 23. Sullivan CA, Konefal SH Jr. Cardiac tamponade in a newborn: a complication of hyperalimentation. *J Parenteral Enteral Nutrition*. May/June 1987;11:319-321.
 24. Thielen JB, Nyquist J. Subclavian catheter removal. *J Intravenous Nurs*. March/April 1991;14:114-118.
- Arrow International, Inc. recommends that the user be acquainted with the reference literature.
- * If you have any questions or would like additional reference information, please contact Arrow International, Inc.

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Arrow International, Inc. rekommenderar att användaren studerar referenslitteraturen.
*Om du har frågor eller önskar ytterligare referensinformation, kontakta Arrow International, Inc.