ARROW®

INTERNATIONAL

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Single-Lumen Central Venous Catheterization Product with Arrow-Johans™ RAECG Adapter

Safety and Efficacy Considerations:

Do not use if package has been previously opened or damaged.

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Do not alter the catheter, spring-wire guide, or any other kit/set component during insertion, use, or removal.

Procedure must be performed by trained personnel well versed in anatomical landmarks, safe technique, and potential complications.

Warning: Do not place the catheter into or allow it to remain in the right atrium or right ventricle. If the catheter is used for aspiration of air emboli, upon completion the tip must be withdrawn to proper placement within the SVC (refer to Fig. 1).

Fig. 1

Cardiac Tamponade: It has been documented by many authors that placement of indwelling catheters in the right atrium is a dangerous practice^{1,2,3,5,7,14,17} that may lead to cardiac perforation and tamponade^{1,2,3,5,14,17}. Although cardiac tamponade secondary to pericardial effusion is uncommon, there is a high mortality rate associated with it¹⁹. Practitioners placing central venous catheters must be aware of this potentially fatal complication before advancing the catheter too far relative to patient size.

No particular route or catheter type is exempt from this potentially fatal complication¹⁷. The actual position of the tip of the indwelling catheter should be confirmed by x-ray after insertion^{1,2,5,14,17,20}. Central venous catheters should be placed in the superior vena cava^{1,2,3,5,7,14,22} above its junction with the right atrium and parallel to the vessel wall^{19,22} and its distal tip positioned at a level above either the azygos vein or the carina of the trachea, whichever is better visualized.

Central venous catheters should not be placed in the right atrium unless specifically required for special relatively short term procedures, such as aspiration of air emboli during neurosurgery. Such procedures are nevertheless risk prone and should be closely monitored and controlled.

Indications for Use:

The single-lumen catheter permits venous access to the central circulation.

Contraindications:

None known.

Warnings and Precautions:*

- Warning: Sterile, Single use:** Do not reuse, reprocess or sterilize. Reuse of device creates a potential risk of serious injury and/or infection which may lead to death.
- Warning: For general central venous use, do not place the catheter into or allow it to remain in the right atrium or right ventricle.** Central vein catheters should be positioned so that the distal tip of the catheter is in the superior vena cava (SVC) above the junction of the SVC and the right atrium and lies parallel to the vessel wall. For femoral vein approach, the catheter should be advanced into the vessel so that the catheter tip lies parallel to the vessel wall and does not enter the right atrium.
- Warning: Practitioners must be aware of complications associated with central vein catheters including cardiac tamponade secondary to vessel wall, atrial or ventricular perforation, pleural and mediastinal injuries, air embolism, catheter embolism, thoracic duct laceration, bacteremia, septicemia, thrombosis, inadvertent arterial puncture, nerve damage, hematoma, hemorrhage, and dysrhythmias.** If catheter is to be used for aspiration of air emboli, practitioner must be aware of complications associated with right atrial and inadvertent right ventricular catheterization before intentionally passing the catheter tip into the right atrium. These complications include dysrhythmia and cardiac tamponade secondary to pericardial effusion of fluids or blood. Do not allow the catheter to remain in the right atrium or right ventricle after it is no longer needed for aspiration of air emboli.
- Warning: Do not apply excessive force in removing guide wire or catheters.** If withdrawal cannot be easily accomplished, a chest x-ray should be obtained and further consultation requested.
- Warning: The practitioner must be aware of potential air embolism associated with leaving open needles or catheters in central venous puncture sites or as a consequence of inadvertent disconnects.** To lessen the risk of disconnects, only securely tightened Luer-Lock connections should be used with this device. Follow hospital protocol to guard against air embolism for all catheter maintenance.
- Warning: Passage of the guide wire into the right heart can cause dysrhythmias, right bundle branch block,⁸ and vessel wall, atrial or ventricular perforation.**
- Warning: Due to the risk of exposure to HIV (Human Immunodeficiency Virus) or other blood borne pathogens, health care workers should routinely use universal blood and body-fluid precautions in the care of all patients.**

8. **Precaution:** Indwelling catheters should be routinely inspected for desired flow rate, security of dressing, correct catheter position, and for secure Luer-Lock connection. Use centimeter markings to identify if the catheter position has changed.
9. **Precaution:** Only x-ray examination of the catheter placement can ensure that the catheter tip has not entered the heart or no longer lies parallel to the vessel wall in catheters left indwelling. If catheter position has changed, immediately perform chest x-ray examination to confirm catheter tip position.
10. **Precaution:** Alcohol and acetone can weaken the structure of polyurethane materials. Check ingredients of prep sprays and swabs for acetone and alcohol content. Acetone: Do not use acetone on catheter surface. Acetone may be applied to skin but must be allowed to dry completely prior to applying dressing. Alcohol: Do not use alcohol to soak catheter surface or to restore catheter patency. Care should be taken when instilling drugs containing high concentration of alcohol. Always allow alcohol to dry completely prior to applying dressing.
11. **Precaution:** Some disinfectants used at the catheter insertion site contain solvents, which can attack the catheter material. Assure insertion site is dry before dressing.
12. **Precaution:** Use of a syringe smaller than 10 mL to irrigate or declot an occluded catheter may cause intraluminal leakage or catheter rupture.⁶

A Suggested Procedure:

Use sterile technique.

1. **Precaution:** Place patient in slight Trendelenburg position as tolerated to reduce the risk of air embolism. If femoral approach is used, place patient in supine position.
2. Prep and drape puncture site as required.
3. Perform skin wheal with desired needle (25 Ga. or 22 Ga. needle). In kits where provided, a SharpsAway® disposal cup is used for the disposal of needles. Push needles into foam after use. Discard entire cup at completion of procedure. **Precaution: Do not reuse needles after they have been placed into the disposal cup. Particulate matter may adhere to needle tip.**
4. Locate central vein with a 22 Ga. needle and syringe.
5. Insert introducer catheter/needle with attached syringe into vein beside locator needle and aspirate. Remove locator needle. Withdraw needle from introducer catheter. If no free flow of venous blood is observed after needle is removed, attach syringe to the catheter and aspirate until good venous blood flow is established. **Precaution: The color of the blood aspirated is not always a reliable indicator of venous access.¹⁰ Do not reinsert needle into introducer catheter.** Because of the potential for inadvertent arterial placement, verify venous access via a wave form obtained by a calibrated pressure transducer (refer to Fig. 2).


Fig. 2

If hemodynamic monitoring equipment is not available to permit transducing a central venous wave form, disconnect the syringe and check for pulsatile flow. Pulsatile flow is usually an indicator of inadvertent arterial puncture.

Alternate Technique:

Introducer needle may be used in the standard manner as alternative to catheter/needle assembly.

6. Insert desired tip of spring-wire guide through introducer needle or catheter into vein. If the "J"-tip portion of the spring-wire guide is used, prepare for insertion by sliding the plastic tube over the "J" to straighten. The spring-wire guide should then be advanced in the routine fashion to the desired depth. Advancement of "J" tip may require a gentle rotating motion. **Warning: Do not cut spring-wire guide to alter length. Do not withdraw spring-wire guide against needle bevel to minimize the risk of possible severing or damaging of spring-wire guide.**
7. Hold spring-wire guide in place and remove introducer needle or catheter. **Precaution: Maintain firm grip on spring-wire guide at all times.**
8. Enlarge cutaneous puncture site with cutting edge of scalpel positioned away from the spring-wire guide. **Precaution: Do not cut guide wire.** Use vessel dilator to enlarge site as required. **Warning: Do not leave vessel dilator in place as an indwelling catheter to minimize the risk of possible vessel wall perforation.**
9. Thread tip of catheter over spring-wire guide. Sufficient guide wire length must remain exposed at hub end of catheter to maintain a firm grip on guide wire. Grasping near skin, advance catheter into vein with slight twisting motion. **Precaution: If provided, catheter clamp and fastener must not be attached to catheter until spring-wire guide is removed.**
10. Advance the catheter to the final indwelling position using centimeter marks (where provided) as positioning reference points.
11. Hold catheter at desired depth and remove spring-wire guide. The Arrow® catheter included in this product has been designed to freely pass over the spring-wire guide. If resistance is encountered when attempting to remove the spring-wire guide after catheter placement, the spring-wire may be kinked about the tip of the catheter within the vessel (refer to Fig. 3).


Fig. 3

In this circumstance, pulling back on the spring-wire guide may result in undue force being applied resulting in spring-wire guide breakage. If resistance is encountered, withdraw the catheter relative to the spring-wire guide about 2-3 cm and attempt to remove the spring-wire guide. If resistance is again encountered, remove the spring-wire guide and catheter simultaneously. **Warning: Although the incidence of spring-wire guide failure is extremely low, practitioner should be aware of the potential for breakage if undue force is applied to the wire.**

12. Verify that the entire spring-wire guide is intact upon removal.
13. Check lumen placement by attaching a syringe to catheter and aspirate until free flow of venous blood is observed. Catheter tip placement can be achieved with Right Atrial Electrocardiography (RAECG) using the Arrow-Johans™ RAECG adapter.^{11,15}

Arrow-Johans™ RAECG

Adapter Electro-Shock Warnings:

1. **Warning: Strict electrical isolation of monitoring equipment and good grounding techniques are mandatory when using a saline-filled catheter as a transvenous ECG electrode. All equipment should comply with American Heart Association recommendation of being incapable of delivering more than 10 microamp leakage current to the patient. An electrocautery protection filter should be inserted between ECG cable and ECG machine. Rubber cap must be over stainless steel male RAECG button when adapter is in-line but not in use.**
2. **Warning: Physician must be aware of complications associated with right atrial and inadvertent right ventricular catheterization whenever central venous catheter being introduced is intentionally passed beyond junction of superior vena cava and right atrium. The Arrow-Johans™ adapter should be used to monitor the initial position of the CVC.**
3. **Warning: Any premature heart beat might indicate catheter touching heart; withdraw catheter and check position.**
4. **Warning: Final catheter position to be confirmed by chest x-ray.**

Right Atrial Electrocardiography

Instructions for Use:

- Sterile Operator – Pass distal or female end of a sterile length of IV tubing to non-sterile assistant.
- Non-sterile Assistant – Luer-Lock adapter to IV tubing and flush system with an ionic salt solution, e.g. 0.9% saline solution, 4% sodium chloride or 8.4% sodium bicarbonate solution. Meticulously remove all air bubbles and blood clots.

- Sterile Operator – Luer-Lock male or proximal end of sterile IV tubing to catheter.
- Non-sterile Assistant – Flush entire assembly.
- Sterile Operator – Remove rubber protective cap from stainless steel male button on adapter.
- Non-sterile Assistant – For ECG machines without lead selectors, connect negative terminal lead of ECG to adapter, the positive terminal to left leg, and the ground terminal lead to left arm. For ECG machines with lead selectors, choose lead II on ECG machine and place right arm lead to adapter.

- Sterile Operator – Advance or withdraw saline filled catheter until size of the P wave of ECG trace just begins to progressively get larger. At this point, catheter tip should be at cavoatrial junction. If tip is advanced further into right atrium, P wave will either be as large as R wave or become biphasic. If catheter enters right ventricle, P wave will disappear and a wider QRS will appear.

Trouble Shooting Tips:

- Salt bridge solution – The solution in the catheter must be highly ionic to conduct the electrical impulse. Do not use dextrose in water since it is non-ionic.
 - Flush fluids – Discontinue all flush fluids to minimize the risk of electrical deflection.
 - Grounding – Ensure good skin connection with ECG electrodes to optimize ECG tracing.
 - Internal catheter movement – Movement of chest during respiration or movement of heart during contraction can cause catheter to move and produce an undulating trace. Repositioning chest leads may alleviate problem.
14. Connect appropriate Luer-Lock line, as required, or “lock” through injection cap using standard hospital protocol. Where provided, slide clamp is used on extension line to occlude flow through catheter during line and injection cap changes. **Precaution: To minimize the risk of damage to extension line from excessive pressure, clamp must be opened prior to infusing through catheter.**
 15. Secure and dress catheter temporarily.
 16. Verify catheter tip position by chest x-ray immediately after placement. **Precaution: For general central venous use, x-ray exam must show the catheter located in the right side of the mediastinum in the SVC with the distal end of the catheter parallel to the vena cava wall and its distal tip positioned at a level above either the azygos vein or the carina of the trachea, whichever is better visualized. If catheter tip is malpositioned, reposition and re-verify.**
 17. Secure catheter to patient. Use integral suture ring and/or side wings as primary suture site. In kits where provided, the catheter clamp and fastener should be utilized as a secondary suture site as necessary. **Precaution: Do not suture directly to the outside diameter of the catheter to minimize the risk of cutting or damaging the catheter or impeding catheter flow.**

Catheter Clamp and Fastener (where provided)

Instructions for Use:

- After spring-wire guide has been removed and the line has been connected or locked, spread wings of rubber clamp and position on catheter, as required, to ensure proper tip location (refer to Fig. 4).


Fig. 4

- Snap rigid fastener onto catheter clamp (refer to Fig. 5).


Fig. 5

- Secure catheter to patient by suturing the catheter clamp and fastener to the skin, using side wings to minimize the risk of catheter migration (refer to Fig. 6).


Fig. 6

18. Dress puncture site per hospital protocol. **Precaution: Maintain the insertion site with regular meticulous redressing using aseptic technique.**
19. Record on the patient's chart the indwelling catheter length using centimeter markings (where provided) on catheter where it enters the skin. Frequent visual reassessment should be made to ensure that the catheter has not moved.

Catheter Exchange Procedure:

1. Use sterile technique.
2. **Precaution: Prior to attempting a catheter exchange procedure, remove catheter clamp and fastener, where provided.**
3. Proceed per hospital protocol. Cutting the catheter is not recommended due to the potential for catheter embolism.

Catheter Removal Procedure:

1. **Precaution: Place patient in a supine position.**
2. Remove dressing. **Precaution: To minimize the risk of cutting catheter, do not use scissors to remove dressing.**
3. **Warning: Exposure of central vein to atmospheric pressure may result in entry of air into central venous system.** Remove suture(s) from catheter clamp, where provided, and primary suture site. Be careful not to cut catheter. Remove catheter slowly, pulling parallel to skin. As catheter exits site, apply pressure with a dressing impermeable to air, e.g. VASELINE®[†] gauze. Because residual catheter track remains an air entry point until completely sealed, the occlusive dressing should remain in place for at least 24-72 hours dependent upon amount of time catheter was indwelling.^{12,16,18,21}
4. Upon removal of catheter, inspect it to make sure entire length has been withdrawn.
5. Document removal procedure.

References:

1. Bar-Joseph G, Galvis AG. Perforation of the heart by central venous catheters in infants: guidelines to diagnosis and management. *J Pediatr Surg.* 1983;18:284-287.
2. Blitt CD, ed. *Monitoring in Anesthesia and Critical Care Medicine.* Central venous pressure monitoring. New York, NY: Churchill Livingstone; 1985:121-165.
3. Brandt RL, Foley WJ, Fink GH, Regan WJ. Mechanism of perforation of the heart with production of hydropericardium by a venous catheter and its prevention. *Am J Surg.* 1970;119:311-316.
4. Carbone K, Gimenez LF, Rogers WH, Watson AJ. Hemothorax due to vena caval erosion by a subclavian dual-lumen dialysis catheter. *South Med J.* 1987;80:795-796.
5. Collier PE, Ryan JJ, Diamond DL. Cardiac tamponade from central venous catheters – report of a case and review of the English literature. *Angiology.* September 1984;35:595-600.
6. Conn C. The importance of syringe size when using implanted vascular access devices. *J Vasc Access Nurs.* Winter 1993;3:11-18.
7. Curelaru J, Linder LE, Gustavsson B. Displacement of catheters inserted through internal jugular veins with neck flexion and extension. *Intens Care Med.* 1980;6:179-183.
8. Eissa NT, Kvetan V. Guide wire as a cause of complete heart block in patients with preexisting left bundle branch block. *Anesthesiology.* 1990;73:772-774.
9. Iberti TJ, Katz LB, Reiner MA, Brownie T, Kwun KB. Hydrothorax as a late complication of central venous indwelling catheters. *Surgery.* November 1983;842-846.
10. Jobes DR, Schwartz AJ, Greenhow DE, Stephenson LW, Ellison N. Safer jugular vein cannulation: recognition of arterial punctures and preferential use of the external jugular route. *Anesthesiology.* 1983;59:353-355.
11. Johans TG. Multiorificated catheter placement with an intravascular electrocardiographic technique. [letter] *Anesthesiology.* 1986;64:411-413.

12. Kashuk JL, Penn I. Air embolism after central venous catheterization. *Surg Gynecol Obstet.* September 1984;159:249-252.
 13. Kozeny GA, Bansal VK, Vertuno LL, Hano JE. Contralateral hemothorax secondary to chronic subclavian dialysis catheter. *Am J Nephrol.* 1984;4:312-314.
 14. Maschke SP, Rogove HJ. Cardiac tamponade associated with a multilumen central venous catheter. *Crit Care Med.* 1984;12:611-612.
 15. McGee W, Ackerman BL, Rouben LR, Prasad VM, Bandi V, Mallory DL. Accurate placement of central venous catheters: a prospective, randomized, multicenter trial. *Crit Care Med.* 1993;21:1118-1123.
 16. Paskin DL, Hoffman WS, Tuddenham WJ. A new complication of subclavian vein catheterization. *Ann Surg.* March 1974;179:266-268.
 17. Peters JL, ed. *A Manual of Central Venous Catheterization and Parenteral Nutrition.* Boston, MA: John Wright PSG; 1983:58-61, 155-157.
 18. Phifer TJ, Bridges M, Conrad SA. The residual central venous catheter track – an occult source of lethal air embolism: case report. *J Trauma.* 1991;31:1558-1560.
 19. Sheep RE, Guiney WB Jr. Fatal cardiac tamponade. *JAMA.* 1982;248:1632-1635.
 20. Sigurdsson J, Riba P, Sigurdsson S. The wandering central venous catheter. *Intensive Care Med.* 1985;11:263-264.
 21. Thielen JB, Nyquist J. Subclavian catheter removal. *J Intravenous Nurs.* March/April 1991;14:114-118.
 22. Tocino IM, Watanabe A. Impending catheter perforation of superior vena cava: radiographic recognition. *Am J Roentgenology.* March 1986;146:487-490.
 23. Vaziri ND, Maksy M, Lewis M, Martin D, Edwards K. Massive mediastinal hematoma caused by a double-lumen subclavian catheter. *Artif Organs.* 1984;8:223-226.
 24. Wanschler M, Frifelt JJ, Smith-Sivertsen C, et al. Thrombosis caused by polyurethane double-lumen subclavian superior vena cava catheter and hemodialysis. *Crit Care Med.* 1988;16:624-628.
- Arrow International, Inc. recommends that the user be acquainted with the reference literature.
- *If you have any questions or would like additional reference information, please contact Arrow International, Inc.
- † A registered trademark of Unilever Supply Chain, Inc.

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Arrow International, Inc. kullanıcının referans literatüre aşina olmasını önerir.
*Sorumuz varsa veya ek referans bilgisi isterseniz lütfen Arrow International, Inc. ile irtibat kurun.