ARROW®

INTERNATIONAL

ENG	Single-Lumen Central Venous Catheterization Product	1
CZ	Jednolumenové zařízení pro centrální žilní katetrizaci	5
DK	Centralvenekateteriseringsprodukt med enkelt lumen	8
NL	Enkellumen product voor centraal-veneuze katheterisatie	11
F	Produit de cathétérisme veineux central à lumière unique	14
D	Katheterisierungsbesteck mit einlumigem zentralen Venenkatheter	17
GR	Προϊόν κεντρικού φλεβικού καθετηριασμού μονού αυλού	20
H	Egylumenű centrális vénás katéter	23
I	Prodotto a lume singolo per cateterismo della vena centrale	26
PL	Jednoświatłowy cewnik cewnikowania żył centralnych	29
P	Produto para Cateterização Venosa Central de Lúmen Simples	32
RUS	Однопросветное изделие для катетеризации главных вен	35
SK	Jednolúmenový centrálny venózny katetrizačný produkt	38
E	Dispositivo de una sola luz para cateterismo venoso central	41
S	Enkanalig kateter för central ven	44

Single-Lumen Central Venous Catheterization Product

Safety and Efficacy Considerations:

Do not use if package has been previously opened or damaged.

Warning: Prior to use read all package insert warnings, precautions, and instructions. Failure to do so may result in severe patient injury or death.

Do not alter the catheter, spring-wire guide, or any other kit/set component during insertion, use, or removal.

Procedure must be performed by trained personnel well versed in anatomical landmarks, safe technique, and potential complications.

Warning: Do not place the catheter into or allow it to remain in the right atrium or right ventricle (refer to Fig. 1).

No particular route or catheter type is exempt from this potentially fatal complication.¹⁵ The actual position of the tip of the indwelling catheter should be confirmed by x-ray after insertion.^{1,2,3,13,15,18} Central venous catheters should be placed in the superior vena cava^{1,2,3,5,7,13,20} above its junction with the right atrium and parallel to the vessel wall^{9,20} and its distal tip positioned at a level above either the azygos vein or the carina of the trachea, whichever is better visualized.

Central venous catheters should not be placed in the right atrium unless specifically required for special relatively short-term procedures such as aspiration of air emboli during neurosurgery. Such procedures are, nevertheless, risk prone and should be closely monitored and controlled.

Indications for Use:

The Arrow Single-Lumen Catheter permits venous access to the central circulation.

Contraindications:

None Known

Warnings and Precautions:*

- Warning:** Sterile, Single use: Do not reuse, reprocess or resterilize. Reuse of device creates a potential risk of serious injury and/or infection which may lead to death.
- Warning:** Do not place the catheter into or allow it to remain in the right atrium or right ventricle. Central vein catheters should be positioned so that the distal tip of the catheter is in the superior vena cava (SVC) above the junction of the SVC and the right atrium and lies parallel to the vessel wall. For femoral vein approach, the catheter should be advanced into the vessel so that the catheter tip lies parallel to the vessel wall and does not enter the right atrium.
- Warning:** Practitioners must be aware of complications associated with central vein catheters including cardiac tamponade secondary to vessel wall, atrial or ventricular perforation, pleural and mediastinal injuries, air embolism, catheter embolism, catheter occlusion, thoracic duct laceration, bacteremia, septicemia, thrombosis, inadvertent arterial puncture, nerve damage, hematoma, hemorrhage, and dysrhythmias.
- Warning:** Do not apply excessive force in removing guide wire or catheters. If withdrawal cannot be easily accomplished, a chest x-ray should be obtained and further consultation requested.
- Warning:** The practitioner must be aware of potential air embolism associated with leaving open needles or catheters in central venous puncture sites or as a consequence of inadvertent disconnects. To lessen the risk of disconnects, only securely tightened Luer-Lock connections should be used with this device. Follow hospital protocol to guard against air embolism for all catheter maintenance.
- Warning:** Passage of the guide wire into the right heart can cause dysrhythmias, right bundle branch block,⁸ and vessel wall, atrial or ventricular perforation.


Fig. 1

Cardiac Tamponade: It has been documented by many authors that placement of indwelling catheters in the right atrium is a dangerous practice^{1,2,3,5,7,13,15} that may lead to cardiac perforation and tamponade.^{1,2,3,5,13,15} Although cardiac tamponade secondary to pericardial effusion is uncommon, there is a high mortality rate associated with it.¹⁷ Practitioners placing central venous catheters must be aware of this potentially fatal complication before advancing the catheter too far relative to patient size.

7. **Warning:** Practitioners must be aware of the potential for entrapment of guide wire by any implanted device in the circulatory system (i.e., vena cava filters, stents). Review patient's history before catheterization procedure to assess for possible implants. Care should be taken regarding the length of spring-wire guide inserted. It is recommended that if patient has a circulatory system implant, catheter procedure be done under direct visualization to minimize the risk of guidewire entrapment.²
8. **Warning:** Due to the risk of exposure to HIV (Human Immunodeficiency Virus) or other blood borne pathogens, health care workers should routinely use universal blood and body-fluid precautions in the care of all patients.
9. **Precaution:** Indwelling catheters should be routinely inspected for desired flow rate, security of dressing, correct catheter position, and for secure Luer-Lock connection. Use centimeter markings to identify if the catheter position has changed.
10. **Precaution:** Only x-ray examination of the catheter placement can ensure that the catheter tip has not entered the heart or no longer lies parallel to the vessel wall. If catheter position has changed, immediately perform chest x-ray examination to confirm catheter tip position.
11. **Precaution:** Alcohol and acetone can weaken the structure of polyurethane materials. Check ingredients of prep sprays and swabs for acetone and alcohol content. **Acetone:** Do not use acetone on catheter surface. Acetone may be applied to skin but must be allowed to dry completely prior to applying dressing. **Alcohol:** Do not use alcohol to soak catheter surface or to restore catheter patency. Care should be taken when instilling drugs containing high concentration of alcohol. Always allow alcohol to dry completely prior to applying dressing.
12. **Precaution:** Some disinfectants used at the catheter insertion site contain solvents, which can attack the catheter material. Assure insertion site is dry before dressing.
13. **Precaution:** Use of a syringe smaller than 10 mL to irrigate or declot an occluded catheter may cause intraluminal leakage or catheter rupture.⁶

A Suggested Procedure: Use sterile technique.

1. **Precaution:** Place patient in slight Trendelenburg position, as tolerated, to reduce risk of air embolism. If femoral approach is used, place patient in supine position.
2. Prep and drape puncture site as required.
3. Perform skin wheal with desired needle. In kits where provided, a SharpsAway® disposal cup is used for disposal of needles. Push needles into foam after use. Discard entire cup at completion of procedure. **Precaution:** Do not reuse needles after they have been placed into disposal cup. Particulate matter may adhere to needle tip.
4. Locate central vein with introducer needle and syringe. **Precaution:** The color of blood aspirated is not always a reliable indicator of venous access.¹⁰ Because of potential for inadvertent arterial placement, verify venous access via a waveform obtained by a calibrated pressure transducer (refer to Fig. 2).


Fig. 2

If hemodynamic monitoring equipment is not available to permit transducing a central venous waveform, disconnect syringe and check for pulsatile flow. Pulsatile flow is usually an indicator of inadvertent arterial puncture.

5. Insert desired tip of spring-wire guide through introducer needle into vein. **Precaution:** Due to small size of the 0.018" diameter spring-wire guide contained in this product, inspect "J" tip for damage prior to insertion. If "J" tip of spring-wire guide is used, prepare for insertion by sliding plastic tube over "J" to straighten. The spring-wire guide should be advanced a desired depth. Advancement of "J" tip may require a gentle rotating motion. **Warning:** Do not cut spring-wire guide to alter length. Do not withdraw spring-wire guide against needle bevel to minimize the risk of severing or damaging guide wire.
6. Hold spring-wire guide in place and remove introducer needle. **Precaution:** Maintain firm grip on spring-wire guide at all times.
7. Enlarge cutaneous puncture site with cutting edge of scalpel positioned away from spring-wire guide. **Precaution:** Do not cut guide wire.
8. Remove insert tube from catheter. Thread tip of catheter over spring-wire guide. Sufficient guide wire length must remain exposed at hub end of catheter to maintain firm grip on guide wire. Grasping near skin, advance catheter into vein with slight twisting motion. **Precaution:** If provided, catheter clamp must not be attached to catheter until spring-wire guide is removed.
9. Advance catheter to final indwelling position.
10. Hold catheter at desired depth and remove spring-wire guide. The catheter in this product has been designed to freely pass over the spring-wire guide. If resistance is encountered when attempting to remove spring-wire guide after catheter placement, it may be kinked about tip of catheter within vessel (refer to Fig. 3).


Fig. 3

In this circumstance, pulling back on spring-wire guide may result in undue force being applied resulting in spring-wire guide breakage. If resistance is encountered, withdraw catheter relative to spring-wire guide about 2-3 cm and attempt to remove spring-wire guide. If resistance is again encountered, remove spring-wire guide and catheter simultaneously. **Warning:** Although the incidence of spring-wire guide failure is extremely low, practitioner should be aware of potential for breakage if undue force is applied to wire.

11. Verify that entire spring-wire guide is intact on removal.
12. Check lumen placement by attaching a syringe to catheter and aspirate until free flow of venous blood is observed. Connect appropriate Luer-Lock line, as required, or "lock" through injection cap using standard hospital protocol.
13. Secure and dress catheter temporarily.

- Verify catheter tip position by chest x-ray immediately after placement. **Precaution: X-ray exam must show catheter located in right side of mediastinum in SVC with distal end of catheter parallel to vena cava wall and its distal tip positioned at a level above either the azygos vein or carina of the trachea, whichever is better visualized. If catheter tip is malpositioned, reposition and reverify.**
- Secure catheter to patient. Use integral suture ring as primary suture site. In kits where provided, the catheter clamp should be utilized as a secondary suture site as necessary. **Precaution: Do not suture directly to outside diameter of catheter to minimize the risk of cutting or damaging the catheter or impeding flow.**

Catheter Clamp Instructions:

- After spring-wire guide has been removed and line has been connected or locked, spread wings of rubber clamp and position on catheter, as required, to ensure proper tip location (refer to Fig. 4).


Fig. 4

- Secure catheter to patient by suturing catheter clamp to skin using side wings to minimize the risk of catheter migration.
- Dress puncture site per hospital protocol. **Precaution: Maintain insertion site with regular, meticulous redressing using aseptic technique.**
 - Record on patient's chart the indwelling catheter length per centimeter markings on catheter where it enters the skin. Frequent visual reassessment should be made to ensure catheter has not moved.

Catheter Exchange Procedure:

- Use sterile technique.
- Precaution: Prior to attempting a catheter exchange procedure, remove catheter clamp, where provided.**
- Proceed per hospital protocol. Cutting catheter is not recommended due to potential for catheter embolism.

Catheter Removal Procedure:

- Precaution: Place patient in supine position.**
- Remove dressing. **Precaution: To minimize the risk of cutting catheter, do not use scissors to remove dressing.**
- Warning: Exposure of central vein to atmospheric pressure may result in entry of air into central venous system.** Remove suture(s) from catheter clamp, where provided, and primary suture site. Be careful not to cut catheter. Remove catheter slowly, pulling parallel to skin. As catheter exits site, apply pressure with a dressing impermeable to air, e. g. Vaseline® gauze. Because the residual catheter track remains an air entry point until completely sealed, the occlusive dressing should remain in place for at least 24-72

hours dependent upon amount of time catheter was indwelling^{11,14,16,19}

- Upon removal of catheter, inspect it to make sure entire length has been withdrawn.
- Document removal procedure.

References:

- Bar-Joseph G, Galvis AG. Perforation of the heart by central venous catheters in infants: guidelines to diagnosis and management. *J Pediatr Surg.* 1983;18:284-287.
- Blitt CD, ed. *Monitoring in Anesthesia and Critical Care Medicine.* Central venous pressure monitoring. New York, NY: Churchill Livingstone;1985:121-165.
- Brandt RL, Foley WJ, Fink GH, Regan WJ. Mechanism of perforation of the heart with production of hydropericardium by a venous catheter and its prevention. *Am J Surg.* 1970;119:311-316.
- Carbone K, Gimenez LF, Rogers WH, Watson AJ. Hemothorax due to vena caval erosion by a subclavian dual-lumen dialysis catheter. *South Med J.* 1987;80:795-796.
- Collier PE, Ryan JJ, Diamond DL. Cardiac tamponade from central venous catheters –report of a case and review of the English literature. *Angiology.* September 1984;35:595-600.
- Conn C. The importance of syringe size when using implanted vascular access devices. *J Vasc Access Nurs.* Winter 1993;3:11-18.
- Curelaru J, Linder LE, Gustavsson B. Displacement of catheters inserted through internal jugular veins with neck flexion and extension. *Intens Care Med.* 1980;6:179-183.
- Eissa NT, Kvetan V. Guide wire as a cause of complete heart block in patients with preexisting left bundle branch block. *Anesthesiology.* 1990;73:772-774.
- Iberty TJ, Katz LB, Reiner MA, Brownie T, Kwon KB. Hydrothorax as a late complication of central venous indwelling catheters. *Surgery.* November 1983;842-846.
- Jobs DR, Schwartz AJ, Greenhow DE, Stephenson LW, Ellison N. Safer jugular vein cannulation: recognition of arterial punctures and preferential use of the external jugular route. *Anesthesiology.* 1983;59:353-355.
- Kashuk JL, Penn I. Air embolism after central venous catheterization. *Surg Gynecol Obstet.* September 1984;159:249-252.
- Kozeny GA, Bansal VK, Vertuno LL, Hano JE. Contralateral hemothorax secondary to chronic subclavian dialysis catheter. *Am J Nephrol.* 1984;4:312-314.
- Maschke SP, Rogove HJ. Cardiac tamponade associated with a multilumen central venous catheter. *Crit Care Med.* 1984;12:611-612.
- Paskin DL, Hoffman WS, Tudenham WJ. A new complication of subclavian vein catheterization. *Ann Surg.* March 1974;179:266-268.
- Peters JL, ed. *A Manual of Central Venous Catheterization and Parenteral Nutrition.* Boston, MA: John Wright PSG;1983:58-61, 155-157.
- Phifer TJ, Bridges M, Conrad SA. The residual central venous catheter track – an occult source of lethal air embolism: case report. *J Trauma.* 1991;31:1558-1560.
- Sheep RE, Guiney WB Jr. Fatal cardiac tamponade. *JAMA.* 1982;246:1632-1635.
- Sigurdsson J, Riba P, Sigurdsson S. The wandering central venous catheter. *Intens Care Med.* 1985;1:263-264.
- Thielen JB, Nyquist J. Subclavian catheter removal. *J Intravenous Nurs.* March/April 1991;14:114-118.

20. Tocino IM, Watanabe A. Impending catheter perforation of superior vena cava: radiographic recognition. *Am J Roentgenology*. March 1986;146:487-490. vena cava catheter and hemodialysis. *Crit Care Med*. 1988;16:624-628.
21. Vaziri ND, Maksy M, Lewis M, Martin D, Edwards K. Massive mediastinal hematoma caused by a double-lumen subclavian catheter. *Artif Organs*. 1984;8:223-226. Arrow International, Inc. recommends that the user be acquainted with the reference literature.
22. Wanscher M, Frifelt JJ, Smith-Sivertsen C, et al. Thrombosis caused by polyurethane double-lumen subclavian superior * If you have any questions or would like additional reference information, please contact Arrow International, Inc.

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