Inspected Dimensions:  
Folded width: 5-1/2" (14 cm)  
Folded length: 7-1/2" (19 cm)

# ARROW®

## Single Lumen Infusion Catheter (SLIC®)

### Safety and Efficacy Considerations:

Do not use if package has been previously opened or damaged. **Warning: Prior to use read all package insert warnings, precautions, and instructions. Failure to do so may result in severe patient injury or death.**

Do not alter the SLIC® during insertion, use, or removal.

Procedure must be performed by trained personnel well versed in anatomical landmarks, safe technique, and potential complications.

**Warning: Do not place the catheter into or allow it to remain in the right atrium or right ventricle (refer to Fig. 1).**

Central venous catheters should not be placed in the right atrium unless specifically required for special relatively short term procedures, such as aspiration of air emboli during neurosurgery. Such procedures are, nevertheless, risk prone and should be closely monitored and controlled.

### Indications:

The SLIC® is a two-piece assembly consisting of an infusion catheter and an obturator. With SLIC® obturator removed, the Arrow SLIC® permits access to the central venous circulation through an indwelling sheath/hemostasis valve. With the SLIC® obturator in place, the SLIC® occludes the hemostasis valve preventing air entry and blood loss through the valve.

### Contraindications:

None known.

### Warnings and Precautions:\*

- Warning: Sterile, Single use: Do not reuse, reprocess or resterilize.** Reuse of device creates a potential risk of serious injury and/or infection which may lead to death.
- Warning: The practitioner must be aware of potential air embolism problems associated with leaving open needles, sheaths, or catheters in venous puncture sites or as a consequence of inadvertent disconnects.** To lessen the risk of disconnects, only securely tightened Luer-Lock connections should be used with this device. Follow hospital protocol for all sheath and side port maintenance to guard against air embolism.
- Warning: Practitioners must be aware of complications associated with central vein catheters including cardiac tamponade secondary to vessel wall, atrial or ventricular perforation, pleural and mediastinal injuries, air embolism, catheter embolism, catheter occlusion, thoracic duct laceration, bacteremia, septicemia, thrombosis, inadvertent arterial puncture, nerve damage, hematoma, hemorrhage, and dysrhythmias.**
- Warning: Hemostasis valve/side port assembly to SLIC® connection and SLIC® to obturator connection must be secured and routinely examined to minimize the risk of disconnection and possible air embolism, hemorrhage, or exsanguination.**
- Warning: Due to the risk of exposure to HIV (Human Immunodeficiency Virus) or other bloodborne pathogens, health care workers should routinely use universal blood and body fluid precautions in the care of all patients.**
- Precaution: Alcohol and acetone can weaken the structure of polyurethane materials. Check ingredients of prep sprays and swabs for acetone**


Fig. 1

**Cardiac Tamponade:** It has been documented by many authors that placement of indwelling catheters in the right atrium is a dangerous practice<sup>1,2,3,4,5,7,8</sup> that may lead to cardiac perforation and tamponade.<sup>1,2,3,4,7,8</sup> Although cardiac tamponade secondary to pericardial effusion is uncommon, there is a high mortality rate associated with it.<sup>9</sup> Practitioners placing central venous catheters must be aware of this potentially fatal complication before advancing the catheter too far relative to patient size.

No particular route or catheter type is exempt from this potentially fatal complication.<sup>8</sup> The actual position of the tip of the indwelling catheter should be confirmed by x-ray after insertion.<sup>1,2,4,7,8,10</sup> Central venous catheters should be placed in the superior vena cava<sup>1,2,3,4,5,7,11</sup> above its junction with the right atrium and parallel to the vessel wall<sup>6,11</sup> and its distal tip positioned at a level above either the azygos vein or the carina of the trachea, whichever is better visualized.

and alcohol content.

**Acetone:** Do not use acetone on catheter surface. Acetone may be applied to skin but must be allowed to dry completely prior to applying dressing.

**Alcohol:** Do not use alcohol to soak catheter surface or to restore catheter patency. Care should be taken when instilling drugs containing high concentration of alcohol. Always allow alcohol to dry completely prior to applying dressing.

## A Suggested Procedure: Use sterile technique.

1. **Precaution:** Place patient in slight Trendelenburg position, as tolerated, to reduce the risk of air embolism. If femoral approach is used, place patient in supine position.
2. If catheter is removed from sheath or catheter insertion is delayed, introduce the entire length of the SLIC<sup>®</sup> assembly through the hemostasis valve/sheath assembly. Twist to lock (refer to Fig. 2).


Fig. 2

- Orient slot in hub with locking pin on assembly cap.
- Slide hub forward over cap and twist.

**Precaution:** When SLIC<sup>®</sup> assembly is inserted after catheter removal, use an appropriate antiseptic to prep the hemostasis valve housing prior to inserting the SLIC<sup>®</sup>. Include the exposed portion of the valve on the top of the housing. The SLIC<sup>®</sup>, with the obturator in place, occludes the hemostasis valve preventing air entry or blood loss through the valve. **Warning:** Connection between SLIC<sup>®</sup> obturator and SLIC<sup>®</sup> must be tightened securely and routinely examined to minimize the risk of disconnection and possible air embolism, hemorrhage or exsanguination.

3. To use the SLIC<sup>®</sup> for intravenous infusion, remove the blue-capped SLIC<sup>®</sup> obturator by twisting counter clockwise. Hold the infusion port to maintain positive lock to hemostasis valve housing. Pull the SLIC<sup>®</sup> obturator from the infusion catheter (refer to Fig. 3). Immediately attach desired line to Luer-Lock hub. **Warning:** Exposure of the central vein to atmospheric pressure may result in entry of air into the central venous system. Inspect the SLIC<sup>®</sup> obturator to ensure the entire length has been withdrawn. Document the SLIC<sup>®</sup> obturator withdrawal and start of infusion.


Fig. 3

4. If the infusion through the SLIC<sup>®</sup> is discontinued, the hub should be capped with a Luer-Lock injection cap and handled per hospital flushing protocol, or the SLIC<sup>®</sup> should be withdrawn and replaced with a sterile Arrow obturator, sold separately, to ensure that leakage does not occur and inner seal is protected from contamination. **Warning:** Cover the lumen during any manipulation to minimize the risk of blood loss or the introduction of air into the sheath.

\* If you have any questions or would like reference information, please contact Arrow International, Inc.

## Reference:

1. Bar-Joseph G, Galvis AG. Perforation of the heart by central venous catheters in infants: guidelines to diagnosis and management. *J Pediatr Surg.* 1983;18:284-287.
2. Bliitt CD, ed. *Monitoring in Anesthesia and Critical Care Medicine.* Central venous pressure monitoring. New York, NY: Churchill Livingstone; 1985:121-165.
3. Brandt RL, Foley WJ, Fink GH, Regan WJ. Mechanism of perforation of the heart with production of hydropericardium by a venous catheter and its prevention. *Am J Surg.* 1970;119:311-316.
4. Collier PE, Ryan JJ, Diamond DL. Cardiac tamponade from central venous catheters – report of a case and review of the English literature. *Angiology.* September 1984;35:595-600.
5. Curelaru J, Linder LE, Gustavsson B. Displacement of catheters inserted through internal jugular veins with neck flexion and extension. *Intens Care Med.* 1980;6:179-183.
6. Iberti TJ, Katz LB, Reiner MA, Brownie T, Kwun KB. Hydrothorax as a late complication of central venous indwelling catheters. *Surgery.* November 1983:842-846.
7. Maschke SP, Rogove HJ. Cardiac tamponade associated with a multilumen central venous catheter. *Crit Care Med.* 1984;12:611-612.
8. Peters JL, ed. *A Manual of Central Venous Catheterization and Parenteral Nutrition.* Boston, MA: John Wright PSG; 1983:58-61,155-157.
9. Sheep RE, Guiney WB Jr. Fatal cardiac tamponade. *JAMA.* 1982;248:1632-1635.
10. Sigurdsson J, Riba P, Sigurdsson S. The wandering central venous catheter. *Intensive Care Med.* 1985;11:263-264.
11. Tocino IM, Watanabe A. Impending catheter perforation of superior vena cava: radiographic recognition. *Am J Roentgenology.* March 1986;146:487-490.

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\* Om du har frågor eller önskar referensinformation, kontakta Arrow International, Inc.