EN

ARROW®

HANDS-OFF® Central Venous Catheterization Product with Cath-Gard® Catheter Contamination Shield

Recommended for use with Arrow® One-Piece Percutaneous Sheath Introducer

Safety and Efficacy Considerations:

Do not use if package has been previously opened or damaged.

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Do not alter the catheter, catheter contamination shield, or any other kit/set component during insertion, use, or removal.

Central venous catheterization must be performed by trained personnel well versed in anatomical landmarks, safe technique, and potential complications.

Warning: Do not place the catheter into or allow it to remain in the right atrium or right ventricle (refer to Fig. 1).

Fig. 1

Cardiac Tamponade:

It has been documented by many authors that placement of indwelling catheters in the right atrium is a dangerous practice^{1,2,3,4,6,8,11,12} that may lead to cardiac perforation and tamponade.^{1,2,3,5,6,11,12} Although cardiac tamponade secondary to pericardial effusion is uncommon, there is a high mortality rate associated with it.^{1,3} Practitioners placing central venous catheters must be aware of this potentially fatal complication before advancing the catheter too far relative to patient size.

No particular route or catheter type is exempt from this potentially fatal complication.^{1,2} The actual position of the tip of the indwelling catheter should be confirmed by x-ray after insertion.^{1,2,5,6,11,12,14} Central venous catheters should be placed in the superior vena cava^{1,2,3,5,6,8,11,15} above its junction with the right atrium and parallel to the vessel wall^{6,9,15} and its distal tip positioned at a level above either the azygos vein or the carina of the trachea, whichever is better visualized.

Central venous catheters should not be placed in the right atrium unless specifically required for special relatively short term procedures, such as aspiration of air emboli during neurosurgery. Such procedures are nevertheless risk prone and should be closely monitored and controlled.

Indications for Use:

The central venous catheter permits venous access to the central circulation through an indwelling Arrow® percutaneous sheath introducer.

Contraindications:

None known.

Central Venous Catheterization

Warnings and Precautions:*

1. **Warning: Sterile, Single use: Do not reuse, reprocess or resterilize. Reuse of device creates a potential risk of serious injury and/or infection which may lead to death.**
2. **Warning: Do not place the catheter into or allow it to remain in the right atrium or right ventricle. Central vein catheters should be positioned so that the distal tip of the catheter is in the superior vena cava (SVC) above the junction of the SVC and the right atrium and lies parallel to the vessel wall. For femoral vein approach the catheter should be advanced into the vessel so that the catheter tip lies parallel to the vessel wall and does not enter the right atrium.**
3. **Warning: Practitioners must be aware of complications associated with central vein catheters including cardiac tamponade secondary to vessel wall, atrial or ventricular perforation, pleural and mediastinal injuries, air embolism, catheter embolism, bacteremia, septicemia, thrombosis, hemorrhage, and dysrhythmias.**
4. **Warning: The practitioner must be aware of potential air embolism associated with leaving open needles or catheters in central venous puncture sites or as a consequence of inadvertent disconnects. To lessen the risk of disconnects, only securely tightened Luer-Lock connections should be used with this device. Follow hospital protocol to guard against air embolism for all catheter maintenance.**
5. **Warning: Due to the risk of exposure to HIV (Human Immunodeficiency Virus) or other blood borne pathogens, health care workers should routinely use universal blood and body-fluid precautions in the care of all patients.**
6. **Precaution: Indwelling catheters should be routinely inspected for desired flow rate, security of dressing, correct catheter position, and for secure Luer-Lock connection. Use centimeter markings to identify if the catheter position has changed.**
7. **Precaution: Only x-ray examination of the catheter placement can ensure that the catheter tip has not entered the heart or no longer lies parallel to the vessel wall. If catheter position has changed, immediately perform chest x-ray examination to confirm catheter tip position.**
8. **Precaution: For blood sampling, temporarily shut off remaining port(s), if applicable, through which solutions are being infused.**
9. **Precaution: Alcohol and acetone can weaken the structure of polyurethane material. Therefore, care should be taken**

when instilling drugs containing alcohol or when using high concentration of alcohol or acetone when performing routine catheter care and maintenance. Alcohol should not be utilized to decontaminate polyurethane catheters.

10. **Precaution:** Use of a syringe smaller than 10 mL to irrigate or decontaminate an occluded catheter may cause intraluminal leakage or catheter rupture.⁷

Catheter Insertion Procedure:

Use sterile technique.

1. Prepare the catheter for insertion by flushing each lumen and clamping or attaching injection caps to the appropriate extension lines, if applicable. **Warning: Do not cut catheter to alter length.**
2. Prep hemostasis valve/side port assembly cap of the percutaneous sheath introducer with appropriate antiseptic per hospital protocol. Include the exposed portion of the valve on the top of the cap. Occlude hemostasis valve with sterile-gloved finger.
3. Press distal hub of catheter contamination shield over hemostasis valve/side port assembly cap. Twist to lock (refer to Fig. 2).


Fig. 2

- Orient slot in hub with locking pin on assembly cap.
 - Slide hub forward over cap and twist.
4. Grasp catheter through catheter contamination shield and advance into sheath introducer. Continue advancing catheter to final indwelling position. **Warning: The catheter insertion depth must be determined relative to patient size. Centimeter marks on catheter may be used as positioning reference points taking into consideration the amount of catheter external to the body passing through the sheath hemostasis valve/side port assembly (approximately 6 cm).**
 5. While maintaining catheter position, twist the upper half of the distal hub in clockwise direction to lock catheter in place. Test the adapter by gently tugging on the catheter to ensure a secure grip on the catheter (refer to Fig. 3).


Fig. 3

6. Check catheter lumen placement by attaching a syringe to each extension line and aspirate until free flow of venous blood is observed. **Precaution: The proximal lumen on a multiple**

lumen catheter is located approximately 5 cm above the distal tip of the catheter. Depending on the patient's anatomy and site of sheath placement, be aware that the proximal lumen of a multiple lumen catheter may be inside the sheath when the catheter tip is properly placed within the superior vena cava. This condition may affect flow rate through proximal lumen, and incompatible medications should not be administered simultaneously through the sheath side port and the proximal lumen within the sheath.

7. Connect all extension lines to appropriate Luer-Lock line(s), if applicable. Unused ports may be "locked" through injection cap(s) using standard hospital protocol. **Warning: This product contains slide clamps which may be inadvertently removed, and potentially aspirated, by children or confused adults. In such situations, it is recommended that the practitioner remove the clamps when not in use.** Slide clamps are provided on extension line(s) to occlude flow through each lumen during line and injection cap changes. **Precaution: To avoid damage to extension lines from excessive pressure, each clamp must be opened prior to infusing through that lumen.**
8. Verify catheter tip position by chest x-ray immediately after placement. **Precaution: X-ray exam must show the catheter located in the right side of the mediastinum in the SVC with the distal end of the catheter parallel to the vena cava wall and its distal tip positioned at a level above either the azygos vein or the carina of the trachea, whichever is better visualized. If catheter tip is malpositioned, reposition and re-verify.**
9. Record on the patient's chart the indwelling catheter length using centimeter markings on the catheter. Frequent visual reassessment should be made to ensure that the catheter has not moved.

Catheter Removal Procedure:

1. **Precaution: Place the patient in a supine position.**
2. Twist upper half of distal hub of catheter contamination shield in counterclockwise direction to unlock the catheter from position.
3. Withdraw catheter from sheath into catheter contamination shield.
4. Twist distal hub of catheter contamination shield to allow removal from hemostasis valve/side port assembly. Temporarily occlude hemostasis valve with sterile-gloved finger.
5. **Precaution: Insert Arrow® obturator, sold separately, as a dummy catheter into hemostasis valve/side port. This will ensure that leakage does not occur and inner seal is protected from contamination.**

References:

1. Bar-Joseph G, Galvis AG. Perforation of the heart by central venous catheters in infants: guidelines to diagnosis and management. *J Pediatr Surg.* 1983;18:284-287.
2. Blitt CD, ed. Monitoring in anesthesia and critical care medicine. *Central venous pressure monitoring.* New York, NY: Churchill Livingstone, 1985. pp.121-165.
3. Brandt RL, Foley WJ, Fink GH, Regan WJ. Mechanism of perforation of the heart with production of hydropericardium by a venous catheter and its prevention. *Am J Surg.* 1970;119:311-316.
4. Carbone K, Gimenez LF, Rogers WH, Watson AJ. Hemothorax due to vena caval erosion by a subclavian dual-lumen dialysis catheter. *South Med J.* 1987;80:795-796.

5. Collier PE, Goodman GB. Cardiac tamponade caused by central venous catheter perforation of the heart: a preventable complication. *J Am Coll Surg*. 1995;181:459-463.
6. Collier PE, Ryan JJ, Diamond DL. Cardiac tamponade from central venous catheters – report of a case and review of the English literature. *Angiology*. September 1984;35: 595-600.
7. Conn C. The importance of syringe size when using implanted vascular access devices. *J Vasc Access Nurs*. Winter 1993;3:11-18.
8. Curelaru J, Linder LE, Gustavsson B. Displacement of catheters inserted through internal jugular veins with neck flexion and extension. *Intens Care Med*. 1980;6:179-183.
9. Iberti TJ, Katz LB, Reiner MA, Brownie T, Kwun KB. Hydrothorax as a late complication of central venous indwelling catheters. *Surgery*. November 1983:842-846.
10. Kozeny GA, Bansal VK, Vertuno LL, Hano JE. Contralateral hemothorax secondary to chronic subclavian dialysis catheter. *Am J Nephrol*. 1984;4:312-314.
11. Maschke SP, Rogove HJ. Cardiac tamponade associated with a multilumen central venous catheter. *Crit Care Med*. 1984;12:611-612.
12. Peters JL, ed. *A Manual of Central Venous Catheterization and Parenteral Nutrition*. Boston Mass:John Wright PSG; 1983:58-61, 155-157.
13. Sheep RE, Guiney WB Jr. Fatal cardiac tamponade. *JAMA*. 1982;248:1632-1635.
14. Sigurdsson J, Riba P, Sigurdsson S. The wandering central venous catheter. *Intens Care Med*. 1985;11:263-264.
15. Tocino IM, Watanabe A. Impending catheter perforation of superior vena cava: radiographic recognition. *Am J Roentgenology*. March 1986;146:487-490.
16. Vaziri ND, Maksy M, Lewis M, Martin D, Edwards K. Massive mediastinal hematoma caused by a double-lumen subclavian catheter. *Artif Organs*. 1984;8:223-226.
17. Wanscher M, Frifelt JJ, Smith-Sivertsen C, et al. Thrombosis caused by polyurethane double-lumen subclavian superior vena cava catheter and hemodialysis. (*Crit Care Med*) 1988;16:624-628.

Arrow International, Inc. recommends that the user be acquainted with the reference literature.

*If you have any questions or would like additional reference information, please contact Arrow International, Inc.

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Arrow International, Inc. kullanıcının referans literatüre aşına olmasını önerir.
*Sorumuz varsa veya ekle referans bilgisi isterseniz lütfen Arrow International, Inc. ile irtibat kurun.