Inspected Dimensions:
Folded Length: 7-1/2" (19 cm)
Folded Width: 5-1/2" (14 cm)

ARROW®

INTERNATIONAL

EN	Arterial Catheterization Product	1
CS	Výrobek určený k arteriální katetrizaci	3
DA	Arterielt kateteriseringsprodukt	5
NL	Product voor arteriële katheterisatie	7
ET	Arterikanüül	9
FI	Valtimokatetrituote	11
FR	Produit de cathétérisme artériel	13
DE	Besteck zur Katheterisierung von Arterien	15
EL	Προϊόν αρτηριακού καθετηριασμού	17
HU	Artériás katéterező eszköz	19
IT	Prodotto per cateterizzazione arteriosa	21
LT	Arterijų kateterizavimo įtaisai	23
NO	Arteriell kateteriseringsprodukt	25
PL	Przyrząd do cewnikowania tętnic	27
PT	Produto de Cateterização Arterial	29
RO	Dispozitiv de cateterizare arterială	31
RU	Устройство для артериальной катетеризации	33
SK	Arteriálny katetrizačný produkt	35
SL	Izdelek za arterijsko kateterizacijo	37
ES	Producto para cateterización arterial	39
SV	Anordning för artärkatetrisering	41
TR	Arteriyel Kateterizasyon Ürünü	43

Safety and Efficacy Considerations:

Do not use if package has been previously opened or damaged.

Warning: Prior to use read all package insert warnings, precautions, and instructions. Failure to do so may result in severe patient injury or death.

Do not alter catheter or any kit/set components during insertion, use or removal.

Procedure must be performed by trained personnel well versed in anatomical landmarks, safe technique and potential complications.

Indications for Use:

The Arrow® Arterial Catheterization Device permits access to the peripheral arterial circulation.

Contraindications:

Inadequate collateral arterial circulation.

Warnings and Precautions:*

1. **Warning:** Sterile, Single use: Do not reuse, reprocess or resterilize. Reuse of device creates a potential risk of serious injury and/or infection which may lead to death.
2. **Warning:** Practitioners must be aware of complications associated with arterial procedures^{1,7,8,11} including bacteremia,² septicemia, vessel wall perforation, intravascular clotting and embolization, hematoma, arterial spasm, tissue necrosis, hemorrhage, thrombosis, peripheral ischemia and infarction, peripheral nerve damage,⁹ air embolism,³ and occlusion.
3. **Warning:** In brachial procedures, collateral flow cannot be guaranteed, and therefore intravascular clotting can result in tissue necrosis.¹²
4. **Warning:** In radial artery procedures, practitioners must ascertain that definite evidence of collateral ulnar flow exists.
5. **Warning:** Care should be exercised that indwelling catheter is not inadvertently kinked at hub area when securing catheter to patient. Kinking may weaken wall of catheter and cause a fraying or fatigue of material, leading to possible separation of the catheter.
6. **Warning:** Due to the risk of exposure to HIV (Human Immunodeficiency Virus) or other blood borne pathogens, health care workers should routinely use universal blood and body-fluid precautions in the care of all patients.
7. **Precaution:** Do not suture directly to outside diameter of catheter body to minimize the risk of damaging the catheter or adversely affecting monitoring capabilities.
8. **Precaution:** In order to minimize the risk of problems associated with disconnects, it is recommended that only Luer-Lock connecting tubing be used with this device.⁵ It is recommended to regularly check all connections of the system.
9. **Precaution:** Alcohol and acetone can weaken the structure of polyurethane materials. Check ingredients of prep sprays and swabs for acetone and alcohol content.

ARROW®

Arterial Catheterization Product

Acetone: Do not use acetone on catheter surface. Acetone may be applied to skin but must be allowed to dry completely prior to applying dressing.

Alcohol: Do not use alcohol to soak catheter surface or to restore catheter patency. Care should be taken when instilling drugs containing high concentration of alcohol. Always allow alcohol to dry completely prior to applying dressing.

10. **Precaution:** Some disinfectants used at catheter insertion site contain solvents, which can attack catheter material. Assure insertion site is dry before dressing.
11. **Precaution:** Indwelling catheter should be routinely inspected for desired flow rate, security of dressing, and possible migration. Do not use scissors to remove dressing to minimize the risk of cutting catheter.
12. **Precaution:** The entire system (dome of pressure transducer, extension tube, etc.) must be permanently filled with isotonic saline or heparinised isotonic saline solution and must not contain any air bubbles or blood.
13. **Precaution:** It is recommended to check patency of catheter at least every three hours by delivering a bolus of isotonic saline. It is recommended to use heparinised isotonic saline.
14. **Precaution:** It is recommended to check correct blood reflux at least twice every shift.
15. **Precaution:** Never inject drugs via an arterial catheter.
16. **Precaution:** Remove insert tube prior to use, where applicable.

A Suggested Procedure:

Use sterile technique.

1. Assess for adequate collateral arterial circulation.^{5,12}
2. Prep and drape anticipated puncture site per hospital protocol.
3. Utilize local anesthetic as required.
4. **Precaution:** Remove insert tube from catheter, where applicable. Place on sterile field.
5. Locate and puncture artery with desired introducer needle or introducer catheter/needle assembly, where provided. Avoid transfixing both vessel walls. **Precaution:** If both vessel walls are punctured, subsequent advancement of spring-wire guide could result in inadvertent sub-arterial placement. Pulsatile blood flow from needle hub indicates successful entry into artery.
6. Insert desired tip of spring-wire guide through introducer needle into artery (until depth-marking [if provided] on wire enters hub of needle). At this point, the depth-marking entering hub shows that tip of wire leaves tip of introducer needle and enters vessel. If catheter/needle assembly is used, remove needle and insert spring-wire guide through catheter into artery as described above. If "J" tip spring-wire guide is used, prepare for insertion by sliding plastic tube over "J" to straighten. Advance spring-wire guide to required depth. Advancement of "J" tip may require a gentle rotating motion.

Precaution: Maintain firm grip on spring-wire guide at all times.

7. Hold spring-wire guide in place and remove introducer needle. **Warning: Do not withdraw spring-wire guide against needle bevel to minimize the risk of possible severing of spring-wire guide.** Hold spring-wire guide in place and remove catheter if catheter/needle assembly is used.
8. Enlarge puncture site with scalpel if desired. **Warning: Do not cut spring-wire guide.**
9. Thread tip of indwelling catheter over spring-wire guide. **Precaution: Allow sufficient guide wire length to remain exposed at hub end of catheter to maintain firm grip on guide wire.** Grasping near skin, advance catheter into vessel with slight twisting motion.
10. Advance catheter to final indwelling position using centimeter marks (where provided) as positioning reference points. All centimeter marks are referenced from the catheter tip. Marking symbology is as follows: (1) numerical: 5, 15, 25, etc.; (2) bands: each band denotes 10 cm intervals, with one band indicating 10 cm, two bands indicating 20 cm, etc.; (3) each dot denotes a 1 cm interval.
11. Hold catheter at desired depth and remove spring-wire guide. Pulsatile blood flow indicates positive arterial placement. **Precaution: Use care when removing spring-wire guide. If resistance is encountered, remove spring-wire guide and catheter simultaneously. Use of excessive force may damage catheter or spring-wire guide.**
12. For blood containment purposes, in accordance with standard hospital protocol, use slide clamp provided on extension line to occlude flow through the lumen during line and injection cap changes. **Precaution: To minimize the risk of damage to extension line from excessive pressure, each clamp must be opened prior to infusion through that lumen to prevent any possible adverse affects on monitoring capabilities.**
13. Attach desired Luer-Lock stopcock, injection cap or connecting tubing to catheter hub. Secure catheter to patient in preferred manner using suture wings, suture groove or wing clip, where provided. **Warning: Care should be exercised that the indwelling catheter is not inadvertently kinked at hub area when securing catheter to patient. Kinking may weaken wall of catheter and cause a fraying or fatigue of the material, leading to possible separation of catheter. Precaution: Do not suture directly to outside diameter of catheter body to minimize the risk of damaging the catheter or adversely affecting monitoring capabilities.**
14. Dress puncture site according to hospital protocol.

Catheter Removal Procedure:

1. Remove dressing. **Precaution: To minimize the risk of cutting catheter do not use scissors to remove dressing.**

2. **Warning: Exposure of arterial circulation to atmospheric pressure may result in entry of air into circulation.** Remove suture(s). Be careful not to cut catheter. Remove catheter slowly, pulling it parallel to skin. As catheter exits site, apply pressure with a dressing until hemostasis occurs. Apply light pressure dressing to site.

3. Upon removal of catheter, inspect it to make sure that entire length has been withdrawn.

4. Document removal procedure.

References:

1. Abadir AR, Ung KA. Complications of radial artery cannulation. *Anesthesiology Rev.* 1980;7:11-16.
 2. Band JD, Maki DG. Infections caused by arterial catheters used for hemodynamic monitoring. *Am J Med.* 1979;67:735-741.
 3. Chang C, Dughi J, Shitabata P, Johnson G, Coel M, McNamara JJ. Air embolism and the radial arterial line. *Crit Care Med.* 1988;16:141-143.
 4. Clark CA, Harman EM. Hemodynamic monitoring: arterial catheters. In: Taylor RW, Civetta JM, Dirby RR, eds. *Techniques and Procedures in Critical Care.* Philadelphia, PA: JB Lippincott; 1990:218-230.
 5. Daily EK, Schroeder JS. *Techniques in Bedside Hemodynamic Monitoring.* 5th ed. St. Louis, MO: Mosby; 1994:71-77.
 6. Falk PS, Scuderi PE, Sherertz RJ, Motsinger SM. Infected radial artery pseudoaneurysms occurring after percutaneous cannulation. *Chest.* 1992;101:490-495.
 7. Kaye W. Invasive monitoring techniques: arterial cannulation, bedside pulmonary artery catheterization, and arterial puncture. *Heart Lung.* 1983;12:395-427.
 8. Kaye WE, Dubin HG. Vascular cannulation. In: Taylor RW, Civetta JM, Dirby RR, eds. *Techniques and Procedures in Critical Care.* Philadelphia, PA: JB Lippincott; 1990:204-208.
 9. Marshall G, Edelstein G, Hirshman CA. Median nerve compression following radial arterial puncture. *Anesth Analg.* 1980;59:953-954.
 10. Norwood SH, Cormier B, McMahon NG, Moss A, Moore V. Prospective study of catheter-related infection during prolonged arterial catheterization. *Crit Care Med.* 1988;16:836-839.
 11. Puri VK, Carlson RW, Bander JJ, Weil MH. Complications of vascular catheterization in the critically ill. *Crit Care Med.* 1980;8:495-499.
 12. Venus B, Satish P. Vascular Cannulation. In: Civetta JM, Taylor RW, Kirby RR, eds. *Critical Care.* 3rd ed. Philadelphia, PA: Lippincott-Raven Publishers; 1997:539-542.
- Arrow International, Inc. recommends that the user be acquainted with the reference literature.

*If you have any questions or would like additional information, please contact Arrow International, Inc.

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