Inspected Dimensions:
Folded Length: 7-1/2" (19 cm)
Folded Width: 5-1/2" (14 cm)

ARROW®

Central Venous Catheterization Product with Catheter Contamination Guard

for use only with Arrow MAC™ Two-Lumen Central Venous Access Device

Safety and Efficacy Considerations:

Do not use if package has been previously opened or damaged.

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Do not resterilize or reuse. Do not alter the catheter during insertion, use, or removal.

Procedure must be performed by trained personnel well versed in anatomical landmarks, safe technique, and potential complications.

Warning: Do not place the catheter into or allow it to remain in the right atrium or right ventricle (refer to Fig. 1).

Fig. 1

Cardiac Tamponade: It has been documented by many authors that placement of indwelling catheters in the right atrium is a dangerous practice^{1,2,3,4,6,8,9} that may lead to cardiac perforation and tamponade.^{1,2,3,4,8,9} Although cardiac tamponade secondary to pericardial effusion is uncommon, there is a high mortality rate associated with it¹⁰. Practitioners placing central venous catheters must be aware of this potentially fatal complication before advancing the catheter too far relative to patient size.

No particular route or catheter type is exempt from this potentially fatal complication.⁹ The actual position of the tip of the indwelling catheter should be confirmed by x-ray after insertion.^{1,2,4,8,9,11} Central venous catheters should be placed in the superior vena cava^{1,2,3,4,6,8,12} above its junction with the right atrium and parallel to the vessel wall^{7,12} and its distal tip positioned at a level above either the azygos vein or the carina of the trachea, whichever is better visualized.

Central venous catheters should not be placed in the right atrium unless specifically required for special relatively short term procedures, such as aspiration of air emboli during neurosurgery. Such procedures are nevertheless risk prone and should be closely monitored and controlled.

Indications for Use:

The Central Venous Catheter permits venous access to the central venous circulation through a hemostasis valve on an indwelling two-lumen central venous access device.

Contraindications:

None known.

Warnings and Precautions:*

1. **Warning: Sterile, Single use: Do not reuse, reprocess or resterilize. Reuse of device creates a potential risk of serious injury and/or infection which may lead to death.**
2. **Warning: Do not place the catheter into or allow it to remain in the right atrium or right ventricle. Central vein catheters should be positioned so that the distal tip of the catheter is in the superior vena cava (SVC) above the junction of the SVC and the right atrium and lies parallel to the vessel wall. For femoral vein approach the catheter should be advanced into the vessel so that the catheter tip lies parallel to the vessel wall and does not enter the right atrium.**
3. **Warning: Practitioners must be aware of complications associated with central vein catheters including cardiac tamponade secondary to vessel wall, atrial or ventricular perforation, pleural and mediastinal injuries, air embolism, catheter embolism, catheter occlusion, thoracic duct laceration, bacteremia, septicemia, thrombosis, inadvertent arterial puncture, nerve damage, hematoma, hemorrhage, and dysrhythmias.**
4. **Warning: The practitioner must be aware of potential air embolism associated with leaving open needles, sheaths, or catheters in venous puncture sites or as a consequence of inadvertent disconnects. To lessen the risk of disconnects, only securely tightened Luer-Lock connections should be used with this device. Follow hospital protocol for all sheath and side port maintenance to guard against air embolism.**
5. **Warning: Hemostasis valve to two-lumen central venous access device connection must be secured and routinely examined to minimize the risk of disconnection and possible air embolism, hemorrhage, or exsanguination.**
6. **Warning: Due to the risk of exposure to HIV (Human Immunodeficiency Virus) or other bloodborne pathogens, health care workers should routinely use universal blood and body fluid precautions in the care of all patients.**
7. **Precaution: For blood sampling, temporarily shut off remaining port(s), if applicable, through which solutions are being infused.**
8. **Precaution: Indwelling catheters should be routinely inspected for desired flow rate, security of dressing, correct catheter position, and for secure Luer-Lock connection.**
9. **Precaution: Use of a syringe smaller than 10 mL to irrigate or declot an occluded catheter may cause intraluminal leakage or catheter rupture.⁵**

10. **Precaution:** Alcohol and acetone can weaken the structure of polyurethane materials. Check ingredients of prep sprays and swabs for acetone and alcohol content.

Acetone: Do not use acetone on catheter surface. Acetone may be applied to skin but must be allowed to dry completely prior to applying dressing.

Alcohol: Do not use alcohol to soak catheter surface or to restore catheter patency. Care should be taken when instilling drugs containing high concentration of alcohol. Always allow alcohol to dry completely prior to applying dressing.

11. **Precaution:** Some disinfectants used at the catheter insertion site contain solvents, which can attack the catheter material. Assure insertion site is dry before dressing.

A Suggested Procedure:

Use sterile technique.

1. **Precaution:** Place patient in slight Trendelenburg position, as tolerated, to reduce the risk of air embolism. If femoral approach is used, place patient in supine position.
2. If catheter is removed from hemostasis valve of two-lumen central venous access device or catheter insertion is delayed, insert Arrow® obturator, sold separately, as a dummy catheter through hemostasis valve. This will ensure that leakage does not occur and inner seal is protected from contamination.
3. Prepare the central venous catheter for insertion by flushing each lumen and clamping or attaching injection caps to the appropriate extension lines, as required. **Warning:** Do not cut catheter to alter length.
4. Prep hemostasis valve cap with appropriate antiseptic per hospital protocol. Include the exposed portion of the valve on the top of the cap. Occlude hemostasis valve with sterile-gloved finger.
 - With sterile-gloved fingers, peel back tail of peel-away catheter contamination guard to expose tip of catheter (Fig. 2).


Fig. 2

- Grasping catheter body, still covered by the peel-away catheter contamination guard, advance catheter into hemostasis valve of two-lumen central venous access device. Continue advancing catheter and peeling away the contamination guard until catheter hub locks on hemostasis valve (Fig. 3).


Fig. 3

5. Verify catheter tip position by chest x-ray immediately after placement. **Precaution:** X-ray exam must show the catheter located in the right side of the mediastinum in the SVC with the distal end of the catheter parallel to the vena cava wall and its distal tip positioned at a level above either the azygos vein or the carina of the trachea, whichever is better visualized. If catheter tip is malpositioned, reposition and re-verify.
6. Connect all extension lines to appropriate Luer-Lock line(s), as required. Unused ports may be “locked” through injection cap(s) using standard hospital protocol. **Warning:** This product contains slide clamps which may be inadvertently removed, and potentially aspirated, by children or confused adults. In such situations, it is recommended that the practitioner remove the clamps when not in use. Slide clamps are provided on extension lines to occlude flow through each lumen during line and injection cap changes. **Precaution:** To minimize the risk of damage to extension lines from excessive pressure, each clamp must be opened prior to infusing through that lumen.

Catheter Removal Procedure:

1. **Precaution:** Place the patient in a supine position.
2. Withdraw catheter from hemostasis valve.
3. Temporarily occlude hemostasis valve with sterile-gloved finger.
4. **Warning:** Hemostasis valve must be occluded at all times to minimize the risk of air embolism or hemorrhage.
5. **Precaution:** Insert Arrow® obturator, sold separately, as a dummy catheter into hemostasis valve. This will ensure that leakage does not occur and inner seal is protected from contamination.

References:

1. Bar-Joseph G, Galvis AG. Perforation of the heart by central venous catheters in infants: guidelines to diagnosis and management. *J Pediatr Surg.* 1983;18:284-287.
2. Blitt CD, ed. *Monitoring in Anesthesia and Critical Care Medicine.* Central venous pressure monitoring. New York, NY: Churchill Livingstone; 1985:121-165.
3. Brandt RL, Foley WJ, Fink GH, Regan WJ. Mechanism of perforation of the heart with production of hydropericardium by a venous catheter and its prevention. *Am J Surg.* 1970;119:311-316.
4. Collier PE, Ryan JJ, Diamond DL. Cardiac tamponade from central venous catheters – report of a case and review of the English literature. *Angiology.* September 1984;35:595-600.
5. Conn C. The importance of syringe size when using implanted vascular access devices. *J Vasc Access Nurs.* Winter 1993;3:11-18.
6. Curelaru J, Linder LE, Gustavsson B. Displacement of catheters inserted through internal jugular veins with neck flexion and extension. *Intens Care Med.* 1980;6:179-183.
7. Iberti TJ, Katz LB, Reiner MA, Brownie T, Kwun KB. Hydrothorax as a late complication of central venous indwelling catheters. *Surgery.* November 1983:842-846.

8. Maschke SP, Rogove HJ. Cardiac tamponade associated with a multilumen central venous catheter. *Crit Care Med.* 1984;12:611-612.
9. Peters JL, ed. *A Manual of Central Venous Catheterization and Parenteral Nutrition*. Boston, MA: John Wright PSG; 1983:58-61, 155-157.
10. Sheep RE, Guiney WB Jr. Fatal cardiac tamponade. *JAMA.* 1982;248:1632-1635.
11. Sigurdsson J, Riba P, Sigurdsson S. The wandering central venous catheter. *Intensive Care Med.* 1985;11:263-264.
12. Tocino IM, Watanabe A. Impending catheter perforation of superior vena cava: radiographic recognition. *Am J Roentgenology.* March 1986;146:487-490.

Arrow International, Inc. recommends that the user be acquainted with the reference literature.

*If you have any questions or would like additional reference information, please contact Arrow International, Inc.

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Arrow International, Inc. kullanıcının referans literatüre aşına olmasını önerir.
*Sorularınız varsa veya ek referans bilgisi isterseniz lütfen Arrow International, Inc. ile irtibat kurun.