Inspected Dimensions:  
Folded Length: 7-1/2" (19 cm)  
Folded Width: 5-1/2" (14 cm)

# ARROW®

|           |  |    |
|-----------|--|----|
| <b>en</b> | ECG (Electrocardiogram) Extension Cable          | 3  |
| <b>be</b> | Кабель-падаўжальнік для ЭКГ (электракардыёграмы) | 4  |
| <b>hr</b> | Produžni kabel za ECG (elektrokardiogram)        | 5  |
| <b>cs</b> | Prodlužovací kabel pro EKG (elektrokardiografii) | 6  |
| <b>da</b> | EKG (elektrokardiogram) forlænger kabel          | 7  |
| <b>nl</b> | ECG-verlengsnoer (Elektrocardiogram)             | 8  |
| <b>et</b> | EKG (elektrokardiogrammi) pikendusjuhe           | 9  |
| <b>fi</b> | EKG-jatkokaapeli                                 | 10 |
| <b>fr</b> | Rallonge à ECG (Électrocardiogramme)             | 11 |
| <b>de</b> | EKG (Elektrokardiogramm) Verlängerungskabel      | 12 |
| <b>el</b> | Καλώδιο επέκτασης για ΗΚΓ (ηλεκτροκαρδιογράφημα) | 13 |
| <b>hu</b> | EKG (Elektrokardiogram) hosszabbító kábel        | 14 |
| <b>is</b> | EKG (Hjartarafrit) framlengingarleiðsla          | 15 |
| <b>it</b> | Cavo di prolunga per ECG (elettrocardiogramma)   | 16 |
| <b>lt</b> | EKG (elektrokardiogramos) ilginamasis laidas     | 17 |
| <b>no</b> | Forlengelseskabel til EKG (elektrokardiogram)    | 18 |
| <b>pl</b> | Przedłużacz do EKG                               | 19 |

|           |  |    |
|-----------|--|----|
| <b>pt</b> | Cabo de Extensão para ECG (Electrocardiograma)     | 20 |
| <b>ro</b> | Cablu de extensie ECG (Electrocardiogramă)         | 21 |
| <b>ru</b> | Удлинительный кабель для ЭКГ (электрокардиограмма) | 22 |
| <b>sr</b> | Produžni kabl za EKG (elektrokardiogram)           | 23 |
| <b>sk</b> | Predlžovací kábel pre EKG (elektrokardiogram)      | 24 |
| <b>sl</b> | Kabelski podaljšek za EKG (elektrokardiogram)      | 25 |
| <b>es</b> | Cable de extensión para ECG (electrocardiograma)   | 26 |
| <b>sv</b> | EKG-förlängningskabel (Elektrokardiogramkabel)     | 27 |
| <b>tr</b> | EKG (Elektrokardiyografi) Uzatma Kablosu           | 28 |
| <b>uk</b> | Подовжувальний кабель для ЕКГ (електрокардіограма) | 29 |

## ECG (Electrocardiogram) Extension Cable

Rx only.

### Safety and Efficacy Considerations:

Do not use if package has been previously opened or damaged.

**Warning:** Prior to use, read all package insert warnings, precautions, and instructions. Failure to do so may result in severe patient injury or death.

Do not alter cable during use.

ECG Extension Cable use must be performed by trained personnel well versed in anatomical landmarks, safe technique, and potential complications.

### Indications for Use:

The ECG Extension Cable provides confirmation of central venous catheter placement by right atrial electrophysiology.

### Contraindications:

None known.

### Warnings and Precautions:\*

- Warning: Sterile, Single use: Do not reuse, reprocess or resterilize. Reuse of device creates a potential risk of serious injury and/or infection which may lead to death.**
- Warning: Strict electrical isolation and grounding techniques must be used for monitoring equipment.**
- Warning: Do not leave cable attached to spring-wire guide when cable not attached to monitor. An electrical shock hazard is created due to a direct electrical path to the patient. Always detach from spring-wire guide first when disconnecting ECG extension cable.**
- Warning: Due to the risk of exposure to HIV (Human Immunodeficiency Virus) or other blood borne pathogens, health care workers should routinely use universal blood and body-fluid precautions in the care of all patients.**

### A Suggested Procedure:

#### Use sterile technique.

- a) for Arrow catheters 20 and 30 cm:*

Place the catheter into vessel using over-the-wire technique per manufacturer's instructions. When 1 cm-wide mark on spring wire guide exits catheter hub, the wire "J" is just outside of catheter tip.

- b) for Arrow catheters 16 cm:*

Place the catheter into vessel using over-the-wire technique per manufacturer's instructions. When 32 cm mark on spring wire guide exits catheter hub, the wire "J" is just outside of catheter tip.

- c) for non-Arrow catheters:*

Follow instructions for use of the manufacturer of the catheter and spring wire.

- Connect cable to monitor.
- Connect alligator clamp of ECG extension cable to spring-wire guide at 1 cm-wide mark.
- Precaution: Do not allow spring-wire guide to move within the catheter.** Advance catheter and spring-wire guide as a unit until P wave of ECG trace begins to progressively get larger. At this point, the catheter tip should be at the cavo-atrial junction. If the tip is advanced further into the right atrium, the P wave will either be as large as the R wave or become biphasic. For catheters 16 cm or shorter the catheter may not reach the right atrium, and this technique can only be used to determine if the catheter is facing the proper direction. Attach extension cable as described above and advance catheter and spring-wire guide as a unit. When the P wave becomes higher than normal the catheter is facing the proper direction.
- Withdraw catheter and spring-wire guide as a unit 2-3 cm from cavo-atrial junction. At this position the tip should be within the superior vena cava.
- Holding catheter in place remove spring-wire guide.
- Continue procedure per manufacturer's instructions.

\*If you have any questions or would like additional reference information, please contact Arrow International, Inc.

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