ECG (Electrocardiogram) Extension Cable

Rx only.

Indications for Use:

The ECG Extension Cable permits verification of central venous catheter placement by right atrial electrophysiology.

Contraindications:

None known.

General Warnings and Precautions

Warnings:

1. **Sterile, Single use: Do not reuse, reprocess or resterilize.** Reuse of device creates a potential risk of serious injury and/or infection which may lead to death. Reprocessing of medical devices intended for single use only may result in degraded performance or a loss of functionality.
2. **Read all package insert warnings, precautions and instructions prior to use.** Failure to do so may result in severe patient injury or death.
3. **Strict electrical isolation and grounding techniques must be used for monitoring equipment.**
4. **Do not leave cable attached to guidewire when cable not attached to monitor.** An electrical shock hazard is created due to a direct electrical path to the patient. Always detach from guidewire first when disconnecting ECG extension cable.

Precautions:

1. **Do not alter the cable during use.**
2. **Procedure must be performed by trained personnel well versed in anatomical landmarks, safe technique and potential complications.**
3. **Use standard precautions and follow institutional policies for all procedures including safe disposal of devices.**

A Suggested Procedure: Use sterile technique.

1. Place the catheter into vessel using over-the-wire technique per manufacturer's instructions. When using a marked guidewire, refer to Table 1 to determine when the wire is just outside the catheter tip based on Arrow CVC catheter lengths.

16 cm catheter	27 cm mark
20 cm catheter	32 cm mark
30 cm catheter	40 cm mark

Table 1

Non-Arrow catheters:

Follow instructions for use of the manufacturer of the catheter and guidewire.

2. Connect cable to monitor.
3. Connect alligator clamp of ECG extension cable to guidewire at referenced mark per Table 1 based on Arrow CVC catheter length.

Precaution: Do not allow guidewire to move within the catheter.

4. Advance catheter and guidewire as a unit until P wave of ECG trace begins to progressively get larger. At this point, the catheter tip should be at the cavo-atrial junction. If the tip is advanced further into the right atrium, the P wave will either be as large as the R wave or become biphasic.

For catheters 16 cm or shorter the catheter may not reach the right atrium, and this technique can only be used to determine if the catheter is facing the proper direction. Attach extension cable as described above and advance catheter and guidewire as a unit. When the P wave becomes higher than normal the catheter is facing the proper direction.

5. Withdraw catheter and guidewire as a unit 2-3 cm from cavo-atrial junction. At this position the tip should be within the superior vena cava.
6. Holding catheter in place remove guidewire.
7. Continue procedure per manufacturer's instructions.

For reference literature concerning patient assessment, clinician education, insertion technique, and potential complications associated with this procedure, consult standard textbooks, medical literature, and Arrow International LLC website: www.teleflex.com

A pdf copy of this IFU is located at www.teleflex.com/IFU

For a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation 2017/745/EU on Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority. The contacts of national competent authorities (Vigilance Contact Points) and further information can be found on the following European Commission website: https://ec.europa.eu/growth/sectors/medical-devices/contacts_en

en

Symbol Glossary: Symbols are in compliance with ISO 15223-1

								
Caution	Medical device	Consult instructions for use	Do not reuse	Do not resterilize	Sterilized by ethylene oxide	Single sterile barrier system	Keep away from sunlight	Keep dry
								
Do not use if package is damaged	Not made with natural rubber latex	Catalogue number	Lot number	Use by	Manufacturer	Date of manufacture	Importer	

Teleflex, the Teleflex logo, Arrow and the Arrow logo are trademarks or registered trademarks of Teleflex Incorporated or its affiliates, in the U.S. and/or other countries. © 2020 Teleflex Incorporated. All rights reserved.

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