SecondSite™ Adjustable Hub: Catheter Clamp and Fastener Rx only.

Indications for Use:

The catheter clamp and fastener are indicated for use in securing catheter to patient or as a secondary suture site.

Contraindications:

None known.

Warnings and Precautions:

Warnings:

1. **Sterile, Single use: Do not reuse, reprocess or resterilize. Reuse of device creates a potential risk of serious injury and/or infection which may lead to death.**
2. **Read all package insert warnings, precautions and instructions prior to use. Failure to do so may result in severe patient injury or death.**
3. **Do not attach catheter clamp and fastener until guidewire is removed.**

Cautions:

1. **Procedure must be performed by trained personnel well versed in anatomical landmarks, safe technique and potential complications.**
2. **Do not alter the catheter clamp and fastener during use or removal.**

A Suggested Procedure:

Use sterile technique.

1. Follow manufacturer's instructions for placing catheter.
2. After guidewire has been removed and necessary lines have been connected or locked, spread wings of rubber clamp and position on catheter making sure catheter is not moist, as required, to maintain proper tip location (refer to Fig. 1).


Fig. 1

3. Snap rigid fastener onto catheter clamp (refer to Fig. 2).


Fig. 2

4. Secure catheter clamp and fastener as a unit to patient by using either catheter stabilization device, stapling or suturing. Both catheter clamp and fastener need to be secured to reduce risk of catheter migration (refer to Fig. 3).


Fig. 3

5. Continue procedure per manufacturer's instructions.

For reference literature concerning patient assessment, clinician education, insertion technique, and potential complications associated with this procedure, consult standard textbooks, medical literature, and Arrow International, Inc. website: www.teleflex.com

Teleflex and the Teleflex logo are trademarks or registered trademarks of Teleflex Incorporated or its affiliates, in the U.S. and/or other countries. © 2017 Teleflex Incorporated. All rights reserved.

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