SecondSite™ Adjustable Hub: Catheter Clamp and Fastener

Rx only

Intended Purpose:

The catheter clamp and fastener are indicated for use to secure catheter when an additional securement site other than the catheter hub is required for catheter stabilization.

Patient Target Group:

Intended to be used when a vascular catheter requires secondary securement.

Contraindications:

None known.

Clinical Benefits to be Expected:

The Arrow® Catheter clamp and fastener are indicated for use in securing a catheter to patient as a secondary securement site.

Performance Characteristics:

This system provides a method to secure the exposed portion of a partially inserted intravascular catheter to prevent catheter migration.

⚠ General Warnings and Precautions

Warnings:

1. Sterile, Single use: Do not reuse, reprocess or resterilize. Reuse of device creates a potential risk of serious injury and/or infection which may lead to death. Reprocessing of medical devices intended for single use only may result in degraded performance or a loss of functionality.
2. Read all package insert warnings, precautions and instructions prior to use. Failure to do so may result in severe patient injury or death.
3. Do not attach catheter clamp and fastener until either guidewire or placement wire is removed.

Precautions:

1. Do not alter any component during use or removal.
2. Procedure must be performed by trained personnel well versed in anatomical landmarks, safe technique and potential complications.
3. Use standard precautions and follow institutional policies for all procedures including safe disposal of devices.
4. If the package is damaged or unintentionally opened before use do not use the device. Dispose of the device.
5. Storage conditions for these devices require that they are kept dry and out of direct sunlight.

A Suggested Procedure: Use sterile technique.

1. Follow manufacturer's instructions for placing catheter.
2. After guidewire has been removed and necessary lines have been connected or locked, spread wings of rubber clamp and position on catheter body making sure catheter surface is not moist to maintain proper securement (refer to Figure 1).


Figure 1

3. Snap rigid fastener onto catheter clamp (refer to Figure 2).


Figure 2

4. Secure catheter clamp and fastener as a unit to patient by using either catheter stabilization device, stapling or suturing. Catheter juncture hub and both catheter clamp and fastener need to be secured to minimize the risk of catheter migration (refer to Figure 3).


Figure 3

5. Continue procedure per manufacturer's instructions.

For reference literature concerning patient assessment, clinician education, insertion technique, and potential complications associated with this procedure, consult standard textbooks, medical literature, and Arrow International LLC website: www.teleflex.com

A pdf copy of this IFU is located at www.teleflex.com/IFU

For a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation 2017/745/EU on Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority. The contacts of national competent authorities (Vigilance Contact Points) and further information can be found on the following European Commission website: https://ec.europa.eu/growth/sectors/medical-devices/contacts_en

Symbol Glossary: Symbols are in compliance with ISO 15223-1.

							
Caution	Medical device	Consult instructions for use	Do not reuse	Do not resterilize	Sterilized by ethylene oxide	Single sterile barrier system	Keep away from sunlight
							
Keep dry	Do not use if package is damaged	Not made with natural rubber latex	Catalogue number	Lot number	Use by	Manufacturer	Date of manufacture
	Arrow, the Arrow logo, Teleflex, and the Teleflex logo, are trademarks or registered trademarks of Teleflex Incorporated or its affiliates, in the U.S. and/or other countries. © 2023 Teleflex Incorporated. All rights reserved. "Rx only" is used within this labeling to communicate the following statement as presented in the FDA CFR : Caution: Federal law restricts this device to sale by or on the order of a licensed healthcare practitioner.						
Importer							

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