Inspected Dimensions:  
Folded Length: 7-1/2" (19 cm)  
Folded Width: 5-1/2" (14 cm)

# ARROW®

INTERNATIONAL

<b>EN</b>	<b>Pediatric Multiple-Lumen Central Venous Catheterization Product</b>	<b>1</b>
<b>FR</b>	<b>Produit pédiatrique de cathétérisme veineux central à une lumière multiples</b>	<b>5</b>
<b>DE</b>	<b>Pädiatisches multilumige Zentralvenen-Katheterisierungsprodukt</b>	<b>8</b>
<b>IT</b>	<b>Prodotto per cateterizzazione venosa centrale a lume multiplo per uso pediatrico</b>	<b>11</b>
<b>PL</b>	<b>Pediatriczny wieloświatłowy cewnik do cewnikowania żył centralnych</b>	<b>14</b>
<b>PT</b>	<b>Produto para Cateterização Venosa Central com Lúmen Múltiplo Pediátrico</b>	<b>17</b>
<b>RU</b>	<b>Педиатрическое многопросветное изделие для катетеризации главных вен</b>	<b>20</b>
<b>SL</b>	<b>Pripomoček za pediatrično centralno venozno kateterizacijo z več svetlinami</b>	<b>23</b>
<b>ES</b>	<b>Producto para cateterización de vena central de luz múltiple para uso pediátrico</b>	<b>26</b>
<b>SV</b>	<b>Pediatriisk flerkanaliga centralvenkatetreringsprodukt</b>	<b>29</b>
<b>TR</b>	<b>Pediyatrik Çoklu Lümenli Santral Venöz Kateterizasyon Ürünü</b>	<b>32</b>

## Pediatric Multiple-Lumen Central Venous Catheterization Product

### Safety and Efficacy Considerations:

Do not use if package has been previously opened or damaged.

**Warning:** Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Do not alter the catheter, spring-wire guide, or any other kit/set component during insertion, use, or removal.

Procedure must be performed by trained personnel well versed in anatomical landmarks, safe technique, and potential complications.

**Warning:** Do not place the catheter into or allow it to remain in the right atrium or right ventricle (refer to Fig. 1).


Fig. 1

### Cardiac Tamponade:

It has been documented by many authors that placement of indwelling catheters in the right atrium is a dangerous practice<sup>1,2,3,5,6,11,14,15,16,23</sup> that may lead to cardiac perforation and tamponade.<sup>1,2,3,5,6,11,14,15,16,23</sup> Although cardiac tamponade secondary to pericardial effusion is uncommon, there is a high mortality rate associated with it.<sup>22</sup> Practitioners placing central venous catheters must be aware of this potentially fatal complication before advancing the catheter too far relative to patient size.

No particular route or catheter type is exempt from this potentially fatal complication.<sup>18</sup> The actual position of the tip of the indwelling catheter should be confirmed by x-ray after insertion.<sup>1,2,3,6,11,16,23</sup> Central venous catheters should be placed in the superior vena cava<sup>1,2,3,6,11,15,16,20,23</sup> above its junction with the right atrium and parallel to the vessel wall and its distal tip positioned at a level above either the azygos vein or the carina of the trachea, whichever is better visualized.

Central venous catheters should not be placed in the right atrium unless specifically required for special relatively short term procedures, such as aspiration of air emboli during neurosurgery. Such procedures are nevertheless risk prone and should be closely monitored and controlled.

### Indications for Use:

The multiple-lumen catheter permits venous access to the central circulation of pediatric patients.

### Contraindications:

None known.

### Warnings and Precautions:\*

1. **Warning:** Sterile, Single use: Do not reuse, reprocess or resterilize. Reuse of device creates a potential risk of serious injury and/or infection which may lead to death.

2. **Warning:** Do not place the catheter into or allow it to remain in the right atrium or right ventricle. Central vein catheters should be positioned so that the distal tip of the catheter is in the lower 1/2 to 1/3 of the superior vena cava (SVC) above the junction of the SVC and the right atrium and lies parallel to the vessel wall. For femoral vein approach the catheter should be advanced into the vessel so that the catheter tip lies parallel to the vessel wall and does not enter the right atrium.

Due to the variety of pediatric catheter lengths available, patient size must be carefully considered relative to actual length of catheter introduced.

3. **Warning:** Practitioners must be aware of complications associated with central vein catheters including cardiac tamponade secondary to vessel wall, atrial or ventricular perforation, pleural and mediastinal injuries, air embolism, catheter embolism, thoracic duct laceration, bacteremia, septicemia, thrombosis, inadvertent arterial puncture, nerve damage, hematoma, hemorrhage, and dysrhythmias.

4. **Warning:** Do not apply excessive force in removing guide wire or catheters. If withdrawal cannot be easily accomplished, a chest x-ray should be obtained and further consultation requested.

5. **Warning:** The practitioner must be aware of potential air embolism associated with leaving open needles or catheters in central venous puncture sites or as a consequence of inadvertent disconnects. To lessen the risk of disconnects, only securely tightened Luer-Lock connections should be used with this device. Follow hospital protocol to guard against air embolism for all catheter maintenance.

6. **Warning:** Passage of the guide wire into the right heart can cause dysrhythmias, right bundle branch block,<sup>9</sup> and vessel wall, atrial or ventricular perforation.

7. **Warning:** Due to the risk of exposure to HIV (Human Immunodeficiency Virus) or other blood borne pathogens, health care workers should routinely use universal blood and body-fluid precautions in the care of all patients.

8. **Precaution:** Indwelling catheters should be routinely inspected for desired flow rate, security of dressing, correct catheter position and for secure Luer-Lock connection. Use centimeter markings to identify if the catheter position has changed.

9. **Precaution:** Only x-ray examination of the catheter placement can ensure that the catheter tip has not entered the heart or no longer lies parallel to the vessel wall. If catheter position has changed, immediately perform chest x-ray examination to confirm catheter tip position.

10. **Precaution:** For blood sampling, temporarily shut off remaining port(s) through which solutions are being infused.

11. **Precaution:** Alcohol and acetone can weaken the structure of polyurethane materials. Check ingredients of prep sprays and swabs for acetone and alcohol content.

**Acetone:** Do not use acetone on catheter surface. Acetone may be applied to skin but must be allowed to dry completely prior to applying dressing.

**Alcohol:** Do not use alcohol to soak catheter surface or to restore catheter patency. Care should be taken when instilling drugs containing high concentration of alcohol.

Always allow alcohol to dry completely prior to applying dressing.

12. **Precaution:** Some disinfectants used at the catheter insertion site contain solvents, which can attack the catheter material. Assure insertion site is dry before dressing.

13. **Precaution:** Use of a syringe smaller than 10 mL to irrigate or declot an occluded catheter may cause intraluminal leakage or catheter rupture.<sup>7</sup>

#### A Suggested Procedure:

##### Use sterile technique.

1. **Precaution:** Place patient in slight Trendelenburg position as tolerated to reduce the risk of air embolism. If femoral approach is used, place patient in supine position.
2. Measure the patient's external anatomy to estimate the length of catheter required to place the distal tip in the SVC.
3. Prep and drape puncture site as required.
4. Perform skin wheal with desired needle. In kits where provided, a SharpsAway® disposal cup is used for the disposal of needles. Push needles into foam after use. Discard entire cup at completion of procedure. **Precaution: Do not reuse needles after they have been placed into the disposal cup. Particulate matter may adhere to needle tip.**
5. Prepare the catheter for insertion by flushing each lumen and clamping or attaching the injection caps to the appropriate extension lines. Leave the distal extension line uncapped for guide wire passage. **Warning: Do not cut the catheter to alter length.**
6. Locate central vein with a needle and syringe appropriate to the size of the patient.
7. Insert introducer catheter/needle with attached syringe into vein along side of locator needle and aspirate. Remove locator needle. Withdraw needle from introducer catheter. If no free flow of venous blood is observed after needle is removed, attach syringe to the catheter and aspirate until good venous blood flow is established. **Precaution: The color of the blood aspirated is not always a reliable indicator of venous access.<sup>12</sup> Do not reinsert needle into introducer catheter.** Because of the potential for inadvertent arterial placement, verify venous access via a wave form obtained by a calibrated pressure transducer (refer to Fig. 2).


Fig. 2

If hemodynamic monitoring equipment is not available to permit transducing a central venous wave form, disconnect the syringe and check for pulsatile flow. Pulsatile flow is usually an indicator of inadvertent arterial puncture.

##### Alternate Technique:

Introducer needle may be used in the standard manner as alternative to catheter/needle assembly.

8. Insert desired tip of spring-wire guide through the introducer needle or catheter into vein. **Precaution: Due to the fragile nature of the spring-wire guide contained in this product, inspect "J"-tip for damage prior to insertion.** If the "J"-tip portion of the spring-wire guide is used, prepare for insertion by sliding the plastic tube over the "J" to straighten. The spring-wire guide should then be advanced in the routine fashion to a depth appropriate to patient size. Advancement of "J"-tip may require a gentle rotating motion. **Warning: Do not cut spring-wire guide to alter length. Do not withdraw spring-wire guide against needle bevel to minimize the risk of possible severing or damaging of spring-wire guide.**
9. Hold spring-wire guide in place and remove introducer needle or catheter. **Precaution: Maintain firm grip on spring-wire guide at all times.**
10. Enlarge cutaneous puncture site with cutting edge of scalpel positioned away from the spring-wire guide. **Precaution: Do not cut guide wire.** Use vessel dilator to enlarge site as required. **Warning: Do not leave vessel dilator in place as an indwelling catheter to minimize the risk of possible vessel wall perforation.**
11. Thread tip of multiple-lumen catheter over spring-wire guide. Sufficient guide wire length must remain exposed at hub end of catheter to maintain a firm grip on guide wire. Grasping near skin, advance catheter into vein with slight twisting motion. **Precaution: Where provided, catheter clamp and fastener must not be attached to catheter until spring-wire guide is removed.**
12. Using centimeter marks on catheter as positioning reference points, advance catheter to final indwelling position.
13. Hold catheter at desired depth and remove spring-wire guide. The Arrow® catheter included in this product has been designed to freely pass over the spring-wire guide. If resistance is encountered when attempting to remove the spring-wire guide after catheter placement, the spring-wire may be kinked about the tip of the catheter within the vessel (refer to Fig. 3). In this circumstance, pulling back on the spring-wire guide may result in undue force being applied resulting in spring-wire guide breakage. If resistance is encountered, withdraw the catheter relative to the spring-wire guide about 2-3 cm and attempt to remove the spring-wire guide. If resistance is again encountered remove the spring-wire guide and catheter simultaneously. **Warning: Although the incidence of spring-wire guide failure is extremely low, practitioner should**

be aware of the potential for breakage if undue force is applied to the wire.


Fig. 3

14. Verify that the entire spring-wire guide is intact upon removal.
15. Check lumen placement by attaching a syringe to each extension line and aspirate until free flow of venous blood is observed. Connect all extension lines to appropriate Luer-Lock line(s) as required. Unused port(s) may be "locked" through injection cap(s) using standard hospital protocol. Slide clamps are provided on extension lines to occlude flow through each lumen during line and injection cap changes. **Precaution: To minimize the risk of damage to extension lines from excessive pressure, each clamp must be opened prior to infusing through that lumen.**
16. Secure and dress catheter temporarily.
17. Verify catheter tip position by chest x-ray immediately after placement. **Precaution: X-ray exam must show the catheter located in the right side of the mediastinum in the SVC with the distal end of the catheter parallel to the vena cava wall and its distal tip positioned at a level above either the azygos vein or the carina of the trachea, whichever is better visualized.** If catheter tip is malpositioned, reposition and re-verify.
18. Secure catheter to patient. Use juncture hub with integral suture ring and side wings as primary suture site. In kits where provided, the catheter clamp and fastener should be utilized as a secondary suture site as necessary. **Precaution: Do not suture directly to the outside diameter of the catheter to minimize the risk of cutting or damaging the catheter or impeding catheter flow.**

**Catheter Clamp and Fastener Instructions:**

After spring-wire guide has been removed and the necessary lines have been connected or locked, spread wings of rubber clamp and position on catheter as required to ensure proper tip location (refer to Fig. 4).


Fig. 4

- Snap rigid fastener onto catheter clamp (refer to Fig. 5).


Fig. 5

- Secure catheter to patient by suturing the catheter clamp and fastener together to the skin, using side wings to prevent catheter migration (refer to Fig. 6).


Fig. 6

19. Dress puncture site per hospital protocol. **Precaution: Maintain the insertion site with regular meticulous redressing using aseptic technique.**
20. Record on the patient's chart the indwelling catheter length as to centimeter markings on catheter where it enters the skin. Frequent visual reassessment should be made to ensure that the catheter has not moved.

**Catheter Exchange Procedure:**

1. Use sterile technique.
2. **Precaution: Prior to attempting a catheter exchange procedure, remove catheter clamp and fastener.**
3. Proceed per hospital protocol. Cutting the catheter is not recommended due to the potential for catheter embolism.

**Catheter Removal Procedure:**

1. **Precaution: Place the patient in a supine position.**
2. Remove dressing. **Precaution: To minimize the risk of cutting the catheter do not use scissors to remove the dressing.**
3. **Warning: Exposure of the central vein to atmospheric pressure may result in entry of air into the central venous**

**system.** Remove suture(s) from catheter clamp and primary suture site. Be careful not to cut the catheter. Remove catheter slowly, pulling it parallel to the skin. As catheter exits the site, apply pressure with a dressing impermeable to air, e.g. VASELINE<sup>®†</sup> gauze. Because the residual catheter track remains an air entry point until completely sealed, the occlusive dressing should remain in place for at least 24-72 hours dependent upon the amount of time the catheter was indwelling.<sup>3,17,19,24</sup>

4. Upon removal of the catheter, inspect it to make sure that the entire length has been withdrawn.
5. Document removal procedure.

#### References:

1. Agarwal KC, Khan MAA, Falla A, Amato JJ. Cardiac perforation from central venous catheters: survival after cardiac tamponade in an infant. *Pediatrics*. March 1984;73:333-338.
2. Bar-Joseph G, Galvis AG. Perforation of the heart by central venous catheters in infants: guidelines to diagnosis and management. *J Pediatr Surg*. June 1983;18:284-287.
3. Brandt RL, Foley WJ, Fink GH, Regan WJ. Mechanism of perforation of the heart with production of hydropericardium by a venous catheter and its prevention. *Am J Surg*. 1970;119:311-316.
4. Carey BE. Major complications of central lines in neonates. *Neonatal Network*. June 1989;7:17-28.
5. Casado-Flores J, Valdivielso-Serna A, Perez-Jurado L, et al. Subclavian vein catheterization in critically ill children: analysis of 322 cannulations. *Intensive Care Med*. 1991;17:350-354.
6. Collier PE, Ryan JJ, Diamond DL. Cardiac tamponade from central venous catheters – report of a case and review of the English literature. *Angiology*. September 1984;35:595-600.
7. Conn C. The importance of syringe size when using implanted vascular access devices. *J Vasc Access Nurs*. Winter 1993;3:11-18.
8. Daniels SR, Hannon DW, Meyer RA, Kaplan S. Paroxysmal supraventricular tachycardia. *AJDC*. 1984;138:474-475.
9. Eissa NT, Kvetan V. Guide wire as a cause of complete heart block in patients with preexisting left bundle branch block. *Anesthesiology*. 1990;73:772-774.
10. Franciosi RA, Ellefson RD, Uden D, Drake RM. Sudden unexpected death during central hyperalimentation. *Pediatrics*. March 1982;69:305-307.
11. Fuhrman BP. Perforation by vascular catheters. *Perspectives in Critical Care*. 1988;1:147-154.
12. Jobs DR, Schwartz AJ, Greenhow DE, Stephenson LW, Ellison N. Safer jugular vein cannulation: recognition of arterial punctures and preferential use of the external jugular route. *Anesthesiology*. 1983;59:353-355.
13. Kashuk JL, Penn I. Air embolism after central venous catheterization. *Surg Gynecol Obstet*. September 1984;159:249-252.
14. Leibovitz E, Ashkenazi A, Levin S, Nissim F. Fatal pericardial tamponade complicating total parenteral nutrition via a silastic central vein catheter. *J Pediatr Gastroenterol Nutr*. 1988;7:306-307.
15. Mupanemunda RH, Mackanjee HR. A life-threatening complication of percutaneous central venous catheters in neonates. *AJDC*. 1992;146:1414-1415. Letter.
16. Otto CW. Central venous pressure monitoring. In: Blitt CD, ed. *Monitoring in Anesthesia and Critical Care Medicine*. New York, NY: Churchill Livingstone, 1985:139.
17. Paskin DL, Hoffman WS, Tuddenham WJ. A new complication of subclavian vein catheterization. *Ann Surg*. March 1974;179:266-268.
18. Peters JL, ed. *A Manual of Central Venous Catheterization and Parenteral Nutrition*. Boston, MA: John Wright PSG, 1983;58-61, 155-157.
19. Phifer TJ, Bridges M, Conrad SA. The residual central venous catheter track – an occult source of lethal air embolism: case report. *J Trauma*. 1991;31:1558-1560.
20. Rogers BB, Berns SD, Maynard EC, Hansen TWR. Pericardial tamponade secondary to central venous catheterization and hyperalimentation in a very low birth weight infant. *Pediatr Pathology*. 1990;10:819-823.
21. Ross P Jr, Seashore JH. Bilateral hydrothorax complicating central venous catheterization in a child: case report. *J Pediatr Surg*. 1989;24:263-264.
22. Sheep RE, Guiney WB Jr. Fatal cardiac tamponade. *JAMA*. 1982;248:1632-1635.
23. Sullivan CA, Konefal SH Jr. Cardiac tamponade in a newborn: a complication of hyperalimentation. *J Parenteral Enteral Nutrition*. May/June 1987;11:319-321.
24. Thielen JB, Nyquist J. Subclavian catheter removal. *J Intravenous Nurs*. March/April 1991;14:114-118.

Arrow International, Inc. recommends that the user be acquainted with the reference literature.

\*If you have any questions or would like additional reference information, please contact Arrow International, Inc.

†A registered trademark of Unilever Supply Chain, Inc.

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Arrow International, Inc. kullanıcının referans literatüre aşına olmasını önerir.
\*Sorularınız varsa veya ek referans bilgisi isterseniz lütfen Arrow International, Inc. ile irtibat kurun.