Inspected Dimensions:
Folded Length: 7-1/2" (19 cm)
Folded Width: 5-1/2" (14 cm)

ARROW®

INTERNATIONAL

EN	Spring-Wire Guide Product	2
NL	Voerdraad met veer	4
FR	Fil de guidage métallisé spiralé	5
DE	Federführungsdraht	6
EL	Προϊόν Συρμάτινου Οδηγού με Έλασμα	7
IT	Prodotto per guida metallica a molla	8
PL	Prowadnik	9
PT	Produto do fio-guia	10
RU	Пружинный проводник	11
ES	Guía de alambre flexible	12
SV	Fjädertrådsledareprodukt	13
TR	Yaylı Kılavuz Tel Ürünü	14

ARROW®

Spring-Wire Guide Product

Safety and Efficacy Considerations:

Do not use if package has been previously opened or damaged.

Warning: Prior to use read all package insert warnings, precautions, and instructions. Failure to do so may result in severe patient injury or death.

Do not alter spring-wire guide during insertion, use, or removal.

Spring-wire guide use must be performed by trained personnel well versed in anatomical landmarks, safe technique and potential complications.

Indications for Use:

The spring-wire guide permits central venous catheter introduction or exchange.

Contraindications:

None known.

Warnings and Precautions:*

1. **Warning:** Sterile, Single use: Do not reuse, reprocess or resterilize. Reuse of device creates a potential risk of serious injury and/or infection which may lead to death.
2. **Warning:** Practitioners must be aware of complications associated with spring-wire guide use including knotting^{8,10}, separation of the coil and core wires^{1,3,5,7,12,13}, vessel perforation¹¹, malposition¹⁶, displacement of vena cava filters²⁻¹⁴, cardiac tamponade^{4,9,17} and cardiac dysrhythmias^{6,15}.
3. **Warning:** Do not apply excessive force in removing guide wire. If withdrawal cannot be easily accomplished, a chest x-ray should be obtained and further consultation requested.
4. **Warning:** Due to the risk of exposure to HIV (Human Immunodeficiency Virus) or other blood borne pathogens, health care workers should routinely use universal blood and bodyfluid precautions in the care of all patients.
5. **Precaution:** If resistance is encountered when attempting to remove spring-wire guide after catheter placement, spring-wire may be kinked about tip of catheter within vessel (refer to Fig. 1).


Fig. 1

To minimize the risk of spring-wire guide damage withdraw catheter relative to spring-wire guide about 2-3 cm and attempt to remove spring-wire guide. If resistance is again encountered remove spring-wire guide and catheter simultaneously.

6. **Precaution:** Prior to attempting catheter exchange procedure remove catheter clamp and fastener.

A Suggested Procedure:

Use sterile technique.

1. Follow manufacturer's instructions for patient preparation and venous access.
2. Insert desired tip of spring-wire guide through the introducer needle or catheter into vein.
Straightening Tube: If the "J" tip portion of the spring-wire guide is used, prepare for insertion by sliding the plastic tube over the "J" to straighten.
Arrow Advancer™: Straighten spring-wire guide "J" by retracting into Arrow Advancer™ with thumb.
3. The spring-wire guide should then be advanced in the routine fashion to the desired depth. Advancement of "J" tip may require a gentle rotating motion. **Warning:** Do not cut spring-wire guide to alter length. Do not withdraw spring-wire guide against needle bevel to minimize the risk of possible severing or damaging of spring-wire guide.
4. Hold spring-wire guide in place and remove introducer needle or catheter. **Precaution:** Maintain firm grip on spring-wire guide at all times.
5. Continue procedure per manufacturer's instructions.
6. Verify that entire spring-wire guide is intact upon removal.

References:

1. Akazawa S, Nakaigawa Y, Hotta K, Shimizu R, Kashiwagi H, Takahashi K. Unrecognized migration of an entire guidewire on insertion of a central venous catheter into the cardiovascular system. *Anesthesiology*. 1996;84:241-242. Letter.
2. Amesbury S, Vargish T, Hall J. An unusual complication of central venous catheterization. *Chest*. 1994;105:905-907.
3. Bhanumurthy S, McMurray TJ, Bannon C. Early venous catheter blockade: a cause for immediate concern. *Anaesthesia*. 1995;50:473. Letter.
4. Blake PG, Uldall R. Cardiac perforation by a guide wire during subclavian catheter insertion. *Int J Artif Organs*. 1989;12:111-113.
5. Cucchiara RF, Muzzi DA. Guide-wire retention after right atrial catheter insertion. *Anesth Analg*. 1992;74:303-304.
6. Eissa NT, Kvetan V. Guide wire as a cause of complete heart block in patients with preexisting left bundle branch block. *Anesthesiology*. 1990;73:772-774.
7. Hickman JA, McCrerrick A. Failure of a spring wire guide. *Anaesthesia and Intensive Care*. 1990;18:587. Letter.
8. Hoogsteden L, Filshie J, Querci della Rovere G. Knotted guidewire – a complication of Hickman line insertion. *Anaesthesia*. 1996;51:713. Letter.
9. Khan IH. Pericardial tamponade resulting from changing a central venous catheter over a guide wire. *Postgraduate Med J*. 1991;67:591.
10. Kwok SPY, Lau WY, Li AKC. Insertion of subclavian line: don't get knotted. *Aust NZ J Surg*. 1994;64:64-65.
11. Masud T, Tapson JS. Mediastinal hematoma caused by subclavian catheterization for hemodialysis. *Int J Artif Organs*. 1989;12:708-710.
12. Polos PG, Sahn SA. Complication of central venous catheter insertion: fragmentation of a guidewire with pulmonary artery embolism. *Crit Care Med*. 1991;19:438-440.
13. Reynen K. 14-year follow-up of central embolization by a guide wire. *N Engl J Med*. 1993;329:970-971. Letter.
14. Rosenblum JD, Boyle CM. Percutaneous retrieval of a Vena-Tech filter displaced during central line placement. *AJR*. 1996;166:994-995. Letter.
15. Stuart RK, Shikora SA, Akerman P, et al. Incidence of arrhythmia with central venous catheter insertion and exchange. *J Parenteral Enteral Nutrition*. 1990;14:152-155.
16. Trujillo MH, Arai K. Hydrothorax after inadvertent placement of a central venous catheter in the left pericardiophrenic vein. *J Intensive Care Med*. 1994;9:257-260.
17. van Haeflen TW, van Pampus ECM, Boot H, Strack van Schijndel RJM, Thijs LG. Cardiac tamponade from misplaced central venous line in pericardiophrenic vein. *Arch Intern Med*. 1988;148:1649-1650.

Arrow International, Inc. recommends that the user be acquainted with the reference literature.

*If you have any questions or would like additional reference information, please contact Arrow International, Inc.

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Arrow International, Inc. kullanıcının referans literatürüne aşina olmasını tavsiye eder.
* Herhangi bir sorunuz olursa veya ek referans bilgisi isterseniz lütfen Arrow International, Inc. ile irtibata geçin.
Arrow International Inc.