Inspected Dimensions:  
Folded Length: 7-1/2" (19 cm)  
Folded Width: 5-1/2" (14 cm)

# ARROW®

INTERNATIONAL

<b>EN</b>	<b>Pediatric Multiple-Lumen Central Venous Catheterization with ARROWg<sup>ard</sup>® Antimicrobial Surface</b>	<b>1</b>
<b>FR</b>	<b>Produit de cathétérisme veineux central à lumières multiples pédiatrique avec surface antimicrobes ARROWg<sup>ard</sup></b>	<b>7</b>
<b>DE</b>	<b>Einführbesteck für viellumigen, zentralen Venenkatheter für Kinder mit antimikrobieller ARROWg<sup>ard</sup>-Oberflächenbeschichtung</b>	<b>11</b>
<b>IT</b>	<b>Prodotto multilume per cateterismo della vena centrale, pediatrico con superficie antimicrobica ARROWg<sup>ard</sup></b>	<b>15</b>
<b>PL</b>	<b>Wielokanalowy cewnik pediatryczny do żył centralnych z powłoką przeciwbakteryjną ARROWg<sup>ard</sup></b>	<b>19</b>
<b>PT</b>	<b>Produto para Cateterização Venosa Central com Lúmen Múltiplo Pediátrico com Superfície Antimicrobiana ARROWg<sup>ard</sup></b>	<b>23</b>
<b>RU</b>	<b>Педиатрическое многопросветное изделие для катетеризации центральных вен с противомикробной поверхностью ARROWg<sup>ard</sup></b>	<b>27</b>
<b>SL</b>	<b>Pediatrični pripomoček z več svetlinami za katetrizacijo centralne vene s antimikrobno površino ARROWg<sup>ard</sup></b>	<b>31</b>
<b>ES</b>	<b>Dispositivo de cateterización venosa central de luz múltiple para uso pediátrico con superficie antimicrobiana ARROWg<sup>ard</sup></b>	<b>35</b>
<b>SV</b>	<b>Flerkanalig produkt för pediatrik central venkatetrering med ARROWg<sup>ard</sup> antimikrobiell yta</b>	<b>39</b>
<b>TR</b>	<b>Pediyatrik Çoklu Lümen Santral Venöz Kateterizasyon Ürünü, ARROWg<sup>ard</sup> Antimikrobiyel Yüzevli</b>	<b>43</b>

EN

ARROW®

## Pediatric Multiple-Lumen Central Venous Catheterization Product with ARROWg<sup>ard</sup>® Antimicrobial Surface

### Safety and Efficacy Considerations:

Do not use if package has been previously opened or damaged.

**Warning:** Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Do not alter the catheter, spring-wire guide, or any other kit/ set component during insertion, use, or removal.

Procedure must be performed by trained personnel well versed in anatomical landmarks, safe technique, and potential complications.

**Warning:** Do not place the catheter into or allow it to remain in the right atrium or right ventricle. Read instructions (refer to Fig. 1).


Fig. 1

**Cardiac Tamponade:** It has been documented by many authors that placement of indwelling catheters in the right atrium is a dangerous practice<sup>1,3,4,6,8,14,22,26,29,36,41</sup> that may lead to cardiac perforation and tamponade<sup>1,3,4,6,8,14,22,26,29,36</sup>. Although cardiac tamponade secondary to pericardial effusion is uncommon, there is a high mortality rate associated with it<sup>35</sup>. Practitioners placing central venous catheters must be aware of this potentially fatal complication before advancing the catheter too far relative to patient size.

No particular route or catheter type is exempt from this potentially fatal complication<sup>31</sup>. The actual position of the tip of the indwelling catheter should be confirmed by x-ray after insertion<sup>1,3,4,8,14,29,36</sup>. Central venous catheters should be placed in the superior vena cava<sup>1,3,4,8,14,26,29,33,36</sup> above its junction with the right atrium and parallel to the vessel wall and its distal tip positioned at a level above either the azygos vein or the carina of the trachea, whichever is better visualized.

Central venous catheters should not be placed in the right atrium unless specifically required for special relatively short term procedures, such as aspiration of air emboli during neurosurgery. Such procedures are nevertheless risk prone and should be closely monitored and controlled.

**ARROWg<sup>ard</sup>® Antimicrobial Surface:** The Arrow® antimicrobial catheter consists of our standard polyurethane catheter with Blue FlexTip®, plus an exterior antimicrobial surface treatment. Substantial antimicrobial activity associated with the 7 Fr. triple-lumen ARROWg<sup>ard</sup>® catheter has been demonstrated in the following ways:

• Significant antimicrobial activity associated with the Arrow® catheter has been demonstrated using zone of inhibition bioassays against the following organisms:<sup>24</sup>

*Escherichia coli*

*Pseudomonas aeruginosa*

*Staphylococcus epidermidis*

*Staphylococcus aureus*

*Klebsiella pneumoniae*

*Candida albicans*

• Contact inhibition of microbial growth on the surface of the Arrow® catheter has been demonstrated against organisms commonly associated with nosocomial infections; e.g. *Staphylococcus epidermidis* and *Staphylococcus aureus*.<sup>24</sup>

• Antimicrobial activity on the surface of the catheter during handling and placement has been demonstrated *in situ* in limited animal studies.<sup>16</sup>

• The ARROWg<sup>ard</sup>® catheter has demonstrated a significant decrease in the rate of bacterial colonization along the catheter in limited animal studies.<sup>16</sup>

• A prospective randomized clinical trial of 403 catheter insertions in adult patients in a medical-surgical ICU showed that the antimicrobial catheters were 50% less likely to be colonized than control catheters (p=0.003) and 80% less likely to produce catheter related bacteremia (p=0.02).<sup>23</sup>

• Arrow® antimicrobial catheters retained antibacterial activity with zones of inhibition of 4 to 10 mm against *Staphylococcus aureus* and *Escherichia coli* after 10 days of implantation in rats.<sup>24</sup>

• Complete data were obtained for 403 catheters (195 control catheters and 208 antimicrobial catheters) in 158 patients. Control catheters removed from patients who were receiving systemic antibiotic therapy occasionally showed low-level surface activity that was unrelated to the length of time the catheter had been in place (mean zone of inhibition ± SD, 1.7 ± 2.8 mm); in contrast, antimicrobial catheters uniformly showed residual surface activity (mean zone of inhibition, 5.4 ± 2.2 mm; P < 0.002), which declined after prolonged periods *in situ*. Antimicrobial activity was seen with antimicrobial catheters that had been in place for as long as 15 days.<sup>23</sup>

• Arrow® antimicrobial catheters produced large zones of inhibition *in vitro* (range 10 to 18 mm) against the following microbes:

Methicillin-resistant *Staphylococcus aureus*

Gentamicin/methicillin-resistant

*Staphylococcus aureus*

*Staphylococcus aureus*

*Staphylococcus epidermidis*

*Escherichia coli*

*Pseudomonas aeruginosa*

*Klebsiella pneumoniae*

*Candida albicans*

- After 7 days of implantation the catheters retained 6-7 mm zones of inhibition against *Staphylococcus aureus*.<sup>16</sup>
- Antibacterial activity was retained against *Staphylococcus epidermidis* (10<sup>6</sup> bacterial concentration) from subcutaneous segments of ARROWg<sup>ard</sup> antimicrobial surface catheters for at least 120 hours and some up to 520 hours after insertion of the catheters into cardiac surgical patients (both double- and triple-lumen catheters). The zone of inhibition size varied in 7 Fr. triple-lumen catheters from 2.5 to 10 mm at 500 hours.<sup>2</sup>

If the total amount of silver sulfadiazine and chlorhexidine contained in the antimicrobial surface was released from the catheter as a single dose, the blood levels of silver, sulfadiazine, and chlorhexidine that would be found would be less than the blood levels found after clinical usage of these compounds in established safe dosages as administered via mucous membranes and skin.<sup>12</sup>

The potential exposure of patients to the two agents, silver sulfadiazine and chlorhexidine, on the antimicrobial surface is significantly less than that encountered when these compounds are used on burn wounds, on cutaneous wounds, or as mucosal irrigants.<sup>12</sup>

No adverse effects of a toxicologic nature have been associated with the clinical use of this antimicrobial surface in spite of the fact that catheters have been placed in patients sensitive to sulfonamides but who were unaware of their sensitivity.<sup>12</sup> However, the ARROWg<sup>ard</sup> Blue<sup>®</sup> antimicrobial catheter has been reported to cause severe anaphylactic reactions in a limited number of patients in Japan and the UK (first case reported May 1996). Refer to the Contraindications section for additional information.

#### Indications for Use:

The multiple-lumen catheter permits venous access to the central circulation of pediatric patients.

The ARROWg<sup>ard</sup> antimicrobial surface is intended to help provide protection against catheter-related infections. The catheter is not intended to be used as a treatment for existing infections nor as a substitute for a tunneled catheter in those patients requiring long-term therapy. One clinical study indicates that the antimicrobial properties of the catheter may not be effective when it is used to administer TPN.<sup>7</sup>

The use of the 4 Fr. antimicrobial surface catheter in infants weighing less than 2.0 kilograms and the use of the 5 Fr. antimicrobial surface catheter in infants weighing less than 2.5 kilograms have not been evaluated. Therefore, the benefits of these catheters in infants weighing less than the indicated weights should be weighed against any possible risks.

#### Contraindications:

The ARROWg<sup>ard</sup> Blue<sup>®</sup> antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine, and/or sulfa drugs. The ARROWg<sup>ard</sup> Blue<sup>®</sup> antimicrobial catheter has been reported to cause severe anaphylactic reactions in a limited number of patients in Japan and the UK (first case reported May 1996). There have been no reported incidents of hypersensitivity in the United States. Since its introduction in 1990 through October 1999, 20 cases of potential hypersensitivity have been reported worldwide with 17 cases occurring in individuals of Japanese extraction living in Japan. The literature indicates that individuals of Japanese extraction are known to have had similar anaphylactic reactions following topical chlorhexidine administration.<sup>15,17,20,21,27,28,37,39</sup>

**If adverse reactions occur after catheter placement, remove catheter immediately.**

#### Special Patient Populations:

Controlled studies of this product have not been conducted in pregnant women,<sup>25</sup> pediatric or neonatal patients, and patients with known sulfonamide hypersensitivity, erythema multiforme, and Stevens-Johnson syndrome deficiency.<sup>12</sup> The benefits of the use of this catheter should be weighed against any possible risk.

#### Warnings and Precautions:\*

1. **Warning: Sterile, Single use: Do not reuse, reprocess or resterilize. Reuse of device creates a potential risk of serious injury and/or infection which may lead to death.**
2. **Warning: Do not place the catheter into or allow it to remain in the right atrium or right ventricle. Central vein catheters should be positioned so that the distal tip of the catheter is in the lower 1/2 to 1/3 of the superior vena cava (SVC) above the junction of the SVC and the right atrium and lies parallel to the vessel wall. For femoral vein approach, the catheter should be advanced into the vessel so that the catheter tip lies parallel to the vessel wall and does not enter the right atrium.**  
**Due to the variety of pediatric catheter lengths available, patient size must be carefully considered relative to actual length of catheter introduced.**
3. **Warning: Practitioners must be aware of complications associated with central vein catheters including cardiac tamponade secondary to vessel wall, atrial or ventricular perforation, pleural and mediastinal injuries, air embolism, catheter embolism, thoracic duct laceration, bacteremia, septicemia, thrombosis, inadvertent arterial puncture, nerve damage, hematoma, hemorrhage, and dysrhythmias.**
4. **Warning: Do not apply excessive force in removing guide wire or catheters. If withdrawal cannot be easily accomplished, a chest x-ray should be obtained and further consultation requested.**
5. **Warning: The practitioner must be aware of potential air embolism associated with leaving open needles or catheters in central venous puncture sites or as a consequence of inadvertent disconnects. To lessen the risk of disconnects, only securely tightened Luer-Lock connections should be used with this device. Follow hospital protocol to guard against air embolism for all catheter maintenance.**
6. **Warning: Passage of the guide wire into the right heart can cause dysrhythmias, right bundle branch block,<sup>41</sup> and vessel wall, atrial or ventricular perforation.**
7. **Warning: Due to the risk of exposure to HIV (Human Immunodeficiency Virus) or other blood borne pathogens, health care workers should routinely use universal blood and body-fluid precautions in the care of all patients.**
8. **Precaution: Indwelling catheters should be routinely inspected for desired flow rate, security of dressing, correct catheter position, and for secure Luer-Lock connection. Use centimeter markings to identify if the catheter position has changed.**
9. **Precaution: Only x-ray examination of the catheter placement can ensure that the catheter tip has not entered the heart or no longer lies parallel to the vessel wall. If catheter position has changed, immediately perform chest x-ray examination to confirm catheter tip position.**

10. **Precaution:** For blood sampling, temporarily shut off remaining port(s) through which solutions are being infused.

11. **Precaution:** Alcohol and acetone can weaken the structure of polyurethane materials. Check ingredients of prep sprays and swabs for acetone and alcohol content.

**Acetone:** Do not use acetone on catheter surface. Acetone may be applied to skin but must be allowed to dry completely prior to applying dressing.

**Alcohol:** Do not use alcohol to soak catheter surface or to restore catheter patency. Care should be taken when instilling drugs containing high concentration of alcohol. Always allow alcohol to dry completely prior to applying dressing.

12. **Precaution:** Some disinfectants used at the catheter insertion site contain solvents, which can attack the catheter material. Assure insertion site is dry before dressing.

13. **Precaution:** Use of a syringe smaller than 10 mL to irrigate or declot an occluded catheter may cause intraluminal leakage or catheter rupture.<sup>9</sup>

#### A Suggested Procedure:

##### Use sterile technique.

1. **Precaution:** Place patient in slight Trendelenburg position as tolerated to reduce the risk of air embolism. If femoral approach is used, place patient in supine position.
2. Measure the patient's external anatomy to estimate the length of catheter required to place the distal tip in the SVC.
3. Prep and drape puncture site as required.
4. Perform skin wheal with desired needle.
5. Prepare the catheter for insertion by flushing each lumen and clamping or attaching the injection caps to the appropriate extension lines. Leave the distal extension line uncapped for guide wire passage. **Warning:** Do not cut the catheter to alter length.

##### Arrow UserGard® Needle-Free Injection Hub (where provided)

###### Instructions for Use:

- Attach Luer end of UserGard® hub to syringe.
- Prepare injection site with alcohol or betadine per standard hospital protocol.
- Remove red dust cap.
- Press UserGard® hub onto injection site and twist to lock on pin (refer to Fig. 2).


Fig. 2

- Inject or withdraw fluid as required.
- Disengage UserGard® hub from injection site and discard. **Warning:** To minimize the risk of possible air embolism, do not leave UserGard® hub connected to injection site. Single use only.

6. Locate central vein with a needle and syringe appropriate to the size of the patient.

7. Insert introducer catheter/needle with attached syringe into vein along side of locator needle and aspirate. Remove locator needle. Withdraw needle from introducer catheter. If no free flow of venous blood is observed after needle is removed, attach syringe to the catheter and aspirate until good venous blood flow is established. **Precaution:** The color of the blood aspirated is not always a reliable indicator of venous access.<sup>18</sup> Do not reinsert needle into introducer catheter. Because of the potential for inadvertent arterial placement, verify venous access via a wave form obtained by a calibrated pressure transducer (refer to Fig. 3).


Fig. 3

If hemodynamic monitoring equipment is not available to permit transducing a central venous wave form, disconnect the syringe and check for pulsatile flow. Pulsatile flow is usually an indicator of inadvertent arterial puncture.

##### Alternate Technique:

Introducer needle may be used in the standard manner as alternative to catheter/needle assembly.

8. Insert desired tip of spring-wire guide through the introducer needle or catheter into vein. **Precaution:** Due to the fragile nature of the spring-wire guide contained in this product, inspect "J"-tip for damage prior to insertion. If the "J"-tip portion of the spring-wire guide is used, prepare for insertion by sliding the plastic tube over the "J" to straighten. The spring-wire guide should then be advanced in the routine fashion to a depth appropriate to patient size. Advancement of "J"-tip may require a gentle rotating motion. **Warning:** Do not cut spring-wire guide to alter length. Do not withdraw spring-wire guide against needle bevel to minimize the risk of possible severing or damaging of spring-wire guide.
9. Hold spring-wire guide in place and remove introducer needle or catheter. **Precaution:** Maintain firm grip on spring-wire guide at all times.
10. Enlarge cutaneous puncture site with cutting edge of scalpel positioned away from the spring-wire guide. **Precaution:** Do not cut guide wire. Use tissue dilator to enlarge site as required. **Warning:** Do not leave tissue dilator in place as an indwelling catheter to minimize the risk of possible vessel wall perforation.
11. Thread tip of multiple-lumen catheter over spring-wire guide. Sufficient guide wire length must remain exposed at hub end of catheter to maintain a firm grip on guide wire. Grasping near skin, advance catheter into vein with slight twisting motion. **Precaution:** Where provided, catheter clamp and fastener must not be attached to catheter until spring-wire guide is removed.
12. Using centimeter marks on catheter as positioning reference points, advance catheter to final indwelling position. All

centimeter marks are referenced from the catheter tip. Marking symbology is as follows: (1) numerical: 5, 15, 25, etc.; (2) bands: each band denotes 10 cm intervals, with one band indicating 10 cm, two bands indicating 20 cm, etc.; (3) each dot denotes a 1 cm interval.

13. Hold catheter at desired depth and remove spring-wire guide. The Arrow® catheter included in this product has been designed to freely pass over the spring-wire guide. If resistance is encountered when attempting to remove the spring-wire guide after catheter placement, the spring-wire may be kinked about the tip of the catheter within the vessel (refer to Fig. 4).


Fig. 4

In this circumstance, pulling back on the spring-wire guide may result in undue force being applied resulting in spring-wire guide breakage. If resistance is encountered, withdraw the catheter relative to the spring-wire guide about 2-3 cm and attempt to remove the spring-wire guide. If resistance is again encountered, remove the spring-wire guide and catheter simultaneously. **Warning: Although the incidence of spring-wire guide failure is extremely low, practitioner should be aware of the potential for breakage if undue force is applied to the wire.**

14. Verify that the entire spring-wire guide is intact upon removal.
15. Check lumen placement by attaching a syringe to each extension line and aspirate until free flow of venous blood is observed. Connect all extension lines to appropriate Luer-Lock line(s) as required. Unused port(s) may be "locked" through injection cap(s) using standard hospital protocol. Slide clamps are provided on extension lines to occlude flow through each lumen during line and injection cap changes. **Precaution: To minimize the risk of damage to extension lines from excessive pressure, each clamp must be opened prior to infusing through that lumen.**
16. Secure and dress catheter temporarily.
17. Verify catheter tip position by chest x-ray immediately after placement. **Precaution: X-ray exam must show the catheter located in the right side of the mediastinum in the SVC with the distal end of the catheter parallel to the vena cava wall and its distal tip positioned at a level above either the azygos vein or the carina of the trachea, whichever is better visualized.** If catheter tip is malpositioned, reposition and reverify.
18. Secure catheter to patient. Use triangular juncture hub with integral suture ring and side wings as primary suture site. In kits where provided, the catheter clamp and fastener should be utilized as a secondary suture site as necessary. **Precaution: Do not suture directly to the outside diameter of the catheter to minimize the risk of cutting or damaging the catheter or impeding catheter flow.**

**Catheter Clamp and Fastener Instructions:**

- After spring-wire guide has been removed and the necessary lines have been connected or locked, spread wings of rubber

clamp and position on catheter as required to ensure proper tip location (refer to Fig. 5).


Fig. 5

- Snap rigid fastener onto catheter clamp (refer to Fig. 6).


Fig. 6

- Secure catheter to patient by suturing the catheter clamp and fastener together to the skin, using side wings to prevent catheter migration (refer to Fig. 7).


Fig. 7

19. Dress puncture site per hospital protocol. **Precaution: Maintain the insertion site with regular meticulous redressing using aseptic technique.**
20. Record on the patient's chart the indwelling catheter length as to centimeter markings on catheter where it enters the skin.

Frequent visual reassessment should be made to ensure that the catheter has not moved.

#### Catheter Exchange Procedure:

1. Use sterile technique.
2. **Precaution: Prior to attempting a catheter exchange procedure, remove the catheter clamp and fastener.**
3. Proceed per hospital protocol. Cutting the catheter is not recommended due to the potential for catheter embolism.

#### Catheter Removal Procedure:

1. **Precaution: Place the patient in a supine position.**
2. Remove dressing. **Precaution: To minimize the risk of cutting the catheter do not use scissors to remove the dressing.**
3. **Warning: Exposure of the central vein to atmospheric pressure may result in entry of air into the central venous system.** Remove suture(s) from primary suture site. Be careful not to cut the catheter. Remove catheter slowly, pulling it parallel to the skin. As catheter exits the site, apply pressure with a dressing impermeable to air, e.g. VASELINE® gauze. Because the residual catheter track remains an air entry point until completely sealed, the occlusive dressing should remain in place for at least 24-72 hours dependent upon the amount of time the catheter was indwelling.<sup>19,30,32,38</sup>
4. Upon removal of the catheter, inspect it to make sure that the entire length has been withdrawn.
5. Document removal procedure.

#### References:

1. Agarwal KC, Khan MAA, Falla A, Amato JJ. Cardiac perforation from central venous catheters: survival after cardiac tamponade in an infant. *Pediatrics*. March 1984;73:333-338.
2. Bach A, Schmidt H, Bottinger B, et al. Retention of antibacterial activity and bacterial colonization of antiseptic-bonded central venous catheters. *J Antimicrob Chemother*. 1996;37:315-322.
3. Bar-Joseph G, Galvis AG. Perforation of the heart by central venous catheters in infants: guidelines to diagnosis and management. *J Pediatr Surg*. June 1983;18:284-287.
4. Brandt RL, Foley WJ, Fink GH, Regan WJ. Mechanism of perforation of the heart with production of hydropericardium by a venous catheter and its prevention. *Am J Surg*. 1970;119:311-316.
5. Carey BE. Major complications of central lines in neonates. *Neonatal Network*. June 1989;7:17-28.
6. Casado-Flores J, Valdivielso-Serna A, Perez-Jurado L, et al. Subclavian vein catheterization in critically ill children: analysis of 322 cannulations. *Intensive Care Med*. 1991;17:350-354.
7. Ciresi DL, Albrecht RM, Volkens PA, Scholten DJ. Failure of antiseptic bonding to prevent central venous catheter-related infection and sepsis. *Am Surgeon*. 1996;62:641-646.
8. Collier PE, Ryan JJ, Diamond DL. Cardiac tamponade from central venous catheters – report of a case and review of the English literature. *Angiology*. September 1984;35:595-600.
9. Conn C. The importance of syringe size when using implanted vascular access devices. *J Vasc Access Nurs*. Winter 1993;3:11-18.
10. Daniels SR, Hannon DW, Meyer RA, Kaplan S. Paroxysmal supraventricular tachycardia. *AJDC*. 1984;138:474-475.
11. Eissa NT, Kvetan V. Guide wire as a cause of complete heart block in patients with preexisting left bundle branch block. *Anesthesiology*. 1990;73:772-774.

12. Farber T. ARROWgard™ antiseptic surface – toxicology review. Monograph. Published by Arrow International, Inc. April, 1992.
13. Franciosi RA, Ellefson RD, Uden D, Drake RM. Sudden unexpected death during central hyperalimentation. *Pediatrics*. March 1982;69:305-307.
14. Fuhrman BP. Perforation by vascular catheters. *Perspectives in Critical Care*. 1988;1:147-154.
15. Fukui A, Ohsumi A, Takaori M. A case of anaphylactic shock induced by chlorhexidine glucuronate. *J Japan Society Clin Anesthesia*. 1989;9:356-360.
16. Greenfield JI, Sampath L, Popilskis SJ, Brunnert SR, Stylianos S, Modak S. Decreased bacterial adherence and biofilm formation on chlorhexidine and silver sulfadiazine-impregnated central venous catheters implanted in swine. *Crit Care Med*. 1995;23:894-900.
17. Harukuni I, Ishizawa Y, Nishikawa T, Takeshima R, Dohi S, Naito H. Anaphylactic shock with ventricular fibrillation induced by chlorhexidine. *Japanese J Anesthesiology*. 1992;41:455-459.
18. Jobs DR, Schwartz AJ, Greenhow DE, Stephenson LW, Ellison N. Safer jugular vein cannulation: recognition of arterial punctures and preferential use of the external jugular route. *Anesthesiology*. 1983;59:353-355.
19. Kashuk JL, Penn I. Air embolism after central venous catheterization. *Surg Gynecol Obstet*. September 1984;159:249-252.
20. Kubo H, Akiyama Y, Honda K, Nakajo N. Anaphylaxis following oral irrigation with chlorhexidine gluconate. *J Japanese Dental Society Anesthesiology*. 1985;13:659-663.
21. Layton GT, Stanworth DR, Amos HE. The incidence of IgE and IgG antibodies to chlorhexidine. *Clin Experimental Allergy*. 1989;19:307-314.
22. Leibovitz E, Ashkenazi A, Levin S, Nissim F. Fatal pericardial tamponade complicating total parenteral nutrition via a silastic central vein catheter. *J Pediatr Gastroenterol Nutr*. 1988;7:306-307.
23. Maki DG, Wheeler SJ, Stolz SM, Mermel LA. Prevention of central venous catheter-related bloodstream infection by use of an antiseptic-impregnated catheter. *Ann Int Med*. August 15, 1997;127:257-266.
24. Modak SM, Sampath L. Development and evaluation of a new polyurethane central venous antiseptic catheter: reducing central venous catheter infections. *Infections in Medicine*. June 1992;23-29.
25. Modak SM. (Written communication, June 1991).
26. Mupanemunda RH, Mackanjee HR. A life-threatening complication of percutaneous central venous catheters in neonates. *AJDC*. 1992;146:1414-1415. Letter.
27. Okano M, Nomura M, Hata S, et al. Anaphylactic symptoms due to chlorhexidine gluconate. *Arch Dermatol*. 1989;125:50-52.
28. Okano M, Nomura M, Okada N, Sato K, Tashiro M. Four cases presenting anaphylactic reactions due to topical application of Hibitane®. *Skin Research*. 1983;25:587-592.
29. Otto CW. Central venous pressure monitoring. In: Blitt CD, ed. *Monitoring in Anesthesia and Critical Care Medicine*. New York, NY: Churchill Livingstone; 1985:139.
30. Paskin DL, Hoffman WS, Tuddenham WJ. A new complication of subclavian vein catheterization. *Ann Surg*. March 1974;179:266-268.
31. Peters JL, ed. *A Manual of Central Venous Catheterization and Parenteral Nutrition*. Boston, MA: John Wright PSG; 1983:58-61, 155-157.

32. Phifer TJ, Bridges M, Conrad SA. The residual central venous catheter track – an occult source of lethal air embolism: case report. *J Trauma*. 1991;31:1558-1560.
33. Rogers BB, Berns SD, Maynard EC, Hansen TWR. Pericardial tamponade secondary to central venous catheterization and hyperalimentation in a very low birth weight infant. *Pediatr Pathology*. 1990;10:819-823.
34. Ross P Jr, Seashore JH. Bilateral hydrothorax complicating central venous catheterization in a child: case report. *J Pediatr Surg*. 1989;24:263-264.
35. Sheep RE, Guiney WB Jr. Fatal cardiac tamponade. *JAMA*. 1982;248:1632-1635.
36. Sullivan CA, Konefal SH Jr. Cardiac tamponade in a newborn: a complication of hyperalimentation. *J Parenteral Enteral Nutrition*. May/June 1987;11:319-321.
37. Takeda K, Inoue K, Matsuya T, et al. An allergic shock possibly induced by the chlorhexidine: report of a case. *J Osaka Univ Dent Soc*. 1985;30:221-225.
38. Thielen JB, Nyquist J. Subclavian catheter removal. *J Intravenous Nurs*. March/April 1991;14:114-118.
39. Tsuneto S, Watanabe S, Koyama K, Nakayama K, Saito H, Saito K. Anaphylax shock induced by chlorhexidine mixed in the vial of lidocaine. *J Japan Society Clin Anesthesia*. 1987;7:272-277.

Arrow International, Inc. recommends that the user be acquainted with the reference literature.

\*If you have any questions or would like additional reference information, please contact Arrow International, Inc.

†A registered trademark of Unilever Supply Chain, Inc.

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Arrow International, Inc. kullanıcının referans literatüre aşına olmasını önerir.
\*Sorularınız varsa veya ek referans bilgisi isterseniz lütfen Arrow International, Inc. ile irtibat kurun.