Inspected Dimensions:
Folded Length: 7-1/2" (19 cm)
Folded Width: 5-1/2" (14 cm)

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Multiple-Lumen Central Venous Catheterization Product with ARROWg^{ard}® Antimicrobial Surface

Safety and Efficacy Considerations:

Do not use if package has been previously opened or damaged. **Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.**

Do not alter the catheter, spring-wire guide, or any other kit/ set component during insertion, use, or removal.

Procedure must be performed by trained personnel well versed in anatomical landmarks, safe technique, and potential complications.

Warning: Do not place the catheter into or allow it to remain in the right atrium or right ventricle (refer to Fig. 1).

No particular route or catheter type is exempt from this potentially fatal complication.²⁸ The actual position of the tip of the indwelling catheter should be confirmed by x-ray after insertion.^{2,3,7,22,28,31} Central venous catheters should be placed in the superior vena cava^{2,3,4,7,9,22,34} above its junction with the right atrium and parallel to the vessel wall^{15,34} and its distal tip positioned at a level above either the azygos vein or the carina of the trachea, whichever is better visualized.

Central venous catheters should not be placed in the right atrium unless specifically required for special relatively short term procedures such as aspiration of air emboli during neurosurgery. Such procedures are, nevertheless, risk prone and should be closely monitored and controlled.

ARROWg^{ard}® Antimicrobial Surface: The Arrow antimicrobial catheter consists of our standard polyurethane catheter with Blue FlexTip®, plus an exterior antimicrobial surface treatment. Substantial antimicrobial activity associated with the 7 Fr. triple-lumen ARROWg^{ard}® catheter has been demonstrated in the following ways:

- Significant antimicrobial activity associated with the Arrow catheter has been demonstrated using zone of inhibition bioassays against the following organisms²³:

Escherichia coli
Pseudomonas aeruginosa
Staphylococcus epidermidis
Staphylococcus aureus
Klebsiella pneumoniae
Candida albicans

- Contact inhibition of microbial growth on the surface of the Arrow catheter has been demonstrated against organisms commonly associated with nosocomial infections; e.g. *Staphylococcus epidermidis* and *Staphylococcus aureus*.²³
- Antimicrobial activity on the surface of the catheter during handling and placement has been demonstrated *in situ* in limited animal studies.¹³
- The ARROWg^{ard}® catheter has demonstrated a significant decrease in the rate of bacterial colonization along the catheter in limited animal studies.¹³
- A prospective randomized clinical trial of 403 catheter insertions in adult patients in a medical-surgical ICU showed that the antimicrobial catheters were 50% less likely to be colonized than control catheters (p=0.003) and 80% less likely to produce catheter related bacteremia (p=0.02).²¹
- Arrow antimicrobial catheters retained antibacterial activity with zones of inhibition of 4 to 10 mm against *Staphylococcus aureus* and *Escherichia coli* after 10 days of implantation in rats.²³
- Complete data were obtained for 403 catheters (195 control catheters and 208 antimicrobial catheters) in 158 patients. Control catheters removed from patients who were receiving systemic antibiotic therapy occasionally showed low-level


Fig. 1

Cardiac Tamponade: It has been documented by many authors that placement of indwelling catheters in the right atrium is a dangerous practice^{2,3,4,7,9,22,28} that may lead to cardiac perforation and tamponade.^{2,3,4,7,22,28} Although cardiac tamponade secondary to pericardial effusion is uncommon, there is a high mortality rate associated with it.³⁰ Practitioners placing central venous catheters must be aware of this potentially fatal complication before advancing the catheter too far relative to patient size.

surface activity that was unrelated to the length of time the catheter had been in place (mean zone of inhibition \pm SD, 1.7 ± 2.8 mm); in contrast, antimicrobial catheters uniformly showed residual surface activity (mean zone of inhibition, 5.4 ± 2.2 mm; $P < 0.002$), which declined after prolonged periods *in situ*. Antimicrobial activity was seen with antimicrobial catheters that had been in place for as long as 15 days.²¹

- Arrow antimicrobial catheters produced large zones of inhibition *in vitro* (range 10 to 18 mm) against the following microbes:

Methicillin-resistant
Staphylococcus aureus
Gentamicin/methicillin-resistant
Staphylococcus aureus
Staphylococcus aureus
Staphylococcus epidermidis
Escherichia coli
Pseudomonas aeruginosa
Klebsiella pneumoniae
Candida albicans

After 7 days of implantation the catheters retained 6-7 mm zones of inhibition against *Staphylococcus aureus*.¹³

- Antibacterial activity was retained against *Staphylococcus epidermidis* (10^6 bacterial concentration) from subcutaneous segments of ARROWg^{ard}® antimicrobial surface catheters for at least 120 hours and some up to 520 hours after insertion of the catheters into cardiac surgical patients (both double- and triple-lumen catheters). The zone of inhibition size varied in 7 Fr. triple-lumen catheters from 2.5 to 10 mm at 500 hours.¹

If the total amount of silver sulfadiazine and chlorhexidine contained in the antimicrobial surface was released from the catheter as a single dose, the blood levels of silver, sulfadiazine, and chlorhexidine that would be found would be less than the blood levels found after clinical usage of these compounds in established safe dosages as administered via mucous membranes and skin.¹¹

The potential exposure of patients to the two agents, silver sulfadiazine and chlorhexidine, on the antimicrobial surface is significantly less than that encountered when these compounds are used on burn wounds, on cutaneous wounds, or as mucosal irrigants.¹¹

No adverse effects of a toxicologic nature have been associated with the clinical use of this antimicrobial surface in spite of the fact that catheters have been placed in patients sensitive to sulfonamides but who were unaware of their sensitivity.¹¹ However, the ARROWg^{ard} Blue® antimicrobial catheter has been reported to cause severe anaphylactic reactions in a limited number of patients in Japan and the UK (first case reported May 1996). Refer to the Contraindications section for additional information.

Indications for Use:

The Arrow® Multiple-Lumen Catheter permits venous access to the central circulation. The ARROWg^{ard}® antimicrobial surface is intended to provide protection against catheter-related infections.

The catheter is not intended to be used as a treatment for existing infections nor as a substitute for a tunneled catheter in those patients requiring long-term therapy. One clinical study indicates that the antimicrobial properties of the catheter may not be effective when it is used to administer TPN.⁶

Contraindications:

The ARROWg^{ard} Blue® antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine acetate, silver sulfadiazine, and/or sulfa drugs. The ARROWg^{ard} Blue® antimicrobial catheter has been reported to cause severe anaphylactic reactions in a limited number of patients in Japan and the UK (first case reported May 1996). There have been no reported incidents of hypersensitivity in the United States. Since its introduction in 1990 through October 1999, 20 cases of potential hypersensitivity have been reported worldwide with 17 cases occurring in individuals of Japanese extraction living in Japan. The literature indicates that individuals of Japanese extraction are known to have had similar anaphylactic reactions following topical chlorhexidine administration.^{12,14,19,20,25,26,32,35}

If adverse reactions occur after catheter placement, remove catheter immediately.

Special Patient Populations:

Controlled studies of this product have not been conducted in pregnant women,²⁴ pediatric or neonatal patients, and patients with known sulfonamide hypersensitivity, erythema multiforme, and Stevens-Johnson syndrome deficiency.¹¹ The benefits of the use of this catheter should be weighed against any possible risk.

Warnings and Precautions:*

1. **Warning: Sterile, Single use: Do not reuse, reprocess or resterilize. Reuse of device creates a potential risk of serious injury and/or infection which may lead to death.**
2. **Warning: Do not place the catheter into or allow it to remain in the right atrium or right ventricle. Central vein catheters should be positioned so that the distal tip of the catheter is in the superior vena cava (SVC) above the junction of the SVC and the right atrium and lies parallel to the vessel wall. For femoral vein approach, the catheter should be advanced into the vessel so that the catheter tip lies parallel to the vessel wall and does not enter the right atrium.**
3. **Warning: Practitioners must be aware of complications associated with central vein catheters including cardiac tamponade secondary to vessel wall, atrial or ventricular perforation, pleural and mediastinal injuries, air embolism, catheter embolism, catheter occlusion, thoracic duct laceration, bacteremia, septicemia, thrombosis, inadvertent arterial puncture, nerve damage, hematoma, hemorrhage, and dysrhythmias.**
4. **Warning: Do not apply excessive force in removing guide wire or catheters. If withdrawal cannot be easily accomplished, a chest x-ray should be obtained and further consultation requested.**
5. **Warning: The practitioner must be aware of potential air embolism associated with leaving open needles or catheters in central venous puncture sites or as a consequence of inadvertent disconnects. To lessen the risk of disconnects, only securely tightened Luer-Lock connections should be used with this device. Follow hospital protocol to guard against air embolism for all catheter maintenance.**
6. **Warning: Passage of the guide wire into the right heart can cause dysrhythmias, right bundle branch block,¹⁰ and a perforation of the vessel wall, atrial or ventricular.**
7. **Warning: Practitioners must be aware of the potential for entrapment of guide wire by any implanted device in the circulatory system (i.e. vena cava filters, stents). Review patient's history before catheterization procedure to assess for possible implants. Care should be taken regarding the length of spring-wire guide inserted. It is recommended**

that if patient has a circulatory system implant, catheter procedure be done under direct visualization to minimize the risk of guidewire entrapment.¹

8. **Warning:** Due to the risk of exposure to HIV (Human Immunodeficiency Virus) or other blood borne pathogens, health care workers should routinely use universal blood and body-fluid precautions in the care of all patients.
9. **Precaution:** Indwelling catheters should be routinely inspected for desired flow rate, security of dressing, correct catheter position and for secure Luer-Lock connection. Use centimeter markings to identify if the catheter position has changed.
10. **Precaution:** Only x-ray examination of the catheter placement can ensure that the catheter tip has not entered the heart or no longer lies parallel to the vessel wall. If catheter position has changed, immediately perform chest x-ray examination to confirm catheter tip position.
11. **Precaution:** For blood sampling, temporarily shut off remaining port(s) through which solutions are being infused.
12. **Precaution:** Alcohol and acetone can weaken the structure of polyurethane materials. Check ingredients of prep sprays and swabs for acetone and alcohol content. **Acetone:** Do not use acetone on catheter surface. Acetone may be applied to skin but must be allowed to dry completely prior to applying dressing. **Alcohol:** Do not use alcohol to soak catheter surface or to restore catheter patency. Care should be taken when instilling drugs containing high concentration of alcohol. Always allow alcohol to dry completely prior to applying dressing.
13. **Precaution:** Some disinfectants used at the catheter insertion site contain solvents, which can attack the catheter material. Assure insertion site is dry before dressing.
14. **Precaution:** Use of a syringe smaller than 10 mL to irrigate or declot an occluded catheter may cause intraluminal leakage or catheter rupture.⁸

A Suggested Procedure: Use sterile technique.

1. **Precaution:** Place patient in slight Trendelenburg position as tolerated to reduce the risk of air embolism. If femoral approach is used, place patient in supine position.
2. Prep and drape puncture site as required.
3. Perform skin wheal with desired needle (25 Ga. or 22 Ga. needle). A SharpsAway® disposal cup is provided for the disposal of needles. Push needles into foam after use. Discard entire cup at completion of procedure. **Precaution: Do not reuse needles after they have been placed into the disposal cup. Particulate matter may adhere to needle tip.**
4. Prepare catheter for insertion by flushing each lumen and clamping or attaching injection caps to appropriate extension lines. Leave distal extension line uncapped for guide wire passage. **Warning: Do not cut catheter to alter length.**
5. Locate central vein with a 22 Ga. needle and syringe.
6. Insert introducer catheter/needle with attached syringe into vein along side of locator needle and aspirate. Remove locator needle. Withdraw needle from introducer catheter. If no free flow of venous blood is observed after needle is removed, attach syringe to catheter and aspirate until good venous blood flow is established. **Precaution: The color of the**

blood aspirated is not always a reliable indicator of venous access.¹⁶ Do not reinsert needle into introducer catheter. Because of the potential for inadvertent arterial placement, verify venous access via a wave form obtained by a calibrated pressure transducer (refer to Fig. 2).


Fig. 2

If hemodynamic monitoring equipment is not available to permit transducing a central venous wave form, disconnect the syringe and check for pulsatile flow. Pulsatile flow is usually an indicator of inadvertent arterial puncture.

Alternate Technique:

Introducer needle may be used in the standard manner as alternative to catheter/needle assembly.

7. Using the two-piece Arrow Advancer™, advance spring-wire guide through guide wire introducer needle or catheter into vein. Advance spring-wire guide to required depth. Advancement of “J” tip may require a gentle rotating motion. **Warning: Do not cut spring-wire guide to alter length. Do not withdraw spring-wire guide against needle bevel to minimize the risk of possible severing or damaging of spring-wire guide.**

Two-Piece Arrow Advancer™ Instructions:

- Using your thumb, straighten the “J” by retracting spring-wire guide into the Arrow Advancer™ (refer to Figs. 3, 4).


Fig. 3


Fig. 4

When tip is straightened, spring-wire guide is ready for insertion. Centimeter marks on guide wire are referenced from “J” end. One band indicates 10 cm, two bands 20 cm, and three bands 30 cm.

Alternate Technique:

If a simple straightening tube is preferred, the straightening tube portion of the Arrow Advancer™ can be disconnected from the unit and used separately.

Separate the Arrow Advancer™ tip or straightening tube from blue Arrow Advancer™ unit. If “J” tip portion of spring-wire guide is used, prepare for insertion by sliding the plastic tube over the “J” to straighten. The spring-wire guide should then be advanced in routine fashion to desired depth.

8. Hold spring-wire guide in place and remove introducer needle or catheter. **Precaution: Maintain firm grip on spring-wire guide at all times.** Use centimeter markings on spring-wire guide to adjust indwelling length according to desired depth of indwelling catheter placement.
9. Enlarge cutaneous puncture site with cutting edge of scalpel positioned away from the spring-wire guide. **Precaution: Do not cut guide wire.** Use tissue dilator to enlarge site as required. **Warning: Do not leave tissue dilator in place as an indwelling catheter to minimize the risk of possible vessel wall perforation.**
10. Thread tip of multiple-lumen catheter over spring-wire guide. Sufficient guide wire length must remain exposed at hub end of catheter to maintain a firm grip on guide wire. Grasping near skin, advance catheter into vein with slight twisting motion. **Precaution: If provided, catheter clamp and fastener must not be attached to catheter until spring-wire guide is removed.**
11. Using centimeter marks on catheter as positioning reference points, advance catheter to final indwelling position. All centimeter marks are referenced from the catheter tip. Marking symbology is as follows: (1) numerical: 5, 15, 25, etc.; (2) bands: each band denotes 10 cm intervals, with one band indicating 10 cm, two bands indicating 20 cm, etc.; (3) each dot denotes a 1 cm interval.
12. Hold catheter at desired depth and remove spring-wire guide. The Arrow® catheter included in this product has been designed to freely pass over the spring-wire guide. If resistance is encountered when attempting to remove spring-wire guide after catheter placement, the spring-wire may be kinked about the tip of the catheter within the vessel (refer to Fig. 5).


Fig. 5

In this circumstance, pulling back on the spring-wire guide may result in undue force being applied resulting in spring-wire guide breakage. If resistance is encountered, withdraw catheter relative to spring-wire guide about 2-3 cm and attempt to remove spring-wire guide. If resistance is again encountered remove spring-wire guide and catheter simultaneously.

Warning: Although the incidence of spring-wire guide failure is extremely low, practitioner should be aware of the potential for breakage if undue force is applied to the wire.

13. Verify that entire spring-wire guide is intact upon removal.
14. Check lumen placement by attaching a syringe to each extension line and aspirate until free flow of venous blood is observed. Connect all extension lines to appropriate Luer-Lock line(s) as required. Unused port(s) may be “locked” through injection cap(s) using standard hospital protocol. Slide clamps are provided on extension lines to occlude flow through each lumen during line and injection cap changes. **Precaution: To minimize the risk of damage to extension lines from excessive pressure, each clamp must be opened prior to infusing through that lumen.**
15. Secure and dress catheter temporarily.
16. Verify catheter tip position by chest x-ray immediately after placement. **Precaution: X-ray exam must show the catheter located in the right side of the mediastinum in the SVC with the distal end of the catheter parallel to the vena cava wall and its distal tip positioned at a level above either the azygos vein or the carina of the trachea, whichever is better visualized.** If catheter tip is malpositioned, reposition and re-verify.
17. Secure catheter to patient. Use triangular juncture hub with integral suture ring and side wings as primary suture site. In kits where provided, the catheter clamp and fastener should be utilized as a secondary suture site as necessary. **Precaution: Do not suture directly to outside diameter of catheter to minimize the risk of cutting or damaging the catheter or impeding catheter flow.**

Catheter Clamp and Fastener (where provided)

Instructions for Use:

- After spring-wire guide has been removed and the necessary lines have been connected or locked, spread wings of rubber clamp and position on catheter, as required, to ensure proper tip location (refer to Fig. 6).


Fig. 6

- Snap rigid fastener onto catheter clamp (refer to Fig. 7).


Fig. 7

- Secure catheter to patient by suturing catheter clamp and fastener together to skin, using side wings to minimize the risk of catheter migration (refer to Fig. 8).


Fig. 8

18. Dress puncture site per hospital protocol. **Precaution: Maintain the insertion site with regular meticulous redressing using aseptic technique.**
19. Record on patient's chart the indwelling catheter length as to centimeter markings on catheter where it enters the skin. Frequent visual reassessment should be made to ensure that the catheter has not moved.

Catheter Exchange Procedure:

1. Use sterile technique.
2. **Precaution: Prior to attempting a catheter exchange procedure, remove catheter clamp and fastener (where provided).**
3. Proceed per hospital protocol. Cutting the catheter is not recommended due to the potential for catheter embolism.

Catheter Removal Procedure:

1. **Precaution: Place patient in a supine position.**
2. Remove dressing. **Precaution: To minimize the risk of cutting catheter, do not use scissors to remove dressing.**

3. **Warning: Exposure of the central vein to atmospheric pressure may result in entry of air into the central venous system.** Remove suture(s) from catheter clamp (where provided) and primary suture site. Be careful not to cut the catheter. Remove catheter slowly, pulling it parallel to the skin. As catheter exits the site, apply pressure with a dressing impermeable to air, e.g. Vaseline[®] gauze. Because the residual catheter track remains an air entry point until completely sealed, the occlusive dressing should remain in place for at least 24-72 hours dependent upon the amount of time the catheter was indwelling.^{17,27,29,33}
4. Upon removal of catheter, inspect it to make sure that entire length has been withdrawn.
5. Document removal procedure.

References:

1. Bach A, Schmidt H, Bottinger B, et al. Retention of antibacterial activity and bacterial colonization of antiseptic-bonded central venous catheters. *J Antimicrob Chemo.* 1996;37:315-322.
2. Bar-Joseph G, Galvis AG. Perforation of the heart by central venous catheters in infants: guidelines to diagnosis and management. *J Pediatr Surg.* 1983;18:284-287.
3. Blitt CD, ed. *Monitoring in Anesthesia and Critical Care Medicine.* Central venous pressure monitoring. New York, NY: Churchill Livingstone; 1985:121-165.
4. Brandt RL, Foley WJ, Fink GH, Regan WJ. Mechanism of perforation of the heart with production of hydropericardium by a venous catheter and its prevention. *Am J Surg.* 1970;119:311-316.
5. Carbone K, Gimenez LF, Rogers WH, Watson AJ. Hemothorax due to vena caval erosion by a subclavian dual-lumen dialysis catheter. *South Med J.* 1987;80:795-796.
6. Ciresi DL, Albrecht RM, Volkens PA, Scholten DJ. Failure of antiseptic bonding to prevent central venous catheter-related infection and sepsis. *Am Surgeon.* 1996;62:641-646.
7. Collier PE, Ryan JJ, Diamond DL. Cardiac tamponade from central venous catheters – report of a case and review of the English literature. *Angiology.* September 1984;35:595-600.
8. Conn C. The importance of syringe size when using implanted vascular access devices. *J Vasc Access Nurs.* Winter 1993;3:11-18.
9. Curelaru J, Linder LE, Gustavsson B. Displacement of catheters inserted through internal jugular veins with neck flexion and extension. *Intens Care Med.* 1980;6:179-183.
10. Eissa NT, Kvetan V. Guide wire as a cause of complete heart block in patients with preexisting left bundle branch block. *Anesthesiology.* 1990;73:772-774.
11. Farber T. ARROWgard™ antiseptic surface – toxicology review. Monograph. Published by Arrow International, Inc. April 1992.
12. Fukui A, Ohsumi A, Takaori M. A case of anaphylactic shock induced by chlorhexidine glucuronate. *J Japan Society Clin Anesthesia.* 1989;9:356-360.
13. Greenfeld JI, Sampath L, Popilskis SJ, Brunnert SR, Stylianos S, Modak S. Decreased bacterial adherence and biofilm formation on chlorhexidine and silver sulfadiazine-impregnated central venous catheters implanted in swine. *Crit Care Med.* 1995;23:894-900.
14. Harukuni I, Ishizawa Y, Nishikawa T, Takeshima R, Dohi S, Naito H. Anaphylactic shock with ventricular fibrillation induced by chlorhexidine. *Japanese J Anesthesiology.* 1992;41:455-459.
15. Iberti TJ, Katz LB, Reiner MA, Brownie T, Kwun KB. Hydrothorax as a late complication of central venous indwelling catheters. *Surgery.* November 1983:842-846.

16. Jobs DR, Schwartz AJ, Greenhow DE, Stephenson LW, Ellison N. Safer jugular vein cannulation: recognition of arterial punctures and preferential use of the external jugular route. *Anesthesiology*. 1983;59:353-355.
17. Kashuk JL, Penn I. Air embolism after central venous catheterization. *Surg Gynecol Obstet*. September 1984;159:249-252.
18. Kozeny GA, Bansal VK, Vertuno LL, Hano JE. Contralateral hemothorax secondary to chronic subclavian dialysis catheter. *Am J Nephrol*. 1984;4:312-314.
19. Kubo H, Akiyama Y, Honda K, Nakajo N. Anaphylaxis following oral irrigation with chlorhexidine gluconate. *J Japanese Dental Society Anesthesiology*. 1985;13:659-663.
20. Layton GT, Stanworth DR, Amos HE. The incidence of IgE and IgG antibodies to chlorhexidine. *Clin Experimental Allergy*. 1989;19:307-314.
21. Maki DG, Wheeler SJ, Stolz SM, Mermel LA. Prevention of central venous catheter-related bloodstream infection by use of an antiseptic-impregnated catheter. *Ann Int Med*. August 15, 1997;127:257-266.
22. Maschke SP, Rogove HJ. Cardiac tamponade associated with a multilumen central venous catheter. *Crit Care Med*. 1984;12:611-612.
23. Modak SM, Sampath L. Development and evaluation of a new polyurethane central venous antiseptic catheter: reducing central venous catheter infections. *Infections in Medicine*. June 1992;23-29.
24. Modak SM. (Written communication, June 1991).
25. Okano M, Nomura M, Hata S, et al. Anaphylactic symptoms due to chlorhexidine gluconate. *Arch Dermatol*. 1989;125:50-52.
26. Okano M, Nomura M, Okada N, Sato K, Tashiro M. Four cases presenting anaphylactic reactions due to topical application of Hibitane®. *Skin Research*. 1983;25:587-592.
27. Paskin DL, Hoffman WS, Tuddenham WJ. A new complication of subclavian vein catheterization. *Ann Surg*. March 1974;179:266-268.
28. Peters JL, ed. *A Manual of Central Venous Catheterization and Parenteral Nutrition*. Boston, MA: John Wright PSG; 1983:58-61, 155-157.
29. Phifer TJ, Bridges M, Conrad SA. The residual central venous catheter track – an occult source of lethal air embolism: case report. *J Trauma*. 1991;31:1558-1560.
30. Sheep RE, Guiney WB Jr. Fatal cardiac tamponade. *JAMA*. 1982;248:1632-1635.
31. Sigurdsson J, Riba P, Sigurdsson S. The wandering central venous catheter. *Intensive Care Med*. 1985;11:263-264.
32. Takeda K, Inoue K, Matsuya T, et al. An allergic shock possibly induced by the chlorhexidine: report of a case. *J Osaka Univ Dent Soc*. 1985;30:221-225.
33. Thielen JB, Nyquist J. Subclavian catheter removal. *J Intravenous Nurs*. March/April 1991;14:114-118.
34. Tocino IM, Watanabe A. Impending catheter perforation of superior vena cava: radiographic recognition. *Am J Roentgenology*. March 1986;146:487-490.
35. Tsuneto S, Watanabe S, Koyama K, Nakayama K, Saito H, Saito K. Anaphylaxis shock induced by chlorhexidine mixed in the vial of lidocaine. *J Japan Society Clin Anesthesia*. 1987;7:272-277.
36. Vaziri ND, Maksy M, Lewis M, Martin D, Edwards K. Massive mediastinal hematoma caused by a double-lumen subclavian catheter. *Artif Organs*. 1984;8:223-226.
37. Wanscher M, Frifelt JJ, Smith-Sivertsen C, et al. Thrombosis caused by polyurethane double-lumen subclavian superior vena cava catheter and hemodialysis. *Crit Care Med*. 1988;16:624-628.
- Arrow International, Inc. recommends that the user be acquainted with the reference literature.
- * If you have any questions or would like additional reference information, please contact Arrow International, Inc.
- † A registered trademark of Unilever Supply Chain, Inc.

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Arrow International, Inc. kullanıcının referans literatürüne aşina olmasını tavsiye eder.
* Herhangi bir sorunuz olursa veya ek referans bilgisi isterseniz lütfen Arrow International, Inc. ile irtibata geçin.
Arrow International Inc.