# Pressure Injectable CVC Product

Venous Access | Critical Care

EN

LT

NO

PL

PT

RO

RU

SR

SK

SL

ES

SV

TR

UK

ARROW®

# Table of Contents

|   |   |  |   |
|---|---|--|---|
| <i>Product Description</i>                              | 1 | <i>Catheter Insertion Instructions</i>     | 3 |
| Indications / Contraindications . . . . .               | 1 | Preparing for Catheter Insertion . . . . . | 3 |
| <i>Central Venous Catheter Warnings and Precautions</i> | 1 | Gain Initial Venous Access . . . . .       | 3 |
| General . . . . .                                       | 1 | Verify Venous Access . . . . .             | 4 |
| Catheter . . . . .                                      | 2 | Guidewire / SWG Insertion . . . . .        | 4 |
| Guidewire / SWG . . . . .                               | 2 | Advance Catheter . . . . .                 | 4 |
| Tissue Dilator . . . . .                                | 2 | Complete Catheter Insertion . . . . .      | 4 |
| Pressure Injection . . . . .                            | 2 |  |   |
| Possible Complications . . . . .                        | 3 |  |   |
| <i>Accessory Component Instructions</i>                 | 3 | <i>Pressure Injection Instructions</i>     | 5 |
| Arrow Advancer™ . . . . .                               | 3 |  |   |
| Arrow® Raulerson Syringe . . . . .                      | 3 |  |   |
| SharpsAway II™ Locking Disposal Cup . . . . .           | 3 | <i>Catheter Removal Instructions</i>       | 5 |

For convenience, procedural and general Warnings and Precautions are listed at the beginning of the instructions. Please review all content before performing the procedure.

For reference literature concerning patient assessment, clinician education, insertion techniques and potential complications associated with this procedure refer to Arrow International, Inc. website: [www.arrowintl.com](http://www.arrowintl.com)

Caution

Do not reuse

Do not resterilize

Sterilized  
by ethylene oxide

Do not use  
if package is damaged

Does not contain  
natural rubber latex

Consult  
instructions for use

## Arrow International, Inc.

© 2012 Arrow International, Inc. All rights reserved.

An issued or revision date for these instructions is included for user information. If two years have elapsed between this date and product use, the user should contact Arrow International, Inc. to see if additional product information is available.

Issued Date: May 2012

Rx only.

# Pressure Injectable Central Venous Catheter (CVC) Product

## Product Description

The Arrow® Pressure Injectable CVC is a central venous catheter manufactured with medical grade, flexible polyurethane. A Pressure Injectable CVC may vary from two to four noncommunicating lumens. The catheter has a soft Blue FlexTip® that is more pliable than the catheter body. Lumens are connected to separate color-coded extension lines which have hubs on the end that are standard Luer-Lock. Centimeter markings referenced from the tip are placed along length of indwelling catheter body to facilitate proper positioning. The kit components assist the clinician in maintaining maximal sterile barrier precautions (where provided).

## Indications:

The Arrow® CVC is indicated to provide short-term (< 30 days) central venous access for treatment of diseases or conditions requiring central venous access including, but not limited to:

- multiple infusions of fluids, medications, or chemotherapy
- infusion of fluids that are hypertonic, hyperosmolar, or have divergent pH values
- frequent blood sampling or blood/blood component infusions
- infusion of incompatible medications
- central venous pressure monitoring
- lack of usable peripheral IV sites
- replacement of multiple peripheral sites for IV access
- injection of contrast media
  - When used for pressure injection of contrast media, do not exceed the maximum indicated flow rate for each catheter lumen. The maximum pressure of power injector equipment used with pressure injectable CVC may not exceed 400 psi.

## Contraindications:

None known. See additional labeling for product specific contraindications.

## Central Venous Catheter

### *Warnings and Precautions:*

### **WARNING**

DO NOT PLACE THE CATHETER INTO OR ALLOW IT TO REMAIN IN THE RIGHT ATRIUM OR RIGHT VENTRICLE. FAILURE TO FOLLOW THESE INSTRUCTIONS CAN RESULT IN SEVERE PATIENT INJURY OR DEATH.

**READ INSTRUCTIONS**

## General Warnings and Precautions

### Warnings:

1. Sterile, Single use: Do not reuse, reprocess or resterilize. Reuse of device creates a potential risk of serious injury and/or infection which may lead to death.
2. Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.
3. Practitioners must be aware of complications associated with central vein catheters including but not limited to: cardiac tamponade secondary to vessel wall, atrial or ventricular perforation, pleural (i.e., pneumothorax) and mediastinal injuries, air embolism, catheter embolism, catheter occlusion, thoracic duct laceration, bacteremia, septicemia, thrombosis, inadvertent arterial puncture, nerve damage, hematoma, hemorrhage, and dysrhythmias.

4. Do not place Central Venous Catheter (CVC) or Peripherally Inserted Central Catheter (PICC) into or allow them to remain in the right atrium or right ventricle. X-ray exam or other method in compliance with hospital/institutional protocol must show catheter tip located in right side of mediastinum in the Superior Vena Cava (SVC) above its junction with right atrium and parallel to vessel wall and its distal tip positioned at a level above either azygos vein or carina of the trachea, whichever is better visualized. Although cardiac tamponade secondary to pericardial effusion is uncommon, there is a high mortality rate associated with it. Improper advancement of guidewire into the heart has also been implicated in causing cardiac perforation and tamponade. **For femoral vein approach, catheter should be advanced into vessel so catheter tip lies parallel to vessel wall and does not enter right atrium.**
5. Ensure catheter tip has not entered the heart or no longer lies parallel to vessel wall by performing an x-ray exam or other method in compliance with hospital/institutional protocol. If catheter position has changed, immediately re-evaluate.
6. Practitioners must be aware of the potential for entrapment of guidewire by any implanted device in circulatory system (i.e., vena cava filters, stents). Review patient's history before catheterization procedure to assess for possible implants. Care should be taken regarding length of guidewire inserted. It is recommended that if patient has a circulatory system implant, catheter procedure be done under direct visualization to minimize the risk of guidewire entrapment.
7. Choose appropriate sized catheter for size of vessel to be cannulated.
8. Catheter tip must be located in central circulation when administering > 10% glucose solution, total parenteral nutrition, continuous vesicant therapy, infusates with pH less than 5 or greater than 9, and infusates with an osmolality above 600 mOsm/L, or any medication known to be irritating to vessels proximal to the vena cava.
9. Do not leave open needles or uncapped, unclamped catheters in central venous puncture site. Air embolism can occur with these practices.
10. Use only securely tightened Luer-Lock connections with any Venous Access Device (VAD) to guard against inadvertent disconnect.
11. Use Luer-Lock connectors to help guard against air embolism and blood loss.
12. Pulsatile flow is usually an indicator of inadvertent arterial puncture.

### Cautions:

1. Do not use if package has been previously opened or damaged.
2. Do not alter the catheter, guidewire, or any other kit/set component during insertion, use, or removal (except as instructed).
3. Procedure must be performed by trained personnel well versed in anatomical landmarks, safe technique, and potential complications.
4. Assess patient for heparin sensitivity. Heparin-Induced Thrombocytopenia (HIT) has been reported with use of heparin flush solutions.
5. Do not use topical antibiotic ointment or creams on insertion sites (except when using dialysis catheters) because of their potential to promote fungal infections and antimicrobial resistance.
6. Temporarily shut off remaining port(s) through which solutions are being infused before blood sampling.
7. Do not rely on blood aspirate color to indicate venous access.
8. Do not reinsert needle into introducer catheter (where provided) to reduce risk of catheter embolism.
9. Maintain insertion site with regular meticulous redressing using aseptic technique.
10. Engage safety and/or locking feature of scalpel (where provided) when not in use to reduce risk of sharps injury.
11. Perform hand hygiene:
  - before and immediately after all clinical procedures
  - before donning and after removal of gloves

12. Properly handle and dispose of sharps in sharps container in accordance with US OSHA or other governmental standards for blood borne pathogens and/or hospital/institutional policy.
13. Keep hands behind needle at all times during use and disposal.
14. Use universal blood and body-fluid precautions in the care of all patients due to the risk of exposure to Human Immunodeficiency Virus (HIV) or other blood borne pathogens.

## Catheter Warnings and Precautions

### Warnings:

1. Only utilize catheters indicated for high pressure injection applications for such applications. Utilizing catheters not indicated for high pressure applications can result in inter-lumen crossover or rupture with potential for injury.
2. Do not apply excessive force in placing or removing catheter. Excessive force can cause catheter breakage. If placement or withdrawal cannot be easily accomplished, an x-ray should be obtained and further consultation requested.
3. Do not secure, staple, and/or suture directly to outside diameter of catheter body or extension lines to reduce risk of cutting or damaging the catheter or impeding catheter flow. Secure only at indicated stabilization locations.
4. Do not cut catheter to alter catheter length.
5. Do not attach catheter clamp and fastener (where provided) until guidewire is removed.
6. Do not use scissors to remove dressing to reduce risk of cutting catheter.
7. Open catheter clamp prior to infusion through lumen to reduce risk of damage to extension line from excessive pressure.
8. Practitioners should remove slide clamp(s), where provided, when not in use. Slide clamp(s) may be inadvertently removed and aspirated by children or confused adults.
9. Do not routinely replace central venous catheters solely for the purpose of reducing incidence of infection.
10. Do not use guidewire techniques to replace catheters in patient suspected of having catheter-related infection.
11. Residual catheter track remains an air entry point until completely sealed, occlusive dressing should remain in place for at least 24 - 72 hours dependent upon amount of time catheter was indwelling.

### Cautions:

1. Check ingredients of prep sprays and swabs before using. Some disinfectants used at catheter insertion site contain solvents which can attack the catheter material. Alcohol and acetone can weaken the structure of polyurethane materials. These agents may also weaken the adhesive bond between catheter stabilization device and skin.
  - Acetone: Do not use acetone on catheter surface.
  - Alcohol: Do not use alcohol to soak catheter surface or to restore catheter patency.
 Take care when instilling drugs containing high concentration of alcohol. Allow insertion site to dry completely prior to applying dressing.
2. Ensure catheter patency prior to injection. Do not use syringes smaller than 10 mL (a fluid filled 1 mL syringe can exceed 300 psi), to reduce risk of intraluminal leakage or catheter rupture.
3. Remove catheter clamp and fastener (where provided) prior to attempting a catheter exchange procedure.
4. Do not exert excessive force while removing catheter to reduce risk of catheter breakage.
5. Continuously monitor indwelling catheter for:
  - desired flow rate
  - security of dressing
  - adherence of stabilization device to skin and connection to catheter
  - correct catheter position; use centimeter markings to identify if catheter position has changed
  - secure Luer-Lock connection(s)

6. Minimize catheter manipulation throughout procedure to maintain proper catheter tip position.
7. Inject a small amount of radiopaque dye to locate catheter tip if difficulty is encountered in visualizing the catheter tip.

## Guidewire / SWG Warnings and Precautions

### Warnings:

1. Do not cut guidewire to alter length.
2. Do not withdraw guidewire against needle bevel to reduce risk of possible severing or damaging of guidewire.
3. Do not use excessive force when introducing guidewire or tissue dilator as this can lead to vessel perforation and bleeding.
4. Passage of guidewire into the right heart can cause dysrhythmias, right bundle branch block, and a perforation of vessel, atrial, or ventricular wall.
5. Do not apply undue force on guidewire to reduce risk of possible breakage.
6. Do not apply excessive force in removing guidewire or catheter. If withdrawal cannot be easily accomplished, a visual image should be obtained and further consultation requested.
7. Do not cut guidewire with scalpel.
8. Position cutting edge of scalpel away from guidewire.
9. Engage safety and/or locking feature of scalpel (where provided) once cutaneous puncture site is enlarged, to reduce risk of cutting the guidewire.
10. Do not aspirate with guidewire in place or air will enter syringe.

### Cautions:

1. Maintain a firm grip on guidewire at all times. Keep sufficient guidewire length exposed at hub for handling purposes. A non-controlled guidewire can lead to wire embolism.
2. If resistance is encountered when attempting to remove guidewire after catheter placement, guidewire may be kinked about tip of catheter within vessel (refer to Figure 3).

## Tissue Dilator Warnings

### Warnings:

1. Do not leave tissue dilator in place as an indwelling catheter. Leaving tissue dilator in place puts patient at risk for possible vessel wall perforation.
2. Do not use excessive force when introducing guidewire or tissue dilator as this can lead to vessel perforation and bleeding.

## Pressure Injection Warnings and Precautions

### Warnings:

1. Assess each patient for appropriateness of a pressure injection procedure. Pressure injection procedures must be performed by trained personnel well versed in safe technique and potential complications.
2. Obtain a visual image to confirm catheter tip position prior to each pressure injection.
3. Ensure catheter patency prior to pressure injection to reduce risk of catheter failure and/or patient complications.
4. Discontinue pressure injections at first sign of extravasation or catheter deformation. Follow hospital/institutional protocol for appropriate medical intervention.

### Cautions:

1. Do not exceed ten (10) injections or catheter's maximum recommended flow rate located on product labeling and catheter luer hub to reduce risk of catheter failure and/or tip displacement.
2. Warm contrast media to body temperature prior to pressure injection to reduce risk of catheter failure.
3. Pressure limit settings on injector equipment may not prevent over pressurizing an occluded or partially occluded catheter.

4. Use appropriate administration set tubing between catheter and pressure injector equipment to reduce risk of catheter failure.
5. Follow the contrast media manufacturer's specified instructions for use, contraindications, warnings, and precautions.

### Possible Complications (but not limited to):

- ◆ cardiac tamponade secondary to vessel wall, atrial or ventricular perforation
- ◆ pleural injury
- ◆ air embolism
- ◆ catheter embolism
- ◆ bleeding / hemorrhage
- ◆ bacteremia
- ◆ thrombosis
- ◆ hematoma
- ◆ fibrin sheath formation
- ◆ vessel erosion
- ◆ mediastinal injury
- ◆ nerve injury
- ◆ thoracic duct laceration
- ◆ occlusion
- ◆ septicemia
- ◆ inadvertent arterial puncture
- ◆ dysrhythmias
- ◆ exit site infection
- ◆ catheter tip malposition

### Accessory Component Instructions

Review the list of components that will be utilized before beginning the Arrow® Pressure Injectable CVC insertion procedure. Kits / Sets may not contain all accessory components detailed in this section. Become familiar with instructions for individual component(s) before beginning the actual procedure.

#### Arrow Advancer™:

Arrow Advancer™ is used to straighten "J" Tip of Spring-Wire Guide (SWG) for introduction of guidewire into Arrow® Raulerson Syringe or a needle.

- Using thumb, retract "J" (refer to Figure 1).


Figure 1

- Place tip of Arrow Advancer™ – with "J" retracted – into the hole in rear of Arrow® Raulerson Syringe plunger or introducer needle.
- Advance SWG into Arrow® Raulerson Syringe approximately 10 cm until it passes through syringe valves (refer to Figure 2) or into introducer needle.


Figure 2

- Raise thumb and pull Arrow Advancer™ approximately 4 - 8 cm away from Arrow® Raulerson Syringe or introducer needle. Lower thumb onto Arrow Advancer™ and while maintaining a firm grip on SWG, push assembly into syringe barrel to further advance SWG. Continue until SWG reaches desired depth.

#### Arrow® Raulerson Syringe:

Arrow® Raulerson Syringe is used in conjunction with Arrow Advancer™ for SWG insertion.

- |                   |  |
|-------------------|--|
| <b>⚠ Warning:</b> | Do not withdraw guidewire against needle bevel to reduce risk of possible severing or damaging of guidewire.   |
| <b>⚠ Warning:</b> | Do not aspirate with guidewire in place or air may enter syringe.  |
| <b>⚠ Caution:</b> | Do not rely on blood aspirate color to indicate venous access.   |
| <b>⚠ Caution:</b> | Do not reinforce blood to reduce risk of blood leakage from rear (cap) of syringe.   |
| <b>⚠ Caution:</b> | Maintain a firm grip on guidewire at all times. Keep sufficient guidewire length exposed at hub for handling purposes. A non-controlled guidewire can lead to wire embolism. |

- Insert introducer needle attached to Arrow® Raulerson Syringe into vessel and aspirate.
  - Vessel may be prelocated with a 22 Ga. locator needle.
  - For pressure wave form transduction, a separately packaged transduction probe is available.
- Straighten SWG "J".
- Advance SWG through Arrow® Raulerson Syringe into vessel to desired depth.
- Hold SWG in place and remove introducer needle and Arrow® Raulerson Syringe.

#### SharpsAway II™ Locking Disposal Cup:

The SharpsAway II™ Locking Disposal Cup is used for disposal of needles (15 Ga. - 30 Ga.).

- |                   |  |
|-------------------|--|
| <b>⚠ Caution:</b> | Do not attempt to remove needles that have been placed into SharpsAway II™ Locking Disposal Cup. These needles are secured in place. Damage may occur to needles if they are forced out of disposal cup. |
| <b>⚠ Caution:</b> | Do not re-use needles after they have been placed into the foam SharpsAway® system. Particulate matter may adhere to needle tip.   |
| •                 | Discard after single use.  |

### Catheter Insertion Instructions

#### A Suggested Procedure:

- |                   |  |
|-------------------|--|
| <b>⚠ Warning:</b> | Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. |
|-------------------|--|

Clinical assessment of patient must be completed to ensure no contraindications exist e.g. allergies. This device is not recommended for use in the presence of device related infections or previous/current thrombosis.

#### Preparing for Catheter Insertion:

- Use sterile technique.
- 1. Prep and drape puncture site.
- 2. Flush each lumen of catheter with sterile saline solution, to establish patency and prime lumen(s).

- |                   |                                      |
|-------------------|--------------------------------------|
| <b>⚠ Warning:</b> | Do not cut catheter to alter length. |
|-------------------|--------------------------------------|

#### Gain Initial Venous Access:

- |                   |  |
|-------------------|--|
| <b>⚠ Warning:</b> | Do not leave open needles or uncapped, unclamped catheters in central venous puncture site. Air embolism can occur with these practices. |
|-------------------|--|

1. Insert introducer needle or catheter/needle with attached syringe or Arrow® Raulerson Syringe (where provided) into vein and aspirate.

*NOTE: If larger introducer needle is used, vessel may be pre-located with 22 Ga. locator needle and syringe.*

*NOTE: An echogenic needle, if included, is used to provide greater needle visibility under ultrasound. See Arrow® Raulerson Syringe under Accessory Component Instructions section.*

- Remove locator needle.

**Caution:** Do not rely on blood aspirate color to indicate venous access.

**Caution:** Do not reinsert needle into introducer catheter to reduce risk of catheter embolism.

#### Verify Venous Access:

Utilize one of the following techniques to verify venous access, because of the potential for inadvertent arterial placement:

- Central Venous Waveform:
  - Insert fluid primed blunt tip pressure transduction probe into rear plunger and through valves of the Arrow® Raulerson Syringe and observe for central venous pressure waveform.
  - Remove transduction probe if using Arrow® Raulerson Syringe.
- Pulsatile Flow (if hemodynamic monitoring equipment is not available):
  - Use transduction probe to open syringe valving system of Arrow® Raulerson Syringe and observe for pulsatile flow.
  - Disconnect syringe from needle and observe for pulsatile flow.

**Warning:** Pulsatile flow is usually an indicator of inadvertent arterial puncture.

#### Guidewire / SWG Insertion:

- Straighten "J" of SWG by using a straightening tube or Arrow Advancer™ (where provided).

See Arrow Advancer™ under Accessory Component Instructions section.

- Insert tip of SWG into plunger of Arrow® Raulerson Syringe or into introducer needle.
- Advance SWG through Arrow® Raulerson Syringe or introducer needle into vein to desired depth.
  - Advancement of "J" Tip through Arrow® Raulerson Syringe may require a gentle rotating motion.
  - Advance SWG until triple band mark reaches rear of Arrow® Raulerson Syringe plunger.

- Use centimeter markings on SWG as a reference to assist in determining how much guidewire has been inserted.

**NOTE:** When guidewire is used in conjunction with the Arrow® Raulerson Syringe (fully aspirated) and a 2-1/2 inch introducer needle, the following positioning references can be made:

- 20 cm mark (two bands) entering back of plunger = guidewire tip is at end of needle
- 32 cm mark (three bands) entering back of plunger = guidewire tip is approximately 10 cm beyond end of needle

**Caution:** Maintain firm grip on guidewire at all times. Keep sufficient guidewire length exposed at hub for handling purposes. A non-controlled guidewire can lead to wire embolism.

**Warning:** Do not aspirate with guidewire in place or air may enter syringe.

**Caution:** Do not reinfuse blood to reduce risk of blood leakage from rear (cap) of syringe.

**Warning:** Do not withdraw guidewire against needle bevel to reduce risk of possible severing or damaging of guidewire.

- Remove introducer needle and Arrow® Raulerson Syringe (or catheter) while holding SWG in place.

- Use centimeter markings on SWG to adjust indwelling length according to desired depth of indwelling catheter placement.

- Enlarge cutaneous puncture site with cutting edge of scalpel, if necessary, positioned away from guidewire.

**Warning:** Do not cut guidewire to alter length.

**Warning:** Do not cut guidewire with scalpel.

- Position cutting edge of scalpel away from guidewire.
- Engage safety and/or locking feature of scalpel (where provided) once cutaneous puncture site is enlarged, to reduce risk of cutting the guidewire.

- Use tissue dilator to enlarge puncture site as required.

**Warning:** Do not leave tissue dilator in place as an indwelling catheter. Leaving tissue dilator in place puts patient at risk for possible vessel wall perforation.

#### Advance Catheter:

- Thread tip of catheter over SWG. Sufficient SWG length must remain exposed at hub end of catheter to maintain a firm grip on SWG.
  - Grasping near skin, advance catheter into vein with slight twisting motion.
- Warning:** Do not attach catheter clamp and fastener (where provided) until guidewire is removed.
- Using centimeter marks on catheter as positioning reference points, advance catheter to final indwelling position.

**NOTE:** Centimeter marking symbology is referenced from catheter tip.

- numerical: 5, 15, 25, etc.
- bands: each band denotes a 10 cm interval, with one band indicating 10 cm, two bands indicating 20 cm, etc.
- dots: each dot denotes a 1 cm interval

- Hold catheter at desired depth and remove SWG.

**NOTE:** Arrow® catheters are designed to pass freely over SWG.

**Caution:** If resistance is encountered when attempting to remove guidewire after catheter placement, guidewire may be kinked about tip of catheter within vessel (refer to Figure 3).


Figure 3

- In this circumstance, pulling back on SWG may result in undue force being applied resulting in SWG breakage.
- If resistance is encountered, withdraw catheter relative to SWG about 2-3 cm and attempt to remove SWG.
- If resistance is again encountered, remove SWG and catheter simultaneously.

**Warning:** Do not apply undue force on guidewire to reduce risk of possible breakage.

- Verify entire SWG is intact upon removal.

#### Complete Catheter Insertion:

- Check lumen placement by attaching a syringe to each extension line and aspirate until free flow of venous blood is observed.
- Flush lumen(s) to completely clear blood from catheter.
- Connect all extension line(s) to appropriate Luer-Lock line(s) as required. Unused port(s) may be "locked" through luer-activated connector(s) using standard hospital/institutional protocol.
  - Slide clamps are provided on extension lines to occlude flow through each lumen during line and luer-activated connector changes.

**Warning:** Open catheter clamp prior to infusion through lumen to reduce risk of damage to extension line from excessive pressure.

- Secure catheter: Use a catheter clamp and rigid fastener, catheter stabilization device, staples or sutures.
  - Use triangular junction hub with side wings as primary suture site.
  - Use catheter clamp and fastener as a secondary suture site as necessary. Snap rigid fastener onto catheter clamp.

**Caution:** Minimize catheter manipulation throughout procedure to maintain proper catheter tip position.

- Ensure insertion site is dry before applying dressing. Apply skin protectant as needed.

**Caution:** Do not use topical antibiotic ointment or creams on insertion sites (except when using dialysis catheters) because of their potential to promote fungal infections and antimicrobial resistance.

- Assess placement of catheter tip in compliance with hospital/institutional protocol.

7. If catheter tip is malpositioned, reposition, redress, and re-verify.
8. Document procedure.

## Pressure Injection Instructions

**Use sterile technique.**

**⚠ Warning:** Obtain a visual image to confirm catheter tip position prior to each pressure injection.

1. Remove injection cap from appropriate extension line of catheter.
2. Check for catheter patency:
  - Attach 10 mL syringe filled with sterile normal saline.
  - Aspirate catheter for adequate blood return.
  - Vigorously flush catheter.
3. Detach syringe.
4. Attach pressure injection administration set tubing to appropriate extension line of catheter according to manufacturer's recommendations.
5. Inject contrast media in accordance with hospital/institutional protocol.
6. Disconnect catheter from pressure injector equipment.
7. Flush catheter using 10 mL syringe or larger filled with sterile normal saline.
8. Disconnect syringe and replace with sterile injection cap on catheter extension line.

**NOTE: Do not exceed ten (10) pressure injections.**

## Catheter Removal Instructions

**⚠ Warning:** Do not use guidewire techniques to replace catheters in patient suspected of having catheter-related infection.

1. Place patient in supine position, as clinically indicated to reduce risk of potential air embolism.
2. Remove dressing.
3. Remove sutures or staples; or open catheter stabilization device retainer wings and remove catheter from catheter stabilization device posts.
4. Place gauze pad over insertion site and catheter.
5. Remove catheter by slowly pulling it parallel to skin. If resistance is met while removing, catheter should not be forcibly removed and physician should be notified.
6. Upon removal of catheter:
  - inspect for intact Blue FlexTip® or catheter tip
  - ensure entire catheter length has been removed
7. Apply direct pressure to site until hemostasis is achieved.
8. Dress insertion site. Apply sterile air occlusive dressing and assess site every 24 hours until site is epithelialized.

**⚠ Caution:** Do not exert excessive force while removing catheter, to reduce risk of catheter breakage.

**⚠ Warning:** Residual catheter track remains an air entry point until completely sealed, occlusive dressing should remain in place for at least 24 - 72 hours dependent upon amount of time catheter was indwelling.

9. Document catheter removal procedure on patient's chart per hospital/institutional protocol.

For reference literature concerning patient assessment, clinician education, insertion techniques and potential complications associated with this procedure refer to Arrow International, Inc. website: [www.arrowintl.com](http://www.arrowintl.com)

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