Cath-Gard® Contamination Shield

Rx only.

Intended Purpose:

The Cath-Gard® Contamination Shield is intended to cover the exposed portion of a catheter inserted through an introducer to reduce contamination of the indwelling devices during procedural manipulation.

Contraindications:

None known.

⚠ General Warnings and Precautions

Warnings:

1. Sterile, Single use: Do not reuse, reprocess or resterilize. Reuse of device creates a potential risk of serious injury and/or infection which may lead to death. Reprocessing of medical devices intended for single use only may result in degraded performance or a loss of functionality.
2. Read all package insert warnings, precautions and instructions prior to use. Failure to do so may result in severe patient injury or death.

Precautions:

1. Do not alter the Cath-Gard component except as instructed.
2. Procedure must be performed by trained personnel well versed in anatomical landmarks, safe technique and potential complications.
3. Use standard precautions and follow institutional policies for all procedures including safe disposal of devices.

A Suggested Procedure: Use sterile technique.

1. If using a flow-directed catheter prepare according to manufacturer's instructions. Wet balloon with flush solution to facilitate passage through catheter contamination shield.

⚠ **Precaution:** Do not inflate balloon of flow-directed catheter prior to insertion through catheter contamination shield to reduce the risk of balloon damage.

2. Apply Contamination Shield:

- a. If using a catheter contamination shield with Tuohy-Borst adapter, insert tip of desired catheter through Tuohy-Borst adapter end of catheter contamination shield. Advance catheter through tubing and hub at other end (refer to Figure 1).


Figure 1

- b. If using a catheter contamination shield with TwistLock™ adapter, ensure double TwistLock of catheter contamination shield is fully opened (refer to Figure 2).


Figure 2

- 0 Insert tip of desired catheter through proximal end of catheter contamination shield. Advance catheter through tubing and hub at other end (refer to Figure 3).


Figure 3

3. Slide entire catheter contamination shield to proximal end of catheter.
4. If flow directed catheter is used, inflate and deflate balloon with syringe to ensure integrity.

⚠ **Caution:** Do not exceed balloon catheter manufacturer's recommended volume.

Place catheter and catheter contamination shield on sterile field awaiting final placement.

5. After successful positioning of the catheter per manufacturer's instructions, hold catheter in place and reposition catheter contamination shield so that distal hub is approximately five inches (12.7 cm) from hemostasis valve.
6. Hold proximal hub of catheter contamination shield in place. Disengage distal hub from inner feed tube by pulling forward. Advance distal hub forward toward hemostasis valve assembly. Hold assembly in place.
7. Press distal hub of catheter contamination shield over assembly cap. Twist to lock (refer to Figure 4).


Figure 4

- Orient slot in hub with locking pin on assembly cap.
 - Slide hub forward over cap and twist.
8. While maintaining catheter position lock the catheter in place:
 - a. If using a catheter contamination shield with a Tuohy-Borst adapter, grasp catheter through front portion of catheter contamination shield and hold in place while repositioning Tuohy-Borst adapter end as desired.

⚠ **Precaution:** Do not reposition Tuohy-Borst adapter on insertion catheter once moved to this final position.

- Tighten Tuohy-Borst adapter by pressing down on cap and simultaneously turning clockwise to secure hub to catheter. Gently pull catheter to verify securement.

⚠️ Precaution: Do not overtighten Tuohy-Borst adapter to reduce the risk of lumen constriction or insertion catheter damage.

- Tuohy-Borst adapter end of catheter contamination shield should be secured with sterile tape to inhibit catheter movement (refer to Figure 5).

⚠️ Precaution: Do not apply tape to the transparent sheathing on the shield to reduce the risk of tearing material.


Figure 5

- b. If using a catheter contamination shield with a TwistLock adapter, twist the upper half of the distal hub in clockwise direction to lock catheter in place. Reposition proximal end of catheter shield as desired. Twist upper and lower halves in opposite directions to lock in place. Test the adapter by gently tugging on the catheter to ensure a secure grip on the catheter (refer to Figure 6).

⚠️ Precaution: Do not reposition proximal hub once locked in final position.


Figure 6

For reference literature concerning patient assessment, clinician education, insertion technique, and potential complications associated with this procedure, consult standard textbooks, medical literature, and Arrow International LLC website: www.teleflex.com

A pdf copy of this IFU is located at www.teleflex.com/IFU

For a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation 2017/745/EU on Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority. The contacts of national competent authorities (Vigilance Contact Points) and further information can be found on the following European Commission website: https://ec.europa.eu/growth/sectors/medical-devices/contacts_en

en

Symbol Glossary: Symbols are in compliance with ISO 15223-1.

Some symbols may not apply to this product. Refer to product labeling for symbols that apply specifically to this product.

Caution	Medical device	Consult instructions for use	Do not reuse	Do not resterilize	Sterilized by ethylene oxide	Single sterile barrier system	Keep away from sunlight	Keep dry
Do not use if package is damaged	Not made with natural rubber latex	Catalogue number	Lot number	Use by	Manufacturer	Date of manufacture	Importer	

Arrow, the Arrow logo, Cath-Gard, Teleflex, the Teleflex logo, and TwistLock are trademarks or registered trademarks of Teleflex Incorporated or its affiliates, in the U.S. and/or other countries. © 2021 Teleflex Incorporated. All rights reserved.

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