EN

ARROW®

Arterial/Venous Catheterization Product

Safety and Efficacy Considerations:

Do not use if package has been previously opened or damaged.

Warning: Read all package insert warnings, precautions and instructions prior to use. Failure to do so may result in severe patient injury or death.

Do not alter catheter or any other kit/set component during insertion, use or removal.

Procedure must be performed by trained personnel well versed in anatomical landmarks, safe technique and potential complications.

Indications for Use:

The single-lumen catheter permits access to the arterial or peripheral venous circulation.

Contraindications:

Inadequate collateral arterial circulation.

Arterial/Venous Warnings and Precautions:*

- Warning: Sterile, Single use: Do not reuse, reprocess or resterilize. Reuse of device creates a potential risk of serious injury and/or infection which may lead to death.**
- Warning: Care should be exercised that the indwelling catheter is not inadvertently kinked at the hub area when securing catheter to the patient. Kinking may weaken the wall of the catheter and cause a fraying or fatigue of the material, leading to possible separation of the catheter.**
- Warning: Do not apply excessive force in removing guide wire or catheter.**
- Warning: The practitioner must be aware of potential air embolism associated with leaving open needles or catheters in venous or arterial puncture sites or as a consequence of inadvertent disconnects. To lessen the risk of disconnects, only securely tightened Luer-Lock connections should be used with this device. Follow hospital protocol to guard against air embolism for all catheter maintenance.**
- Warning: Due to the risk of exposure to Human Immunodeficiency Virus (HIV) or other blood borne pathogens, health care workers should routinely use universal blood and body-fluid precautions in the care of all patients.**
- Warning: Do not suture directly to outside diameter of catheter body to minimize the risk of damaging the catheter or adversely affecting monitoring capabilities.**
- Precaution: Indwelling catheters should be routinely inspected for desired flow rate, security of dressing, correct catheter position and for secure Luer-Lock connection.**
- Precaution: Alcohol and acetone can weaken the structure of polyurethane materials. Check ingredients of prep sprays and swabs for acetone and alcohol content.**

Acetone: Do not use acetone on catheter surface. Acetone may be applied to skin but must be allowed to dry completely prior to applying dressing.

Alcohol: Do not use alcohol to soak catheter surface or to restore catheter patency. Care should be taken

when instilling drugs containing high concentration of alcohol. Always allow alcohol to dry completely prior to applying dressing.

- 9. Precaution: Some disinfectants used at the catheter insertion site contain solvents which can attack the catheter material. Assure insertion site is dry before dressing.**

Arterial Catheterization

Warnings and Precautions:*

- Warning: Practitioners must be aware of complications associated with arterial procedures^{1,10,11,15} including bacteremia,² septicemia, vessel wall perforation, intravascular clotting and embolization, hematoma, arterial spasm, tissue necrosis, hemorrhage, thrombosis, peripheral ischemia and infarction, peripheral nerve damage¹² and air embolism.³**

A Suggested Procedure:

Use sterile technique.

- Prep and drape anticipated puncture site per hospital protocol.
- Utilize local anesthetic as required.
- Vessel may be accessed using catheter/needle with vented plug only, with spring-wire guide tube or with a syringe. Catheter hub wing clip may be removed, if desired.
 - Vented plug only – Remove needle guard.
 - Spring-wire guide tube – Remove vented plug. Attach spring-wire guide tube assembly (AW-04018) to hub of needle. To ensure proper function, advance and retract spring-wire guide using actuating lever. **Precaution: Prior to insertion actuating lever must be retracted proximally as far as possible or blood flashback may be inhibited.** Remove needle guard.
 - Syringe – Remove vented plug. Attach syringe to hub of needle. Remove needle guard.
- Make a 2 mm skin nick with scalpel at the desired puncture site.
- Puncture vessel using a continuous, controlled, slow, forward motion. Avoid transfixing both vessel walls. **Precaution: If both vessel walls are punctured, subsequent advancement of spring-wire guide could result in inadvertent sub-arterial placement.** Blood flashback indicates successful entry into vessel.
- Advance entire placement device a maximum of 1 to 2 mm further into vessel.
- Stabilize position of introducer needle. Insert spring-wire guide.
 - Vented plug only – Remove needle and plug. Pass 0.018" guide wire through catheter.
 - Spring-wire guide tube – Carefully advance spring-wire guide as far as required into vessel using actuating lever (refer to Fig. 1).


Fig. 1

Precaution: Do not advance guide wire unless there is free blood flashback. When reference mark on clear spring-wire tube coincides with edge of internal cylinder of actuation lever, the tip of the spring-wire guide is located at needle tip. **Warning:** Do not retract spring-wire guide against edge of needle while in vessel to minimize the risk of spring-wire guide damage.

- Syringe – Remove needle and syringe. Pass 0.018" spring-wire guide through catheter. **Precaution:** If resistance is encountered while advancing spring-wire guide, do not force feed. If resistance is encountered during spring-wire guide advancement, withdraw entire unit and attempt new puncture.
- 8. Firmly hold introducer needle hub and/or spring-wire guide in position and advance catheter forward with a slight rotating motion over spring-wire guide into vessel.
- 9. Hold catheter in place and remove spring-wire guide or introducer needle/spring-wire guide assembly. Pulsatile blood flow indicates positive arterial placement. **Precaution:** Do not reinsert needle into catheter to minimize risk of catheter damage.
- 10. Attach desired stopcock, injection cap or connecting tubing to catheter hub. Secure catheter to patient in preferred manner using suture wings, suture groove or wing clip, where provided. **Warning:** Care should be exercised that the indwelling catheter is not inadvertently kinked at the hub area when securing catheter to the patient. Kinking may weaken the wall of the catheter and cause a fraying or fatigue of the material, leading to possible separation of the catheter. **Warning:** Do not suture directly to outside diameter of catheter body to minimize the risk of damaging the catheter or adversely affecting monitoring capabilities.
- 11. Dress according to hospital protocol. **Precaution:** Indwelling catheters should be routinely inspected for desired flow rate, security of dressing, correct catheter position and for secure Luer-Lock connection.

Peripheral Venous Catheterization

Warnings and Precautions:*

1. **Warning:** Practitioners must be aware of complications associated with peripheral intravenous catheters^{14,18} including infiltration, catheter embolus, bacteremia, septicemia, thrombosis, inadvertent arterial puncture, nerve damage, hematoma and hemorrhage.
2. **Precaution:** Use of a syringe smaller than 10 mL to irrigate or declot an occluded catheter may cause catheter rupture.⁵

A Suggested Procedure:

Use sterile technique.

1. Prep puncture site as required.
2. Puncture vessel using a continuous, controlled, slow, forward motion. Avoid transfixing both vessel walls. Blood flashback in hub of introducer needle indicates successful entry into vessel. **Precaution:** The color of the blood aspirated is not always a reliable indicator of venous access.⁹ Aspiration may be required.

3. Advance catheter and needle, as a unit, approximately 1-2 mm further into vessel.
4. Firmly hold needle in position and advance catheter forward with a slight rotating motion into vessel. Remove introducer needle. **Precaution:** Do not reinsert needle into catheter to minimize the risk of catheter damage.
5. Attach stopcock, injection cap or connecting tubing to catheter hub. **Precaution:** To lessen the risk of disconnects, only securely tightened Luer-Lock connections should be used with this device.
6. Secure catheter to patient. Suture groove or wing clip may be used. **Warning:** Care should be exercised that the indwelling catheter is not inadvertently kinked at the hub area when securing catheter to the patient. Kinking may weaken the wall of the catheter and cause a fraying or fatigue of the material, leading to possible separation of the catheter. **Warning:** Do not suture directly to outside diameter of catheter body to minimize the risk of damaging the catheter or adversely affecting monitoring capabilities. **Precaution:** If using suture ring, suture should be snug but not tight.
7. Dress puncture site per hospital protocol. **Precaution:** Indwelling catheters should be routinely inspected for desired flow rate, security of dressing, correct catheter position and for secure Luer-Lock connection.

Catheter Removal Procedure:

1. Remove dressing. **Precaution:** To minimize the risk of cutting the catheter, do not use scissors to remove dressing.
2. **Warning:** Exposure of vessel to atmospheric pressure may result in entry of air into the circulation. Remove suture(s). Be careful not to cut catheter. Remove catheter slowly, pulling it parallel to skin. As catheter exits site, apply pressure with a dressing until hemostasis occurs. Apply light pressure dressing to site.
3. Upon removal of catheter, inspect it to make sure entire length has been withdrawn.
4. Verify catheter was intact upon removal. Document removal procedure.

References:

1. Abadir AR, Ung KA. Complications of radial artery cannulation. *Anesthesiology Rev.* 1980;7:11-16.
2. Band JD, Maki DG. Infections caused by arterial catheters used for hemodynamic monitoring. *Am J Med.* 1979;67:735-741.
3. Chang C, Dughi J, Shitabata P, Johnson G, Coel M, McNamara JJ. Air embolism and the radial arterial line. *Crit Care Med.* 1988;16:141-143.
4. Clark CA, Harman EM. Hemodynamic monitoring: arterial catheters. In: Taylor RW, Civetta JM, Dirby RR, eds. *Techniques and Procedures in Critical Care.* Philadelphia, PA: JB Lippincott; 1990:218-230.
5. Conn C. The importance of syringe size when using implanted vascular access devices. *J Vasc Access Nurs.* Winter 1993;3:11-18.
6. Daily EK, Schroeder JS. *Techniques in Bedside Hemodynamic Monitoring.* 5th ed. St. Louis, MO: Mosby; 1994:71-77.
7. Eissa NT, Kvetan V. Guide wire as a cause of complete heart block in patients with preexisting bundle branch block. *Anesthesiology.* 1990;73:772-774.
8. Falk PS, Scuderi PE, Sherertz RJ, Motsinger SM. Infected radial artery pseudoaneurysms occurring after percutaneous cannulation. *Chest.* 1992;101:490-495.

9. Jobs DR, Schwartz AJ, Greenhow DE, Stephenson LW, Ellison N. Safer jugular vein cannulation: recognition of arterial punctures and preferential use of the external jugular route. *Anesthesiology*. 1983;59:353-355.
10. Kaye W. Invasive monitoring techniques: arterial cannulation, bedside pulmonary artery catheterization, and arterial puncture. *Heart Lung*. 1983;12:395-427.
11. Kaye WE, Dubin HG. Vascular cannulation. In: Taylor RW, Civetta JM, Dirby RR, eds. *Techniques and Procedures in Critical Care*. Philadelphia, PA: JB Lippincott; 1990:204-208.
12. Marshall G, Edelstein G, Hirschman CA. Median nerve compression following radial arterial puncture. *Anesth Analg*. 1980;59:953-954.
13. Norwood SH, Cormier B, McMahon NG, Moss A, Moore V. Prospective study of catheter related infection during prolonged arterial catheterization. *Crit Care Med*. 1988;16:836-839.
14. Perdue M. Intravenous complications. In: Terry J, Baranowski L, Lonsway RA, Hedrick C, eds. *Intravenous Therapy*. Philadelphia PA: WB Saunders;1995:419-436.
15. Puri VK, Carlson RW, Bander JJ, Weil MH. Complications of vascular catheterization in the critically ill. *Crit Care Med*. 1980;8:495-499.
16. Sigurdsson J, Riba P, Sigurdsson S. The wandering central venous catheter. *Intensive Care Med*. 1985;11:263-264.
17. Venus B, Satish P. Vascular Cannulation. In: Civetta JM, Taylor RW, Kirby RR, eds. *Critical Care*. 3rd ed. Philadelphia, PA: Lippincott-Raven Publishers; 1997:539-542.
18. Weinstein SM. *Plumer's Principles & Practice of Intravenous Therapy*. 5th ed. Philadelphia PA: JB Lippincott Co; 1993: 83-91.

Arrow International, Inc. recommends that the user be acquainted with the reference literature.

* If you have any questions or would like additional reference information, please contact Arrow International, Inc.

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Arrow International, Inc. rekommenderar att användaren studerar referenslitteraturen.
* Om du har frågor eller önskar ytterligare referensinformation, kontakta Arrow International, Inc.