Inspected Dimensions:
Folded width: 5-1/2" (14 cm)
Folded length: 7-1/2" (19 cm)

ARROW®

Spring-Wire Guide Product

Safety and Efficacy Considerations:

Do not use if package has been previously opened or damaged. **Warning:** Prior to use read all package insert warnings, precautions, and instructions. Failure to do so may result in severe patient injury or death.

Do not alter spring-wire guide during insertion, use, or removal.

Assembly procedure must be performed by trained personnel well versed in anatomical landmarks, safe technique and potential complications.

Indications for Use:

The spring-wire guide with arterial/venous (AV) catheter introduction.

Contraindications:

None known.

Warnings and Precautions:*

1. **Warning: Sterile, Single use: Do not reuse, reprocess or resterilize. Reuse of device creates a potential risk of serious injury and/or infection which may lead to death.**
2. **Warning: Practitioners must be aware of complications associated with spring-wire guide use including knotting^{5,6} separation of the coil and core wires^{1,2,3,4,8,9} vessel perforation⁷ and malposition.¹⁰**
3. **Warning: Do not apply excessive force in removing guide wire. If withdrawal cannot be easily accomplished, a radiograph should be obtained and further consultation requested.**
4. **Warning: Due to the risk of exposure to HIV (Human Immunodeficiency Virus) or other blood borne pathogens, health care workers should routinely use universal blood and body-fluid precautions in the care of all patients.**

AV Catheter Introduction

1. Attach spring-wire tube assembly to hub of needle.
2. Remove protective shield from catheter (if applicable). Trial advance and retract spring-wire

guide through needle using actuating lever to ensure proper function. **Precaution: Prior to insertion actuating lever must be retracted proximally as far as possible or blood flashback may be inhibited.**

3. Puncture vessel using a continuous, controlled, slow, forward motion. Avoid transfixing both vessel walls. **Precaution: If both vessel walls are punctured, subsequent advancement of spring-wire guide could result in inadvertent sub-vascular placement.** Blood flashback in hub of introducer needle indicates successful entry into vessel (refer to Fig. 1).


Fig. 1

4. Stabilize position of introducer needle and carefully advance spring-wire guide as far as required into vessel using actuating lever (refer to Fig. 2). **Precaution: Do not advance guide wire unless there is free blood flashback.**


Fig. 2

When reference mark of clear feed tube coincides with edge of internal cylinder of actuating lever, the tip of spring-wire guide is located at the needle tip. **Precaution: If resistance is encountered while advancing spring-wire guide, do not force feed. Warning: Do not retract spring-wire guide against edge of needle while in vessel to minimize the risk of spring-wire guide damage.**

If resistance is encountered during spring-wire guide advancement, withdraw entire unit and attempt new puncture.

5. Advance entire placement device a maximum of 1 to 2 mm further into vessel.
6. Firmly hold introducer needle hub in position and advance catheter forward with a slight rotating motion over spring-wire guide into vessel (refer to Fig. 3).


Fig. 3

7. Hold catheter in place and remove needle and spring-wire guide assembly. Pulsatile blood flow indicates positive arterial placement. **Precaution: Do not reinsert needle into catheter to minimize risk of catheter damage.**
8. Continue procedure per manufacturer's instructions.
9. Verify that entire spring-wire guide is intact upon removal.

References:

1. Akazawa S, Nakaigawa Y, Hotta K, Shimizu R, Kashiwagi H, Takahashi K. Unrecognized migration of an entire guidewire on insertion of a central venous catheter into the cardiovascular system. *Anesthesiology*. 1996;84:241-242. Letter.
2. Bhanumurthy S, McMurray TJ, Bannon C. Early venous catheter blockade: a cause for immediate concern. *Anaesthesia*. 1995;50:473. Letter.

3. Cucchiara RF, Muzzi DA. Guide-wire retention after right atrial catheter insertion. *Anesth Analg*. 1992;74:303-304.
4. Hickman JA, McCrerrick A. Failure of a spring wire guide. *Anaesthesia and Intens Care*. 1990;18:587. Letter.
5. Hoogsteden L, Filshie J, Querci della Rovere G. Knotted guidewire – a complication of Hickman line insertion. *Anaesthesia*. 1996;51:713. Letter.
6. Kwok SPY, Lau WY, Li AKC. Insertion of subclavian line: don't get knotted. *Aust NZ J Surg*. 1994;64:64-65.
7. Masud T, Tapson JS. Mediastinal hematoma caused by subclavian catheterization for hemodialysis. *Int J Artif Organs*. 1989;12:708-710.
8. Polos PG, Sahn SA. Complication of central venous catheter insertion: fragmentation of a guidewire with pulmonary artery embolism. *Crit Care Med*. 1991;19:438-440.
9. Reynen K. 14-year follow-up of central embolization by a guide wire. *N Engl J Med*. 1993;329:970-971. Letter.
10. Trujillo MH, Arai K. Hydrothorax after inadvertent placement of a central venous catheter in the left pericardiophrenic vein. *J Intensive Care Med*. 1994;9:257-260.

Arrow International, Inc. recommends that the user be acquainted with the reference literature.

*If you have any questions or would like additional reference information, please contact Arrow International, Inc.

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Arrow International, Inc. rekommenderar att användaren studerar referenslitteraturen.
*Om du har frågor eller önskar ytterligare referensinformation, kontakta Arrow International, Inc.