# ARROW®

## Spring-Wire Guide Product

### Safety and Efficacy Considerations:

Do not use if package has been previously opened or damaged. **Warning:** Prior to use read all package insert warnings, precautions, and instructions. Failure to do so may result in severe patient injury or death.

Do not alter spring-wire guide during insertion, use, or removal.

Spring-wire guide use must be performed by trained personnel well versed in anatomical landmarks, safe technique and potential complications.

### Indications for Use:

The spring-wire guide permits central venous catheter introduction or exchange.

### Contraindications:

None known.

### Warnings and Precautions:\*

1. **Warning: Sterile, Single use: Do not reuse, reprocess or resterilize. Reuse of device creates a potential risk of serious injury and/or infection which may lead to death.**
2. **Warning: Practitioners must be aware of complications associated with spring-wire guide use including knotting<sup>8,10</sup>, separation of the coil and core wires<sup>1,3,5,7,12,13</sup>, vessel perforation<sup>11</sup>, malposition<sup>16</sup>, displacement of vena cava filters<sup>2,14</sup>, cardiac tamponade<sup>4,9,17</sup> and cardiac dysrhythmias<sup>6,15</sup>.**
3. **Warning: Do not apply excessive force in removing guide wire. If withdrawal cannot be easily accomplished, a chest x-ray should be obtained and further consultation requested.**
4. **Warning: Due to the risk of exposure to HIV (Human Immunodeficiency Virus) or other blood borne pathogens, health care workers should routinely use universal blood and body-fluid precautions in the care of all patients.**
5. **Precaution: If resistance is encountered when attempting to remove spring-wire guide after catheter placement, spring-wire may be kinked about tip of catheter within vessel (refer to Fig. 1).**


Fig. 1

To minimize the risk of spring-wire guide damage withdraw catheter relative to spring-wire guide about 2-3 cm and attempt to remove spring-wire guide. If resistance is again encountered remove spring-wire guide and catheter simultaneously.

6. **Precaution: Prior to attempting catheter exchange procedure remove catheter clamp and fastener.**

### A Suggested Procedure: Use sterile technique.

1. Follow manufacturer's instructions for patient preparation and venous access.
2. Insert desired tip of spring-wire guide through the introducer needle or catheter into vein.

**Straightening Tube:** If the "J" tip portion of the spring-wire guide is used, prepare for insertion by sliding the plastic tube over the "J" to straighten.

**Arrow Advancer™:** Straighten spring-wire guide "J" by retracting into Arrow Advancer™ with thumb.

3. The spring-wire guide should then be advanced in the routine fashion to the desired depth. Advancement of "J" tip may require a gentle rotating motion. **Warning: Do not cut spring-wire guide to alter length. Do not withdraw spring-wire guide against needle bevel to minimize the risk of possible severing or damaging of spring-wire guide.**
4. Hold spring-wire guide in place and remove introducer needle or catheter. **Precaution: Maintain firm grip on spring-wire guide at all times.**
5. Continue procedure per manufacturer's instructions.
6. Verify that entire spring-wire guide is intact upon removal.

## References:

1. Akazawa S, Nakaigawa Y, Hotta K, Shimizu R, Kashiwagi H, Takahashi K. Unrecognized migration of an entire guidewire on insertion of a central venous catheter into the cardiovascular system. *Anesthesiology*. 1996;84:241-242. Letter.
2. Amesbury S, Vargish T, Hall J. An unusual complication of central venous catheterization. *Chest*. 1994;105:905-907.
3. Bhanumurthy S, McMurray TJ, Bannon C. Early venous catheter blockade: a cause for immediate concern. *Anaesthesia*. 1995;50:473. Letter.
4. Blake PG, Uldall R. Cardiac perforation by a guide wire during subclavian catheter insertion. *Int J Artif Organs*. 1989;12:111-113.
5. Cucchiara RF, Muzzi DA. Guide-wire retention after right atrial catheter insertion. *Anesth Analg*. 1992;74:303-304.
6. Eissa NT, Kvetan V. Guide wire as a cause of complete heart block in patients with preexisting left bundle branch block. *Anesthesiology*. 1990;73:772-774.
7. Hickman JA, McCrirrick A. Failure of a spring wire guide. *Anaesthesia and Intens Care*. 1990;18:587. Letter.
8. Hoogsteden L, Filshie J, Querci della Rovere G. Knotted guidewire – a complication of Hickman line insertion. *Anaesthesia*. 1996;51:713. Letter.
9. Khan IH. Pericardial tamponade resulting from changing a central venous catheter over a guide wire. *Postgraduate Med J*. 1991;67:591.
10. Kwok SPY, Lau WY, Li AKC. Insertion of subclavian line: don't get knotted. *Aust NZ J Surg*. 1994;64:64-65.
11. Masud T, Tapson JS. Mediastinal hematoma caused by subclavian catheterization for hemodialysis. *Int J Artif Organs*. 1989;12:708-710.
12. Polos PG, Sahn SA. Complication of central venous catheter insertion: fragmentation of a guidewire with pulmonary artery embolism. *Crit Care Med*. 1991;19:438-440.
13. Reynen K. 14-year follow-up of central embolization by a guide wire. *N Engl J Med*. 1993;329:970-971. Letter.
14. Rosenblum JD, Boyle CM. Percutaneous retrieval of a Vena-Tech filter displaced during central line placement. *AJR*. 1996;166:994-995. Letter.
15. Stuart RK, Shikora SA, Akerman P, et al. Incidence of arrhythmia with central venous catheter insertion and exchange. *J Parenteral Enteral Nutrition*. 1990;14:152-155.
16. Trujillo MH, Arai K. Hydrothorax after inadvertent placement of a central venous catheter in the left pericardiophrenic vein. *J Intensive Care Med*. 1994;9:257-260.
17. van Haeften TW, van Pampus ECM, Boot H, Strack van Schijndel RJM, Thijs LG. Cardiac tamponade from misplaced central venous line in pericardiophrenic vein. *Arch Intern Med*. 1988;148:1649-1650.

Arrow International, Inc. recommends that the user be acquainted with the reference literature.

\*If you have any questions or would like additional reference information, please contact Arrow International, Inc.

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