Guidewire / SWG Product

FR	Fil guide / Fil guide spirale	4
DE	Führungsdraht / SWG	9
IT	Guida / Prodotto SWG	14
PT	Fio-guia	19
ES	Guía / producto SWG	23
SV	Ledare / SWG-produkt	27

ARROW[®]

Arrow International, Inc.

© 2009 Arrow International, Inc. All rights reserved. Printed in the USA.

An issued or revision date for these instructions is included for user information. In the event two years have elapsed between this date and product use, the user should contact Arrow International, Inc. to see if additional product information is available.

Issued Date: September 2009

Rx only.

Guidewire / SWG Product

Table of Contents

<i>Product Description</i>	1
Indications / Contraindications	1
<hr/>	
<i>Guidewire Warnings and Precautions</i>	1
General	1
Positive Placement	2
<hr/>	
<i>Guidewire Product Instructions</i>	2
General Guidewire Usage and Removal	2
Positive Placement Guidewire Usage and Removal	3
<hr/>	
<i>Documentation</i>	3

For convenience, procedural and general Warnings and Precautions are listed at the beginning of the instructions. Please review all content before performing the procedure.

For reference literature concerning patient assessment, clinician education, insertion techniques and potential complications associated with this procedure refer to Arrow International, Inc. website: www.arrowintl.com

Guidewire / SWG Product

Product Description

A guidewire is used to aid in the placement of a catheter into the vasculature or into a body cavity. Once placed, the guidewire provides a pathway for the catheter to be advanced over the guidewire and into correct position. After the catheter is inserted the guidewire is removed and is appropriately discarded.

Indications:

The guidewire permits vessel or cavity access for catheter introduction or exchange.

Contraindications:

None known.

Guidewire

Warnings and Precautions:

General Warnings and Precautions

Warnings:

1. Sterile, Single use: Do not reuse, reprocess or resterilize. Reuse of device creates a potential risk of serious injury and/or infection which may lead to death.
2. Prior to use read all package insert warnings, precautions, and instructions. Failure to do so may result in severe patient injury or death.
3. Practitioners must be aware of complications associated with guidewire use including knotting, separation of the coil and core wires, vessel perforation, malposition, displacement of vena cava filters, cardiac tamponade and cardiac dysrhythmias.
4. Do not apply excessive force in removing guidewire to minimize the risk of breakage. If withdrawal cannot be easily accomplished, a visual image should be obtained and further consultation requested.
5. Do not insert the solid proximal end (without coil) of guidewire, where provided, as it may lead to vessel damage.
6. Do not cut guidewire to alter length.

7. Do not cut guidewire with scalpel.
 - Position cutting edge of scalpel away from guidewire.
 - Retract blade of safety scalpel to protected position once cutaneous puncture site is enlarged, to minimize the risk of cutting the guidewire.
8. Do not withdraw guidewire against needle bevel to minimize the risk of possible severing or damaging of guidewire.
9. Do not use excessive force when introducing guidewire as this can lead to vessel perforation and bleeding.
10. Passage of guidewire into the right heart can cause dysrhythmias, right bundle branch block, and a perforation of vessel wall, atrial, or ventricular.
11. Practitioners must be aware of the potential for entrapment of guidewire by any implanted device in the circulatory system (i.e., vena cava filters, stents). Review patient's history before catheterization procedure to assess for possible implants. Care should be taken regarding length of guidewire inserted. It is recommended that if patient has a circulatory implant, catheter procedure to be done under direct visualization to minimize the risk of guidewire entrapment.

Cautions:

1. Do not use if package has been previously opened or damaged.
2. Do not alter the guidewire or any other kit/set component during insertion, use, or removal (except as instructed).
3. Prior to attempting catheter exchange procedure remove catheter clamp and fastener so that a guidewire can be passed freely.
4. Guidewire use must be performed by trained personnel well versed in anatomical landmarks, safe technique and potential complications.
5. If resistance is encountered when attempting to remove guidewire after catheter placement, guidewire may be

kinked about tip of catheter within vessel (refer to Figure 1).

Figure 1

To minimize the risk of guidewire damage withdraw catheter relative to guidewire about 2-3 cm and attempt to remove guidewire. If resistance is again encountered remove guidewire and catheter simultaneously.

6. Maintain a firm grip on guidewire at all times. Keep sufficient guidewire length exposed at hub for handling purposes. A non-controlled guidewire can lead to wire embolism.
7. Perform hand hygiene before and immediately after all clinical procedures and before and after donning and removal of gloves.
8. Properly handle and dispose of sharps in sharps container in accordance with US OSHA or other governmental standards for blood borne pathogens and/or hospital/institutional policy.
9. Hands must remain behind the needle at all times during use and disposal.
10. Use universal blood and body-fluid precautions in the care of all patients due to the risk of exposure to HIV (Human Immunodeficiency Virus) or other blood borne pathogens.

Positive Placement Warnings and Precautions

Warnings:

1. Do not use excessive force when introducing guidewire as this can lead to vessel perforation and bleeding.
2. Do not retract guidewire against edge of needle while in vessel to minimize the risk of guidewire damage.
3. Do not apply excessive force in removing guidewire to minimize the risk of breakage.

Caution:

1. Prior to insertion actuating lever must be retracted proximally as far as possible or blood flashback may be inhibited.
2. If both vessel walls are punctured, subsequent advancement of guidewire could result in inadvertent subvascular placement.
3. Do not advance guidewire unless there is free blood flashback.

Guidewire Product Instructions

Clinicians should use sterile technique throughout a procedure.

General Guidewire Usage and Removal

1. Follow manufacturer's instructions for patient preparation and vessel or body cavity access.
2. Insert desired tip of guidewire through the introducer needle or catheter into vessel or body cavity.

Warning: Do not insert the solid proximal end (without coil) of guidewire, where provided, as it may lead to vessel damage.

- "J" Tip Usage:
 - Straightening Tube (where provided): If the "J" Tip portion of the guidewire is used, prepare for insertion by sliding the plastic tube over the "J" to straighten.
 - Arrow Advancer™ (where provided): Straighten guidewire "J" by retracting into Arrow Advancer™ with thumb.
- Use centimeter markings on guidewire to assist in determining how much guidewire has been inserted.

NOTE: If guidewire has three sets of markings, they will be located as follows:

- One band - 10 cm from "J" Tip
- Two bands - 20 cm from "J" Tip
- Three bands - 32 cm from "J" Tip

NOTE: When guidewire is used in conjunction with the Arrow® Raulerson Syringe (fully aspirated) and a 2-1/2 inch introducer needle, the following positioning references can be made:

- 20 cm mark entering back of plunger = guidewire tip is at end of needle
- 32 cm mark entering back of plunger = guidewire tip is approximately 10 cm beyond end of needle

3. The guidewire should then be advanced in the routine fashion to the desired depth.
 - Advancement of “J” Tip may require a gentle rotating motion.

Warning: Do not cut guidewire to alter length.

Warning: Do not withdraw guidewire against needle bevel to minimize the risk of possible severing or damaging of guidewire.

4. Hold guidewire in place and remove introducer needle or catheter.

Caution: Maintain firm grip on guidewire at all times.

5. Continue procedure per manufacturer’s instructions.

Warning: Do not apply excessive force in removing guidewire to minimize the risk of breakage. If withdrawal cannot be easily accomplished, a visual image should be obtained and further consultation requested.

6. Verify that entire guidewire is intact upon removal.

Positive Placement Guidewire Usage and Removal

1. Follow manufacturer’s instructions for patient preparation and vessel access.
2. Attach guidewire tube assembly to hub of needle.
3. Remove protective shield from catheter (if applicable). Trial advance and retract guidewire through needle using actuating lever to ensure proper function.

Caution: Prior to insertion actuating lever must be retracted proximally as far as possible or blood flashback may be inhibited.

Caution: If both vessel walls are punctured, subsequent advancement of guidewire could result in inadvertent subvascular placement.

4. Use activating lever to advance guidewire into vessel (refer to Figure 2).


Figure 2

- If resistance is encountered during guidewire advancement, withdraw entire unit and attempt new puncture.

Caution: Do not advance guidewire unless there is free blood flashback.

Warning: Do not use excessive force when introducing guidewire as this can lead to vessel perforation and bleeding.

Warning: Do not retract guidewire against edge of needle while in vessel to minimize the risk of guidewire damage.

5. Hold catheter in place and remove needle and guidewire assembly.

Warning: Do not apply excessive force in removing guidewire to minimize the risk of breakage.

6. Verify that entire guidewire is intact upon removal.

Documentation: Hospitals/institutions must establish a permanent medical record that documents the entire procedure, based upon their policy, procedures, and Best Practices. The actual format can differ from institution to institution. Report any product defects/failures to organization risk management, manufacturers, and appropriate regulatory agencies.

For reference literature concerning patient assessment, clinician education, insertion techniques and potential complications associated with this procedure refer to Arrow International, Inc. website: www.arrowintl.com

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