



REF CDC-44063-HPK1A (IPNXXXXXX)

6 Fr. 3 Lumen 40 cm catheter length .018 inch dia. spring-wire guide

ARROW®

ErgoPack®

Pressure Injectable Arrowg+ard Blue Advance® Three-Lumen PICC



Contents:

- | | |
|---|---|
| 1: Three Lumen TaperFree® Catheter with Arrowg+ard Blue Advance®
Antimicrobial/Antithrombogenic Protection ¹ : 6 Fr. (2.11 mm OD) x 40 cm,
Pressure Injectable, T-Port Connector, Blue FlexTip® and Placement Wire | 1: Towel |
| 1: GlideThru™ Peel-Away Sheath: 6 Fr. x 2-3/4" (7 cm) Radiopaque over 6 Fr.
Dilator | 1: Safety Scalpel: #11 |
| 1: Spring-Wire Guide, Nitinol, Marked: .018" (0.46 mm) dia. x 17-3/4" (45 cm)
(Straight Soft Tip on One End - Straight Stiff Tip on Other) with Arrow Advancer | 1: Sterile Procedure Sign |
| 1: Injection Needle: SafetyGlide™ ² 25 Ga. x 5/8" (1.60 cm) | 1: Patient ID Card |
| 1: Safety Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW | 1: Chart Sticker |
| 1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm) | 1: Patient Information Booklet |
| 1: Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW | 1: Medication Label: 1% Lidocaine |
| 1: Syringe: 10 mL Luer-Lock | 1: Checklist/CLIP Sheet |
| 1: Syringe: 3 mL Luer-Lock | 2: Paper Tape Measure |
| 2: 10 mL Luer-Lock Pre-Filled Saline Syringe | 1: Tourniquet |
| 1: Skin Protectant Prep | 10: Gauze Pad: 2" x 2" (5 cm x 5 cm) |
| 1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep | 10: Gauze Pad: 4" x 4" (10 cm x 10 cm) |
| 1: 3 mL Applicator 2% CHG and 70% IPA ChloraPrep® ³ One-Step Solution | 1: Dressing: Tegaderm® ⁴ 15.5 cm x 10 cm |
| 1: 3 mL Applicator 2% CHG and 70% IPA ChloraPrep® ³ One-Step Solution with
Hi-Lite Orange™ Tint | 1: Surgical Apparel: Impervious Gown |
| 3: Dust Cap: Non-Vented | 1: Dressing: BIOPATCH® ⁵ |
| 1: SecondSite™ Adjustable Hub: Fastener | 1: Dressing: STATLOCK® ⁶ Catheter Stabilization Device |
| 1: SecondSite™ Adjustable Hub: Catheter Clamp | 1: Surgical Apparel: Mask with Eye Shield |
| 1: SharpsAway® II Locking Disposal Cup | 1: Surgical Apparel: Bouffant Cap |
| 1: SharpsAway® Disposal Cup | 1: Surgical Apparel: Mask |
| 1: Catheter Trimmer | 1: Tape: Steri-Strip® ⁴ |
| 1: Towel: 24" x 36" (61 cm x 91 cm) | 1: HemoHopper® Fluid Receptacle |
| 1: Maximal Barrier Drape™ with two 13 cm (5") x 15 cm (6") fenestrations | |

¹Licensed under US Patent No. 7,329,412.²A trademark of Becton, Dickinson and Company.³A registered trademark of CareFusion or one of its subsidiaries.⁴A registered trademark of 3M Company.⁵A registered trademark of Johnson & Johnson Corporation.⁶A registered trademark of C. R. Bard, Inc.

Not made with natural rubber latex.

Store between 20 – 25°C (68 – 77°F).

Arrow International provides the enclosed medication label(s) for your convenience. Please ensure that these labels are applied to the correct syringe and corresponding medication.

Contraindications: The Pressure Injectable Arrowg+ard Blue Advance antimicrobial/antithrombogenic catheter is contraindicated:

- for patients with known hypersensitivity to chlorhexidine
- in the presence of device related infection in the intended insertion vessel or catheter pathway
- in the presence of thrombosis in the intended insertion vessel or catheter pathway



Lumen

	Priming Volume* (mL)	Gravity Flow Rate† (mL/hr)	Pump Flow Rate†† (mL/hr)	MAX Pressure Injection Flow Rate** (mL/sec)
Distal (18 Ga.)	0.5	1137	6970	6
Medial (19 Ga.)	0.39	323	2810	not rated
Proximal (19 Ga.)	0.41	337	2880	not rated

* Priming volumes are approximate and are done without accessories.

† Flow rate values are approximate and are determined using deionized water at 100 cm head height.

†† Pump flow rates are determined at maximum pump pressure of 10 psig and represent approximate flow capabilities.

** Pressure injection flow rates are determined at the injector pressure setting of 300 psi maximum using media of 11.8 centipoise viscosity, with 152 cm pressure tubing.

