



# ARROWgard Blue<sup>®</sup> You-Bend<sup>™</sup> Two-Lumen Hemodialysis Catheterization Kit for High Volume Infusions

# ARROW

REF CDA-25122-U1A

**12** Fr. | **2** Lumen | **20** cm catheter length | **.035** inch dia. spring-wire guide

**ErgoPack**<sup>®</sup>

**Contents:**

- |   |  |  |
|---|--|--|
| <ul style="list-style-type: none"> <li>1: Two-Lumen Indwelling Catheter<sup>1</sup>: 12 Fr. x 20 cm Radiopaque Polyurethane with Blue FlexTip<sup>®</sup>, ARROWgard<sup>®</sup> Antimicrobial Surface Treatment<sup>1</sup></li> <li>1: Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 26-3/4" (68 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer</li> <li>1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle</li> <li>1: Pressure Transduction Probe</li> <li>1: Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow<sup>®</sup> Raulerson Spring-Wire Introduction Syringe</li> <li>1: Injection Needle: Eclipse<sup>™2</sup> 22 Ga. x 1-1/2" (3.81 cm) RW and 5 mL Luer-Slip Syringe</li> <li>1: Injection Needle: Eclipse<sup>™2</sup> 25 Ga. x 1" (2.54 cm) RW and 3 mL Luer-Lock Syringe</li> <li>1: Syringe: 10 mL Luer-Lock</li> </ul> | <ul style="list-style-type: none"> <li>1: Tissue Dilator: 10 Fr. (3.5 mm) x 10.2 cm</li> <li>1: Tissue Dilator: 12 Fr. (4.0 mm) x 14 cm</li> <li>1: SharpsAway<sup>®</sup> Disposal Cup</li> <li>1: SharpsAway II<sup>™</sup> Locking Disposal Cup</li> <li>1: Maximal Barrier Drape<sup>™</sup> with 4" fenestration</li> <li>1: Towel</li> <li>1: Needle Holder</li> <li>1: Forceps</li> <li>1: Safety Scalpel: #11</li> <li>1: Checklist/CLIP Sheet</li> <li>1: Sterile Procedure Sign</li> <li>2: Gauze Pad: 2" x 2" (5 cm x 5 cm)</li> <li>5: Gauze Pad: 4" x 4" (10 cm x 10 cm)</li> <li>1: Dressing: Tegaderm<sup>®3</sup> 10 cm x 12 cm</li> <li>1: Surgical Apparel: Impervious Gown</li> </ul> | <ul style="list-style-type: none"> <li>1: Surgical Apparel: Bouffant Cap</li> <li>1: Surgical Apparel: Mask with Eye Shield</li> <li>1: Suture: 3-0 Braided Silk with Straight Needle</li> <li>1: Suture: 3-0 Silk with Curved Needle</li> <li>1: HemoHopper<sup>®</sup> Fluid Receptacle</li> </ul> |
|---|--|--|

<sup>1</sup>Licensed under US Patent No. 6,706,024.  
<sup>2</sup>A trademark of Becton, Dickinson and Company.  
<sup>3</sup>A registered trademark of 3M Company.

**Contraindications:** The Arrow large-bore two-lumen catheter is not designed for long-term hemodialysis or for use in patients with thrombosed vessels. The ARROWgard Blue antimicrobial surface catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine, and/or sulfa drugs (refer to product instructions for references).




Not made with natural rubber latex



‡ MR Conditional



Store below 25°C (77°F).  
Avoid excessive heat above 40°C (104°F).



\*Priming volumes are approximate and are done without injection caps.

<sup>†</sup> In vitro tests performed using 3:1 glycerine/saline at 37°C.