



Arrowgard Blue[®] PSI Kit

Cath-Gard[®] for use with 7.5 - 8 Fr. Catheters

REF CDA-29903-1A

9^{Fr.} | 10^{cm} sheath length | .035^{inch dia.} spring-wire guide

ErgoPack[®]

Contents:

- 1: Sheath: 9 Fr. x 10 cm Arrow-Flex[®] Radiopaque Polyurethane with Arrowgard[®] Antimicrobial Surface Treatment¹, Integral Hemostasis Valve/Side Port and Tissue Dilator
- 1: Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- 1: Pressure Transduction Probe
- 1: Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow[®] Raulerson Spring-Wire Introduction Syringe
- 1: Injection Needle: Eclipse^{™2} 22 Ga. x 1-1/2" (3.81 cm) RW and 5 mL Luer-Slip Syringe
- 1: Injection Needle: Eclipse^{™2} 25 Ga. x 1" (2.54 cm) RW and 3 mL Luer-Lock Syringe
- 1: Syringe: 10 mL Luer-Lock
- 1: Obturator: 8 Fr.

- 1: Cath-Gard[®] Catheter Contamination Shield: 80 cm with TwistLock[™] Adapter
- 1: Dust Cap: Non-Vented
- 1: SharpsAway[®] Disposal Cup
- 1: SharpsAway II[™] Locking Disposal Cup
- 1: Maximal Barrier Drape[™] with 4" fenestration
- 1: Towel
- 1: Needle Holder
- 1: Forceps
- 1: Safety Scalpel: #11
- 1: Stopcock: 4-way High-Flow
- 1: Checklist/CLIP Sheet
- 1: Sterile Procedure Sign
- 2: Gauze Pad: 2" x 2" (5 cm x 5 cm)
- 5: Gauze Pad: 4" x 4" (10 cm x 10 cm)
- 1: Dressing: Tegaderm^{®3} 10 cm x 12 cm

- 1: Surgical Apparel: Impervious Gown
- 1: Surgical Apparel: Bouffant Cap
- 1: Surgical Apparel: Mask with Eye Shield
- 1: Suture: 3-0 Silk with Curved Needle
- 1: Tubing²: 20" (50 cm) Extension
- 1: HemoHopper[®] Fluid Receptacle

¹Licensed under US Patent No. 6,706,024.

²A trademark of Becton, Dickinson and Company.

³A registered trademark of 3M Company.

Sheath OD 3.9 mm

Contraindications: The Arrowgard Blue antimicrobial sheath introducer is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs (refer to product instructions for references).

Warning: The results of certain animal experiments have shown phthalates to be potentially toxic to reproduction. Proceeding from the present state of scientific knowledge, risks for male premature infants cannot be excluded in the case of long-term exposure or application. Medical products containing phthalates should be used only temporarily with pregnant women, nursing mothers, babies and infants.



Not made with natural rubber latex



± Precaution: Contains Phthalate: DEHP



Store below 25°C (77°F). Avoid excessive heat above 40°C (104°F).