



Pressure Injectable Multi-Lumen CVC Kit

REF CDA-42703-P1A

7 Fr. | 3 Lumen | 16 cm catheter length | .032 inch dia. spring-wire guide

ErgoPack®

Contents:

- 1: Multi-Lumen Indwelling Catheter: 7 Fr. x 16 cm Pressure Injectable Radiopaque Polyurethane with Blue FlexTip®, ARROWgãrd Blue PLUS® Antimicrobial Surface Treatment¹, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- 1: Pressure Transduction Probe
- 1: Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe
- 1: Injection Needle: Eclipse™² 22 Ga. x 1-1/2" (3.81 cm) RW and 5 mL Luer-Slip Syringe
- 1: Injection Needle: Eclipse™² 25 Ga. x 1" (2.54 cm) RW and 3 mL Luer-Lock Syringe
- 1: Syringe: 10 mL Luer-Lock
- 1: Tissue Dilator: 8.5 Fr. (2.8 mm) x 10.2 cm
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp
- 1: SharpsAway® Disposal Cup
- 1: SharpsAway II™ Locking Disposal Cup

- 1: Maximal Barrier Drape™ with 4" fenestration
- 1: Towel
- 1: Needle Holder
- 1: Forceps
- 1: Safety Scalpel: #11
- 1: Checklist/CLIP Sheet
- 1: Sterile Procedure Sign
- 1: Flow Rate Information Card with Injection Log
- 2: Gauze Pad: 2" x 2" (5 cm x 5 cm)
- 5: Gauze Pad: 4" x 4" (10 cm x 10 cm)
- 1: Dressing: Tegaderm®³ 10 cm x 12 cm
- 1: Surgical Apparel: Impervious Gown
- 1: Surgical Apparel: Bouffant Cap
- 1: Surgical Apparel: Mask with Eye Shield
- 1: Suture: 3-0 Silk with Curved Needle
- 1: HemoHopper® Fluid Receptacle

¹Licensed under US Patent Nos. 6,706,024; 6,872,195 and 7,329,412.

²A trademark of Becton, Dickinson and Company.

³A registered trademark of 3M Company.

Contraindications: The ARROWgãrd Blue PLUS® antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine, and/or sulfa drugs (refer to product instructions for references).

LATEX
Not made with natural rubber latex

Store below 25°C (77°F).
Avoid excessive heat above 40°C (104°F).

pt

Kit CVC ARROWgãrd BLUE PLUS Multi-Lúmen Para Infusão Sobre Pressão

Advertência: Leia todas as advertências, precauções e instruções do folheto informativo antes da utilização. Se não o fizer, poderá provocar lesões graves ou mesmo a morte do doente.

Contra-indicações: A utilização do cateter antimicrobiano ARROWgãrd Blue PLUS está contra-indicada em doentes com hipersensibilidade conhecida à cloro-hexidina, à sulfadiazina de prata e/ou a fármacos sulfa (consulte as instruções do produto para obter referências).

Esterilizado. Cuidado. Consultar as instruções de utilização. Não reutilizar. Não reesterilizar. Esterilizado por óxido de etileno. Manter afastado da luz solar. Manter seco. Não utilizar se a embalagem estiver danificada. Não fabricado em látex de borracha natural. Os componentes do trajecto do fluido são apirógenos.

Armazenar abaixo de 25°C (77°F). Evitar calor excessivo acima de 40°C (104°F).

es


Kit CVC ARROWgãrd BLUE PLUS Multi-Lumen de Inyección a Presión

Advertencia: Lea todas las advertencias, precauciones e instrucciones del prospecto antes de utilizar el dispositivo. El no hacerlo puede ocasionar lesiones graves o la muerte del paciente.

Contraindicaciones: El catéter antimicrobiano ARROWgãrd Blue PLUS está contraindicado en pacientes con hipersensibilidad conocida a la clorhexidina, la sulfadiazina de plata o las sulfamidas (consulte las instrucciones del producto).

Aséptico. Aviso. Consulte las instrucciones de uso. No reutilizar. No reesterilizar. Esterilizado mediante óxido de etileno. Mantener alejado de la luz del sol. Mantener seco. No utilizar si el envase está dañado. No fabricado con látex de caucho natural. Los componentes de la vía de líquido son apirógenos.

Almacenar a menos de 25 °C (77 °F). Evite las temperaturas superiores a los 40 °C (104 °F).



Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)	MAX Pressure Injection Flow Rate‡ (mL/sec)
Distal (16 Ga.)	0.38	3000	10
Medial (18 Ga.)	0.37	1400	5
Proximal (18 Ga.)	0.39	1600	5

* Priming volumes are approximate and are done without injection caps.

† Based on 100 cm head height.

‡ Pressure injection flow rates are determined using media of 11.8 centipoise (cP) viscosity, with 152 cm administration set tubing.