



MAC™ Two-Lumen Central Venous Access Kit Cath-Gard® for use with 7.5 - 8 Fr. Catheters

ARROW®

REF **CDC-11242-1A**

9 Fr. distal lumen | **2** Lumen | **11.5** cm device length | **.035** inch dia. spring-wire guide

ErgoPack®

Contents:

- 1: Two-Lumen Central Venous Access Device¹: 115 mm Radiopaque Polyurethane with 9 Fr. Distal Lumen, Integral Hemostasis Valve, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- 1: Injection Needle: SafetyGlide™² 25 Ga. x 1" (2.54 cm)
- 1: Local Anesthetic Component
- 1: Pressure Transduction Probe
- 1: Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe
- 1: Injection Needle: SafetyGlide™² 22 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe
- 1: Syringe: 3 mL Luer-Lock

- 1: Tissue Dilator
- 1: 5 mL 1% Lidocaine HCl Solution
- 1: 3 mL Applicator 2% CHG and 70% IPA ChloraPrep®³ One-Step Solution with Hi-Lite Orange™ Tint (non-sterile solution)
- 1: Obturator: 8 Fr.
- 1: Cath-Gard® Catheter Contamination Shield: 80 cm with TwistLock™ Adapter
- 1: Dust Cap: Vented
- 1: Dust Cap: Non-Vented
- 1: SharpsAway® Disposal Cup
- 1: SharpsAway II™ Locking Disposal Cup
- 1: Maximal Barrier Drape™ with 4" fenestration
- 1: Needle Holder
- 1: Safety Scalpel: #11
- 1: Stopcock: 4-way High-Flow

- 1: Checklist/CLIP Sheet
 - 1: Sterile Procedure Sign
 - 1: Medication Label: 1% Lidocaine
 - 2: Gauze Pad: 2" x 2" (5 cm x 5 cm)
 - 5: Gauze Pad: 4" x 4" (10 cm x 10 cm)
 - 1: Surgical Apparel: Impervious Gown
 - 1: Surgical Apparel: Bouffant Cap
 - 1: Surgical Apparel: Mask with Eye Shield
 - 1: Suture: 3-0 Silk with Curved Needle
 - 1: Tubing[†]: 20" (50 cm) Extension
 - 1: HemoHopper® Fluid Receptacle
- ¹Licensed under US Patent No. 5,250,038.
²A trademark of Becton, Dickinson and Company.
³A registered trademark of CareFusion or one of its subsidiaries.

DuPont™ Tyvek®
polypropylene




STERILE EO



~~LATEX~~

PHT DEHP



Rx only. **Warning: Read all package insert warnings, cautions and instructions prior to use. Failure to do so may result in severe patient injury or death.**

Sterilized using ethylene oxide unless otherwise indicated in the contents list.

California Prop. 65

Warning: Cancer and Reproductive Harm www.P65Warnings.ca.gov.

Fluid path components are non-pyrogenic.

Arrow International provides the enclosed medication label(s) for your convenience. Please ensure that these labels are applied to the correct syringe and corresponding medication.


Not made with natural rubber latex

± Precaution: Contains Phthalate: DEHP.

Store between 20 - 25°C (68 - 77°F).

Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)
Distal (9 Fr.) no catheter	1.80	30450
Distal (9 Fr.) with 8 Fr. catheter	1.15	12390
Distal (9 Fr.) with 7 Fr. catheter	1.45	19660
Proximal (12 Ga.) no catheter	0.60	11950

* Priming volumes are approximate and are done without accessories.
 † Based on 100 cm head height.



Manufacturer:
 Arrow International LLC
 Subsidiary of Teleflex Incorporated
 3015 Carrington Mill Blvd.
 Morrisville, NC 27560 USA

Teleflex®
www.teleflex.com

C-11242-200E, Rev. 01 (2020-12)



MAC™ Two-Lumen Central Venous Access Kit Cath-Gard® for use with 7.5 - 8 Fr. Catheters

ARROW®

REF **CDC-11242-1A**

ErgoPack®

9 Fr. distal lumen | **2** Lumen | **11.5** cm device length | **.035** inch dia. spring-wire guide