



ARROW®

ARROWgard Blue® Two-Lumen CVC Kit

REF CDC-22702-XP1A

7 Fr. | 2 Lumen | 16 cm catheter length | .025 inch dia. spring-wire guide

ErgoPack®

Contents:

- | | | |
|--|---|--|
| <ul style="list-style-type: none">1: Two-Lumen Indwelling Catheter: 7 Fr. x 16 cm Pressure Injectable Radiopaque Polyurethane with Blue FlexTip®, ARROWgard® Antimicrobial Surface Treatment¹, Extension Line Clamps1: Spring-Wire Guide, Marked: .025" (0.64 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer1: Catheter: 20 Ga. x 1-3/4" (4.45 cm) Radiopaque over 22 Ga. TW Introducer Needle1: Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) RW1: Injection Needle: SafetyGlide™² 25 Ga. x 1" (2.54 cm)1: Local Anesthetic Component1: Pressure Transduction Probe1: Introducer Needle: Echogenic 20 Ga. x 1-1/2" (3.81 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe1: Injection Needle: SafetyGlide™² 22 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe1: Syringe: 3 mL Luer-Lock | <ul style="list-style-type: none">1: Tissue Dilator: 8.5 Fr. (2.8 mm) x 10.2 cm2: 10 mL Luer-Lock Pre-Filled Saline Syringe1: 5 mL 1% Lidocaine HCl Solution1: 3 mL Applicator 2% CHG and 70% IPA Chloraprep®³ One-Step Solution with Hi-Lite Orange™ Tint (non-sterile solution)1: SecondSite™ Adjustable Hub: Fastener1: SecondSite™ Adjustable Hub: Catheter Clamp1: SharpsAway® Disposal Cup1: SharpsAway II™ Locking Disposal Cup1: Maximal Barrier Drape™ with 4" fenestration1: Towel1: Needle Holder1: Safety Scalpel: #111: Checklist/CLIP Sheet1: Flow Rate Information Card with Injection Log1: Sterile Procedure Sign | <ul style="list-style-type: none">1: Medication Label: 1% Lidocaine2: Gauze Pad: 2" x 2" (5 cm x 5 cm)5: Gauze Pad: 4" x 4" (10 cm x 10 cm)1: Dressing: Tegaderm®⁴ 10 cm x 12 cm1: Surgical Apparel: Impervious Gown1: Dressing: BIOPATCH®⁵1: Surgical Apparel: Bouffant Cap1: Surgical Apparel: Mask with Eye Shield1: Suture: 3-0 Silk with Curved Needle1: Tubing[‡]: 20" (50 cm) Extension1: HemoHopper® Fluid Receptacle |
|--|---|--|

¹Licensed under US Patent No. 6,706,024.
²A trademark of Becton, Dickinson and Company.
³A registered trademark of CareFusion or one of its subsidiaries.
⁴A registered trademark of 3M Company.
⁵A registered trademark of Johnson & Johnson Corporation.

Sterilized using ethylene oxide unless otherwise indicated in the contents list.

Contraindications: The ARROWgard Blue antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine, and/or sulfa drugs (refer to product instructions for references).

Arrow International provides the enclosed medication label(s) for your convenience. Please ensure that these labels are applied to the correct syringe and corresponding medication.




Not made with natural rubber latex



‡ Precaution:
Contains Phthalate: DEHP



Store between
20-25°C (68-77°F).



Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)	MAX Pressure Injection Flow Rate‡ (mL/sec)
Distal (18 Ga.)	0.35	1300	5
Proximal (14 Ga.)	0.50	5600	10

* Priming volumes are approximate and are done without injection caps.
† Based on 100 cm head height.
‡ Pressure injection flow rates are determined using media of 11.8 centipoise (cP) viscosity, with 152 cm administration set tubing.