



ARROW⁺gard Blue[®] Pediatric Multi-Lumen CVC Kit



REF CDC-26553-X1A

5.5 Fr. | **3** Lumen | **13** cm catheter length | **.018** inch dia. spring-wire guide



Contents:

- 1: Multi-Lumen Indwelling Catheter: 5.5 Fr. x 13 cm Radiopaque Polyurethane with Blue FlexTip[®], ARROW⁺gard[®] Antimicrobial Surface Treatment¹, Extension Line Clamps
- 1: Spring-Wire Guide: .018" (0.46 mm) dia. x 17-3/4" (45 cm) PTFE-Coated (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Catheter: 20 Ga. x 1.16" (30 mm) Insite™ Autoguard™ Radiopaque over 22 Ga. Introducer Needle
- 1: Injection Needle: SafetyGlide™ 25 Ga x 1" (2.54 cm)
- 1: Local Anesthetic Component
- 1: Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW
- 1: Introducer Needle: 21 Ga. x 1-1/2" (3.81 cm) TW and 5 mL Luer-Slip Syringe
- 1: Injection Needle: SafetyGlide™ 25 Ga. x 1" (2.54 cm) and 5 mL Luer-Slip Syringe
- 1: Syringe: 3 mL Luer-Lock
- 1: Tissue Dilator: 6 Fr. (2.0 mm) x 7.6 cm
- 1: Tissue Dilator: 6 Fr. (2.0 mm) x 4.1 cm
- 2: 10 mL Luer-Lock Pre-Filled Saline Syringe
- 1: 5 mL 1% Lidocaine HCl Solution
- 1: 3 mL Applicator 2% CHG and 70% IPA Chloraprep^{®3} One-Step Solution with Hi-Lite Orange™ Tint (non-sterile solution)
- 3: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp
- 1: SharpsAway[®] Disposal Cup
- 1: SharpsAway II™ Locking Disposal Cup
- 1: Drape: 41" x 55" with 4" fenestration
- 1: Towel
- 1: Needle Holder
- 1: Safety Scalpel: #11
- 1: Checklist/CLIP Sheet
- 1: Sterile Procedure Sign
- 1: Medication Label: 1% Lidocaine
- 2: Gauze Pad: 2" x 2" (5 cm x 5 cm)
- 5: Gauze Pad: 4" x 4" (10 cm x 10 cm)
- 1: Dressing: Tegaderm^{®4} 10 cm x 12 cm
- 1: Surgical Apparel: Impervious Gown
- 1: Dressing: BIOPATCH^{®5}
- 1: Dressing: STATLOCK^{®6} CV Plus Catheter Stabilization Device (prep pad solution non-sterile)
- 1: Surgical Apparel: Bouffant Cap
- 1: Surgical Apparel: Mask with Eye Shield
- 1: Suture: 3-0 Silk with Curved Needle
- 1: HemoHopper[®] Fluid Receptacle

¹Licensed under US Patent No. 6,706,024.
²A trademark of Becton, Dickinson and Company.
³A registered trademark of CareFusion or one of its subsidiaries.
⁴A registered trademark of 3M Company.
⁵A registered trademark of Johnson & Johnson Corporation.
⁶A registered trademark of C. R. Bard, Inc.

Sterilized using ethylene oxide unless otherwise indicated in the contents list.

Contraindications: The ARROW⁺gard Blue antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine, and/or sulfa drugs (refer to product instructions for references).


Arrow International provides the enclosed medication label(s) for your convenience. Please ensure that these labels are applied to the correct syringe and corresponding medication.



Not made with natural rubber latex



Store between 20-25°C (68-77°F).



Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr.)
Distal (20 Ga.)	0.30	1220
Medial (22 Ga.)	0.30	560
Proximal (22 Ga.)	0.40	670

* Priming volumes are approximate and are done without accessories.

† Based on 100 cm head height.