



Arrowg+ard Blue® PSI Kit Cath-Gard® for use with 7.5 - 8 Fr. Catheters

ARROW®

REF **CDC-29803-XN1A**

8.5 Fr. | **10** cm sheath length | **.035** inch dia. spring-wire guide

ErgoPack[®] Complete

Contents:

- | | | |
|--|--|---|
| 1: Sheath: 8.5 Fr. x 4" (10 cm) Arrow-Flex® Radiopaque Polyurethane with Arrowg+ard® Antimicrobial Surface Treatment ¹ , Integral Hemostasis Valve/Side Port and Tissue Dilator | 1: Syringe: 3 mL Luer-Lock | 1: Stopcock: 4-way High-Flow |
| 1: Spring-Wire Guide, Nitinol, Marked: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow® GlideWheel™ Wire Advancer | 2: 10 mL Luer-Lock Pre-Filled Saline Syringe | 1: Checklist/CLIP Sheet |
| 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle | 1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep | 1: Sterile Procedure Sign |
| 1: Injection Needle: SafetyGlide™ ² 25 Ga. x 1" (2.54 cm) | 2: 3 mL Applicator 2% CHG and 70% IPA Chloraprep® ³ One-Step Solution with Hi-Lite Orange™ Tint | 1: Medication Label: 1% Lidocaine |
| 1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm) | 1: Obturator: 8 Fr. | 2: Gauze Pad: 2" x 2" (5 cm x 5 cm) |
| 1: Pressure Transduction Probe | 1: Cath-Gard® Catheter Contamination Shield: 80 cm with TwistLock™ Adapter | 5: Gauze Pad: 4" x 4" (10 cm x 10 cm) |
| 1: Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe | 1: Dust Cap: Non-Vented | 1: Surgical Apparel: Impervious Gown |
| 1: Injection Needle: SafetyGlide™ ² 22 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe | 1: SharpsAway® Disposal Cup | 1: Dressing: I.V. Advanced 8.5 cm x 11.5 cm |
| | 1: SharpsAway® II Locking Disposal Cup | 1: Dressing: BIOPATCH® ⁴ |
| | 1: Maximal Barrier Drape™ with 4" fenestration | 1: Transducer Cover: 14 cm x 147 cm |
| | 1: Towel | 1: Surgical Apparel: Mask with Eye Shield |
| | 1: Needle Holder | 1: Surgical Apparel: Bouffant Cap |
| | 1: Safety Scalpel: #11 | 1: Suture: 3-0 Silk with Curved Needle |
| | | 1: Tubing±: 20" (50 cm) Extension |
| | | 1: HemoHopper® Fluid Receptacle |

Sheath OD 3.6 mm

Sterilized using ethylene oxide unless otherwise indicated in the contents list.

Contraindications: The Arrowg+ard Blue antimicrobial sheath introducer is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs (refer to product instructions for references).

Warning: The results of certain animal experiments have shown phthalates to be potentially toxic to reproduction. Proceeding from the present state of scientific knowledge, risks for male premature infants cannot be excluded in the case of long-term exposure or application. Medical products containing phthalates should be used only temporarily with pregnant women, nursing mothers, babies and infants.

±**Precaution: Contains Phthalate: DEHP.**

Not made with natural rubber latex.

Store between 20 - 25°C (68 - 77°F).

Arrow International provides the enclosed medication label(s) for your convenience. Please ensure that these labels are applied to the correct syringe and corresponding medication.

¹Licensed under US Patent No. 6,706,024.

²A trademark of Becton, Dickinson and Company.

³A registered trademark of CareFusion or one of its subsidiaries.

⁴A registered trademark of Johnson & Johnson Corporation.

