



ARROW®

Pressure Injectable Arrowg[†]ard Blue Plus® Two-Lumen CVC

REF CDC-42802-XPCN1A

8^{Fr.} | 2^{Lumen} | 16 cm catheter length | .032 inch dia. spring-wire guide

Ergo[®]Pack[®] Complete

Contents:

- | | | |
|---|--|--|
| 1: Arrowg [†] ard Blue Plus® Two-Lumen Catheter: 8 Fr. (2.8 mm OD) x 16 cm, Pressure Injectable | 2: 10 mL Luer-Lock Pre-Filled Saline Syringe | 2: Gauze Pad: 2" x 2" (5 cm x 5 cm) |
| 1: Spring-Wire Guide, Nitinol, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow® GlideWheel™ Wire Advancer | 1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep | 5: Gauze Pad: 4" x 4" (10 cm x 10 cm) |
| 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle | 2: 3 mL Applicator 2% CHG and 70% IPA ChloraPrep® ³ One-Step Solution with Hi-Lite Orange™ Tint | 1: Surgical Apparel: Impervious Gown |
| 1: Injection Needle: SafetyGlide™ ² 25 Ga. x 1" (2.54 cm) | 1: SecondSite™ Adjustable Hub: Fastener | 1: Dressing: Tegaderm® ⁴ CHG 3-1/2" x 4-1/2" (8.5 cm x 11.5 cm) |
| 1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm) | 1: SecondSite™ Adjustable Hub: Catheter Clamp | 1: Transducer Cover: 14 cm x 147 cm |
| 1: Pressure Transduction Probe | 1: SharpsAway® Disposal Cup | 1: Surgical Apparel: Mask with Eye Shield |
| 1: Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe | 1: SharpsAway II™ Locking Disposal Cup | 1: Surgical Apparel: Bouffant Cap |
| 1: Injection Needle: SafetyGlide™ ² 22 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe | 1: Maximal Barrier Drape™ with 4" fenestration | 1: Suture: 3-0 Silk with Curved Needle |
| 1: Syringe: 3 mL Luer-Lock | 1: Towel | 1: Tubing [‡] : 20" (50 cm) Extension |
| 1: Tissue Dilator: 9 Fr. (3.0 mm) x 10.2 cm | 1: Needle Holder | 1: HemoHopper® Fluid Receptacle |
| | 1: Safety Scalpel: #11 | |
| | 1: Checklist/CLIP Sheet | |
| | 1: Flow Rate Information Card with Injection Log | |
| | 1: Sterile Procedure Sign | |
| | 1: Medication Label: 1% Lidocaine | |

¹Licensed under US Patent Nos. 6,706,024; 6,872,195 and 7,329,412.

²A trademark of Becton, Dickinson and Company.

³A registered trademark of CareFusion or one of its subsidiaries.


⁴A registered trademark of 3M Company.

Contraindications: The Arrowg[†]ard Blue Plus antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine, and/or sulfa drugs (refer to product instructions for references).

Warning: The results of certain animal experiments have shown phthalates to be potentially toxic to reproduction.

Proceeding from the present state of scientific knowledge, risks for male premature infants cannot be excluded in the case of long-term exposure or application. Medical products containing phthalates should be used only temporarily with pregnant women, nursing mothers, babies and infants.

Arrow International provides the enclosed medication label(s) for your convenience. Please ensure that these labels are applied to the correct syringe and corresponding medication.



Lumen	Priming Volume* (mL)	Gravity Flow Rate† (mL/hr)	MAX Pressure Injection Flow Rate‡ (mL/sec)
Distal (14 Ga.)	0.68	7650	10
Proximal (14 Ga.)	0.70	6000	10

* Priming volumes are approximate and are done without accessories.

† Based on 100 cm head height.

‡ Pressure injection flow rates are determined using media of 11.8 centipoise (cP) viscosity, with 152 cm administration set tubing.



Not made with natural rubber latex



± Precaution: Contains Phthalate: DEHP



Store between 20-25°C (68-77°F).