



ARROWgard Blue PLUS[®] Pressure Injectable Quad-Lumen CVC Kit

REF CDC-42854-P1A **8.5** Fr. | **4** Lumen | **16** cm catheter length | **.032** inch dia. spring-wire guide

Contents:

- 1: Quad-Lumen Indwelling Catheter: 8.5 Fr. x 16 cm Pressure Injectable Radiopaque Polyurethane with Blue FlexTip[®], ARROWgard Blue PLUS[®] Antimicrobial Surface Treatment¹, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- 1: Injection Needle: SafetyGlide[™] 25 Ga. x 1" (2.54 cm)
- 1: Local Anesthetic Component
- 1: Pressure Transduction Probe
- 1: Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow[®] Raulerson Spring-Wire Introduction Syringe
- 1: Injection Needle: SafetyGlide[™] 22 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe
- 1: Syringe: 3 mL Luer-Lock
- 1: Tissue Dilator: 10 Fr. (3.5 mm) x 10.2 cm
- 1: 5 mL 1% Lidocaine HCl Solution
- 1: 3 mL Applicator 2% CHG and 70% IPA Chloraprep[®] One-Step Solution with Hi-Lite Orange[™] Tint (non-sterile solution)
- 1: SecondSite[™] Adjustable Hub: Fastener
- 1: SecondSite[™] Adjustable Hub: Catheter Clamp
- 1: SharpsAway II[™] Locking Disposal Cup

- 1: SharpsAway[®] Disposal Cup
- 1: Maximal Barrier Drape[™] with 4" fenestration
- 1: Needle Holder
- 1: Safety Scalpel: #11
- 1: Flow Rate Information Card with Injection Log
- 1: Sterile Procedure Sign
- 1: Medication Label: 1% Lidocaine
- 2: Gauze Pad: 2" x 2" (5 cm x 5 cm)
- 5: Gauze Pad: 4" x 4" (10 cm x 10 cm)
- 1: Surgical Apparel: Impervious Gown
- 1: Surgical Apparel: Bouffant Cap
- 1: Surgical Apparel: Mask with Eye Shield
- 1: Suture: 3-0 Silk with Curved Needle

¹Licensed under US Patent Nos. 6,706,024; 6,872,195 and 7,329,412.

²A trademark of Becton, Dickinson and Company.

³A registered trademark of CareFusion or one of its subsidiaries.



STERILE EO



Not made with natural rubber latex



Store between 20 - 25°C (68 - 77°F).

Rx only. **Warning: Read all package insert warnings, cautions and instructions prior to use. Failure to do so may result in severe patient injury or death.**

Sterilized using ethylene oxide unless otherwise indicated in the contents list.


California Prop. 65

Warning: Cancer and Reproductive Harm www.P65Warnings.ca.gov.

Fluid path components are non-pyrogenic.

Arrow International provides the enclosed medication label(s) for your convenience. Please ensure that these labels are applied to the correct syringe and corresponding medication.

Contraindications: The ARROWgard Blue Plus antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine, and/or sulfa drugs (refer to product instructions for references).



Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)	MAX Pressure Injection Flow Rate‡ (mL/sec)
Distal (16 Ga.)	0.39	3350	5
Medial 1 (14 Ga.)	0.47	6000	10
Medial 2 (18 Ga.)	0.36	1750	not rated
Proximal (18 Ga.)	0.38	1950	not rated

* Priming volumes are approximate and are done without the injection cap.

† Based on 100 cm head height.

‡ Pressure injection flow rates are determined using media of 11.8 centipoise viscosity, with 60" administration set tubing.

Manufacturer:
Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA

Teleflex[®]
www.teleflex.com

C-42854-101E (10/17)



ARROWgard Blue PLUS[®] Pressure Injectable Quad-Lumen CVC Kit

REF CDC-42854-P1A

8.5 Fr. | **4** Lumen | **16** cm catheter length | **.032** inch dia. spring-wire guide