

ARROW



Quad-Lumen Central Venous Catheterization Kit

REF **CDC-45854-XP1A**

ErgoPack

Contents:

- 1: Quad-Lumen Indwelling Catheter: 8.5 Fr. x 8" (20 cm) Pressure Injectable Radiopaque Polyurethane with Blue FlexTip®, ARROWgard Blue PLUS® Antimicrobial Surface Treatment¹, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- 1: Injection Needle: SafetyGlide™² 25 Ga. x 1" (2.54 cm)
- 1: Local Anesthetic Components
- 1: Pressure Transduction Probe
- 1: Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe
- 1: Injection Needle: SafetyGlide™² 22 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe

- 1: Syringe: 3 mL Luer-Lock
- 1: Tissue Dilator: 10 Fr. (3.5 mm) x 10.2 cm
- 1: 5 mL 1% Lidocaine HCl Solution
- 2: 10 mL Luer-Lock Pre-Filled Saline Syringe
- 1: 3 mL Applicator 2% CHG and 70% IPA ChloroPrep®³ One-Step Solution with Hi-Lite Orange™ Tint
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp
- 1: SharpsAway® Disposal Cup
- 1: SharpsAway II™ Locking Disposal Cup
- 1: Maximal Barrier Drape™ with 4" fenestration
- 1: Towel
- 1: Needle Holder
- 1: Safety Scalpel: #11
- 1: Checklist/CLIP Sheet
- 1: Sterile Procedure Sign
- 1: Flow Rate Information Card with Injection Log

- 1: Medication Label: 1% Lidocaine
- 2: Gauze Pad: 2" x 2" (5 cm x 5 cm)
- 5: Gauze Pad: 4" x 4" (10 cm x 10 cm)
- 1: Dressing: Tegaderm®⁴ 10 cm x 12 cm
- 1: Surgical Apparel: Impervious Gown
- 1: Dressing: BIOPATCH®⁵
- 1: Surgical Apparel: Bouffant Cap
- 1: Surgical Apparel: Mask with Eye Shield
- 1: Suture: 3-0 Silk with Curved Needle
- 1: Tubing±: 20" (50 cm) Extension
- 1: HemoHopper® Fluid Receptacle

¹ Licensed under US Patent Nos. 6,706,024; 6,872,195 and 7,329,412.
² A trademark of Becton, Dickinson and Company.
³ A registered trademark of CareFusion or one of its subsidiaries.
⁴ A registered trademark of 3M Company.
⁵ A registered trademark of Johnson & Johnson Corporation.

Caution	Consult instructions for use	Do not reuse	Do not resterilize	Sterilized by ethylene oxide	Keep away from sunlight	Keep dry	Do not use if package is damaged	Not made with natural rubber latex	± Precaution: Contains Phthalate: DEHP.	Store between 20 - 25°C (68 - 77°F).

Rx only. Warning: Read all package insert warnings, precautions and instructions prior to use. Failure to do so may result in severe patient injury or death.

California Prop. 65
Warning: Cancer and Reproductive Harm www.P65Warnings.ca.gov.

Contraindications: The ARROWgard Blue PLUS® antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs (refer to product instructions for references).

All contents sterile unless otherwise indicated.

Arrow International provides the enclosed medication label(s) for your convenience. Please ensure that these labels are applied to the correct syringe and corresponding medication.

Fluid path components are non-pyrogenic.

DuPont™ Tyvek®
package

Manufacturer:
 Arrow International, Inc.
 Subsidiary of Teleflex Incorporated
 2400 Bernville Road
 Reading, PA 19605 USA

Teleflex
www.teleflex.com

C-45854-107F (5/17)



Quad-Lumen Central Venous Catheterization Kit

REF **CDC-45854-XP1A**

8.5 Fr.

4 Lumen

20 cm catheter length

.032 inch dia. spring-wire guide

ARROW