



ARROW ENDURANCE™ Extended Dwell Peripheral Catheter System

ARROW

REF **EDC-00622**
Contents:

22 Ga. catheter | **6** cm catheter length | **.010** inch dia. wire guide

- 1: Catheter System: 22 Ga. (0.91 mm) x 60 mm (2.36")
ARROW ENDURANCE™ over 24 Ga. (0.60 mm) TW
Introducer Needle with Integral 0.010" (0.26 mm) Straight
Soft Tip Wire Guide
- 1: Label: Venous/Arterial Placement

Rx only. Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Warning: This product contains a chemical known to the State of California to cause cancer, birth defects or other reproductive harm.

Sterile. Fluid path components are non-pyrogenic.

Manufacturer:
Arrow International, Inc.
Subsidiary of
Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA

Teleflex®
www.teleflex.com
EDC-00622-100A (1/17)



Caution



Consult
instructions for use



Do not reuse



Do not resterilize



Sterilized
by ethylene oxide



Keep away
from sunlight



Keep dry



Do not use
if package is damaged



Not made with
natural rubber latex



ARROW ENDURANCE™ Extended Dwell Peripheral Catheter System

ARROW

REF **EDC-00622**
Contents:

22 Ga. catheter | **6** cm catheter length | **.010** inch dia. wire guide

- 1: Catheter System: 22 Ga. (0.91 mm) x 60 mm (2.36")
ARROW ENDURANCE™ over 24 Ga. (0.60 mm) TW
Introducer Needle with Integral 0.010" (0.26 mm) Straight
Soft Tip Wire Guide
- 1: Label: Venous/Arterial Placement

Rx only. Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Warning: This product contains a chemical known to the State of California to cause cancer, birth defects or other reproductive harm.

Sterile. Fluid path components are non-pyrogenic.

Manufacturer:
Arrow International, Inc.
Subsidiary of
Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA

Teleflex®
www.teleflex.com
EDC-00622-100A (1/17)



Caution



Consult
instructions for use



Do not reuse



Do not resterilize



Sterilized
by ethylene oxide



Keep away
from sunlight



Keep dry



Do not use
if package is damaged



Not made with
natural rubber latex

DuPont™ Tyvek®
packaging

DuPont™ Tyvek®
packaging