



ARROW ENDURANCE™ Extended Dwell Peripheral Catheter System



REF **EDC-00820**

20 Ga. catheter | **8** cm catheter length | **.010** inch dia. wire guide

Contents:

- 1: Catheter System: 20 Ga. (1.10 mm) x 80 mm (3.15") ARROW ENDURANCE™ over 22 Ga. (0.80 mm) RW Introducer Needle with Integral 0.010" (0.26 mm) Straight Soft Tip Wire Guide
- 1: Label: Venous/Arterial Placement

Rx only. Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.
Warning: This product contains a chemical known to the State of California to cause cancer, birth defects or other reproductive harm.
Sterile. Fluid path components are non-pyrogenic.

Manufacturer:
Arrow International, Inc.
Subsidiary of
Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA



www.teleflex.com
EDC-00820-100A (1/17)



Caution



Consult instructions for use



Do not reuse



Do not resterilize



Sterilized by ethylene oxide



Keep away from sunlight



Keep dry



Do not use if package is damaged



Not made with natural rubber latex



ARROW ENDURANCE™ Extended Dwell Peripheral Catheter System



REF **EDC-00820**

20 Ga. catheter | **8** cm catheter length | **.010** inch dia. wire guide

Contents:

- 1: Catheter System: 20 Ga. (1.10 mm) x 80 mm (3.15") ARROW ENDURANCE™ over 22 Ga. (0.80 mm) RW Introducer Needle with Integral 0.010" (0.26 mm) Straight Soft Tip Wire Guide
- 1: Label: Venous/Arterial Placement

Rx only. Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.
Warning: This product contains a chemical known to the State of California to cause cancer, birth defects or other reproductive harm.
Sterile. Fluid path components are non-pyrogenic.

Manufacturer:
Arrow International, Inc.
Subsidiary of
Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA



www.teleflex.com
EDC-00820-100A (1/17)



Caution



Consult instructions for use



Do not reuse



Do not resterilize



Sterilized by ethylene oxide



Keep away from sunlight



Keep dry



Do not use if package is damaged



Not made with natural rubber latex

DuPont™ Tyvek®
packaging

DuPont™ Tyvek®
packaging