



ARROWgard Blue⁺ Two-Lumen Hemodialysis Catheterization Kit for High Volume Infusions

ARROW

REF CDC-22142-1A

14^{Fr.} | 2^{Lumen} | 15^{cm} catheter length | .035^{inch dia.} spring-wire guide

ErgoPack[®]

Contents:

- | | | |
|---|--|---|
| <ul style="list-style-type: none">1: Two-Lumen Indwelling Catheter: 14 Fr. x 15 cm Radiopaque Polyurethane with Blue FlexTip[®], ARROWgard[®] Antimicrobial Surface Treatment¹, Extension Line Clamps1: Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle1: Injection Needle: SafetyGlide[™] 25 Ga. x 1" (2.54 cm)1: Local Anesthetic Component1: Pressure Transduction Probe1: Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow[®] Raulerson Spring-Wire Introduction Syringe1: Injection Needle: SafetyGlide[™] 22 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe1: Syringe: 3 mL Luer-Lock | <ul style="list-style-type: none">1: Tissue Dilator: 14 Fr. (4.7 mm) x 14 cm1: Tissue Dilator: 12 Fr. (4.0 mm) x 14 cm1: 5 mL 1% Lidocaine HCl Solution1: 3 mL Applicator 2% CHG and 70% IPA Chloraprep^{®3} One-Step Solution with Hi-Lite Orange[™] Tint (non-sterile solution)1: SharpsAway[®] Disposal Cup1: SharpsAway II[™] Locking Disposal Cup1: Maximal Barrier Drape[™] with 4" fenestration1: Needle Holder1: Safety Scalpel: #111: Checklist/CLIP Sheet1: Sterile Procedure Sign1: Medication Label: 1% Lidocaine2: Gauze Pad: 2" x 2" (5 cm x 5 cm)5: Gauze Pad: 4" x 4" (10 cm x 10 cm) | <ul style="list-style-type: none">1: Surgical Apparel: Impervious Gown1: Surgical Apparel: Bouffant Cap1: Surgical Apparel: Mask with Eye Shield1: Suture: 3-0 Braided Silk with Straight Needle1: Suture: 3-0 Silk with Curved Needle1: HemoHopper[®] Fluid Receptacle |
|---|--|---|

¹Licensed under US Patent No. 6,706,024.

²A trademark of Becton, Dickinson and Company.

³A registered trademark of CareFusion or one of its subsidiaries.


Sterilized using ethylene oxide unless otherwise indicated in the contents list.

Contraindications: The Arrow large-bore two-lumen catheter is not designed for long-term hemodialysis or for use in patients with thrombosed vessels. The ARROWgard Blue antimicrobial surface catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine, and/or sulfa drugs (refer to product instructions for references).

Arrow International provides the enclosed medication label(s) for your convenience. Please ensure that these labels are applied to the correct syringe and corresponding medication.

LATEX
Not made with natural rubber latex

Store between 20-25°C (68-77°F)



H-22142-400A (03/18)