



# ARROW<sup>+</sup>gard Blue<sup>®</sup> You-Bend<sup>™</sup> Two-Lumen Hemodialysis Catheterization Kit for High Volume Infusions



REF CDC-25122-XU1A

12 Fr. | 2 Lumen | 20 cm catheter length | .035 inch dia. spring-wire guide

ErgoPack<sup>®</sup>

**Contents:**

- 1: Two-Lumen Indwelling Catheter<sup>†</sup>: 12 Fr. x 20 cm Radiopaque Polyurethane with Blue FlexTip<sup>®</sup>, ARROW<sup>+</sup>gard<sup>®</sup> Antimicrobial Surface Treatment<sup>†</sup>
- 1: Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 26-3/4" (68 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- 1: Injection Needle: SafetyGlide<sup>™</sup> 25 Ga. x 1" (2.54 cm)
- 1: Local Anesthetic Component
- 1: Pressure Transduction Probe
- 1: Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow<sup>®</sup> Raulerson Spring-Wire Introduction Syringe
- 1: Injection Needle: SafetyGlide<sup>™</sup> 22 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe
- 1: Syringe: 3 mL Luer-Lock

- 1: Tissue Dilator: 10 Fr. (3.5 mm) x 10.2 cm
- 1: Tissue Dilator: 12 Fr. (4.0 mm) x 14 cm
- 2: 10 mL Luer-Lock Pre-Filled Saline Syringe
- 1: 5 mL 1% Lidocaine HCl Solution
- 1: 3 mL Applicator 2% CHG and 70% IPA Chloraprep<sup>®</sup> One-Step Solution with Hi-Lite Orange<sup>™</sup> Tint (non-sterile solution)
- 1: SharpsAway<sup>®</sup> Disposal Cup
- 1: SharpsAway II<sup>™</sup> Locking Disposal Cup
- 1: Maximal Barrier Drape<sup>™</sup> with 4" fenestration
- 1: Towel
- 1: Needle Holder
- 1: Safety Scalpel: #11
- 1: Checklist/CLIP Sheet
- 1: Sterile Procedure Sign
- 1: Medication Label: 1% Lidocaine

- 2: Gauze Pad: 2" x 2" (5 cm x 5 cm)
- 5: Gauze Pad: 4" x 4" (10 cm x 10 cm)
- 1: Dressing: Tegaderm<sup>®</sup> 10 cm x 12 cm
- 1: Surgical Apparel: Impervious Gown
- 1: Dressing: BIOPATCH<sup>®</sup>
- 1: Surgical Apparel: Bouffant Cap
- 1: Surgical Apparel: Mask with Eye Shield
- 1: Suture: 3-0 Braided Silk with Straight Needle
- 1: Suture: 3-0 Silk with Curved Needle
- 1: HemoHopper<sup>®</sup> Fluid Receptacle

<sup>†</sup>Licensed under US Patent No. 6,706,024.  
<sup>‡</sup>A trademark of Becton, Dickinson and Company.  
<sup>§</sup>A registered trademark of CareFusion or one of its subsidiaries.  
<sup>¶</sup>A registered trademark of 3M Company.  
<sup>§§</sup>A registered trademark of Johnson & Johnson Corporation.

Sterilized using ethylene oxide unless otherwise indicated in the contents list.

**Contraindications: The Arrow large-bore two-lumen catheter is not designed for long-term hemodialysis or for use in patients with thrombosed vessels.**

**The ARROW<sup>+</sup>gard Blue antimicrobial surface catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine, and/or sulfa drugs (refer to product instructions for references).**

Arrow International provides the enclosed medication label(s) for your convenience. Please ensure that these labels are applied to the correct syringe and corresponding medication.




Not made with natural rubber latex



<sup>†</sup> MR Conditional




Store between 20 - 25°C (68 - 77°F).



Priming Volume*	
Distal	1.80 mL
Proximal	1.80 mL

\*Priming volumes are approximate and are done without injection caps.



<sup>†</sup> In vitro tests performed using 3:1 glycerine/saline at 37°C.