

Femoral Arterial Line Catheterization Kit



Arrow Select Kits

ASK-04018-HCA

Contents:

- 1: Indwelling Catheter: 18 Ga. x 6" (16 cm) with Integral Suture Wing
- 1: Spring-Wire Guide: .025" (.64 mm) dia. x 17-13/16" (45 cm) (Straight Soft Tip on One End – "J" Tip on Other)
- 1: Catheter: 20 Ga. x 1-3/4" (4.45 cm) Radiopaque over 22 Ga. TW Introducer Needle
- 1: Injection Needle: Eclipse™¹ 22 Ga. x 1-1/2" (3.81 cm) RW
- 1: Injection Needle: Eclipse™¹ 25 Ga. x 1" (2.54 cm) RW
- 1: Introducer Needle: 18 Ga. x 3" (7.62 cm) RW
- 1: Introducer Needle: 20 Ga. x 1-1/2" (3.81 cm) XTW
- 2: Syringes: 5 mL Luer-Slip
- 1: **10.5 mL Applicator Pouch, 2% CHG and 70% IPA Chloraprep®² One-Step Solution with Hi-Lite Orange™ Tint**
- 1: 5 mL HCl, 1% Lidocaine Solution
- 1: Local Anesthetic Components
- 1: Drape: 24" x 36" with 4" fenestration
- 1: **Drape: 41" x 55" with 4" fenestration**
- 1: Needle Holder
- 1: Prep Razor
- 1: Safety Scalpel: #11
- 4: Gauze Pads: 2" x 2"
- 8: Gauze Pads: 4" x 4"
- 2: Prep Sponge Swabs
- 1: Suture: 2-0 Silk with Curved Needle
- 1: Suture: 2-0 Silk with Straight Needle
- 1: Tray: Prep

¹ A trademark of Becton, Dickinson and Company.

² A registered trademark of Cardinal Health, Inc., or one of its subsidiaries.

Infection Protection Features:

- > **Chloraprep®** One-Step skin antiseptic, 2% chlorhexidine gluconate in 70% isopropyl alcohol with Hi-Lite Orange™ Tint
- > Extra-Large Drape size

Rx only.

Warning: Prior to use read all package insert warnings, precautions, and instructions. Failure to do so may result in severe patient injury or death.

All contents sterile unless otherwise indicated.

Single use – do not resterilize. Sterilized by ethylene oxide. Do not use if package has been previously opened or damaged. Fluid path components are non-pyrogenic. Contains medication.

Store between 20 - 25°C (68 - 77°F). Avoid freezing and excessive heat above 40°C (104°F).

Not made with natural rubber latex.

Manufacturer:
Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA

Teleflex®

www.teleflex.com

IA-04018-109D (9/18)

Femoral Arterial Line Catheterization Kit

ASK-04018-HCA

Arrow Select Kits



18 Ga.

16 cm catheter length

.025 inch diameter spring-wire guide