



18<sup>Ga.</sup> 16<sup>cm</sup> catheter length 025 inch dia. spring-w guide

## **Contents:**

- 1: Indwelling Catheter: 18 Ga. x 16 cm with Integral Suture Wing
  1: Spring-Wire Guide: .025" (0.64 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End "J"
- 1: Catheter: 20 Ga. x 1-3/4" (4.45 cm) Radiopaque over 22 Ga. TW Introducer Needle
  1: Introducer Needle: 20 Ga. x 1-1/2" (3.81 cm) XTW
  1: Injection Needle: Eclipse<sup>™1</sup> 22 Ga. x 1-1/2" (3.81 cm) RW
  1: Injection Needle: Eclipse<sup>™1</sup> 25 Ga. x 1" (2.54 cm) RW
  1: Introducer Needle: 18 Ga. x 3" (7.62 cm) RW

- Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)
- 2: Syringe: 5 mL Luer-Slip1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep
- 10.5 mL Applicator 2% CHG and 70% IPA Chlora Prep®2 One-Step Solution with Hi-Lite Orange™Tint
  1: Drape: 24" x 36" (61 cm x 91 cm) with 4" (10 cm) fenestration
  1: Drape: 41" x 55" with 4" fenestration

- 1: Needle Holder Safety Scalpel: #11
- Prep Razor
- Prep Sponge Swab Gauze Pad: 2" x 2" (5 cm x 5 cm) Gauze Pad: 4" x 4" (10 cm x 10 cm)
- Suture: 2-0 Silk with Curved Needle Suture: 2-0 Silk with Straight Needle
- Tray: Prep

<sup>1</sup>A trademark of Becton, Dickinson and Company.

<sup>2</sup>A registered trademark of CareFusion or one of its subsidiaries.

Rx only. Warning: Read all package insert warnings, cautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

⚠ California Prop. 65
Warning: Cancer and Reproductive Harm www.P65Warnings.ca.gov.

Fluid path components are non-pyrogenic.





















Manufacturer: Arrow International LLC Subsidiary of Teleflex Incorporated 3015 Carrington Mill Blvd. Morrisville, NC 27560 USA



IA-04018-115A, Rev. 01 (2020-12)





